





AUXILOCK® TITANIUM INTERFERENCE SCREW



The Titanium interference screw is made up of Titanium alloy. The screw has a fully threaded design which provides strong mechanical fixation for both Bone-Tendon Bone (BTB) and soft tissue grafts. The rounded edge of the threads protects tissue grafts. Titanium Interference Screw is a cannulated screw to be used with guide wire and cannulated screw driver.

ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION
6-010-01	AUXILOCK® 7mm x 20mm Titanium Interference Screw
6-010-02	AUXILOCK® 7mm x 25mm Titanium Interference Screw
6-010-03	AUXILOCK® 7mm x 30mm Titanium Interference Screw
6-010-04	AUXILOCK® 7mm x 35mm Titanium Interference Screw
6-010-05	AUXILOCK® 8mm x 20mm Titanium Interference Screw
6-010-06	AUXILOCK® 8mm x 25mm Titanium Interference Screw
6-010-07	AUXILOCK® 8mm x 30mm Titanium Interference Screw
6-010-08	AUXILOCK® 8mm x 35mm Titanium Interference Screw
6-010-09	AUXILOCK® 9mm x 20mm Titanium Interference Screw
6-010-10	AUXILOCK® 9mm x 25mm Titanium Interference Screw
6-010-11	AUXILOCK® 9mm x 30mm Titanium Interference Screw
6-010-12	AUXILOCK® 9mm x 35mm Titanium Interference Screw
6-010-13	AUXILOCK® 10mm x 25mm Titanium Interference Screw
6-010-14	AUXILOCK® 10mm x 30mm Titanium Interference Screw
6-010-15	AUXILOCK® 10mm x 35mm Titanium Interference Screw
6-010-16	AUXILOCK® 11mm x 25mm Titanium Interference Screw
6-010-17	AUXILOCK® 11mm x 30mm Titanium Interference Screw
6-010-18	AUXILOCK® 11mm x 35mm Titanium Interference Screw



ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION
6-010-19	AUXILOCK® 12mm x 25mm Titanium Interference Screw
6-010-20	AUXILOCK® 12mm x 30mm Titanium Interference Screw
6-010-21	AUXILOCK® 12mm x 35mm Titanium Interference Screw

RECOMMENDED INSTRUMENTATION

CODE	PRODUCT DESCRIPTION	COMPATIBLE IMPLANTS
7-014-43	Nitinol Guide Wire Ø1.2mm x Length 300mm (For Interference Screws)	6-010-01, 6-010-02, 6-010-03, 6-010-04, 6-010-05, 6-010-06, 6-010-07, 6-010-08,
7-014-05	Screwdriver For Ti Interference Screw, T25 Star Driver	6-010-09, 6-010-10, 6-010-11, 6-010-12,
7-014-95	Screwdriver Shaft For Titanium Interference Screw, T25 Star Driver, Hudson Coupling	6-010-13, 6-010-14, 6-010-15, 6-010-16, 6-010-17, 6-010-18, 6-010-19, 6-010-20,
7-014-42	Ratchet Handle For Hudson Coupling Screwdriver Shaft	6-010-21