




AUXILOCK® TITANIUM SCREW-IN SUTURE ANCHOR



ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION	SUTURE/TAPE COLOUR	NEEDLE
6-005-16	AUXILOCK® 2.0mm Titanium Screw-In Suture Anchor With One #0 BioBraid: White/Blue		
6-005-11	AUXILOCK® 2.0mm Titanium Screw-In Suture Anchor With One #0 BioBraid: White/Blue, With Needles: MO-7		
6-005-15	AUXILOCK® 2.5mm Titanium Screw-In Suture Anchor With One #1 BioBraid: White/Blue		
6-005-12	AUXILOCK® 2.5mm Titanium Screw-In Suture Anchor With One #1 BioBraid: White/Blue, With Needles: MO-7		

RECOMMENDED INSTRUMENTATION

CODE	PRODUCT DESCRIPTION	COMPATIBLE IMPLANTS
7-139-01	Ø1.2mm X 100mm Drill Bit	6-005-11, 6-005-16
7-139-02	8 Point Slotted Drill Guide, Length 50mm	6-005-11, 6-005-12, 6-005-16, 6-005-15
7-139-66	Ø1.8mm X 100mm Drill Bit	6-005-15, 6-005-12