

Auxein


**ATHOS NAILING
SYSTEM**



NEXT-LEVEL CARE

ABOUT US

Auxein Medical is an integrated, research based, orthopaedic implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavour to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in India and serves customers in over 75 countries worldwide.



IMPACTING HEALTHCARE GLOBALLY

Our Reach & Influence

500+

Employees Worldwide

75+

Operations in Countries

150+

Distributors Worldwide

20^{M+}

Patients Cured

200+

FDA Approved Products

250+

Machines Worldwide

3000+

Innovative Orthopedic
Products & Solutions


Auxin is dedicated to advancing medical technology through innovative orthopedic solutions, including **Trauma, Arthroscopy, and Reconstructive Procedures**, with a commitment to improving patient outcomes worldwide. Our approach integrates advanced research, precision engineering, and high-quality manufacturing to develop reliable implants and surgical instruments that address the evolving needs of modern orthopedic practice.

OUR ACHIEVEMENTS



OVER THE YEARS






At Auxein, it is our goal to provide the best in class quality and value to every individual. Consistency is a virtue we pursue with a passion. From the procurement of raw materials to the manufacturing of our range of 1200+ orthopaedic implant products, we employ world-class systems and technology at par with the highest standards of globally recognized evaluation systems. The Athos Nailing System is part of our orthopedic trauma portfolio designed for the intramedullary fixation of femoral fractures. Developed with input from global experts, the system provides reliable stabilization, helping achieve stable fixation, proper load sharing, and accurate anatomical alignment during fracture management.



INDEX



Introduction	09-10
Athos Nail	11-20
Screws	14-21
Implant Set	22-29
Athos Nail Anti-Rotation	30-42
Screws	43-44
Implant Set	45-51
Instruments	52-65
Instruments (PEEK)	66-80
Instruments (Distal Targeting Device)	81-90



ATHOS NAILING SYSTEM

ATHOS NAILING SYSTEM

Introduction:

The ATHOS Nailing System by Auxein is an intramedullary fixation system developed for the comprehensive management of proximal femoral fractures. Designed to provide biomechanical stability while supporting minimally invasive surgical techniques, the system offers two configurations to address varying clinical requirements:

- ATHOS Nailing System (Type-II)
- ATHOS Nailing System (With Anti-Rotation)






Athos Nailing System

Key Features:

- Anatomically contoured proximal geometry for natural load transfer.
- Ø10.5mm cephalic screw profile designed for angular stability and load sharing.
- Multiple distal locking options (static/dynamic) to optimize fixation.
- Radiolucent targeting device enabling precision screw placement.
- Left and right specific nails to match anatomy.
- Stainless steel and titanium availability, compliant with ISO 5832-1 & ISO 5832-3 standards.

Indications:

- Stable and unstable pertrochanteric fractures (AO/OTA 31-A1, 31-A2).
- Intertrochanteric fractures with subtrochanteric extension (AO/OTA 31-A3).
- Pathological fractures of the proximal femur.
- Complex trochanteric fractures involving the femoral shaft.
- Revision procedures after failed fixation.



125° Athos Nail Type-II, Short



ØDia	Length	Titanium	Stainless Steel
9mm	180mm	TI-1500S-09-180	1500S-09-180
9mm	200mm	TI-1500S-09-200	1500S-09-200
9mm	220mm	TI-1500S-09-220	1500S-09-220
9mm	240mm	TI-1500S-09-240	1500S-09-240
10mm	180mm	TI-1500S-10-180	1500S-10-180
10mm	200mm	TI-1500S-10-200	1500S-10-200
10mm	220mm	TI-1500S-10-220	1500S-10-220
10mm	240mm	TI-1500S-10-240	1500S-10-240
11mm	180mm	TI-1500S-11-180	1500S-11-180
11mm	200mm	TI-1500S-11-200	1500S-11-200
11mm	220mm	TI-1500S-11-220	1500S-11-220
11mm	240mm	TI-1500S-11-240	1500S-11-240
12mm	180mm	TI-1500S-12-180	1500S-12-180
12mm	200mm	TI-1500S-12-200	1500S-12-200
12mm	220mm	TI-1500S-12-220	1500S-12-220
12mm	240mm	TI-1500S-12-240	1500S-12-240
13mm	180mm	TI-1500S-13-180	1500S-13-180
13mm	200mm	TI-1500S-13-200	1500S-13-200
13mm	220mm	TI-1500S-13-220	1500S-13-220
13mm	240mm	TI-1500S-13-240	1500S-13-240



125° Athos Nail Type-II, Short

STERILE


ØDia	Length	Titanium	Stainless Steel
9mm	180mm	TI-1500S-09-180-S	1500S-09-180-S
9mm	200mm	TI-1500S-09-200-S	1500S-09-200-S
9mm	220mm	TI-1500S-09-220-S	1500S-09-220-S
9mm	240mm	TI-1500S-09-240-S	1500S-09-240-S
10mm	180mm	TI-1500S-10-180-S	1500S-10-180-S
10mm	200mm	TI-1500S-10-200-S	1500S-10-200-S
10mm	220mm	TI-1500S-10-220-S	1500S-10-220-S
10mm	240mm	TI-1500S-10-240-S	1500S-10-240-S
11mm	180mm	TI-1500S-11-180-S	1500S-11-180-S
11mm	200mm	TI-1500S-11-200-S	1500S-11-200-S
11mm	220mm	TI-1500S-11-220-S	1500S-11-220-S
11mm	240mm	TI-1500S-11-240-S	1500S-11-240-S
12mm	180mm	TI-1500S-12-180-S	1500S-12-180-S
12mm	200mm	TI-1500S-12-200-S	1500S-12-200-S
12mm	220mm	TI-1500S-12-220-S	1500S-12-220-S
12mm	240mm	TI-1500S-12-240-S	1500S-12-240-S
13mm	180mm	TI-1500S-13-180-S	1500S-13-180-S
13mm	200mm	TI-1500S-13-200-S	1500S-13-200-S
13mm	220mm	TI-1500S-13-220-S	1500S-13-220-S
13mm	240mm	TI-1500S-13-240-S	1500S-13-240-S



125° Athos Nail Type-II, Long



ØDia	Length	Titanium(Left)	Stainless Steel(Left)	Titanium(Right)	Stainless Steel(Right)
9mm	320mm	TI-1500L-09-320L	1500L-09-320L	TI-1500L-09-320R	1500L-09-320R
9mm	340mm	TI-1500L-09-340L	1500L-09-340L	TI-1500L-09-340R	1500L-09-340R
9mm	360mm	TI-1500L-09-360L	1500L-09-360L	TI-1500L-09-360R	1500L-09-360R
9mm	380mm	TI-1500L-09-380L	1500L-09-380L	TI-1500L-09-380R	1500L-09-380R
9mm	400mm	TI-1500L-09-400L	1500L-09-400L	TI-1500L-09-400R	1500L-09-400R
9mm	420mm	TI-1500L-09-420L	1500L-09-420L	TI-1500L-09-420R	1500L-09-420R
9mm	440mm	TI-1500L-09-440L	1500L-09-440L	TI-1500L-09-440R	1500L-09-440R
10mm	320mm	TI-1500L-10-320L	1500L-10-320L	TI-1500L-10-320R	1500L-10-320R
10mm	340mm	TI-1500L-10-340L	1500L-10-340L	TI-1500L-10-340R	1500L-10-340R
10mm	360mm	TI-1500L-10-360L	1500L-10-360L	TI-1500L-10-360R	1500L-10-360R
10mm	380mm	TI-1500L-10-380L	1500L-10-380L	TI-1500L-10-380R	1500L-10-380R
10mm	400mm	TI-1500L-10-400L	1500L-10-400L	TI-1500L-10-400R	1500L-10-400R
10mm	420mm	TI-1500L-10-420L	1500L-10-420L	TI-1500L-10-420R	1500L-10-420R
10mm	440mm	TI-1500L-10-440L	1500L-10-440L	TI-1500L-10-440R	1500L-10-440R
11mm	320mm	TI-1500L-11-320L	1500L-11-320L	TI-1500L-11-320R	1500L-11-320R
11mm	340mm	TI-1500L-11-340L	1500L-11-340L	TI-1500L-11-340R	1500L-11-340R
11mm	360mm	TI-1500L-11-360L	1500L-11-360L	TI-1500L-11-360R	1500L-11-360R
11mm	380mm	TI-1500L-11-380L	1500L-11-380L	TI-1500L-11-380R	1500L-11-380R
11mm	400mm	TI-1500L-11-400L	1500L-11-400L	TI-1500L-11-400R	1500L-11-400R
11mm	420mm	TI-1500L-11-420L	1500L-11-420L	TI-1500L-11-420R	1500L-11-420R
11mm	440mm	TI-1500L-11-440L	1500L-11-440L	TI-1500L-11-440R	1500L-11-440R
12mm	320mm	TI-1500L-12-320L	1500L-12-320L	TI-1500L-12-320R	1500L-12-320R
12mm	340mm	TI-1500L-12-340L	1500L-12-340L	TI-1500L-12-340R	1500L-12-340R
12mm	360mm	TI-1500L-12-360L	1500L-12-360L	TI-1500L-12-360R	1500L-12-360R
12mm	380mm	TI-1500L-12-380L	1500L-12-380L	TI-1500L-12-380R	1500L-12-380R
12mm	400mm	TI-1500L-12-400L	1500L-12-400L	TI-1500L-12-400R	1500L-12-400R
12mm	420mm	TI-1500L-12-420L	1500L-12-420L	TI-1500L-12-420R	1500L-12-420R
12mm	440mm	TI-1500L-12-440L	1500L-12-440L	TI-1500L-12-440R	1500L-12-440R
13mm	320mm	TI-1500L-13-320L	1500L-13-320L	TI-1500L-13-320R	1500L-13-320R
13mm	340mm	TI-1500L-13-340L	1500L-13-340L	TI-1500L-13-340R	1500L-13-340R
13mm	360mm	TI-1500L-13-360L	1500L-13-360L	TI-1500L-13-360R	1500L-13-360R
13mm	380mm	TI-1500L-13-380L	1500L-13-380L	TI-1500L-13-380R	1500L-13-380R
13mm	400mm	TI-1500L-13-400L	1500L-13-400L	TI-1500L-13-400R	1500L-13-400R
13mm	420mm	TI-1500L-13-420L	1500L-13-420L	TI-1500L-13-420R	1500L-13-420R
13mm	440mm	TI-1500L-13-440L	1500L-13-440L	TI-1500L-13-440R	1500L-13-440R




Left Right



125° Athos Nail Type-II, Long

STERILE



ØDia	Length	Titanium(Left)	Stainless Steel(Left)	Titanium(Right)	Stainless Steel(Right)
9mm	320mm	TI-1500L-09-320L-S	1500L-09-320L-S	TI-1500L-09-320R-S	1500L-09-320R-S
9mm	340mm	TI-1500L-09-340L-S	1500L-09-340L-S	TI-1500L-09-340R-S	1500L-09-340R-S
9mm	360mm	TI-1500L-09-360L-S	1500L-09-360L-S	TI-1500L-09-360R-S	1500L-09-360R-S
9mm	380mm	TI-1500L-09-380L-S	1500L-09-380L-S	TI-1500L-09-380R-S	1500L-09-380R-S
9mm	400mm	TI-1500L-09-400L-S	1500L-09-400L-S	TI-1500L-09-400R-S	1500L-09-400R-S
9mm	420mm	TI-1500L-09-420L-S	1500L-09-420L-S	TI-1500L-09-420R-S	1500L-09-420R-S
9mm	440mm	TI-1500L-09-440L-S	1500L-09-440L-S	TI-1500L-09-440R-S	1500L-09-440R-S
10mm	320mm	TI-1500L-10-320L-S	1500L-10-320L-S	TI-1500L-10-320R-S	1500L-10-320R-S
10mm	340mm	TI-1500L-10-340L-S	1500L-10-340L-S	TI-1500L-10-340R-S	1500L-10-340R-S
10mm	360mm	TI-1500L-10-360L-S	1500L-10-360L-S	TI-1500L-10-360R-S	1500L-10-360R-S
10mm	380mm	TI-1500L-10-380L-S	1500L-10-380L-S	TI-1500L-10-380R-S	1500L-10-380R-S
10mm	400mm	TI-1500L-10-400L-S	1500L-10-400L-S	TI-1500L-10-400R-S	1500L-10-400R-S
10mm	420mm	TI-1500L-10-420L-S	1500L-10-420L-S	TI-1500L-10-420R-S	1500L-10-420R-S
10mm	440mm	TI-1500L-10-440L-S	1500L-10-440L-S	TI-1500L-10-440R-S	1500L-10-440R-S
11mm	320mm	TI-1500L-11-320L-S	1500L-11-320L-S	TI-1500L-11-320R-S	1500L-11-320R-S
11mm	340mm	TI-1500L-11-340L-S	1500L-11-340L-S	TI-1500L-11-340R-S	1500L-11-340R-S
11mm	360mm	TI-1500L-11-360L-S	1500L-11-360L-S	TI-1500L-11-360R-S	1500L-11-360R-S
11mm	380mm	TI-1500L-11-380L-S	1500L-11-380L-S	TI-1500L-11-380R-S	1500L-11-380R-S
11mm	400mm	TI-1500L-11-400L-S	1500L-11-400L-S	TI-1500L-11-400R-S	1500L-11-400R-S
11mm	420mm	TI-1500L-11-420L-S	1500L-11-420L-S	TI-1500L-11-420R-S	1500L-11-420R-S
11mm	440mm	TI-1500L-11-440L-S	1500L-11-440L-S	TI-1500L-11-440R-S	1500L-11-440R-S
12mm	320mm	TI-1500L-12-320L-S	1500L-12-320L-S	TI-1500L-12-320R-S	1500L-12-320R-S
12mm	340mm	TI-1500L-12-340L-S	1500L-12-340L-S	TI-1500L-12-340R-S	1500L-12-340R-S
12mm	360mm	TI-1500L-12-360L-S	1500L-12-360L-S	TI-1500L-12-360R-S	1500L-12-360R-S
12mm	380mm	TI-1500L-12-380L-S	1500L-12-380L-S	TI-1500L-12-380R-S	1500L-12-380R-S
12mm	400mm	TI-1500L-12-400L-S	1500L-12-400L-S	TI-1500L-12-400R-S	1500L-12-400R-S
12mm	420mm	TI-1500L-12-420L-S	1500L-12-420L-S	TI-1500L-12-420R-S	1500L-12-420R-S
12mm	440mm	TI-1500L-12-440L-S	1500L-12-440L-S	TI-1500L-12-440R-S	1500L-12-440R-S
13mm	320mm	TI-1500L-13-320L-S	1500L-13-320L-S	TI-1500L-13-320R-S	1500L-13-320R-S
13mm	340mm	TI-1500L-13-340L-S	1500L-13-340L-S	TI-1500L-13-340R-S	1500L-13-340R-S
13mm	360mm	TI-1500L-13-360L-S	1500L-13-360L-S	TI-1500L-13-360R-S	1500L-13-360R-S
13mm	380mm	TI-1500L-13-380L-S	1500L-13-380L-S	TI-1500L-13-380R-S	1500L-13-380R-S
13mm	400mm	TI-1500L-13-400L-S	1500L-13-400L-S	TI-1500L-13-400R-S	1500L-13-400R-S
13mm	420mm	TI-1500L-13-420L-S	1500L-13-420L-S	TI-1500L-13-420R-S	1500L-13-420R-S
13mm	440mm	TI-1500L-13-440L-S	1500L-13-440L-S	TI-1500L-13-440R-S	1500L-13-440R-S

Left Right





130° Athos Nail Type-II, Short



ØDia	Length	Titanium	Stainless Steel
9mm	180mm	TI-1457S-09-180	1457S-09-180
9mm	200mm	TI-1457S-09-200	1457S-09-200
9mm	220mm	TI-1457S-09-220	1457S-09-220
9mm	240mm	TI-1457S-09-240	1457S-09-240
10mm	180mm	TI-1457S-10-180	1457S-10-180
10mm	200mm	TI-1457S-10-200	1457S-10-200
10mm	220mm	TI-1457S-10-220	1457S-10-220
10mm	240mm	TI-1457S-10-240	1457S-10-240
11mm	180mm	TI-1457S-11-180	1457S-11-180
11mm	200mm	TI-1457S-11-200	1457S-11-200
11mm	220mm	TI-1457S-11-220	1457S-11-220
11mm	240mm	TI-1457S-11-240	1457S-11-240
12mm	180mm	TI-1457S-12-180	1457S-12-180
12mm	200mm	TI-1457S-12-200	1457S-12-200
12mm	220mm	TI-1457S-12-220	1457S-12-220
12mm	240mm	TI-1457S-12-240	1457S-12-240
13mm	180mm	TI-1457S-13-180	1457S-13-180
13mm	200mm	TI-1457S-13-200	1457S-13-200
13mm	220mm	TI-1457S-13-220	1457S-13-220
13mm	240mm	TI-1457S-13-240	1457S-13-240






130° Athos Nail Type-II, Short

STERILE

ØDia	Length	Titanium	Stainless Steel
9mm	180mm	TI-1457S-09-180-S	1457S-09-180-S
9mm	200mm	TI-1457S-09-200-S	1457S-09-200-S
9mm	220mm	TI-1457S-09-220-S	1457S-09-220-S
9mm	240mm	TI-1457S-09-240-S	1457S-09-240-S
10mm	180mm	TI-1457S-10-180-S	1457S-10-180-S
10mm	200mm	TI-1457S-10-200-S	1457S-10-200-S
10mm	220mm	TI-1457S-10-220-S	1457S-10-220-S
10mm	240mm	TI-1457S-10-240-S	1457S-10-240-S
11mm	180mm	TI-1457S-11-180-S	1457S-11-180-S
11mm	200mm	TI-1457S-11-200-S	1457S-11-200-S
11mm	220mm	TI-1457S-11-220-S	1457S-11-220-S
11mm	240mm	TI-1457S-11-240-S	1457S-11-240-S
12mm	180mm	TI-1457S-12-180-S	1457S-12-180-S
12mm	200mm	TI-1457S-12-200-S	1457S-12-200-S
12mm	220mm	TI-1457S-12-220-S	1457S-12-220-S
12mm	240mm	TI-1457S-12-240-S	1457S-12-240-S
13mm	180mm	TI-1457S-13-180-S	1457S-13-180-S
13mm	200mm	TI-1457S-13-200-S	1457S-13-200-S
13mm	220mm	TI-1457S-13-220-S	1457S-13-220-S
13mm	240mm	TI-1457S-13-240-S	1457S-13-240-S



130° Athos Nail Type-II, LONG




ØDia	Length	Titanium(Left)	Stainless Steel(Left)	Titanium(Right)	Stainless Steel(Right)
9mm	320mm	TI-1457L-09-320L	1457L-09-320L	TI-1457L-09-320R	1457L-09-320R
9mm	340mm	TI-1457L-09-340L	1457L-09-340L	TI-1457L-09-340R	1457L-09-340R
9mm	360mm	TI-1457L-09-360L	1457L-09-360L	TI-1457L-09-360R	1457L-09-360R
9mm	380mm	TI-1457L-09-380L	1457L-09-380L	TI-1457L-09-380R	1457L-09-380R
9mm	400mm	TI-1457L-09-400L	1457L-09-400L	TI-1457L-09-400R	1457L-09-400R
9mm	420mm	TI-1457L-09-420L	1457L-09-420L	TI-1457L-09-420R	1457L-09-420R
9mm	440mm	TI-1457L-09-440L	1457L-09-440L	TI-1457L-09-440R	1457L-09-440R
10mm	320mm	TI-1457L-10-320L	1457L-10-320L	TI-1457L-10-320R	1457L-10-320R
10mm	340mm	TI-1457L-10-340L	1457L-10-340L	TI-1457L-10-340R	1457L-10-340R
10mm	360mm	TI-1457L-10-360L	1457L-10-360L	TI-1457L-10-360R	1457L-10-360R
10mm	380mm	TI-1457L-10-380L	1457L-10-380L	TI-1457L-10-380R	1457L-10-380R
10mm	400mm	TI-1457L-10-400L	1457L-10-400L	TI-1457L-10-400R	1457L-10-400R
10mm	420mm	TI-1457L-10-420L	1457L-10-420L	TI-1457L-10-420R	1457L-10-420R
10mm	440mm	TI-1457L-10-440L	1457L-10-440L	TI-1457L-10-440R	1457L-10-440R
11mm	320mm	TI-1457L-11-320L	1457L-11-320L	TI-1457L-11-320R	1457L-11-320R
11mm	340mm	TI-1457L-11-340L	1457L-11-340L	TI-1457L-11-340R	1457L-11-340R
11mm	360mm	TI-1457L-11-360L	1457L-11-360L	TI-1457L-11-360R	1457L-11-360R
11mm	380mm	TI-1457L-11-380L	1457L-11-380L	TI-1457L-11-380R	1457L-11-380R
11mm	400mm	TI-1457L-11-400L	1457L-11-400L	TI-1457L-11-400R	1457L-11-400R
11mm	420mm	TI-1457L-11-420L	1457L-11-420L	TI-1457L-11-420R	1457L-11-420R
11mm	440mm	TI-1457L-11-440L	1457L-11-440L	TI-1457L-11-440R	1457L-11-440R
12mm	320mm	TI-1457L-12-320L	1457L-12-320L	TI-1457L-12-320R	1457L-12-320R
12mm	340mm	TI-1457L-12-340L	1457L-12-340L	TI-1457L-12-340R	1457L-12-340R
12mm	360mm	TI-1457L-12-360L	1457L-12-360L	TI-1457L-12-360R	1457L-12-360R
12mm	380mm	TI-1457L-12-380L	1457L-12-380L	TI-1457L-12-380R	1457L-12-380R
12mm	400mm	TI-1457L-12-400L	1457L-12-400L	TI-1457L-12-400R	1457L-12-400R
12mm	420mm	TI-1457L-12-420L	1457L-12-420L	TI-1457L-12-420R	1457L-12-420R
12mm	440mm	TI-1457L-12-440L	1457L-12-440L	TI-1457L-12-440R	1457L-12-440R
13mm	320mm	TI-1457L-13-320L	1457L-13-320L	TI-1457L-13-320R	1457L-13-320R
13mm	340mm	TI-1457L-13-340L	1457L-13-340L	TI-1457L-13-340R	1457L-13-340R
13mm	360mm	TI-1457L-13-360L	1457L-13-360L	TI-1457L-13-360R	1457L-13-360R
13mm	380mm	TI-1457L-13-380L	1457L-13-380L	TI-1457L-13-380R	1457L-13-380R
13mm	400mm	TI-1457L-13-400L	1457L-13-400L	TI-1457L-13-400R	1457L-13-400R
13mm	420mm	TI-1457L-13-420L	1457L-13-420L	TI-1457L-13-420R	1457L-13-420R
13mm	440mm	TI-1457L-13-440L	1457L-13-440L	TI-1457L-13-440R	1457L-13-440R

Left Right



130° Athos Nail Type-II, LONG



Left Right

STERILE

ØDia	Length	Titanium (Left)	Stainless Steel (Left)	Titanium (Right)	Stainless Steel (Right)
9mm	320mm	TI-1457L-09-320L-S	1457L-09-320L-S	TI-1457L-09-320R-S	1457L-09-320R-S
9mm	340mm	TI-1457L-09-340L-S	1457L-09-340L-S	TI-1457L-09-340R-S	1457L-09-340R-S
9mm	360mm	TI-1457L-09-360L-S	1457L-09-360L-S	TI-1457L-09-360R-S	1457L-09-360R-S
9mm	380mm	TI-1457L-09-380L-S	1457L-09-380L-S	TI-1457L-09-380R-S	1457L-09-380R-S
9mm	400mm	TI-1457L-09-400L-S	1457L-09-400L-S	TI-1457L-09-400R-S	1457L-09-400R-S
9mm	420mm	TI-1457L-09-420L-S	1457L-09-420L-S	TI-1457L-09-420R-S	1457L-09-420R-S
9mm	440mm	TI-1457L-09-440L-S	1457L-09-440L-S	TI-1457L-09-440R-S	1457L-09-440R-S
10mm	320mm	TI-1457L-10-320L-S	1457L-10-320L-S	TI-1457L-10-320R-S	1457L-10-320R-S
10mm	340mm	TI-1457L-10-340L-S	1457L-10-340L-S	TI-1457L-10-340R-S	1457L-10-340R-S
10mm	360mm	TI-1457L-10-360L-S	1457L-10-360L-S	TI-1457L-10-360R-S	1457L-10-360R-S
10mm	380mm	TI-1457L-10-380L-S	1457L-10-380L-S	TI-1457L-10-380R-S	1457L-10-380R-S
10mm	400mm	TI-1457L-10-400L-S	1457L-10-400L-S	TI-1457L-10-400R-S	1457L-10-400R-S
10mm	420mm	TI-1457L-10-420L-S	1457L-10-420L-S	TI-1457L-10-420R-S	1457L-10-420R-S
10mm	440mm	TI-1457L-10-440L-S	1457L-10-440L-S	TI-1457L-10-440R-S	1457L-10-440R-S
11mm	320mm	TI-1457L-11-320L-S	1457L-11-320L-S	TI-1457L-11-320R-S	1457L-11-320R-S
11mm	340mm	TI-1457L-11-340L-S	1457L-11-340L-S	TI-1457L-11-340R-S	1457L-11-340R-S
11mm	360mm	TI-1457L-11-360L-S	1457L-11-360L-S	TI-1457L-11-360R-S	1457L-11-360R-S
11mm	380mm	TI-1457L-11-380L-S	1457L-11-380L-S	TI-1457L-11-380R-S	1457L-11-380R-S
11mm	400mm	TI-1457L-11-400L-S	1457L-11-400L-S	TI-1457L-11-400R-S	1457L-11-400R-S
11mm	420mm	TI-1457L-11-420L-S	1457L-11-420L-S	TI-1457L-11-420R-S	1457L-11-420R-S
11mm	440mm	TI-1457L-11-440L-S	1457L-11-440L-S	TI-1457L-11-440R-S	1457L-11-440R-S
12mm	320mm	TI-1457L-12-320L-S	1457L-12-320L-S	TI-1457L-12-320R-S	1457L-12-320R-S
12mm	340mm	TI-1457L-12-340L-S	1457L-12-340L-S	TI-1457L-12-340R-S	1457L-12-340R-S
12mm	360mm	TI-1457L-12-360L-S	1457L-12-360L-S	TI-1457L-12-360R-S	1457L-12-360R-S
12mm	380mm	TI-1457L-12-380L-S	1457L-12-380L-S	TI-1457L-12-380R-S	1457L-12-380R-S
12mm	400mm	TI-1457L-12-400L-S	1457L-12-400L-S	TI-1457L-12-400R-S	1457L-12-400R-S
12mm	420mm	TI-1457L-12-420L-S	1457L-12-420L-S	TI-1457L-12-420R-S	1457L-12-420R-S
12mm	440mm	TI-1457L-12-440L-S	1457L-12-440L-S	TI-1457L-12-440R-S	1457L-12-440R-S
13mm	320mm	TI-1457L-13-320L-S	1457L-13-320L-S	TI-1457L-13-320R-S	1457L-13-320R-S
13mm	340mm	TI-1457L-13-340L-S	1457L-13-340L-S	TI-1457L-13-340R-S	1457L-13-340R-S
13mm	360mm	TI-1457L-13-360L-S	1457L-13-360L-S	TI-1457L-13-360R-S	1457L-13-360R-S
13mm	380mm	TI-1457L-13-380L-S	1457L-13-380L-S	TI-1457L-13-380R-S	1457L-13-380R-S
13mm	400mm	TI-1457L-13-400L-S	1457L-13-400L-S	TI-1457L-13-400R-S	1457L-13-400R-S
13mm	420mm	TI-1457L-13-420L-S	1457L-13-420L-S	TI-1457L-13-420R-S	1457L-13-420R-S
13mm	440mm	TI-1457L-13-440L-S	1457L-13-440L-S	TI-1457L-13-440R-S	1457L-13-440R-S



End Cap For Athos Nail Type-II



Length	Stainless Steel	Titanium
0mm	1457-001E	TI-1457-001E
5mm	1457-005E	TI-1457-005E
10mm	1457-010E	TI-1457-010E
15mm	1457-015E	TI-1457-015E



STERILE

Length	Stainless Steel	Titanium
0mm	1457-001E-S	TI-1457-001E-S
5mm	1457-005E-S	TI-1457-005E-S
10mm	1457-010E-S	TI-1457-010E-S
15mm	1457-015E-S	TI-1457-015E-S

Inner Screw For Athos Nail Type-II



Stainless Steel	Titanium
1457-0S	TI-1457-0S

STERILE

Stainless Steel	Titanium
1457-0S-S	TI-1457-0S-S



End Cap For Ø10.5mm Cephalic Screw For Athos Nail



Stainless Steel	Titanium
GM01.06	TI-GM01.06

STERILE

Stainless Steel	Titanium
GM01.06-S	TI-GM01.06-S





Ø10.5mm Cephalic Screw, For Athos Nail



Length	Titanium	Stainless Steel
70mm	TI-600C-10.5-070	600C-10.5-070
75mm	TI-600C-10.5-075	600C-10.5-075
80mm	TI-600C-10.5-080	600C-10.5-080
85mm	TI-600C-10.5-085	600C-10.5-085
90mm	TI-600C-10.5-090	600C-10.5-090
95mm	TI-600C-10.5-095	600C-10.5-095
100mm	TI-600C-10.5-100	600C-10.5-100
105mm	TI-600C-10.5-105	600C-10.5-105
110mm	TI-600C-10.5-110	600C-10.5-110
115mm	TI-600C-10.5-115	600C-10.5-115
120mm	TI-600C-10.5-120	600C-10.5-120



STERILE

Length	Titanium	Stainless Steel
70mm	TI-600C-10.5-070-S	600C-10.5-070-S
75mm	TI-600C-10.5-075-S	600C-10.5-075-S
80mm	TI-600C-10.5-080-S	600C-10.5-080-S
85mm	TI-600C-10.5-085-S	600C-10.5-085-S
90mm	TI-600C-10.5-090-S	600C-10.5-090-S
95mm	TI-600C-10.5-095-S	600C-10.5-095-S
100mm	TI-600C-10.5-100-S	600C-10.5-100-S
105mm	TI-600C-10.5-105-S	600C-10.5-105-S
110mm	TI-600C-10.5-110-S	600C-10.5-110-S
115mm	TI-600C-10.5-115-S	600C-10.5-115-S
120mm	TI-600C-10.5-120-S	600C-10.5-120-S



SCREWS





Ø4.8mm Locking Bolt, Self-Tapping, For Athos Nail



Length	Titanium	Stainless Steel
26mm	TI-600-026	600-026
28mm	TI-600-028	600-028
30mm	TI-600-030	600-030
32mm	TI-600-032	600-032
34mm	TI-600-034	600-034
36mm	TI-600-036	600-036
38mm	TI-600-038	600-038
40mm	TI-600-040	600-040
42mm	TI-600-042	600-042
44mm	TI-600-044	600-044
46mm	TI-600-046	600-046
48mm	TI-600-048	600-048
50mm	TI-600-050	600-050
52mm	TI-600-052	600-052
54mm	TI-600-054	600-054
56mm	TI-600-056	600-056
58mm	TI-600-058	600-058
60mm	TI-600-060	600-060
65mm	TI-600-065	600-065
70mm	TI-600-070	600-070
75mm	TI-600-075	600-075
80mm	TI-600-080	600-080
85mm	TI-600-085	600-085
90mm	TI-600-090	600-090

STERILE

Length	Titanium	Stainless Steel
26mm	TI-600-026-S	600-026-S
28mm	TI-600-028-S	600-028-S
30mm	TI-600-030-S	600-030-S
32mm	TI-600-032-S	600-032-S
34mm	TI-600-034-S	600-034-S
36mm	TI-600-036-S	600-036-S
38mm	TI-600-038-S	600-038-S
40mm	TI-600-040-S	600-040-S
42mm	TI-600-042-S	600-042-S
44mm	TI-600-044-S	600-044-S
46mm	TI-600-046-S	600-046-S
48mm	TI-600-048-S	600-048-S
50mm	TI-600-050-S	600-050-S
52mm	TI-600-052-S	600-052-S
54mm	TI-600-054-S	600-054-S
56mm	TI-600-056-S	600-056-S
58mm	TI-600-058-S	600-058-S
60mm	TI-600-060-S	600-060-S
65mm	TI-600-065-S	600-065-S
70mm	TI-600-070-S	600-070-S
75mm	TI-600-075-S	600-075-S
80mm	TI-600-080-S	600-080-S
85mm	TI-600-085-S	600-085-S
90mm	TI-600-090-S	600-090-S




IMPLANTS SET



1-003-02L


Implant Box for 125° or 130° Athos Nail Type-II(Left)



Note: 125° or 130° Athos Nail Type-II will be placed in trays as per order



Tray 1



Type	Angle	Dia	Length	Qty
Short	125° or 130°	Ø9mm	180mm to 240mm	1pc
Long	125° or 130°	Ø9mm	320mm to 440mm	1pc

Tray 2



Type	Angle	Dia	Length	Qty
Short	125° or 130°	Ø10mm	180mm to 240mm	1pc
Long	125° or 130°	Ø10mm	320mm to 440mm	1pc

Tray 3



Type	Angle	Dia	Length	Qty
Short	125° or 130°	Ø11mm	180mm to 240mm	1pc
Long	125° or 130°	Ø11mm	320mm to 440mm	1pc

Tray 4



Type	Angle	Dia	Length	Qty
Short	125° or 130°	Ø12mm	180mm to 240mm	1pc
Long	125° or 130°	Ø12mm	320mm to 440mm	1pc

Tray 5



Type	Angle	Dia	Length	Qty
Short	125° or 130°	Ø13mm	180mm to 240mm	1pc
Long	125° or 130°	Ø13mm	320mm to 440mm	1pc

Stainless Steel or Titanium Implants can be placed as per requirement


Note: 125° or 130° Athos Nail Type-II will be placed in trays as per order



**Ø4.8mm Locking Bolt,
Self-Tapping, For Athos Nail**

Dia	Length	Qty
Ø4.8mm	26mm to 90mm	4pcs

Stainless Steel or Titanium Implants can be placed as per requirement



**Cephalic Screw
For Athos Nail**

Dia	Length	Qty
Ø10.5mm	70mm to 120mm	1pc

**Inner Screw For
Athos Nail Type-II**

Qty
5pcs


**End Cap For
Athos Nail Type-II**

Length	Qty
0mm to 15mm	2pcs

Stainless Steel or Titanium Implants can be placed as per requirement

1-003-02R


Implant Box for 125° or 130° Athos Nail Type-II (Right)



Note: 125° or 130° Athos Nail Type-II will be placed in trays as per order



Tray 1



Type	Angle	Dia	Length	Qty
Short	125° or 130°	Ø9mm	180mm to 240mm	1pc
Long	125° or 130°	Ø9mm	320mm to 440mm	1pc

Tray 2



Type	Angle	Dia	Length	Qty
Short	125° or 130°	Ø10mm	180mm to 240mm	1pc
Long	125° or 130°	Ø10mm	320mm to 440mm	1pc

Tray 3



Type	Angle	Dia	Length	Qty
Short	125° or 130°	Ø11mm	180mm to 240mm	1pc
Long	125° or 130°	Ø11mm	320mm to 440mm	1pc

Tray 4



Type	Angle	Dia	Length	Qty
Short	125° or 130°	Ø12mm	180mm to 240mm	1pc
Long	125° or 130°	Ø12mm	320mm to 440mm	1pc

Tray 5



Type	Angle	Dia	Length	Qty
Short	125° or 130°	Ø13mm	180mm to 240mm	1pc
Long	125° or 130°	Ø13mm	320mm to 440mm	1pc

Stainless Steel or Titanium Implants can be placed as per requirement

Note: 125° or 130° Athos Nail Type-II will be placed in trays as per order



**Ø4.8mm Locking Bolt,
Self-Tapping, For Athos Nail**

Dia	Length	Qty
Ø4.8mm	26mm to 90mm	4pcs

Stainless Steel or Titanium Implants can be placed as per requirement



**Cephalic Screw
For Athos Nail**

Dia	Length	Qty
Ø10.5mm	70mm to 120mm	1pc

**Inner Screw For
Athos Nail Type-II**

Qty
5pcs

**End Cap For
Athos Nail Type-II**

Length	Qty
0mm to 15mm	2pcs

Stainless Steel or Titanium Implants can be placed as per requirement



ATHOS

Nailing System (With Anti-Rotation)

ATHOS Nailing System (Anti-Rotation) Type-II

The Type-II design provides robust fixation for trochanteric and subtrochanteric fractures. It offers a balance of strength, anatomical conformity, and versatility. Available in both short (180–240mm) and long (320–440mm) versions, with diameters ranging from 9mm to 13mm, Type-II nails are suited for both simple and complex fractures.

Key Features:

- Incorporates an Ø6.4mm cannulated anti-rotation screw for enhanced fixation.
- Length options available: 60–130mm for tailored stability.
- Provides angular and rotational control in unstable fracture configurations.
- Multiple distal locking choices for rigidity or dynamization.
- Radiolucent targeting device reduces intraoperative radiation exposure.
- Stainless steel and titanium options for durability and biocompatibility.

Indications:

- Unstable pertrochanteric fractures requiring rotational control.
- Intertrochanteric fractures with comminution or subtrochanteric extension.
- Pathological fractures needing additional angular stability.
- Osteoporotic bone with risk of implant migration.
- Complex revision cases requiring enhanced fixation.



125° Anti-Rotation Athos Nail, Short



ØDia	Length	Titanium	Stainless Steel
9mm	180mm	TI-1501S-09-180	1501S-09-180
9mm	200mm	TI-1501S-09-200	1501S-09-200
9mm	220mm	TI-1501S-09-220	1501S-09-220
9mm	240mm	TI-1501S-09-240	1501S-09-240
10mm	180mm	TI-1501S-10-180	1501S-10-180
10mm	200mm	TI-1501S-10-200	1501S-10-200
10mm	220mm	TI-1501S-10-220	1501S-10-220
10mm	240mm	TI-1501S-10-240	1501S-10-240
11mm	180mm	TI-1501S-11-180	1501S-11-180
11mm	200mm	TI-1501S-11-200	1501S-11-200
11mm	220mm	TI-1501S-11-220	1501S-11-220
11mm	240mm	TI-1501S-11-240	1501S-11-240
12mm	180mm	TI-1501S-12-180	1501S-12-180
12mm	200mm	TI-1501S-12-200	1501S-12-200
12mm	220mm	TI-1501S-12-220	1501S-12-220
12mm	240mm	TI-1501S-12-240	1501S-12-240
13mm	180mm	TI-1501S-13-180	1501S-13-180
13mm	200mm	TI-1501S-13-200	1501S-13-200
13mm	220mm	TI-1501S-13-220	1501S-13-220
13mm	240mm	TI-1501S-13-240	1501S-13-240





125° Anti-Rotation Athos Nail, Short

STERILE

ØDia	Length	Titanium	Stainless Steel
9mm	180mm	TI-1501S-09-180-S	1501S-09-180-S
9mm	200mm	TI-1501S-09-200-S	1501S-09-200-S
9mm	220mm	TI-1501S-09-220-S	1501S-09-220-S
9mm	240mm	TI-1501S-09-240-S	1501S-09-240-S
10mm	180mm	TI-1501S-10-180-S	1501S-10-180-S
10mm	200mm	TI-1501S-10-200-S	1501S-10-200-S
10mm	220mm	TI-1501S-10-220-S	1501S-10-220-S
10mm	240mm	TI-1501S-10-240-S	1501S-10-240-S
11mm	180mm	TI-1501S-11-180-S	1501S-11-180-S
11mm	200mm	TI-1501S-11-200-S	1501S-11-200-S
11mm	220mm	TI-1501S-11-220-S	1501S-11-220-S
11mm	240mm	TI-1501S-11-240-S	1501S-11-240-S
12mm	180mm	TI-1501S-12-180-S	1501S-12-180-S
12mm	200mm	TI-1501S-12-200-S	1501S-12-200-S
12mm	220mm	TI-1501S-12-220-S	1501S-12-220-S
12mm	240mm	TI-1501S-12-240-S	1501S-12-240-S
13mm	180mm	TI-1501S-13-180-S	1501S-13-180-S
13mm	200mm	TI-1501S-13-200-S	1501S-13-200-S
13mm	220mm	TI-1501S-13-220-S	1501S-13-220-S
13mm	240mm	TI-1501S-13-240-S	1501S-13-240-S



125° Anti-Rotation Athos Nail, Long



ØDia	Length	Titanium(Left)	Stainless Steel(Left)	Titanium(Right)	Stainless Steel(Right)
9mm	320mm	TI-1501L-09-320L	1501L-09-320L	TI-1501L-09-320R	1501L-09-320R
9mm	340mm	TI-1501L-09-340L	1501L-09-340L	TI-1501L-09-340R	1501L-09-340R
9mm	360mm	TI-1501L-09-360L	1501L-09-360L	TI-1501L-09-360R	1501L-09-360R
9mm	380mm	TI-1501L-09-380L	1501L-09-380L	TI-1501L-09-380R	1501L-09-380R
9mm	400mm	TI-1501L-09-400L	1501L-09-400L	TI-1501L-09-400R	1501L-09-400R
9mm	420mm	TI-1501L-09-420L	1501L-09-420L	TI-1501L-09-420R	1501L-09-420R
9mm	440mm	TI-1501L-09-440L	1501L-09-440L	TI-1501L-09-440R	1501L-09-440R
10mm	320mm	TI-1501L-10-320L	1501L-10-320L	TI-1501L-10-320R	1501L-10-320R
10mm	340mm	TI-1501L-10-340L	1501L-10-340L	TI-1501L-10-340R	1501L-10-340R
10mm	360mm	TI-1501L-10-360L	1501L-10-360L	TI-1501L-10-360R	1501L-10-360R
10mm	380mm	TI-1501L-10-380L	1501L-10-380L	TI-1501L-10-380R	1501L-10-380R
10mm	400mm	TI-1501L-10-400L	1501L-10-400L	TI-1501L-10-400R	1501L-10-400R
10mm	420mm	TI-1501L-10-420L	1501L-10-420L	TI-1501L-10-420R	1501L-10-420R
10mm	440mm	TI-1501L-10-440L	1501L-10-440L	TI-1501L-10-440R	1501L-10-440R
11mm	320mm	TI-1501L-11-320L	1501L-11-320L	TI-1501L-11-320R	1501L-11-320R
11mm	340mm	TI-1501L-11-340L	1501L-11-340L	TI-1501L-11-340R	1501L-11-340R
11mm	360mm	TI-1501L-11-360L	1501L-11-360L	TI-1501L-11-360R	1501L-11-360R
11mm	380mm	TI-1501L-11-380L	1501L-11-380L	TI-1501L-11-380R	1501L-11-380R
11mm	400mm	TI-1501L-11-400L	1501L-11-400L	TI-1501L-11-400R	1501L-11-400R
11mm	420mm	TI-1501L-11-420L	1501L-11-420L	TI-1501L-11-420R	1501L-11-420R
11mm	440mm	TI-1501L-11-440L	1501L-11-440L	TI-1501L-11-440R	1501L-11-440R
12mm	320mm	TI-1501L-12-320L	1501L-12-320L	TI-1501L-12-320R	1501L-12-320R
12mm	340mm	TI-1501L-12-340L	1501L-12-340L	TI-1501L-12-340R	1501L-12-340R
12mm	360mm	TI-1501L-12-360L	1501L-12-360L	TI-1501L-12-360R	1501L-12-360R
12mm	380mm	TI-1501L-12-380L	1501L-12-380L	TI-1501L-12-380R	1501L-12-380R
12mm	400mm	TI-1501L-12-400L	1501L-12-400L	TI-1501L-12-400R	1501L-12-400R
12mm	420mm	TI-1501L-12-420L	1501L-12-420L	TI-1501L-12-420R	1501L-12-420R
12mm	440mm	TI-1501L-12-440L	1501L-12-440L	TI-1501L-12-440R	1501L-12-440R
13mm	320mm	TI-1501L-13-320L	1501L-13-320L	TI-1501L-13-320R	1501L-13-320R
13mm	340mm	TI-1501L-13-340L	1501L-13-340L	TI-1501L-13-340R	1501L-13-340R
13mm	360mm	TI-1501L-13-360L	1501L-13-360L	TI-1501L-13-360R	1501L-13-360R
13mm	380mm	TI-1501L-13-380L	1501L-13-380L	TI-1501L-13-380R	1501L-13-380R
13mm	400mm	TI-1501L-13-400L	1501L-13-400L	TI-1501L-13-400R	1501L-13-400R
13mm	420mm	TI-1501L-13-420L	1501L-13-420L	TI-1501L-13-420R	1501L-13-420R
13mm	440mm	TI-1501L-13-440L	1501L-13-440L	TI-1501L-13-440R	1501L-13-440R



125° Anti-Rotation Athos Nail, Long



STERILE

ØDia	Length	Titanium(Left)	Stainless Steel(Left)	Titanium(Right)	Stainless Steel(Right)
9mm	320mm	TI-1501L-09-320L-S	1501L-09-320L-S	TI-1501L-09-320R-S	1501L-09-320R-S
9mm	340mm	TI-1501L-09-340L-S	1501L-09-340L-S	TI-1501L-09-340R-S	1501L-09-340R-S
9mm	360mm	TI-1501L-09-360L-S	1501L-09-360L-S	TI-1501L-09-360R-S	1501L-09-360R-S
9mm	380mm	TI-1501L-09-380L-S	1501L-09-380L-S	TI-1501L-09-380R-S	1501L-09-380R-S
9mm	400mm	TI-1501L-09-400L-S	1501L-09-400L-S	TI-1501L-09-400R-S	1501L-09-400R-S
9mm	420mm	TI-1501L-09-420L-S	1501L-09-420L-S	TI-1501L-09-420R-S	1501L-09-420R-S
9mm	440mm	TI-1501L-09-440L-S	1501L-09-440L-S	TI-1501L-09-440R-S	1501L-09-440R-S
10mm	320mm	TI-1501L-10-320L-S	1501L-10-320L-S	TI-1501L-10-320R-S	1501L-10-320R-S
10mm	340mm	TI-1501L-10-340L-S	1501L-10-340L-S	TI-1501L-10-340R-S	1501L-10-340R-S
10mm	360mm	TI-1501L-10-360L-S	1501L-10-360L-S	TI-1501L-10-360R-S	1501L-10-360R-S
10mm	380mm	TI-1501L-10-380L-S	1501L-10-380L-S	TI-1501L-10-380R-S	1501L-10-380R-S
10mm	400mm	TI-1501L-10-400L-S	1501L-10-400L-S	TI-1501L-10-400R-S	1501L-10-400R-S
10mm	420mm	TI-1501L-10-420L-S	1501L-10-420L-S	TI-1501L-10-420R-S	1501L-10-420R-S
10mm	440mm	TI-1501L-10-440L-S	1501L-10-440L-S	TI-1501L-10-440R-S	1501L-10-440R-S
11mm	320mm	TI-1501L-11-320L-S	1501L-11-320L-S	TI-1501L-11-320R-S	1501L-11-320R-S
11mm	340mm	TI-1501L-11-340L-S	1501L-11-340L-S	TI-1501L-11-340R-S	1501L-11-340R-S
11mm	360mm	TI-1501L-11-360L-S	1501L-11-360L-S	TI-1501L-11-360R-S	1501L-11-360R-S
11mm	380mm	TI-1501L-11-380L-S	1501L-11-380L-S	TI-1501L-11-380R-S	1501L-11-380R-S
11mm	400mm	TI-1501L-11-400L-S	1501L-11-400L-S	TI-1501L-11-400R-S	1501L-11-400R-S
11mm	420mm	TI-1501L-11-420L-S	1501L-11-420L-S	TI-1501L-11-420R-S	1501L-11-420R-S
11mm	440mm	TI-1501L-11-440L-S	1501L-11-440L-S	TI-1501L-11-440R-S	1501L-11-440R-S
12mm	320mm	TI-1501L-12-320L-S	1501L-12-320L-S	TI-1501L-12-320R-S	1501L-12-320R-S
12mm	340mm	TI-1501L-12-340L-S	1501L-12-340L-S	TI-1501L-12-340R-S	1501L-12-340R-S
12mm	360mm	TI-1501L-12-360L-S	1501L-12-360L-S	TI-1501L-12-360R-S	1501L-12-360R-S
12mm	380mm	TI-1501L-12-380L-S	1501L-12-380L-S	TI-1501L-12-380R-S	1501L-12-380R-S
12mm	400mm	TI-1501L-12-400L-S	1501L-12-400L-S	TI-1501L-12-400R-S	1501L-12-400R-S
12mm	420mm	TI-1501L-12-420L-S	1501L-12-420L-S	TI-1501L-12-420R-S	1501L-12-420R-S
12mm	440mm	TI-1501L-12-440L-S	1501L-12-440L-S	TI-1501L-12-440R-S	1501L-12-440R-S
13mm	320mm	TI-1501L-13-320L-S	1501L-13-320L-S	TI-1501L-13-320R-S	1501L-13-320R-S
13mm	340mm	TI-1501L-13-340L-S	1501L-13-340L-S	TI-1501L-13-340R-S	1501L-13-340R-S
13mm	360mm	TI-1501L-13-360L-S	1501L-13-360L-S	TI-1501L-13-360R-S	1501L-13-360R-S
13mm	380mm	TI-1501L-13-380L-S	1501L-13-380L-S	TI-1501L-13-380R-S	1501L-13-380R-S
13mm	400mm	TI-1501L-13-400L-S	1501L-13-400L-S	TI-1501L-13-400R-S	1501L-13-400R-S
13mm	420mm	TI-1501L-13-420L-S	1501L-13-420L-S	TI-1501L-13-420R-S	1501L-13-420R-S
13mm	440mm	TI-1501L-13-440L-S	1501L-13-440L-S	TI-1501L-13-440R-S	1501L-13-440R-S



130° Anti-Rotation Athos Nail, Short



ØDia	Length	Titanium	Stainless Steel
9mm	180mm	TI-600S-09-180	600S-09-180
9mm	200mm	TI-600S-09-200	600S-09-200
9mm	220mm	TI-600S-09-220	600S-09-220
9mm	240mm	TI-600S-09-240	600S-09-240
10mm	180mm	TI-600S-10-180	600S-10-180
10mm	200mm	TI-600S-10-200	600S-10-200
10mm	220mm	TI-600S-10-220	600S-10-220
10mm	240mm	TI-600S-10-240	600S-10-240
11mm	180mm	TI-600S-11-180	600S-11-180
11mm	200mm	TI-600S-11-200	600S-11-200
11mm	220mm	TI-600S-11-220	600S-11-220
11mm	240mm	TI-600S-11-240	600S-11-240
12mm	180mm	TI-600S-12-180	600S-12-180
12mm	200mm	TI-600S-12-200	600S-12-200
12mm	220mm	TI-600S-12-220	600S-12-220
12mm	240mm	TI-600S-12-240	600S-12-240
13mm	180mm	TI-600S-13-180	600S-13-180
13mm	200mm	TI-600S-13-200	600S-13-200
13mm	220mm	TI-600S-13-220	600S-13-220
13mm	240mm	TI-600S-13-240	600S-13-240





130° Anti-Rotation Athos Nail, Short

STERILE

ØDia	Length	Titanium	Stainless Steel
9mm	180mm	TI-600S-09-180-S	600S-09-180-S
9mm	200mm	TI-600S-09-200-S	600S-09-200-S
9mm	220mm	TI-600S-09-220-S	600S-09-220-S
9mm	240mm	TI-600S-09-240-S	600S-09-240-S
10mm	180mm	TI-600S-10-180-S	600S-10-180-S
10mm	200mm	TI-600S-10-200-S	600S-10-200-S
10mm	220mm	TI-600S-10-220-S	600S-10-220-S
10mm	240mm	TI-600S-10-240-S	600S-10-240-S
11mm	180mm	TI-600S-11-180-S	600S-11-180-S
11mm	200mm	TI-600S-11-200-S	600S-11-200-S
11mm	220mm	TI-600S-11-220-S	600S-11-220-S
11mm	240mm	TI-600S-11-240-S	600S-11-240-S
12mm	180mm	TI-600S-12-180-S	600S-12-180-S
12mm	200mm	TI-600S-12-200-S	600S-12-200-S
12mm	220mm	TI-600S-12-220-S	600S-12-220-S
12mm	240mm	TI-600S-12-240-S	600S-12-240-S
13mm	180mm	TI-600S-13-180-S	600S-13-180-S
13mm	200mm	TI-600S-13-200-S	600S-13-200-S
13mm	220mm	TI-600S-13-220-S	600S-13-220-S
13mm	240mm	TI-600S-13-240-S	600S-13-240-S





130° Anti-Rotation Athos Nail, Long



ØDia	Length	Titanium(Left)	Stainless Steel(Left)	Titanium(Right)	Stainless Steel(Right)
9mm	320mm	TI-600L-09-320L	600L-09-320L	TI-600L-09-320R	600L-09-320R
9mm	340mm	TI-600L-09-340L	600L-09-340L	TI-600L-09-340R	600L-09-340R
9mm	360mm	TI-600L-09-360L	600L-09-360L	TI-600L-09-360R	600L-09-360R
9mm	380mm	TI-600L-09-380L	600L-09-380L	TI-600L-09-380R	600L-09-380R
9mm	400mm	TI-600L-09-400L	600L-09-400L	TI-600L-09-400R	600L-09-400R
9mm	420mm	TI-600L-09-420L	600L-09-420L	TI-600L-09-420R	600L-09-420R
9mm	440mm	TI-600L-09-440L	600L-09-440L	TI-600L-09-440R	600L-09-440R
10mm	320mm	TI-600L-10-320L	600L-10-320L	TI-600L-10-320R	600L-10-320R
10mm	340mm	TI-600L-10-340L	600L-10-340L	TI-600L-10-340R	600L-10-340R
10mm	360mm	TI-600L-10-360L	600L-10-360L	TI-600L-10-360R	600L-10-360R
10mm	380mm	TI-600L-10-380L	600L-10-380L	TI-600L-10-380R	600L-10-380R
10mm	400mm	TI-600L-10-400L	600L-10-400L	TI-600L-10-400R	600L-10-400R
10mm	420mm	TI-600L-10-420L	600L-10-420L	TI-600L-10-420R	600L-10-420R
10mm	440mm	TI-600L-10-440L	600L-10-440L	TI-600L-10-440R	600L-10-440R
11mm	320mm	TI-600L-11-320L	600L-11-320L	TI-600L-11-320R	600L-11-320R
11mm	340mm	TI-600L-11-340L	600L-11-340L	TI-600L-11-340R	600L-11-340R
11mm	360mm	TI-600L-11-360L	600L-11-360L	TI-600L-11-360R	600L-11-360R
11mm	380mm	TI-600L-11-380L	600L-11-380L	TI-600L-11-380R	600L-11-380R
11mm	400mm	TI-600L-11-400L	600L-11-400L	TI-600L-11-400R	600L-11-400R
11mm	420mm	TI-600L-11-420L	600L-11-420L	TI-600L-11-420R	600L-11-420R
11mm	440mm	TI-600L-11-440L	600L-11-440L	TI-600L-11-440R	600L-11-440R
12mm	320mm	TI-600L-12-320L	600L-12-320L	TI-600L-12-320R	600L-12-320R
12mm	340mm	TI-600L-12-340L	600L-12-340L	TI-600L-12-340R	600L-12-340R
12mm	360mm	TI-600L-12-360L	600L-12-360L	TI-600L-12-360R	600L-12-360R
12mm	380mm	TI-600L-12-380L	600L-12-380L	TI-600L-12-380R	600L-12-380R
12mm	400mm	TI-600L-12-400L	600L-12-400L	TI-600L-12-400R	600L-12-400R
12mm	420mm	TI-600L-12-420L	600L-12-420L	TI-600L-12-420R	600L-12-420R
12mm	440mm	TI-600L-12-440L	600L-12-440L	TI-600L-12-440R	600L-12-440R
13mm	320mm	TI-600L-13-320L	600L-13-320L	TI-600L-13-320R	600L-13-320R
13mm	340mm	TI-600L-13-340L	600L-13-340L	TI-600L-13-340R	600L-13-340R
13mm	360mm	TI-600L-13-360L	600L-13-360L	TI-600L-13-360R	600L-13-360R
13mm	380mm	TI-600L-13-380L	600L-13-380L	TI-600L-13-380R	600L-13-380R
13mm	400mm	TI-600L-13-400L	600L-13-400L	TI-600L-13-400R	600L-13-400R
13mm	420mm	TI-600L-13-420L	600L-13-420L	TI-600L-13-420R	600L-13-420R
13mm	440mm	TI-600L-13-440L	600L-13-440L	TI-600L-13-440R	600L-13-440R



130° Anti-Rotation Athos Nail, Long



Left

Right

STERILE

ØDia	Length	Titanium(Left)	Stainless Steel(Left)	Titanium(Right)	Stainless Steel(Right)
9mm	320mm	TI-600L-09-320L-S	600L-09-320L-S	TI-600L-09-320R-S	600L-09-320R-S
9mm	340mm	TI-600L-09-340L-S	600L-09-340L-S	TI-600L-09-340R-S	600L-09-340R-S
9mm	360mm	TI-600L-09-360L-S	600L-09-360L-S	TI-600L-09-360R-S	600L-09-360R-S
9mm	380mm	TI-600L-09-380L-S	600L-09-380L-S	TI-600L-09-380R-S	600L-09-380R-S
9mm	400mm	TI-600L-09-400L-S	600L-09-400L-S	TI-600L-09-400R-S	600L-09-400R-S
9mm	420mm	TI-600L-09-420L-S	600L-09-420L-S	TI-600L-09-420R-S	600L-09-420R-S
9mm	440mm	TI-600L-09-440L-S	600L-09-440L-S	TI-600L-09-440R-S	600L-09-440R-S
10mm	320mm	TI-600L-10-320L-S	600L-10-320L-S	TI-600L-10-320R-S	600L-10-320R-S
10mm	340mm	TI-600L-10-340L-S	600L-10-340L-S	TI-600L-10-340R-S	600L-10-340R-S
10mm	360mm	TI-600L-10-360L-S	600L-10-360L-S	TI-600L-10-360R-S	600L-10-360R-S
10mm	380mm	TI-600L-10-380L-S	600L-10-380L-S	TI-600L-10-380R-S	600L-10-380R-S
10mm	400mm	TI-600L-10-400L-S	600L-10-400L-S	TI-600L-10-400R-S	600L-10-400R-S
10mm	420mm	TI-600L-10-420L-S	600L-10-420L-S	TI-600L-10-420R-S	600L-10-420R-S
10mm	440mm	TI-600L-10-440L-S	600L-10-440L-S	TI-600L-10-440R-S	600L-10-440R-S
11mm	320mm	TI-600L-11-320L-S	600L-11-320L-S	TI-600L-11-320R-S	600L-11-320R-S
11mm	340mm	TI-600L-11-340L-S	600L-11-340L-S	TI-600L-11-340R-S	600L-11-340R-S
11mm	360mm	TI-600L-11-360L-S	600L-11-360L-S	TI-600L-11-360R-S	600L-11-360R-S
11mm	380mm	TI-600L-11-380L-S	600L-11-380L-S	TI-600L-11-380R-S	600L-11-380R-S
11mm	400mm	TI-600L-11-400L-S	600L-11-400L-S	TI-600L-11-400R-S	600L-11-400R-S
11mm	420mm	TI-600L-11-420L-S	600L-11-420L-S	TI-600L-11-420R-S	600L-11-420R-S
11mm	440mm	TI-600L-11-440L-S	600L-11-440L-S	TI-600L-11-440R-S	600L-11-440R-S
12mm	320mm	TI-600L-12-320L-S	600L-12-320L-S	TI-600L-12-320R-S	600L-12-320R-S
12mm	340mm	TI-600L-12-340L-S	600L-12-340L-S	TI-600L-12-340R-S	600L-12-340R-S
12mm	360mm	TI-600L-12-360L-S	600L-12-360L-S	TI-600L-12-360R-S	600L-12-360R-S
12mm	380mm	TI-600L-12-380L-S	600L-12-380L-S	TI-600L-12-380R-S	600L-12-380R-S
12mm	400mm	TI-600L-12-400L-S	600L-12-400L-S	TI-600L-12-400R-S	600L-12-400R-S
12mm	420mm	TI-600L-12-420L-S	600L-12-420L-S	TI-600L-12-420R-S	600L-12-420R-S
12mm	440mm	TI-600L-12-440L-S	600L-12-440L-S	TI-600L-12-440R-S	600L-12-440R-S
13mm	320mm	TI-600L-13-320L-S	600L-13-320L-S	TI-600L-13-320R-S	600L-13-320R-S
13mm	340mm	TI-600L-13-340L-S	600L-13-340L-S	TI-600L-13-340R-S	600L-13-340R-S
13mm	360mm	TI-600L-13-360L-S	600L-13-360L-S	TI-600L-13-360R-S	600L-13-360R-S
13mm	380mm	TI-600L-13-380L-S	600L-13-380L-S	TI-600L-13-380R-S	600L-13-380R-S
13mm	400mm	TI-600L-13-400L-S	600L-13-400L-S	TI-600L-13-400R-S	600L-13-400R-S
13mm	420mm	TI-600L-13-420L-S	600L-13-420L-S	TI-600L-13-420R-S	600L-13-420R-S
13mm	440mm	TI-600L-13-440L-S	600L-13-440L-S	TI-600L-13-440R-S	600L-13-440R-S





End Cap For Anti-Rotation Athos Nail



Length	Stainless Steel	Titanium
0	600-001E	TI-600-001E
5	600-005E	TI-600-005E
10	600-010E	TI-600-010E
15	600-015E	TI-600-015E



STERILE

Length	Stainless Steel	Titanium
0	600-001E-S	TI-600-001E-S
5	600-005E-S	TI-600-005E-S
10	600-010E-S	TI-600-010E-S
15	600-015E-S	TI-600-015E-S

Inner Screw For Anti-Rotation Athos Nail-Long (For Cephalic Screw)



Stainless Steel	Titanium
600-0S	TI-600-0S

STERILE

Stainless Steel	Titanium
600-0S-S	TI-600-0S-S



Inner Screw For Anti-Rotation Athos Nail-Short (For Anti-Rotation Screw)



Stainless Steel	Titanium
601-0S	TI-601-0S

STERILE

Stainless Steel	Titanium
601-0S-S	TI-601-0S-S



Ø6.4mm Cannulated Anti-Rotation Screw for Athos Nail



STERILE

Length	Titanium	Stainless Steel
60	TI-600A-6.4-060	600A-6.4-060
65	TI-600A-6.4-065	600A-6.4-065
70	TI-600A-6.4-070	600A-6.4-070
75	TI-600A-6.4-075	600A-6.4-075
80	TI-600A-6.4-080	600A-6.4-080
85	TI-600A-6.4-085	600A-6.4-085
90	TI-600A-6.4-090	600A-6.4-090
95	TI-600A-6.4-095	600A-6.4-095
100	TI-600A-6.4-100	600A-6.4-100
105	TI-600A-6.4-105	600A-6.4-105
110	TI-600A-6.4-110	600A-6.4-110
115	TI-600A-6.4-115	600A-6.4-115
120	TI-600A-6.4-120	600A-6.4-120
125	TI-600A-6.4-125	600A-6.4-125
130	TI-600A-6.4-130	600A-6.4-130

Length	Titanium	Stainless Steel
60	TI-600A-6.4-060-S	600A-6.4-060-S
65	TI-600A-6.4-065-S	600A-6.4-065-S
70	TI-600A-6.4-070-S	600A-6.4-070-S
75	TI-600A-6.4-075-S	600A-6.4-075-S
80	TI-600A-6.4-080-S	600A-6.4-080-S
85	TI-600A-6.4-085-S	600A-6.4-085-S
90	TI-600A-6.4-090-S	600A-6.4-090-S
95	TI-600A-6.4-095-S	600A-6.4-095-S
100	TI-600A-6.4-100-S	600A-6.4-100-S
105	TI-600A-6.4-105-S	600A-6.4-105-S
110	TI-600A-6.4-110-S	600A-6.4-110-S
115	TI-600A-6.4-115-S	600A-6.4-115-S
120	TI-600A-6.4-120-S	600A-6.4-120-S
125	TI-600A-6.4-125-S	600A-6.4-125-S
130	TI-600A-6.4-130-S	600A-6.4-130-S



Ø10.5mm Cephalic Screw, For Athos Nail



Length	Titanium	Stainless Steel
70	TI-600C-10.5-070	600C-10.5-070
75	TI-600C-10.5-075	600C-10.5-075
80	TI-600C-10.5-080	600C-10.5-080
85	TI-600C-10.5-085	600C-10.5-085
90	TI-600C-10.5-090	600C-10.5-090
95	TI-600C-10.5-095	600C-10.5-095
100	TI-600C-10.5-100	600C-10.5-100
105	TI-600C-10.5-105	600C-10.5-105
110	TI-600C-10.5-110	600C-10.5-110
115	TI-600C-10.5-115	600C-10.5-115
120	TI-600C-10.5-120	600C-10.5-120

STERILE

Length	Titanium	Stainless Steel
70	TI-600C-10.5-070-S	600C-10.5-070-S
75	TI-600C-10.5-075-S	600C-10.5-075-S
80	TI-600C-10.5-080-S	600C-10.5-080-S
85	TI-600C-10.5-085-S	600C-10.5-085-S
90	TI-600C-10.5-090-S	600C-10.5-090-S
95	TI-600C-10.5-095-S	600C-10.5-095-S
100	TI-600C-10.5-100-S	600C-10.5-100-S
105	TI-600C-10.5-105-S	600C-10.5-105-S
110	TI-600C-10.5-110-S	600C-10.5-110-S
115	TI-600C-10.5-115-S	600C-10.5-115-S
120	TI-600C-10.5-120-S	600C-10.5-120-S

End Cap For Ø10.5mm Cephalic Screw For Athos Nail



Stainless Steel	Titanium
GM01.06	TI-GM01.06

STERILE

Stainless Steel	Titanium
GM01.06-S	TI-GM01.06-S



SCREWS





Ø4.8mm Locking Bolt, Self-Tapping, For Athos Nail



STERILE

Length	Titanium	Stainless Steel
26	TI-600-026	600-026
28	TI-600-028	600-028
30	TI-600-030	600-030
32	TI-600-032	600-032
34	TI-600-034	600-034
36	TI-600-036	600-036
38	TI-600-038	600-038
40	TI-600-040	600-040
42	TI-600-042	600-042
44	TI-600-044	600-044
46	TI-600-046	600-046
48	TI-600-048	600-048
50	TI-600-050	600-050
52	TI-600-052	600-052
54	TI-600-054	600-054
56	TI-600-056	600-056
58	TI-600-058	600-058
60	TI-600-060	600-060
65	TI-600-065	600-065
70	TI-600-070	600-070
75	TI-600-075	600-075
80	TI-600-080	600-080
85	TI-600-085	600-085
90	TI-600-090	600-090

Length	Titanium	Stainless Steel
26	TI-600-026-S	600-026-S
28	TI-600-028-S	600-028-S
30	TI-600-030-S	600-030-S
32	TI-600-032-S	600-032-S
34	TI-600-034-S	600-034-S
36	TI-600-036-S	600-036-S
38	TI-600-038-S	600-038-S
40	TI-600-040-S	600-040-S
42	TI-600-042-S	600-042-S
44	TI-600-044-S	600-044-S
46	TI-600-046-S	600-046-S
48	TI-600-048-S	600-048-S
50	TI-600-050-S	600-050-S
52	TI-600-052-S	600-052-S
54	TI-600-054-S	600-054-S
56	TI-600-056-S	600-056-S
58	TI-600-058-S	600-058-S
60	TI-600-060-S	600-060-S
65	TI-600-065-S	600-065-S
70	TI-600-070-S	600-070-S
75	TI-600-075-S	600-075-S
80	TI-600-080-S	600-080-S
85	TI-600-085-S	600-085-S
90	TI-600-090-S	600-090-S




IMPLANT SET





1-003-01L


Implant Box for Athos Nail Anti-Rotation (Left)



Note: 125° or 130° Anti-Rotation Athos Nail will be placed in trays as per order




Tray 1




Type	Angle	Dia	Length	Qty
Short	125° or 130°	Ø9mm	180mm to 240mm	1pc
Long	125° or 130°	Ø9mm	320mm to 440mm	1pc

Tray 2




Type	Angle	Dia	Length	Qty
Short	125° or 130°	Ø10mm	180mm to 240mm	1pc
Long	125° or 130°	Ø10mm	320mm to 440mm	1pc

Tray 3




Type	Angle	Dia	Length	Qty
Short	125° or 130°	Ø11mm	180mm to 240mm	1pc
Long	125° or 130°	Ø11mm	320mm to 440mm	1pc

Tray 4



Type	Angle	Dia	Length	Qty
Short	125° or 130°	Ø12mm	180mm to 240mm	1pc
Long	125° or 130°	Ø12mm	320mm to 440mm	1pc

Tray 5



Type	Angle	Dia	Length	Qty
Short	125° or 130°	Ø13mm	180mm to 240mm	1pc
Long	125° or 130°	Ø13mm	320mm to 440mm	1pc

Stainless Steel or Titanium Implants can be placed as per requirement

Note: 125° or 130° Anti-Rotation Athos Nail will be placed in trays as per order



**Ø4.8mm Locking Bolt,
Self-Tapping, For Athos Nail**

Dia	Length	Qty
Ø4.8mm	26mm to 90mm	3pcs

Stainless Steel or Titanium Implants can be placed as per requirement

**Ø4.8mm Locking Bolt,
Self-Tapping, For Athos Nail**

Dia	Length	Qty
Ø6.4mm	60mm to 130mm	2pcs

Stainless Steel or Titanium Implants can be placed as per requirement



**Cephalic Screw
For Athos Nail**

Dia	Length	Qty
Ø10.5mm	70mm to 120mm	1pc

**Inner Screw For Anti-Rotation
Athos Nail (For Cephalic Screw)**

For Long Nail	For Short Nail
5pcs	5pcs


**End Cap For
Athos Nail Type-II**

Length	Qty
0mm to 15mm	2pcs

Stainless Steel or Titanium Implants can be placed as per requirement

1-003-01R


Implant Box for Athos Nail Anti-Rotation (Right)



Note: 125° or 130° Anti-Rotation Athos Nail will be placed in trays as per order




Tray 1




Type	Angle	Dia	Length	Qty
Short	125° or 130°	Ø9mm	180mm to 240mm	1pc
Long	125° or 130°	Ø9mm	320mm to 440mm	1pc

Tray 2




Type	Angle	Dia	Length	Qty
Short	125° or 130°	Ø10mm	180mm to 240mm	1pc
Long	125° or 130°	Ø10mm	320mm to 440mm	1pc

Tray 3




Type	Angle	Dia	Length	Qty
Short	125° or 130°	Ø11mm	180mm to 240mm	1pc
Long	125° or 130°	Ø11mm	320mm to 440mm	1pc

Tray 4



Type	Angle	Dia	Length	Qty
Short	125° or 130°	Ø12mm	180mm to 240mm	1pc
Long	125° or 130°	Ø12mm	320mm to 440mm	1pc

Tray 5



Type	Angle	Dia	Length	Qty
Short	125° or 130°	Ø13mm	180mm to 240mm	1pc
Long	125° or 130°	Ø13mm	320mm to 440mm	1pc

Stainless Steel or Titanium Implants can be placed as per requirement
 Note: 125° or 130° Anti-Rotation Athos Nail will be placed in trays as per order



**Ø4.8mm Locking Bolt,
Self-Tapping, For Athos Nail**

Dia	Length	Qty
Ø4.8mm	26mm to 90mm	3pcs

Stainless Steel or Titanium Implants can be placed as per requirement

**Ø4.8mm Locking Bolt,
Self-Tapping, For Athos Nail**

Dia	Length	Qty
Ø6.4mm	60mm to 130mm	2pcs

Stainless Steel or Titanium Implants can be placed as per requirement



Cephalic Screw For Athos Nail

Dia	Length	Qty
Ø10.5mm	70mm to 120mm	1pc

**Inner Screw For Anti-Rotation
Athos Nail (For Cephalic Screw)**

For Long Nail	For Short Nail
5pcs	5pcs

End Cap For Athos Nail Type-II

Length	Qty
0mm to 15mm	2pcs

Stainless Steel or Titanium Implants can be placed as per requirement



INSTRUMENTS

ALUMINIUM TARGETING DEVICE



1457-001

Radiographic Ruler for Short Athos Nail




1457-055

Entry Reamer Protection Sleeve Ø20/17mm for Athos Nail



1457-056

Entry Reamer Guide Wire Sleeve Ø17/3.2mm for Athos Nail




1457-004

Guide Wire With Ball Tip Ø2.5/3.65mm x Length 680mm - Nitinol, for Athos Nail



1457-053

Guide Wire Without Ball Tip Ø2.5mm x Length 680mm - Nitinol, for Athos Nail



1457-005-10

Flexible Shaft with Reamer Ø10mm for Athos Nail



1457-006-11

Flexible Shaft with Reamer \varnothing 11mm for Athos Nail



1457-007-12

Flexible Shaft with Reamer \varnothing 12mm for Athos Nail



1457-008-13

Flexible Shaft with Reamer \varnothing 13mm for Athos Nail



1457-009

T Handle with Quick Coupling for Entry Reamer - Athos Nail



1457-010

Entry Reamer \varnothing 16.5mm with Quick Coupling End for Athos Nail



1457-011

Cannulated Awl for Athos Nail



1457-012

Tissue Protector for Athos Nail



1457-013

Insertion Handle, Aluminium, for Athos Nail



1457-014

Ø4.4mm Cannulated Nail Holding Bolt for Aluminium Targeting Device for Athos Nail



1457-015

Guide Wire Holder for Athos Nail



1457-016

Insertion Connecting Rod for Athos Nail- M10 for Athos Nail



1457-054

Block for Insertion Connecting Rod - Athos Nail



1457-017

Combination Wrench 11mm for Athos Nail



1457-018

Insertion Hammer, Solid, for Athos Nail



1457-019-130

130° Proximal Targeting Device, Aluminium, for Athos Nail



1457-019-125

125° Proximal Targeting Device, Aluminium, for Athos Nail



7-062-01

Locking Bolt for Aluminium Targeting Device - Athos Nail



1457-020

Protection Sleeve $\varnothing 10/8.2\text{mm}$ for Athos Nail



1457-021

Trocar Ø8mm for Athos Nail



1457-022-4.3

Drill Sleeve Ø8/4.3mm for Drill Bit Ø4.3mm - Athos Nail



1457-023

Short Aluminium Distal Targeting Device for Length upto 240mm - Athos Nail



1457-024

Sliding Hammer, Slotted, for Athos Nail



1457-025

Flexible Screwdriver, Hex 3.5mm for Inner Screw - Athos Nail



1457-026

Allen Key, Hex 3mm for Athos Nail



1457-027

Universal Screwdriver, Hex 7mm for Athos Nail



1457-028

Protection Sleeve Ø14/11.2mm for Cephalic Screw - Athos Nail



1457-029

Trocar Sleeve Ø11.2/3.2mm for Cephalic Screw - Athos Nail



1457-030

Trocar Ø3.2mm for Cephalic Screw Ø10.5mm - Athos Nail



1457-031

Measuring Device upto 130mm for Cephalic Screw Ø10.5mm - Athos Nail



1457-032

Insertor for Cephalic Screw Ø10.5mm & End Cap - Athos Nail



1457-033

Screwdriver with T-Handle, Hex 3.5mm, for Athos Nail



1457-034

Universal Screwdriver with T-Handle, Hex 8mm, for Athos Nail



1457-035

Proximal Cortex Drill Ø10mm with Quick Coupling End for Cephalic Screw Ø10.5mm - Athos Nail



1457-036

Cannulated Reamer, Ø10.5mm for Cephalic Screw Ø10.5mm - Athos Nail



7-007-08

Stopper for Cannulated Reamer - Athos Nail



1457-037

Measuring Device upto 130mm for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail



1457-038

Protection Sleeve 10/8.2mm for Ø6.4mm Cannulated Anti - Rotation Screw, for Athos Nail



1457-039

Trocar Sleeve 8.1/2.5mm for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail



1457-040

Trocar Ø2.5mm for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail



1457-041

Insertor for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail



1457-042

Tapered Drill 6.4/4.5mm for Ø6.4mm Cannulated Anti -Rotation Screw for Athos Nail



7-007-09

Stopper for Tapered Drill 6.4/4.5mm - Athos Nail



1457-043

Holding Driver, Hex 3.5 for Inner Screw, for Athos Nail



1457-044-2.5

Guide Wire Ø2.5mm x Length 380mm for Ø6.4mm Cannulated Anti -Rotation Screw



1457-045-3.0

Guide Wire with Threaded Tip Ø3.0mm x Thread Length 12mm x Length 400mm for Cephalic Screw Ø10.5mm



1457-046

Extraction Hammering Rod for Athos Nail



1457-047

Depth Gauge measuring upto 100mm for Athos Nail



1457-048

Depth Gauge measuring upto 115mm for Athos Nail



1457-049

Styler, Ø3mm, for Athos Nail




1457-050

Compression Screw for Athos Nail



1457-051-4.3

Drill Bit Plain Shank/Jacob Chuck End Ø4.3mm x Length 300mm for Athos Nail



1457-052

Stopper for Drill Bit Ø4.3mm - Athos Nail



1457-057

Curved Awl- Ø9.5-Ø5.6mm for Athos Nail



1457-00T

Instrument Trays for Athos Nailing Instrument Set,
Aluminium Targeting Device



1457-00B

Container with PPSU Lid for Athos Nailing Instrument Set,
Aluminium Targeting Device



OPTIONAL INSTRUMENT


1457-004S

Guide Wire With Ball Tip $\varnothing 2.5/3.65\text{mm}$ x Length 680mm -
Stainless Steel, for Athos Nail



1457-053S



Guide Wire Without Ball Tip $\varnothing 2.5\text{mm}$ x Length 680mm -
Stainless Steel, for Athos Nail





1457-000S

Athos Nailing Instrument Set, Aluminium Targeting Device



1457-000S**Athos Nailing Instrument Set, Aluminium Targeting Device**

Code	Set Consisting of	Units
1457-001	Radiographic Ruler for Short Athos Nail	1
1457-055	Entry Reamer Protection Sleeve Ø20/17mm for Athos Nail	1
1457-056	Entry Reamer Guide Wire Sleeve Ø17/3.2mm for Athos Nail	1
1457-004	Guide Wire With Ball Tip Ø2.5/3.65mm x Length 680mm - Nitinol, for Athos Nail	1
1457-053	Guide Wire Without Ball Tip Ø2.5mm x Length 680mm - Nitinol, for Athos Nail	1
1457-005-10	Flexible Shaft with Reamer Ø10mm for Athos Nail	1
1457-006-11	Flexible Shaft with Reamer Ø11mm for Athos Nail	1
1457-007-12	Flexible Shaft with Reamer Ø12mm for Athos Nail	1
1457-008-13	Flexible Shaft with Reamer Ø13mm for Athos Nail	1
1457-009	T Handle with Quick Coupling for Entry Reamer - Athos Nail	1
1457-010	Entry Reamer Ø16.5mm with Quick Coupling End for Athos Nail	1
1457-011	Cannulated Awl for Athos Nail	1
1457-012	Tissue Protector for Athos Nail	1
1457-013	Insertion Handle, Aluminium, for Athos Nail	1
1457-014	Nail Holding Bolt for Aluminium Targeting Device for Athos Nail	2
1457-015	Guide Wire Holder for Athos Nail	1
1457-016	Insertion Connecting Rod for Athos Nail	1
1457-054	Block for Insertion Connecting Rod - Athos Nail	1
1457-017	Combination Wrench 11mm for Athos Nail	1
1457-018	Insertion Hammer, Solid, for Athos Nail	1
1457-019-130	130° Proximal Targeting Device, Aluminium, for Athos Nail	1
1457-019-125	125° Proximal Targeting Device, Aluminium, for Athos Nail	1
7-062-01	Locking Bolt for Aluminium Targeting Device - Athos Nail	3
1457-020	Protection Sleeve Ø10/8.2mm for Athos Nail	2
1457-021	Trocar Ø8mm for Athos Nail	1
1457-022-4.3	Drill Sleeve Ø8/4.3mm for Drill Bit Ø4.3mm - Athos Nail	2
1457-023	Short Aluminium Distal Targeting Device for Length upto 240mm - Athos Nail	1
1457-024	Sliding Hammer, Slotted, for Athos Nail	1
1457-025	Flexible Screwdriver, Hex 3.5mm for Inner Screw - Athos Nail	1
1457-026	Allen Key, Hex 3mm for Athos Nail	1
1457-027	Universal Screwdriver, Hex 7mm for Athos Nail	1
1457-028	Protection Sleeve Ø14/11.2mm for Cephalic Screw - Athos Nail	1
1457-029	Trocar Sleeve Ø11.2/3.2mm for Cephalic Screw - Athos Nail	1
1457-030	Trocar Ø3.2mm for Cephalic Screw Ø10.5mm - Athos Nail	1
1457-031	Measuring Device upto 130mm for Cephalic Screw Ø10.5mm - Athos Nail	1
1457-032	Insertion for Cephalic Screw Ø10.5mm & End Cap - Athos Nail	1

**1457-000S****Athos Nailing Instrument Set, Aluminium Targeting Device**

Code	Set Consisting of	Units
1457-033	Screwdriver with T-Handle, Hex 3.5mm, for Athos Nail	1
1457-034	Universal Screwdriver with T-Handle, Hex 8mm, for Athos Nail	1
1457-035	Proximal Cortex Drill Ø10mm with Quick Coupling End for Cephalic Screw Ø10.5mm - Athos Nail	1
1457-036	Cannulated Reamer, Ø10.5mm for Cephalic Screw Ø10.5mm - Athos Nail	1
7-007-08	Stopper for Cannulated Reamer - Athos Nail	1
1457-037	Measuring Device upto 130mm for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail	1
1457-038	Protection Sleeve 10/8.2mm for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail	1
1457-039	Trocar Sleeve 8.1/2.5mm for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail	1
1457-040	Trocar Ø2.5mm for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail	1
1457-041	Insertor for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail	1
1457-042	Tapered Drill 6.4/4.5mm for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail	1
7-007-09	Stopper for Tapered Drill 6.4/4.5mm - Athos Nail	1
1457-043	Holding Driver, Hex 3.5 for Inner Screw, for Athos Nail	1
1457-044-2.5	Guide Wire Ø2.5mm x Length 380mm for Ø6.4mm Cannulated Anti -Rotation Screw	5
1457-045-3.0	Guide Wire with Threaded Tip Ø3.0mm x Thread Length 12mm x Length 400mm for Cephalic Screw Ø10.5mm	5
1457-046	Extraction Hammering Rod for Athos Nail	1
1457-047	Depth Gauge measuring upto 100mm for Athos Nail	1
1457-048	Depth Gauge measuring upto 115mm for Athos Nail	1
1457-049	Stylet, Ø3mm, for Athos Nail	1
1457-050	Compression Screw for Athos Nail	1
1457-051-4.3	Drill Bit Plain Shank/Jacob Chuck End Ø4.3mm x Length 300mm for Athos Nail	3
1457-052	Stopper for Drill Bit Ø4.3mm - Athos Nail	1
1457-057	Curved Awl- Ø9.5-Ø5.6mm for Athos Nail	1
1457-00T	Instrument Trays for Athos Nailing Instrument Set, Aluminium Targeting Device	3
1457-00B	Container with PPSU Lid for Athos Nailing Instrument Set, Aluminium Targeting Device	1



INSTRUMENTS

PEEK



1457-001

Radiographic Ruler for Short Athos Nail



1457-055

Entry Reamer Protection Sleeve Ø20/17mm for Athos Nail



1457-056

Entry Reamer Guide Wire Sleeve Ø17/3.2mm for Athos Nail



1457-004

Guide Wire With Ball Tip Ø2.5/3.65mm x Length 680mm - Nitinol, for Athos Nail



1457-053

Guide Wire Without Ball Tip Ø2.5mm x Length 680mm - Nitinol, for Athos Nail



1457-005-10

Flexible Shaft with Reamer Ø10mm for Athos Nail



1457-006-11

Flexible Shaft with Reamer \varnothing 11mm for Athos Nail



1457-007-12

Flexible Shaft with Reamer \varnothing 12mm for Athos Nail



1457-008-13

Flexible Shaft with Reamer \varnothing 13mm for Athos Nail



1457-009

T Handle with Quick Coupling for Entry Reamer - Athos Nail



1457-010

Entry Reamer \varnothing 16.5mm with Quick Coupling End for Athos Nail



1457-011

Cannulated Awl for Athos Nail



1457-012

Tissue Protector for Athos Nail



7-007-01

Insertion Handle, Carbon PEEK, for Athos Nail



7-007-10

Nail Holding Bolt for Insertion Handle, Carbon Peek, for Athos Nail



1457-015

Guide Wire Holder for Athos Nail



7-007-12

Insertion Connecting Rod for Athos Nail, Length 163mm



1457-017

Combination Wrench 11mm for Athos Nai



1457-018

Insertion Hammer, Solid, for Athos Nail



7-007-02

130° Proximal Targeting Device, Carbon PEEK, for Athos Nail



7-007-07

Locking Bolt for Carbon PEEK Targeting Device - Athos Nail



7-007-03

125° Proximal Targeting Device, Carbon PEEK, for Athos Nail



7-007-13

Protection Sleeve 10/8.2mm for Athos Nail, Length 175mm



7-007-14

Trocar Ø8mm for Athos Nail, Length 190mm



7-007-15

Drill Sleeve 8/4.3mm for Drill Bit Ø4.3mm - Athos Nail, Length 195mm



7-007-16

Short Carbon PEEK Distal Targeting Device for Length upto 180mm - Athos Nail



1457-024

Sliding Hammer, Slotted, for Athos Nail



1457-025

Flexible Screwdriver, Hex 3.5mm for Inner Screw - Athos Nail



1457-026

Allen Key, Hex 3mm for Athos Nail



1457-027

Universal Screwdriver, Hex 7mm for Athos Nail



7-007-17

Protection Sleeve 14/11.2mm for Cephalic Screw - Athos Nail, Length 235mm



7-007-18

Trocar Sleeve 11.2/3.2mm for Cephalic Screw - Athos Nail, Length 250mm



7-007-19

Trocar Ø3.2mm for Cephalic Screw Ø10.5mm - Athos Nail, Length 265mm



1457-031

Measuring Device upto 130mm for Cephalic Screw Ø10.5mm - Athos Nail



7-007-20

Insertor for Cephalic Screw Ø10.5mm & End Cap - Athos Nail, Length 355mm



7-007-20-A

Threaded bar for locking Cephalic screw & End cap, Length 355mm



7-007-21

Screwdriver with T-Handle, Hex 3.5mm, for Athos Nail, Length 270mm



7-007-22

Proximal Cortex Drill Ø10mm with Quick Coupling End for Cephalic Screw Ø10.5mm - Athos Nail, Length 370mm



7-007-23

Cannulated Reamer, Ø10.5mm for Cephalic Screw Ø10.5mm - Athos Nail, Length 450mm



7-007-08

Stopper for Cannulated Reamer - Athos Nail



1457-037

Measuring Device upto 130mm for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail



7-007-24

Protection Sleeve 10/8.2mm for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail, Length 230mm



7-007-25

Trocar Sleeve $\varnothing 8.1/2.5\text{mm}$ for $\varnothing 6.4\text{mm}$ Cannulated Anti - Rotation Screw, for Athos Nail, Length 250mm



7-007-26

Trocar $\varnothing 2.5\text{mm}$ for $\varnothing 6.4\text{mm}$ Cannulated Anti - Rotation Screw, for Athos Nail, Length 260mm



7-007-27

Insertor for $\varnothing 6.4\text{mm}$ Cannulated Anti -Rotation Screw, for Athos Nail, Length 325mm



7-007-27-A

Threaded bar for Anti-rotation Screw, Length 355mm



7-007-28

Tapered Drill 6.4/4.5mm for $\varnothing 6.4\text{mm}$ Cannulated Anti -Rotation Screw, for Athos Nail, Length 415mm



7-007-09

Stopper for Tapered Drill 6.4/4.5mm - Athos Nail



1457-043

Holding Driver, Hex 3.5 for Inner Screw, for Athos Nail



7-007-29

Guide Wire Ø2.5mm x Length 405mm for Ø6.4mm Cannulated Anti - Rotation Screw



7-007-30

Guide Wire with Threaded Tip Ø3.0mm x Thread Length 12mm x Length 420mm for Cephalic Screw Ø10.5mm



1457-046

Extraction Hammering Rod for Athos Nail



1457-047

Depth Gauge measuring upto 100mm for Athos Nail




7-007-31

Depth Gauge measuring upto 90mm for Athos Nail




1457-049

Styilet, Ø3mm, for Athos Nail



7-007-32

Compression Screw for Athos Nail, Length 95mm



7-007-33

Drill Bit Plain Shank/Jacob Chuck End Ø4.3mm x Length 320mm for Athos Nail



1457-057

Curved Awl- Ø9.5-Ø5.6mm for Athos Nail



1457-052

Stopper for Drill Bit Ø4.3mm - Athos Nail




7-007-11

Radiographic Target Block




7-007-05

Instrument Trays for Athos Nailing Instrument Set (PEEK)



7-007-06

Container with PPSU Lid for Athos Nailing Instrument Set (PEEK)





OPTIONAL INSTRUMENT

1457-004S

Guide Wire With Ball Tip $\varnothing 2.5/3.65\text{mm}$ x Length 680mm - Stainless Steel, for Athos Nail






1457-053S

Guide Wire Without Ball Tip $\varnothing 2.5\text{mm}$ x Length 680mm - Stainless Steel, for Athos Nail



7-007	Athos Nailing Instrument Set (PEEK)
-------	-------------------------------------





ORDER INFORMATION

7-007

Athos Nailing Instrument Set (PEEK)

Code	Set Consisting Of	Units
1457-001	Radiographic Ruler for Short Athos Nail	1
1457-004	Guide Wire With Ball Tip Ø2.5/3.65mm x Length 680mm - Nitinol, for Athos Nail	1
1457-053	Guide Wire Without Ball Tip Ø2.5mm x Length 680mm - Nitinol, for Athos Nail	1
1457-005-10	Flexible Shaft with Reamer Ø10mm for Athos Nail	1
1457-006-11	Flexible Shaft with Reamer Ø11mm for Athos Nail	1
1457-007-12	Flexible Shaft with Reamer Ø12mm for Athos Nail	1
1457-008-13	Flexible Shaft with Reamer Ø13mm for Athos Nail	1
1457-009	T Handle with Quick Coupling for Entry Reamer - Athos Nail	1
1457-010	Entry Reamer Ø16.5mm with Quick Coupling End for Athos Nail	1
1457-011	Cannulated Awl for Athos Nail	1
1457-012	Tissue Protector for Athos Nail	1
7-007-01	Insertion Handle, Carbon PEEK, for Athos Nail	1
7-007-10	Nail Holding Bolt for Insertion Handle, Carbon Peek, for Athos Nail	2
1457-015	Guide Wire Holder for Athos Nail	1
7-007-12	Insertion Connecting Rod for Athos Nail, Length 163mm	1
1457-017	Combination Wrench 11mm for Athos Nail	1
1457-018	Insertion Hammer, Solid, for Athos Nail	1
7-007-02	130° Proximal Targeting Device, Carbon PEEK, for Athos Nail	1
7-007-07	Locking Bolt for Carbon PEEK Targeting Device - Athos Nail	2
7-007-03	125° Proximal Targeting Device, Carbon PEEK, for Athos Nail	1
7-007-13	Protection Sleeve 10/8.2mm for Athos Nail, Length 175mm	2
7-007-14	Trocar Ø8mm for Athos Nail, Length 190mm	1
7-007-15	Drill Sleeve 8/4.3mm for Drill Bit Ø4.3mm - Athos Nail, Length 195mm	2
7-007-16	Short Carbon PEEK Distal Targeting Device for Length upto 180mm - Athos Nail	1
1457-024	Sliding Hammer, Slotted, for Athos Nail	1
1457-025	Flexible Screwdriver, Hex 3.5mm for Inner Screw - Athos Nail	1
1457-026	Allen Key, Hex 3mm for Athos Nail	1
1457-027	Universal Screwdriver, Hex 7mm for Athos Nail	1
7-007-17	Protection Sleeve 14/11.2mm for Cephalic Screw - Athos Nail, Length 235mm	1
7-007-18	Trocar Sleeve 11.2/3.2mm for Cephalic Screw - Athos Nail, Length 250mm	1
7-007-19	Trocar Ø3.2mm for Cephalic Screw Ø10.5mm - Athos Nail, Length 265mm	1
1457-031	Measuring Device upto 130mm for Cephalic Screw Ø10.5mm - Athos Nail	1
1457-032-A	Threaded bar for locking Cephalic screw & End cap	1
7-007-20	Insertor for Cephalic Screw Ø6.4mm Cannulated Anti-Rotation Screw , for Athos Nail, Length 325mm	1
7-007-20-A	Threaded bar for locking Cephalic screw & End cap, Length 355mm	1
7-007-21	Screwdriver with T-Handle, Hex 3.5mm, for Athos Nail, Length 270mm	1
1457-034	Universal Screwdriver with T-Handle, Hex 8mm, for Athos Nail	1
7-007-22	Proximal Cortex Drill Ø10mm with Quick Coupling End for Cephalic Screw Ø10.5mm - Athos Nail, Length 370mm	1
7-007-23	Cannulated Reamer, Ø10.5mm for Cephalic Screw Ø10.5mm - Athos Nail, Length 450mm	1
7-007-08	Stopper for Cannulated Reamer - Athos Nail	1
1457-037	Measuring Device upto 130mm for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail	1



7-007

Athos Nailing Instrument Set (PEEK)

Code	Set Consisting Of	Units
7-007-24	Protection Sleeve 10/8.2mm for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail, Length 230mm	1
7-007-25	Trocar Sleeve 8.1/2.5mm for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail, Length 250mm	1
7-007-26	Trocar Ø2.5mm for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail, Length 260mm	1
7-007-27	Insertor for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail, Length 325mm	1
7-007-27-A	Threaded bar for Anti-Rotation Screw, Length 355mm	1
7-007-28	Tapered Drill 6.4/4.5mm for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail, Length 415mm	1
7-007-09	Stopper for Tapered Drill 6.4/4.5mm - Athos Nail	1
1457-043	Holding Driver, Hex 3.5 for Inner Screw, for Athos Nail	1
7-007-29	Guide Wire Ø2.5mm x Length 405mm for Ø6.4mm Cannulated Anti -Rotation Screw	5
7-007-30	Guide Wire with Threaded Tip Ø3.0mm x Thread Length 12mm x Length 420mm for Cephalic Screw Ø10.5mm	5
1457-046	Extraction Hammering Rod for Athos Nail	1
1457-047	Depth Gauge measuring upto 100mm for Athos Nail	1
7-007-31	Depth Gauge measuring upto 90mm for Athos Nail	1
1457-049	Stylet, Ø3mm, for Athos Nail	1
7-007-32	Compression Screw for Athos Nail, Length 95mm	1
7-007-33	Drill Bit Plain Shank/Jacob Chuck End Ø4.3mm x Length 320mm for Athos Nail	3
1457-055	Entry Reamer Protection Sleeve 20/17mm for Athos Nail	1
1457-056	Entry Reamer Guide Wire Sleeve 17/3.2mm for Athos Nail	1
1457-057	Curved Awl- Ø9.5-Ø5.6mm for Athos Nail	1
1457-052	Stopper for Drill Bit Ø4.3mm - Athos Nail	1
7-007-11	Radiographic Target Block	1
7-007-35	Flexible Shaft with Reamer Ø14mm for Athos Nail	1
7-007-05	Instrument Trays for Athos Nailing Instrument Set (PEEK)	3
7-007-06	Container with PPSU Lid for Athos Nailing Instrument Set (PEEK)	1
	Optional Instruments	
1457-004S	Guide Wire With Ball Tip Ø2.5/3.65mm x Length 680mm - Stainless Steel, for Athos Nail	
1457-053S	Guide Wire Without Ball Tip Ø2.5mm x Length 680mm - Stainless Steel, for Athos Nail	



INSTRUMENTS

DISTAL TARGETING DEVICE



1457-001

Radiographic Ruler for Short Athos Nail



1457-055

Entry Reamer Protection Sleeve 20/17mm for Athos Nail



1457-056

Entry Reamer Guide Wire Sleeve 17/3.2mm for Athos Nail



1457-005-10

Flexible Shaft with Reamer $\phi 10$ mm for Athos Nail



1457-006-11

Flexible Shaft with Reamer $\phi 11$ mm for Athos Nail



1457-007-12

Flexible Shaft with Reamer $\phi 12$ mm for Athos Nail



1457-008-13

Flexible Shaft with Reamer Ø13mm for Athos Nail



1457-009

T Handle with Quick Coupling for Entry Reamer - Athos Nail



1457-010

Entry Reamer Ø16.5mm with Quick Coupling End for Athos Nail



1457-011

Cannulated Awl for Athos Nail



1457-012

Tissue Protector for Athos Nail



1457-013

Insertion Handle, Aluminium, for Athos Nail



1457-014

Ø4.4mm Cannulated Nail Holding Bolt for Aluminium Targeting Device for Athos Nail



1457-015

Guide Wire Holder for Athos Nail



1457-016

Insertion Connecting Rod for Athos Nail- M10.



1457-054

Block for Insertion Connecting Rod - Athos Nail



1457-017

Combination Wrench 11mm for Athos Nail



1457-018

Insertion Hammer, Solid, for Athos Nail





1457-019-130

130° Proximal Targeting Device, Aluminium, for Athos Nail



1457-019-125

125° Proximal Targeting Device, Aluminium, for Athos Nail



7-062-01

Locking Bolt for Aluminium Targeting Device - Athos Nail



1457-020

Protection Sleeve 10/8.2mm for Athos Nail



1457-021

Trocar Ø8mm for Athos Nail



1457-022-4.3

Drill Sleeve 8/4.3mm for Drill Bit Ø4.3mm - Athos Nail



1457-023

Short Aluminium Distal Targeting Device for Length upto 240mm - Athos Nail



1457-024

Sliding Hammer, Slotted, for Athos Nail



1457-025

Flexible Screwdriver, Hex 3.5mm for Inner Screw - Athos Nail



1457-026

Allen Key, Hex 3mm for Athos Nail



1457-027

Universal Screwdriver, Hex 7mm for Athos Nail



1457-028

Protection Sleeve 14/11.2mm for Cephalic Screw - Athos Nail



1457-029

Trocar Sleeve 11.2/3.2mm for Cephalic Screw - Athos Nail



1457-030

Trocar Ø3.2mm for Cephalic Screw Ø10.5mm - Athos Nail



1457-031

Measuring Device upto 130mm for Cephalic Screw Ø10.5mm - Athos Nail



1457-032

Inserter for Cephalic Screw Ø10.5mm & End Cap - Athos Nail



1457-033

Screwdriver with T-Handle, Hex 3.5mm, for Athos Nail



1457-034

Universal Screwdriver with T-Handle, Hex 8mm, for Athos Nail



1457-035

Proximal Cortex Drill Ø10mm with Quick Coupling End for Cephalic Screw
Ø10.5mm for Athos Nail



1457-036

Cannulated Reamer, Ø10.5mm for Cephalic Screw Ø10.5mm - Athos Nail



7-007-08

Stopper for Cannulated Reamer - Athos Nail



1457-037

Measuring Device upto 130mm for Ø6.4mm Cannulated Anti -Rotation Screw,
for Athos Nail



1457-038

Protection Sleeve 10/8.2mm for Ø6.4mm Cannulated Anti -Rotation Screw,
for Athos Nail



1457-039

Trocar Sleeve 8.1/2.5mm for Ø6.4mm Cannulated Anti -Rotation Screw,
for Athos Nail



1457-040

Trocar Ø2.5mm for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail



1457-041

Insertor for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail



1457-042

Tapered Drill 6.4/4.5mm for Ø6.4mm Cannulated Anti - Rotation Screw, for Athos Nail



7-007-09

Stopper for Tapered Drill 6.4/4.5mm - Athos Nail



1457-043

Holding Driver, Hex 3.5 for Inner Screw, for Athos Nail



1457-044-2.5

Guide Wire Ø2.5mm x Length 380mm for Ø6.4mm Cannulated Anti -Rotation Screw



1457-045-3.0

Guide Wire with Threaded Tip $\varnothing 3.0\text{mm}$ x Thread Length 12mm x Length 400mm for Cephalic Screw $\varnothing 10.5\text{mm}$



1457-046

Extraction Hammering Rod for Athos Nail



1457-047

Depth Gauge measuring upto 100mm for Athos Nail



1457-048

Depth Gauge measuring upto 115mm for Athos Nail



1457-049

Stylet, $\varnothing 3\text{mm}$, for Athos Nail



1457-050

Compression Screw for Athos Nail



1457-051-4.3

Drill Bit Plain Shank/Jacob Chuck End Ø4.3mm x Length 300mm for Athos Nail



1457-052

Stopper for Drill Bit Ø4.3mm - Athos Nail



7-008-01

Ø5.2mm T-Handle Drill Bit with Flat End for Athos Nail



7-008-02

Drill Bit Plain Shank/Jacob Chuck End Ø5.2mm x Length 250mm for Athos Nail



7-008-03

Radiographic Ruler For Long Athos Nail



7-008-04

Distal Targeting Device for Athos Nail



7-008-05

Allen Key, Hex 5mm for Athos Nail



7-008-06

Stabilizing Device (Ø9, Ø10, Ø11) for Athos Nail



7-008-34

Stabilizing Device (Ø12, Ø13) for Athos Nail



7-008-07

Flexible Reamer Shaft (AO Connection) for Athos Nail



7-008-08

Outer Drill Sleeve 10/8.1mm for Distal Targeting Device, for Athos Nail




7-008-09

Trocar Ø8mm for Long Athos Nail



7-008-10

Drill Sleeve 8.1/5.2mm for Positioning Rod, for Athos Nail



7-008-11

Positioning Rod 8.1/5.2mm for Athos Nail



7-008-12

Reduction Rod for Long Athos Nail



7-008-13

Proximal Aiming Bar for Long Athos Nail



7-008-14

Distal Aiming Bar for Long Athos Nail



Reamer Head

Code	Dia	Code	Dia
7-008-15	Ø8.5mm	7-008-21	Ø11.5mm
7-008-16	Ø9mm	7-008-22	Ø12mm
7-008-17	Ø9.5mm	7-008-23	Ø12.5mm
7-008-18	Ø10mm	7-008-24	Ø13mm
7-008-19	Ø10.5mm	7-008-25	Ø13.5mm
7-008-20	Ø11mm	7-008-26	Ø14mm



7-008-27

Locking Bolt for Distal Aiming Bar - Athos Nail



7-008-28

Locking Bolt for Proximal Aiming Bar - Athos Nail



7-008-29

Distal Bolt, Hex 5mm for Distal Targeting Device - Athos Nail



7-008-32

Guide Wire without Ball Tip Ø2.5mm x Length 1000mm - Nitinol, for Athos Nail



7-008-33

Guide Wire with Ball Tip Ø2.5/3.65mm x Length 1000mm - Nitinol, for Athos Nail



1457-057

Curved Awl- Ø9.5-Ø5.6mm for Athos Nail




7-008-30

Instrument Trays for Athos Nailing instrument set, Distal Targeting Device



7-008-31


Container with PPSU Lid for Athos Nailing instrument set, Distal Targeting Device



Optional Instruments

7-008-32S

Guide Wire without Ball Tip $\varnothing 2.5\text{mm}$ x Length 1000mm - Nitinol, for Athos Nail



7-008-33S

Guide Wire with Ball Tip $\varnothing 2.5/3.65\text{mm}$ x Length 1000mm - Nitinol, for Athos Nail



7-008

Athos Nailing Instrument Set, Distal Targeting Device





ORDER INFORMATION

7-008
Athos Nailing Instrument Set, Distal Targeting Device

Code	Set Consisting Of	Units
1457-001	Radiographic Ruler for Short Athos Nail	1
1457-005-10	Flexible Shaft with Reamer Ø10mm for Athos Nail	1
1457-006-11	Flexible Shaft with Reamer Ø11mm for Athos Nail	1
1457-007-12	Flexible Shaft with Reamer Ø12mm for Athos Nail	1
1457-008-13	Flexible Shaft with Reamer Ø13mm for Athos Nail	1
1457-009	T Handle with Quick Coupling for Entry Reamer - Athos Nail	1
1457-010	Entry Reamer Ø16.5mm with Quick Coupling End for Athos Nail	1
1457-011	Cannulated Awl for Athos Nail	1
1457-012	Tissue Protector for Athos Nail	1
1457-013	Insertion Handle, Aluminium, for Athos Nail	1
1457-014	Nail Holding Bolt for Aluminium Targeting Device for Athos Nail	1
1457-015	Guide Wire Holder for Athos Nail	1
1457-016	Insertion Connecting Rod for Athos Nail-M10	2
1457-054	Block for Insertion Connecting Rod - Athos Nail	1
1457-017	Combination Wrench 11mm for Athos Nail	1
1457-018	Insertion Hammer, Solid, for Athos Nail	1
1457-019-130	130° Proximal Targeting Device, Aluminium, for Athos Nail	1
1457-019-125	125° Proximal Targeting Device, Aluminium, for Athos Nail	1
7-062-01	Locking Bolt for Aluminium Targeting Device - Athos Nail	1
1457-020	Protection Sleeve 10/8.2mm for Athos Nail	1
1457-021	Trocar Ø8mm for Athos Nail	3
1457-022-4.3	Drill Sleeve 8/4.3mm for Drill Bit Ø4.3mm - Athos Nail	2
1457-023	Short Aluminium Distal Targeting Device for Length upto 240mm - Athos Nail	1
1457-024	Sliding Hammer, Slotted, for Athos Nail	2
1457-025	Flexible Screwdriver, Hex 3.5mm for Inner Screw - Athos Nail	1
1457-026	Allen Key, Hex 3mm for Athos Nail	1
1457-027	Universal Screwdriver, Hex 7mm for Athos Nail	1
1457-028	Protection Sleeve 14/11.2mm for Cephalic Screw - Athos Nail	1
1457-029	Trocar Sleeve 11.2/3.2mm for Cephalic Screw - Athos Nail	1
1457-030	Trocar Ø3.2mm for Cephalic Screw Ø10.5mm - Athos Nail	1
1457-031	Measuring Device upto 130mm for Cephalic Screw Ø10.5mm - Athos Nail	1
1457-032	Inserter for Cephalic Screw Ø10.5mm & End Cap - Athos Nail	1
1457-032-A	Threaded bar for locking Cephalic screw & End cap	1
1457-033	Screwdriver with T-Handle, Hex 3.5mm, for Athos Nail	1
1457-034	Universal Screwdriver with T-Handle, Hex 8mm, for Athos Nail	1
1457-035	Proximal Cortex Drill Ø10mm with Quick Coupling End for Cephalic Screw Ø10.5mm - Athos Nail	1



ORDER INFORMATION

7-008	Athos Nailing Instrument Set, Distal Targeting Device
--------------	--

Code	Set Consisting Of	Units
1457-036	Cannulated Reamer, Ø10.5mm for Cephalic Screw Ø10.5mm - Athos Nail	1
7-007-08	Stopper for Cannulated Reamer - Athos Nail	1
1457-037	Measuring Device upto 130mm for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail	1
1457-038	Protection Sleeve 10/8.2mm for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail	1
1457-039	Trocar Sleeve 8.1/2.5mm for Ø6.4mm Cannulated Anti-Rotation Screw, for Athos Nail	1
1457-040	Trocar Ø2.5mm for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail	1
1457-041	Insertor for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail	1
1457-041-A	Threaded bar for locking Cannulated Antirotation screw, for Athos Nail	1
1457-042	Tapered Drill 6.4/4.5mm for Ø6.4mm Cannulated Anti -Rotation Screw, for Athos Nail	1
7-007-09	Stopper for Tapered Drill 6.4/4.5mm - Athos Nail	1
1457-043	Holding Driver, Hex 3.5 for Inner Screw, for Athos Nail	1
1457-044-2.5	Guide Wire Ø2.5mm x Length 380mm for Ø6.4mm Cannulated Anti -Rotation Screw	5
1457-045-3.0	Guide Wire with Threaded Tip Ø3.0mm x Thread Length 12mm x Length 400mm for Cephalic Screw Ø10.5mm	5
1457-046	Extraction Hammering Rod for Athos Nail	1
1457-047	Depth Gauge measuring upto 100mm for Athos Nail	1
1457-048	Depth Gauge measuring upto 115mm for Athos Nail	1
1457-049	Stylet, Ø3mm, for Athos Nail	1
1457-050	Compression Screw for Athos Nail	1
1457-051-4.3	Drill Bit Plain Shank/Jacob Chuck End Ø4.3mm x Length 300mm for Athos Nail	3
1457-055	Entry Reamer Protection Sleeve 20/17mm for Athos Nail	1
1457-056	Entry Reamer Guide Wire Sleeve 17/3.2mm for Athos Nail	1
1457-057	Curved Awl- Ø9.5-Ø5.6mm for Athos Nail	1
1457-052	Stopper for Drill Bit Ø4.3mm - Athos Nail	1
7-008-01	Ø5.2mm T-Handle Drill Bit with Flat End for Athos Nail	1
7-008-02	Drill Bit Plain Shank/Jacob Chuck End Ø5.2mm x Length 250mm for Athos Nail	1
7-008-03	Radiographic Ruler For Long Athos Nail	1
7-008-04	Distal Targeting Device for Athos Nail	1
7-008-05	Allen Key, Hex 5mm for Athos Nail	1
7-008-06	Stabilizing Device (Ø9, Ø10, Ø11) for Athos Nail	1
7-008-34	Stabilizing Device (Ø12, Ø13) for Athos Nail	1
7-008-07	Flexible Reamer Shaft (AO Connection) for Athos Nail	2
7-008-08	Outer Drill Sleeve 10/8.1mm for Distal Targeting Device, for Athos Nail	1
7-008-09	Trocar Ø8mm for Long Athos Nail	1
7-008-10	Drill Sleeve 8.1/5.2mm for Positioning Rod, for Athos Nail	1
7-008-11	Positioning Rod 8.1/5.2mm for Athos Nail	1
7-008-12	Reduction Rod for Long Athos Nail	1
7-008-13	Proximal Aiming Bar for Long Athos Nail	1
7-008-14	Distal Aiming Bar for Long Athos Nail	1
7-008-15	Reamer Head, Ø8.5mm for Athos Nail	1



ORDER INFORMATION

7-008	Athos Nailing Instrument Set, Distal Targeting Device
--------------	--

Code	Set Consisting Of	Units
7-008-16	Reamer Head, Ø9.0mm for Athos Nail	1
7-008-17	Reamer Head, Ø9.5mm for Athos Nail	1
7-008-18	Reamer Head, Ø10.0mm for Athos Nail	1
7-008-19	Reamer Head, Ø10.5mm for Athos Nail	1
7-008-20	Reamer Head, Ø11.0mm for Athos Nail	1
7-008-21	Reamer Head, Ø11.5mm for Athos Nail	1
7-008-22	Reamer Head, Ø12.0mm for Athos Nail	1
7-008-23	Reamer Head, Ø12.5mm for Athos Nail	1
7-008-24	Reamer Head, Ø13.0mm for Athos Nail	1
7-008-25	Reamer Head, Ø13.5mm for Athos Nail	1
7-008-26	Reamer Head, Ø14.0mm for Athos Nail	2
7-008-27	Locking Bolt for Distal Aiming Bar - Athos Nail	1
7-008-28	Locking Bolt for Proximal Aiming Bar - Athos Nail	1
7-008-29	Distal Bolt, Hex 5mm for Distal Targeting Device - Athos Nail	1
7-008-32	Guide Wire without Ball Tip Ø2.5mm x Length 1000mm - Nitinol, for Athos Nail	1
7-008-33	Guide Wire with Ball Tip Ø2.5/3.65mm x Length 1000mm - Nitinol, for Athos Nail	1
7-008-30	Instrument Trays for Athos Nailing instrument set, Distal Targeting Device	4
7-008-31	Container with PPSU Lid for Athos Nailing instrument set, Distal Targeting Device	1
	Optional Instruments	
7-008-32S	Guide Wire without Ball Tip Ø2.5mm x Length 1000mm - Stainless Steel, for Athos Nail	
7-008-33S	Guide Wire with Ball Tip Ø2.5/3.65mm x Length 1000mm - Stainless Steel, for Athos Nail	



ATHOS NAILING SYSTEM

The content of this catalogue is intended for healthcare professionals only. All product information, illustrations, descriptions, and specifications provided here are for general informational purposes and are subject to change without notice. Auxein Medical Pvt. Ltd. reserves the right to improve or discontinue products at any time without prior notification.

This document does not constitute medical advice. Surgeons and medical practitioners must exercise their own professional judgment when using any medical devices or instruments, including those provided by Auxein.

INTENDED USE

The implants and instruments in this catalogue are intended for use by qualified orthopedic surgeons. Improper selection, misuse, or lack of experience may result in complications.

REGULATORY COMPLIANCE

US FDA 510(k)

CE Mark (Class IIb)

MDSAP (Medical Device Single Audit Program)

ISO 13485:2016 Quality Management System

CDSCO – Indian Medical Device Rules, 2017

EU MDR- European Union Medical Device Regulation, 2017

INTELLECTUAL PROPERTY

All trademarks, trade names, brand names, product names, and logos mentioned in this catalogue are the property of Auxein Medical Pvt. Ltd. Unauthorized use or reproduction is strictly prohibited.

V1-2025



SCAN TO KNOW MORE

INDIA
AUXEIN MEDICAL PVT.LTD.

PLOT NO. 168-169-170, PHASE-4,
KUNDLI INDUSTRIAL AREA,
HSIIDC, SECTOR-57,
SONEPAT-131028, HARYANA
TEL: +91 99106 43638
FAX: +91 86077 70197
EMAIL: INFO@AUXEIN.COM

MÉXICO
AUXEIN MÉXICO, S.A DE CV

TEPIC 39 INT 801, COLONIA
ROMA SUR,
ALCADÍA CUAUHTÉMOC, CDMX
MÉXICO, C.P 06760
TEL: +521 55 7261 0318
EMAIL: INFO@AUXEIN.MX

USA
AUXEIN INC.

1500 NW 89TH COURT,
SUITE 107-108
DORAL, FLORIDA 33172
TEL. +1 305 395 6062
E FAX: +1 305 395 62 62
EMAIL: USOFFICE@AUXEIN.COM

DUBAI UAE
AUXEIN MEDICAL

TRADING FZCO OUD METHA
OFFICES BUILDING B BLOCK,
1ST FLOOR, OFFICE NO.26,
NEAR HEALTHCARE CITY,
PO BOX-25572
TEL. +971 438 71333
EMAIL: INFO@AUXEIN.COM

WWW.AUXEIN.COM