

# Auxein

## BROKEN SCREW REMOVAL




NEXT-LEVEL CARE



## ABOUT US

Auxein Medical is an integrated, research based, orthopaedic implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavour to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in India and serves customers in over 75 countries worldwide.



# IMPACTING HEALTHCARE GLOBALLY

Our Reach & Influence

**500+**

Employees Worldwide

**75+**

Operations in Countries

**150+**

Distributors Worldwide

**20**<sup>M+</sup>

Patients Cured

**200+**

FDA Approved Products

**250+**

Machines Worldwide

**3000+**

Innovative Orthopedic  
Products & Solutions


Auxin is dedicated to advancing medical technology through innovative orthopedic solutions, including **Trauma, Arthroscopy, and Reconstructive Procedures**, with a commitment to improving patient outcomes worldwide. Our approach integrates advanced research, precision engineering, and high-quality manufacturing to develop reliable implants and surgical instruments that address the evolving needs of modern orthopedic practice.

# OUR ACHIEVEMENTS



OVER THE YEARS






At Auxein, it is our goal to provide the best in class quality and value to every individual. Consistency is a virtue we pursue with a passion. From the procurement of raw materials to the manufacturing of our range of 1200+ orthopaedic implant products, we employ world-class systems and technology at par with the highest standards of globally recognized evaluation systems. The Broken Screw is part of our advanced orthopedic portfolio designed to support versatile fixation techniques in trauma and reconstructive procedures. Developed with input from global experts, the system provides reliable solutions for fixation of for intramedullary fixation of femoral fractures, helping achieve stable fixation and accurate anatomical alignment.



# INDEX



Instrument .....	11-20
Instrument Set .....	21
Order Information .....	22-23




# **BROKEN SCREW REMOVAL**

# BROKEN SCREW REMOVAL



7-043-01

T-Shaped Handle Quick Coupling



7-043-02

Straight Handle Quick Coupling



7-043-03

Trephine Reamer Ø4mm



7-043-04

Trephine Reamer Ø5mm



7-043-05

Trephine Reamer Ø6.5mm



7-043-06

Trephine Reamer Ø8mm



7-043-07

Extraction Bolt Ø6.5mm



7-043-08

Extraction Bolt Ø5.0mm



7-043-09

Extraction Bolt Ø4.0mm



7-043-10

Trephine Reamer Ø2.7mm



7-043-11

Trephine Reamer Ø2.0mm



7-043-12

Extraction Screw Hexagonal 4.0mm Conical



7-043-13

Extraction Screw Hexagonal 3.5mm Conical



7-043-14

Extraction Screw Hexagonal 2.5mm Conical



7-043-15

Extraction Screw Hexagonal 2.7mm Conical



7-043-16

Extraction Screw Hexagonal 2.0mm Conical



7-043-17

Slotted Screwdriver Small



7-043-18

Slotted Screwdriver Large



7-043-19

Crosshead Screwdriver



7-043-20

Square Screwdriver, Hex 1.2mm



7-043-21

Square Screwdriver, Hex 1.5mm



7-043-22

Star Screwdriver Shaft, T30



7-043-23

Star Screwdriver Shaft, T25



7-043-24

Star Screwdriver Shaft, T20



7-043-25

Star Screwdriver Shaft, T18



7-043-26

Star Screwdriver Shaft, T15



7-043-27

Star Screwdriver Shaft, T10



7-043-28

Star Screwdriver Shaft, T9



7-043-29

Star Screwdriver Shaft, T8



7-043-30

Star Screwdriver Shaft, T7



7-043-31

Star Screwdriver Shaft, T6.2



7-043-32

Star Screwdriver Shaft, T6



7-043-33

Triangle Screwdriver 1.7mm



7-043-34

Hexagonal Screwdriver, Hex 5.5mm



7-043-35

Hexagonal Screwdriver, Hex 5.0mm



7-043-36

Hexagonal Screwdriver, Hex 4.5mm



7-043-37

Hexagonal Screwdriver, Hex 4.0mm



7-043-38

Hexagonal Screwdriver, Hex 3.5mm



7-043-39

Hexagonal Screwdriver, Hex 3.0mm



7-043-40

Hexagonal Screwdriver, Hex 2.5mm



7-043-41

Hexagonal Screwdriver, Hex 2.0mm



7-043-42

Hexagonal Screwdriver, Hex 1.5mm



7-043-43

Universal Screw Extractor



7-043-44

High Speed Drill Bit Ø4 x Length 70mm



7-043-45

High Speed Drill Bit Ø5 x Length 70mm



7-043-46

High Speed Drill Bit Ø6 x Length 70mm



7-043-47

Sharp Hook for Screw Removal



7-043-48

Plier for Screw Removal



7-043-49

Spinal Monoaxial Pedicle Screwdriver



7-043-50

Hollow Gouge



7-043-51

Tray for Broken Screw Removal Instruments



7-043-52

Container for Broken Screw Removal Instruments



**7-043 : Broken Screw Removal Instrument Set**



# ORDER INFORMATION

<b>7-043</b>	<b>Instrument Set for Broken Screw Removal</b>
--------------	--

Code	Set Consisting Of	Unit
7-043-01	T-Shaped Handle Quick Coupling	1
7-043-02	Straight Handle Quick Coupling	1
7-043-03	Trephine Reamer Ø4mm	1
7-043-04	Trephine Reamer Ø5mm	1
7-043-05	Trephine Reamer Ø6.5mm	1
7-043-06	Trephine Reamer Ø8mm	1
7-043-07	Extraction Bolt Ø6.5mm	1
7-043-08	Extraction Bolt Ø5.0mm	1
7-043-09	Extraction Bolt Ø4.0mm	1
7-043-10	Trephine Reamer Ø2.0mm	1
7-043-11	Extraction Bolt Ø2.0mm	1
7-043-12	Extraction Screw Hexagonal 4.0mm Conical	1
7-043-13	Extraction Screw Hexagonal 3.5mm Conical	1
7-043-14	Extraction Screw Hexagonal 2.5mm Conical	1
7-043-15	Extraction Screw Hexagonal 2.7mm Conical	1
7-043-16	Extraction Screw Hexagonal 2.0mm Conical	1
7-043-17	Slotted Screwdriver Small	1
7-043-18	Slotted Screwdriver Large	1
7-043-19	Crosshead Screwdriver	1
7-043-20	Square Screwdriver, Hex 1.2mm	1
7-043-21	Square Screwdriver, Hex 1.5mm	1
7-043-22	Star Screwdriver Shaft, T30	1
7-043-23	Star Screwdriver Shaft, T25	1
7-043-24	Star Screwdriver Shaft, T20	1
7-043-25	Star Screwdriver Shaft, T18	1
7-043-26	Star Screwdriver Shaft, T15	1
7-043-27	Star Screwdriver Shaft, T10	1
7-043-28	Star Screwdriver Shaft, T9	1
7-043-29	Star Screwdriver Shaft, T8	1
7-043-30	Star Screwdriver Shaft, T7	1
7-043-31	Star Screwdriver Shaft, T6.2	1
7-043-32	Star Screwdriver Shaft, T6	1
7-043-33	Triangle Screwdriver 1.7mm	1
7-043-34	Hexagonal Screwdriver, Hex 5.5mm	1
7-043-35	Hexagonal Screwdriver, Hex 5.0mm	1

Code	Set Consisting Of	Unit
7-043-36	Hexagonal Screwdriver, Hex 4.5mm	1
7-043-37	Hexagonal Screwdriver, Hex 4.0mm	1
7-043-38	Hexagonal Screwdriver, Hex 3.5mm	1
7-043-39	Hexagonal Screwdriver, Hex 3.0mm	1
7-043-40	Hexagonal Screwdriver, Hex 2.5mm	1
7-043-41	Hexagonal Screwdriver, Hex 2.0mm	1
7-043-42	Hexagonal Screwdriver, Hex 1.5mm	1
7-043-43	Universal Screw Extractor	2
7-043-44	High Speed Drill Bit Ø4 x Length 70mm	1
7-043-45	High Speed Drill Bit Ø5 x Length 70mm	1
7-043-46	High Speed Drill Bit Ø6 x Length 70mm	1
7-043-47	Sharp Hook for Screw Removal	1
7-043-48	Plier for Screw Removal	1
7-043-49	Spinal Monoaxial Pedicle Screwdriver	1
7-043-50	Hollow Gouge	1
7-043-51	Tray for Broken Screw Removal Instruments	2
7-043-52	Container for Broken Screw Removal Instruments	1

## **BROKEN SCREW REMOVAL**

The content of this catalogue is intended for healthcare professionals only. All product information, illustrations, descriptions, and specifications provided here are for general informational purposes and are subject to change without notice. Auxein Medical Pvt. Ltd. reserves the right to improve or discontinue products at any time without prior notification.

This document does not constitute medical advice. Surgeons and medical practitioners must exercise their own professional judgment when using any medical devices or instruments, including those provided by Auxein.

## **INTENDED USE**

The implants and instruments in this catalogue are intended for use by qualified orthopedic trauma surgeons trained in fracture fixation and intramedullary nailing procedures. Improper selection of the implant, incorrect surgical technique, misuse, or insufficient surgical experience may result in implant failure or other complications.

## **REGULATORY COMPLIANCE CERTIFIED**

CE Approved (EUMDR-2017/745)

USFDA 510(k)

ISO 13485:2016

CDSCO

MDSAP

## **INTELLECTUAL PROPERTY**

All trademarks, trade names, brand names, product names, and logos mentioned in this catalogue are the property of Auxein Medical Pvt. Ltd. Unauthorized use or reproduction is strictly prohibited.

V1-2026





SCAN TO KNOW MORE

**INDIA**  
**AUXEIN MEDICAL PVT.LTD.**

PLOT NO. 168-169-170, PHASE-4,  
KUNDLI INDUSTRIAL AREA,  
HSIIDC, SECTOR-57,  
SONEPAT-131028, HARYANA  
TEL: +91 99106 43638  
FAX: +91 86077 70197  
EMAIL: INFO@AUXEIN.COM

**MÉXICO**  
**AUXEIN MÉXICO, S.A DE CV**

TEPIC 39 INT 801, COLONIA  
ROMA SUR,  
ALCADÍA CUAUHTÉMOC, CDMX  
MÉXICO, C.P 06760  
TEL: +521 55 7261 0318  
EMAIL: INFO@AUXEIN.MX

**USA**  
**AUXEIN INC.**

1500 NW 89TH COURT,  
SUITE 107-108  
DORAL, FLORIDA 33172  
TEL. +1 305 395 6062  
E FAX: +1 305 395 62 62  
EMAIL: USOFFICE@AUXEIN.COM

**DUBAI UAE**  
**AUXEIN MEDICAL**

TRADING FZCO OUD METHA  
OFFICES BUILDING B BLOCK,  
1ST-FLOOR, OFFICE NO.26,  
NEAR HEALTHCARE CITY,  
PO BOX-25572  
TEL. +971 438 71333  
EMAIL: INFO@AUXEIN.COM

[WWW.AUXEIN.COM](http://WWW.AUXEIN.COM)