

Auxein


CABLE PLATE SYSTEM



NEXT-LEVEL CARE

ABOUT US

Auxein Medical is an integrated, research based, orthopaedic implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavour to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in India and serves customers in over 75 countries worldwide.



IMPACTING HEALTHCARE GLOBALLY

Our Reach & Influence

500+

Employees Worldwide

75+

Operations in Countries

150+

Distributors Worldwide

20^{M+}

Patients Cured

200+

FDA Approved Products

250+

Machines Worldwide

3000+

Innovative Orthopedic
Products & Solutions


Auxein is dedicated to advancing medical technology through innovative orthopedic solutions, including **Trauma, Arthroscopy, and Reconstructive Procedures**, with a commitment to improving patient outcomes worldwide. Our approach integrates advanced research, precision engineering, and high-quality manufacturing to develop reliable implants and surgical instruments that address the evolving needs of modern orthopedic practice.

OUR ACHIEVEMENTS



OVER THE YEARS






At Auxein, it is our goal to provide the best in class quality and value to every individual. Consistency is a virtue we pursue with a passion. From the procurement of raw materials to the manufacturing of our range of 1200+ orthopaedic implant products, we employ world-class systems and technology at par with the highest standards of globally recognized evaluation systems. The Wise-Lock Cable Plate System is part of the advanced orthopedic trauma portfolio, designed for fixation of complex fractures and periprosthetic fractures. Developed with clinical input, the system combines cerclage cable fixation with plate-based stabilization to provide enhanced construct stability, controlled compression, and maintenance of anatomical alignment.



INDEX





Introduction	09
Cable plates	10-14
Screws	15-16
Implant Set	17-19
Instruments	20-25




CABLE PLATE SYSTEM

Cable Plate System

Introduction

The Cable Plate System is part of the advanced orthopedic trauma portfolio, designed for fixation of complex and periprosthetic fractures where conventional screw fixation may be limited. The system integrates cerclage cable technology with plate-based stabilization to enhance construct stability, provide controlled circumferential fixation, and maintain anatomical alignment.

Comprising cable-compatible plates, cerclage cables, connectors, bone needles, and locking screws, the system enables versatile application in challenging clinical scenarios, particularly in the presence of implants, poor bone quality, or multifragmentary fractures.




Indication:

- Periprosthetic fractures (e.g., around hip or knee prostheses).
- Fractures where screw fixation is limited due to existing implants.
- Comminuted long bone fractures requiring supplementary cerclage fixation.
- Diaphyseal and metaphyseal fractures of the femur.
- Fractures in osteoporotic bone requiring enhanced fixation.
- Revision fracture fixation procedures.
- Stabilization of bone fragments in complex reconstructive procedures.
- Nonunions and malunions requiring augmentation with cerclage fixation.



5.0mm Wise-Lock Cable Prosthesis and Revision Femur Plate



STERILE

Holes	Left Direction	Right Direction
2mm	10-018-02LTI	10-018-02RTI
4mm	10-018-04LTI	10-018-04RTI
6mm	10-018-06LTI	10-018-06RTI
8mm	10-018-08LTI	10-018-08RTI
10mm	10-018-010LTI	10-018-010RTI
12mm	10-018-012LTI	10-018-012RTI


Holes	Left Direction	Right Direction
2mm	10-018-02LTI-S	10-018-02RTI-S
4mm	10-018-04LTI-S	10-018-04RTI-S
6mm	10-018-06LTI-S	10-018-06RTI-S
8mm	10-018-08LTI-S	10-018-08RTI-S
10mm	10-018-010LTI-S	10-018-010RTI-S
12mm	10-018-012LTI-S	10-018-012RTI-S

5.0mm Wise-Lock Cable Straight Plate, Titanium



- Left
- Right

5.0mm Wise-Lock Cable Straight Plate



Compression & Locking Hole



Holes for Cable



Holes	Left Direction
6mm	10-043-06TI
8mm	10-043-08TI
10mm	10-043-010TI

STERILE

Holes	Left Direction
6mm	10-043-06TI-S
8mm	10-043-08TI-S
10mm	10-043-010TI-S

Cable



STERILE

Dia	Length	Right Direction
1.0mm	500mm	10-019-500TI
1.7mm	500mm	10-020-500TI

Dia	Length	Right Direction
1.0mm	500mm	10-019-500TI-S
1.7mm	500mm	10-020-500TI-S

Cable Connector Curved



STERILE

Dia	Right Direction
1.0mm	10-062-1.0TI
1.7mm	10-062-1.7TI

Dia	Right Direction
1.0mm	10-062-1.0TI-S
1.7mm	10-062-1.7TI-S

Bone Needle



STERILE

Dia	Length	Codes
1.5mm	150mm	10-022-150SS
1.5mm	250mm	10-022-250SS
1.8mm	400mm	10-023-400SS
2.0mm	150mm	10-024-150SS
2.0mm	400mm	10-024-400SS
2.0mm	450mm	10-024-450SS

Dia	Length	Codes
1.5mm	150mm	10-022-150SS-S
1.5mm	250mm	10-022-250SS-S
1.8mm	400mm	10-023-400SS-S
2.0mm	150mm	10-024-150SS-S
2.0mm	400mm	10-024-400SS-S
2.0mm	450mm	10-024-450SS-S

Bone Needle with Hole



Dia	Length	Right Direction
1.6mm	100mm	10-025-100SS
1.6mm	140mm	10-025-140SS
2.0mm	100mm	10-026-100SS
2.0mm	120mm	10-026-120SS
2.0mm	140mm	10-026-140SS
2.5mm	120mm	10-027-120SS
2.5mm	140mm	10-027-140SS



Dia	Length	Right Direction
1.6mm	100mm	10-025-100SS-S
1.6mm	140mm	10-025-140SS-S
2.0mm	100mm	10-026-100SS-S
2.0mm	120mm	10-026-120SS-S
2.0mm	140mm	10-026-140SS-S
2.5mm	120mm	10-027-120SS-S
2.5mm	140mm	10-027-140SS-S



SCREWS



5.0mm Wise-Lock Screw, Self-Tapping, (Hex Head)



Length	Codes
12mm	TI-119.012
14mm	TI-119.014
16mm	TI-119.016
18mm	TI-119.018
20mm	TI-119.020
22mm	TI-119.022
24mm	TI-119.024
26mm	TI-119.026
28mm	TI-119.028
30mm	TI-119.030
32mm	TI-119.032
34mm	TI-119.034
36mm	TI-119.036
38mm	TI-119.038
40mm	TI-119.040
42mm	TI-119.042
44mm	TI-119.044
46mm	TI-119.046
48mm	TI-119.048
50mm	TI-119.050
52mm	TI-119.052
54mm	TI-119.054
56mm	TI-119.056
58mm	TI-119.058
60mm	TI-119.060
65mm	TI-119.065
70mm	TI-119.070
75mm	TI-119.075
80mm	TI-119.080
85mm	TI-119.085
90mm	TI-119.090

STERILE

Length	Codes
12mm	TI-119.012-S
14mm	TI-119.014-S
16mm	TI-119.016-S
18mm	TI-119.018-S
20mm	TI-119.020-S
22mm	TI-119.022-S
24mm	TI-119.024-S
26mm	TI-119.026-S
28mm	TI-119.028-S
30mm	TI-119.030-S
32mm	TI-119.032-S
34mm	TI-119.034-S
36mm	TI-119.036-S
38mm	TI-119.038-S
40mm	TI-119.040-S
42mm	TI-119.042-S
44mm	TI-119.044-S
46mm	TI-119.046-S
48mm	TI-119.048-S
50mm	TI-119.050-S
52mm	TI-119.052-S
54mm	TI-119.054-S
56mm	TI-119.056-S
58mm	TI-119.058-S
60mm	TI-119.060-S
65mm	TI-119.065-S
70mm	TI-119.070-S
75mm	TI-119.075-S
80mm	TI-119.080-S
85mm	TI-119.085-S
90mm	TI-119.090-S







IMPLANT BOX



10-049 : Implant Box for 5.0mm Wise-Lock Cable Plate Implant Set



ORDER INFORMATION

10-049	5.0mm Wise-Lock Cable Plate Implant Set
--------	---

Code	Set Consisting Of	Unit
10-018-02LTI	5.0mm Wise-Lock Cable Prosthesis and Revision Femur Plate, Left, 2 Holes, Titanium	1
10-018-04LTI	5.0mm Wise-Lock Cable Prosthesis and Revision Femur Plate, Left, 4 Holes, Titanium	1
10-018-06LTI	5.0mm Wise-Lock Cable Prosthesis and Revision Femur Plate, Left, 6 Holes, Titanium	1
10-018-08LTI	5.0mm Wise-Lock Cable Prosthesis and Revision Femur Plate, Left, 8 Holes, Titanium	1
10-018-10LTI	5.0mm Wise-Lock Cable Prosthesis and Revision Femur Plate, Left, 10 Holes, Titanium	1
10-018-12LTI	5.0mm Wise-Lock Cable Prosthesis and Revision Femur Plate, Left, 12 Holes, Titanium	1
10-018-02RTI	5.0mm Wise-Lock Cable Prosthesis and Revision Femur Plate, Right, 2 Holes, Titanium	1
10-018-04RTI	5.0mm Wise-Lock Cable Prosthesis and Revision Femur Plate, Right, 4 Holes, Titanium	1
10-018-06RTI	5.0mm Wise-Lock Cable Prosthesis and Revision Femur Plate, Right, 6 Holes, Titanium	1
10-018-08RTI	5.0mm Wise-Lock Cable Prosthesis and Revision Femur Plate, Right, 8 Holes, Titanium	1
10-018-10RTI	5.0mm Wise-Lock Cable Prosthesis and Revision Femur Plate, Right, 10 Holes, Titanium	1
10-018-12RTI	5.0mm Wise-Lock Cable Prosthesis and Revision Femur Plate, Right, 12 Holes, Titanium	1
10-043-06TI	5.0mm Wise-Lock Cable Straight Plate, 6 Holes, Titanium	1
10-043-08TI	5.0mm Wise-Lock Cable Straight Plate, 8 Holes, Titanium	1
10-043-10TI	5.0mm Wise-Lock Cable Straight Plate, 10 Holes, Titanium	1
10-062-1.0TI	Cable Connector, Curved, Ø1.0mm, Titanium	1
10-062-1.7TI	Cable Connector, Curved, Ø1.7mm, Titanium	1
10-022-150SS	Bone Needle, Ø1.5mm x Length 150mm, Stainless Steel	2
10-022-250SS	Bone Needle, Ø1.5mm x Length 250mm, Stainless Steel	2
10-023-400SS	Bone Needle, Ø1.8mm x Length 400mm, Stainless Steel	2
10-024-150SS	Bone Needle, Ø2.0mm x Length 150mm, Stainless Steel	2
10-024-400SS	Bone Needle, Ø2.0mm x Length 400mm, Stainless Steel	2
10-024-450SS	Bone Needle, Ø2.0mm x Length 450mm, Stainless Steel	2
10-025-100SS	Bone Needle with Hole, Ø1.6mm x Length 100mm, Stainless Steel	2
10-025-140SS	Bone Needle with Hole, Ø1.6mm x Length 140mm, Stainless Steel	2
10-026-100SS	Bone Needle with Hole, Ø2.0mm x Length 100mm, Stainless Steel	2
10-026-120SS	Bone Needle with Hole, Ø2.0mm x Length 120mm, Stainless Steel	2
10-026-140SS	Bone Needle with Hole, Ø2.0mm x Length 140mm, Stainless Steel	2
10-027-120SS	Bone Needle with Hole, Ø2.5mm x Length 120mm, Stainless Steel	2
10-027-140SS	Bone Needle with Hole, Ø2.5mm x Length 140mm, Stainless Steel	2



INSTRUMENTS



7-035-01

Flat Connector Crimper for Cable Plate




7-035-02

Spicule Cutter for Cable Plate



7-035-03

Cable Cutter for Cable Plate



7-035-04

Front Pressure Cable Lock, ID: Ø1.0mm, for Cable Plate



7-035-05

Back Pressure Cable Lock, ID: Ø1.0mm, for Cable Plate



7-035-06

Front Pressure Cable Lock, ID: Ø1.7mm, for Cable Plate



7-035-07

Back Pressure Cable Lock, ID: Ø1.7mm, for Cable Plate



7-035-08

Cable Tensioner for Cable Plate



7-035-09

Cable Passer for Cable Plate



7-035-10

Spicule Bender for Cable Plate



7-035-11

Curved Needle Guide for Cable Plate



7-035-12

Straight Needle Guide for Cable Plate



7-035-13

Impactor for Cable Plate



7-035-16

Wire Tensioner for Cable Plate



7-035-14

Trays for Cable Plate Instrument Set





7-035-15

Container for Cable Plate Instrument Set



7-035	Cable Plate Instrument Set
-------	----------------------------



ORDER INFORMATION

7-035	Cable plate instrument Set
-------	----------------------------

Code	Set Consisting Of	Unit
7-035-01	Flat Connector Crimper for Cable Plate	1
7-035-02	Spicule Cutter for Cable Plate	1
7-035-03	Cable Cutter for Cable Plate	1
7-035-04	Front Pressure Cable Lock, ID: Ø1.0mm, for Cable Plate	1
7-035-05	Back Pressure Cable Lock, ID: Ø1.0mm, for Cable Plate	1
7-035-06	Front Pressure Cable Lock, ID: Ø1.7mm, for Cable Plate	3
7-035-07	Back Pressure Cable Lock, ID: Ø1.7mm, for Cable Plate	1
7-035-08	Cable Tensioner for Cable Plate	1
7-035-09	Cable Passer for Cable Plate	1
7-035-10	Spicule Bender for Cable Plate	1
7-035-11	Curved Needle Guide for Cable Plate	1
7-035-12	Straight Needle Guide for Cable Plate	1
7-035-13	Impactor for Cable Plate	1
7-035-16	Wire Tensioner for Cable Plate	1
7-035-14	Trays for Cable Plate Instrument Set	2
7-035-15	Container for Cable Plate Instrument Set	1

CABLE PLATE SYSTEM

The content of this catalogue is intended for healthcare professionals only. All product information, illustrations, descriptions, and specifications provided here are for general informational purposes and are subject to change without notice. Auxein Medical Pvt. Ltd. reserves the right to improve or discontinue products at any time without prior notification.

This document does not constitute medical advice. Surgeons and medical practitioners must exercise their own professional judgment when using any medical devices or instruments, including those provided by Auxein.

INTENDED USE

The implants and instruments in this catalogue are intended for use by qualified orthopedic surgeons trained in trauma and reconstructive procedures. Improper selection, misuse, or lack of experience may result in complications.

REGULATORY COMPLIANCE

US FDA 510(k)

CE Mark (Class IIb)

MDSAP (Medical Device Single Audit Program)

ISO 13485:2016 Quality Management System

CDSCO - Indian Medical Device Rules, 2017

EU MDR- European Union Medical Device Regulation, 2017

INTELLECTUAL PROPERTY

All trademarks, trade names, brand names, product names, and logos mentioned in this catalogue are the property of Auxein Medical Pvt. Ltd. Unauthorized use or reproduction is strictly prohibited.

V1-2026



SCAN TO KNOW MORE

INDIA
AUXEIN MEDICAL PVT.LTD.

PLOT NO. 168-169-170, PHASE-4,
KUNDLI INDUSTRIAL AREA,
HSIIDC, SECTOR-57,
SONEPAT-131028, HARYANA
TEL: +91 99106 43638
FAX: +91 86077 70197
EMAIL: INFO@AUXEIN.COM

MÉXICO
AUXEIN MÉXICO, S.A DE CV

TEPIC 39 INT 801, COLONIA
ROMA SUR,
ALCADÍA CUAUHTÉMOC, CDMX
MÉXICO, C.P 06760
TEL: +521 55 7261 0318
EMAIL: INFO@AUXEIN.MX

USA
AUXEIN INC.

1500 NW 89TH COURT,
SUITE 107-108
DORAL, FLORIDA 33172
TEL. +1 305 395 6062
E FAX: +1 305 395 62 62
EMAIL: USOFFICE@AUXEIN.COM

DUBAI UAE
AUXEIN MEDICAL

TRADING FZCO OUD METHA
OFFICES BUILDING B BLOCK,
1ST FLOOR, OFFICE NO.26,
NEAR HEALTHCARE CITY,
PO BOX-25572
TEL. +971 438 71333
EMAIL: INFO@AUXEIN.COM

WWW.AUXEIN.COM