



AV-WISELOCKTM


3.5MM AV-Wiselock Foot System

AV-WISELOCK™

3.5mm AV-Wiselock Foot System



Plate Features:

- Low-Profile plates with better anatomical Fit
- AV-Wiselock-variable angle of $\pm 15^\circ$ in all Directions in each screw hole
- Rounded edges and a smooth surface for soft tissue protection
- K-wire Holes and slot for temporary plate fixation and compression



Plates Available in:
Titanium as per ISO 5832-3 and Stainless Steel as per ISO 5832-1

3.5mm AV-Wiselock Medial Column Fusion Plate





STERILE

Holes	Direction	Titanium	Stainless Steel
8	Left	11-093-08TIL	11-093-08SSL
8	Right	11-093-08TIR	11-093-08SSR

Holes	Direction	Titanium	Stainless Steel
8	Left	11-093-08TIL-S	11-093-08SSL-S
8	Right	11-093-08TIR-S	11-093-08SSR-S

3.5mm AV-Wiselock Medial Column Fusion Plate, Talus Extension



STERILE

Holes	Direction	Titanium	Stainless Steel
9	Left	11-094-09TIL	11-094-09SSL
9	Right	11-094-09TIR	11-094-09SSR

Holes	Direction	Titanium	Stainless Steel
9	Left	11-094-09TIL-S	11-094-09SSL-S
9	Right	11-094-09TIR-S	11-094-09SSR-S

3.5mm AV-Wiselock Medial Column Fusion Plantar Plate




Plate Features

- Available in multiple variants, including talus extension and plantar fusion designs to address diverse medial column pathologies
- Versatile plate configurations enabling fixation and compression across single or multiple medial column joints
- Multidirectional locking capability with strategic hole placement allowing intramedullary screw orientation for enhanced medial column stability
- Optimized screw distribution, permitting placement of two to three screws per bone to achieve high construct stability



STERILE

Holes	Direction	Titanium	Stainless Steel
7	Left	11-095-07TIL	11-095-07SSL
7	Right	11-095-07TIR	11-095-07SSR

Holes	Direction	Titanium	Stainless Steel
7	Left	11-095-07TIL-S	11-095-07SSL-S
7	Right	11-095-07TIR-S	11-095-07SSR-S

3.5mm AV-Wiselock Calcaneus M Plate



Postoperative X-ray showing calcaneal fracture fixation with AV-Wiselock™ Calcaneus M Plate and locking screws, achieving stable alignment.

Plate Features:

- Engineered frame design promotes uniform stress distribution throughout the plate.
- Extensive calcaneal coverage allows screws to engage dense bone structures with 3.5mm AV-Wiselock and cortical screw for enhanced stability.
- The subtalar joint reduction can be stabilized with up to five screws targeting the sustentaculum tali.



Postoperative X-ray



Small

Holes	Direction	Titanium	Stainless Steel
12	Left	11-060-01TIL	11-060-01SSL
12	Right	11-060-01TIR	11-060-01SSR

STERILE

Holes	Direction	Titanium	Stainless Steel
7	Left	11-060-01TIL-S	11-060-01SSL-S
7	Right	11-060-01TIR-S	11-060-01SSR-S

Medium

Holes	Direction	Titanium	Stainless Steel
13	Left	11-060-02TIL	11-060-02SSL
13	Right	11-060-02TIR	11-060-02SSR

Holes	Direction	Titanium	Stainless Steel
7	Left	11-060-02TIL-S	11-060-02SSL-S
7	Right	11-060-02TIR-S	11-060-02SSR-S

Large

Holes	Direction	Titanium	Stainless Steel
12	Left	11-060-03TIL	11-060-03SSL
12	Right	11-060-03TIR	11-060-03SSR

Holes	Direction	Titanium	Stainless Steel
7	Left	11-060-03TIL-S	11-060-03SSL-S
7	Right	11-060-03TIR-S	11-060-03SSR-S

3.5mm AV-Wiselock MIS Calcaneal Plate

- Left
- Right




Plate Features:

Designed specifically for the MIS sinus tarsi approach, allowing for smaller incisions. The Plates are low-profile & anatomically contoured to reduce irritation to the peroneal tendons and sural nerve. It utilizes 3.5 mm AV-Wiselock & Cortical screws for rigid fixation. The design enables surgeons to address both posterior facet depression and restore Böhler's angle. Smaller incisions with the radiolucent Jig minimize damage to the local blood supply, reducing the risk of wound healing issues and infection.




STERILE

Holes	Direction	Titanium	Stainless Steel
6	Left	11-100-06LTI	11-100-06LSS
6	Right	11-100-06RTI	11-100-06RSS
7	Left	11-100-07LTI	11-100-07LSS
7	Right	11-100-07RTI	11-100-07RSS

Holes	Direction	Titanium	Stainless Steel
6	Left	11-100-06LTI-S	11-100-06LSS-S
6	Right	11-100-06RTI-S	11-100-06RSS-S
7	Left	11-100-07LTI-S	11-100-07LSS-S
7	Right	11-100-07RTI-S	11-100-07RSS-S

3.5mm AV-Wiselock Calcaneal SLIDE Plate



Key Features:

- **Variable-Angle Locking Technology ($\pm 15^\circ$):**
Allows multi-Directional screw placement for optimal bone purchase.
- **Central Hooded Sustentacular Screw Tunnel:**
Directs the cortical screw toward the sustentaculum tali for strongest calcaneal purchase.
- **Dual Anti-Rotation Arms:**
Maintain posterior calcaneal fragment alignment after slide osteotomy.
- **Ultra-Low Profile Geometry:**
Smooth, contoured design minimizes irritation to soft tissues around the heel.
- Smaller incisions with the radiolucent Jig minimize damage to the local blood supply, reducing the risk of wound healing issues and infection

Indications:

- Medial Calcaneal Slide Osteotomy
- Lateral Calcaneal Slide Osteotomy
- Fixation of Osteotomies with High Rotational Instability



Holes	Titanium
3	11-116-01TI



Holes	Titanium
3	11-116-01TI-S

3.5mm AV-Wiselock Plantar Lapidus Plate




Plate Features:

- Positioned on the plantar aspect of the 1st TMT joint, this plate resists compression forces (tension band) rather than shear forces.
- The design, including 3.5 mm AV-Wiselock, locking and Cortical screw an integrated cancellous screw (4.0 mm) enhances stability and bone healing.
- Precontoured for left/right, it protects the tibialis anterior tendon and limits dorsal soft tissue dissection.

TMT-I (Lapidus) fusion stabilized with a 3.5mm AV-Wiselock™ Plantar Lapidus Plate and screws, achieving stable alignment.




Postoperative X-ray



	Direction	Titanium	Stainless Steel
SHORT	Left	11-102-01LTI	11-102-01LSS
	Right	11-102-01RTI	11-102-01RSS
STANDARD	Left	11-101-01LTI	11-102-01LSS
	Right	11-101-01RTI	11-102-01RSS

STERILE

	Direction	Titanium	Stainless Steel
SHORT	Left	11-102-01LTI-S	11-102-01LSS-S
	Right	11-102-01RTI-S	11-102-01RSS-S
STANDARD	Left	11-101-01LTI-S	11-101-01LSS-S
	Right	11-101-01RTI-S	11-101-01RSS-S

The background of the page features a series of thin, light gray, wavy lines that flow from the top left towards the right and then curve downwards towards the bottom right, creating a sense of movement and depth. These lines are more densely packed in some areas, creating a subtle gradient effect.

SCREWS

3.5mm AV-Wiselock Screw, Self-tapping




Dia	Length	Titanium	Stainless Steel
Ø3.5	10mm	11-024-010TI	11-024-010SS
Ø3.5	12mm	11-024-012TI	11-024-012SS
Ø3.5	14mm	11-024-014TI	11-024-014SS
Ø3.5	16mm	11-024-016TI	11-024-016SS
Ø3.5	18mm	11-024-018TI	11-024-018SS
Ø3.5	20mm	11-024-020TI	11-024-020SS
Ø3.5	22mm	11-024-022TI	11-024-022SS
Ø3.5	24mm	11-024-024TI	11-024-024SS
Ø3.5	26mm	11-024-026TI	11-024-026SS
Ø3.5	28mm	11-024-028TI	11-024-028SS
Ø3.5	30mm	11-024-030TI	11-024-030SS
Ø3.5	32mm	11-024-032TI	11-024-032SS
Ø3.5	34mm	11-024-034TI	11-024-034SS
Ø3.5	36mm	11-024-036TI	11-024-036SS
Ø3.5	38mm	11-024-038TI	11-024-038SS
Ø3.5	40mm	11-024-040TI	11-024-040SS
Ø3.5	42mm	11-024-042TI	11-024-042SS
Ø3.5	44mm	11-024-044TI	11-024-044SS
Ø3.5	46mm	11-024-046TI	11-024-046SS
Ø3.5	48mm	11-024-048TI	11-024-048SS
Ø3.5	50mm	11-024-050TI	11-024-050SS
Ø3.5	52mm	11-024-052TI	11-024-052SS
Ø3.5	54mm	11-024-054TI	11-024-054SS
Ø3.5	56mm	11-024-056TI	11-024-056SS
Ø3.5	58mm	11-024-058TI	11-024-058SS
Ø3.5	60mm	11-024-060TI	11-024-060SS
Ø3.5	65mm	11-024-065TI	11-024-065SS
Ø3.5	70mm	11-024-070TI	11-024-070SS
Ø3.5	75mm	11-024-075TI	11-024-075SS
Ø3.5	80mm	11-024-080TI	11-024-080SS
Ø3.5	85mm	11-024-085TI	11-024-085SS
Ø3.5	90mm	11-024-090TI	11-024-090SS
Ø3.5	95mm	11-024-095TI	11-024-095SS

STERILE



Dia	Length	Titanium	Stainless Steel
Ø3.5	10mm	11-024-010TI-S	11-024-010SS-S
Ø3.5	12mm	11-024-012TI-S	11-024-012SS-S
Ø3.5	14mm	11-024-014TI-S	11-024-014SS-S
Ø3.5	16mm	11-024-016TI-S	11-024-016SS-S
Ø3.5	18mm	11-024-018TI-S	11-024-018SS-S
Ø3.5	20mm	11-024-020TI-S	11-024-020SS-S
Ø3.5	22mm	11-024-022TI-S	11-024-022SS-S
Ø3.5	24mm	11-024-024TI-S	11-024-024SS-S
Ø3.5	26mm	11-024-026TI-S	11-024-026SS-S
Ø3.5	28mm	11-024-028TI-S	11-024-028SS-S
Ø3.5	30mm	11-024-030TI-S	11-024-030SS-S
Ø3.5	32mm	11-024-032TI-S	11-024-032SS-S
Ø3.5	34mm	11-024-034TI-S	11-024-034SS-S
Ø3.5	36mm	11-024-036TI-S	11-024-036SS-S
Ø3.5	38mm	11-024-038TI-S	11-024-038SS-S
Ø3.5	40mm	11-024-040TI-S	11-024-040SS-S
Ø3.5	42mm	11-024-042TI-S	11-024-042SS-S
Ø3.5	44mm	11-024-044TI-S	11-024-044SS-S
Ø3.5	46mm	11-024-046TI-S	11-024-046SS-S
Ø3.5	48mm	11-024-048TI-S	11-024-048SS-S
Ø3.5	50mm	11-024-050TI-S	11-024-050SS-S
Ø3.5	52mm	11-024-052TI-S	11-024-052SS-S
Ø3.5	54mm	11-024-054TI-S	11-024-054SS-S
Ø3.5	56mm	11-024-056TI-S	11-024-056SS-S
Ø3.5	58mm	11-024-058TI-S	11-024-058SS-S
Ø3.5	60mm	11-024-060TI-S	11-024-060SS-S
Ø3.5	65mm	11-024-065TI-S	11-024-065SS-S
Ø3.5	70mm	11-024-070TI-S	11-024-070SS-S
Ø3.5	75mm	11-024-075TI-S	11-024-075SS-S
Ø3.5	80mm	11-024-080TI-S	11-024-080SS-S
Ø3.5	85mm	11-024-085TI-S	11-024-085SS-S
Ø3.5	90mm	11-024-090TI-S	11-024-090SS-S
Ø3.5	95mm	11-024-095TI-S	11-024-095SS-S





OPTIONAL SCREW

3.5mm WISE-LOCK Screw, Self-Tapping, (Star Head)




Dia	Length	Titanium	Stainless Steel
Ø3.5	10mm	TI-1473-010	1473-010
Ø3.5	12mm	TI-1473-012	1473-012
Ø3.5	14mm	TI-1473-014	1473-014
Ø3.5	16mm	TI-1473-016	1473-016
Ø3.5	18mm	TI-1473-018	1473-018
Ø3.5	20mm	TI-1473-020	1473-020
Ø3.5	22mm	TI-1473-022	1473-022
Ø3.5	24mm	TI-1473-024	1473-024
Ø3.5	26mm	TI-1473-026	1473-026
Ø3.5	28mm	TI-1473-028	1473-028
Ø3.5	30mm	TI-1473-030	1473-030
Ø3.5	32mm	TI-1473-032	1473-032
Ø3.5	34mm	TI-1473-034	1473-034
Ø3.5	36mm	TI-1473-036	1473-036
Ø3.5	38mm	TI-1473-038	1473-038
Ø3.5	40mm	TI-1473-040	1473-040
Ø3.5	42mm	TI-1473-042	1473-042
Ø3.5	44mm	TI-1473-044	1473-044
Ø3.5	46mm	TI-1473-046	1473-046
Ø3.5	48mm	TI-1473-048	1473-048
Ø3.5	50mm	TI-1473-050	1473-050
Ø3.5	52mm	TI-1473-052	1473-052
Ø3.5	54mm	TI-1473-054	1473-054
Ø3.5	56mm	TI-1473-056	1473-056
Ø3.5	58mm	TI-1473-058	1473-058
Ø3.5	60mm	TI-1473-060	1473-060
Ø3.5	65mm	TI-1473-065	1473-065
Ø3.5	70mm	TI-1473-070	1473-070
Ø3.5	75mm	TI-1473-075	1473-075
Ø3.5	80mm	TI-1473-080	1473-080
Ø3.5	85mm	TI-1473-085	1473-085
Ø3.5	90mm	TI-1473-090	1473-090
Ø3.5	95mm	TI-1473-095	1473-095

STERILE



Dia	Length	Titanium	Stainless Steel
Ø3.5	10mm	TI-1473-010-S	1473-010-S
Ø3.5	12mm	TI-1473-012-S	1473-012-S
Ø3.5	14mm	TI-1473-014-S	1473-014-S
Ø3.5	16mm	TI-1473-016-S	1473-016-S
Ø3.5	18mm	TI-1473-018-S	1473-018-S
Ø3.5	20mm	TI-1473-020-S	1473-020-S
Ø3.5	22mm	TI-1473-022-S	1473-022-S
Ø3.5	24mm	TI-1473-024-S	1473-024-S
Ø3.5	26mm	TI-1473-026-S	1473-026-S
Ø3.5	28mm	TI-1473-028-S	1473-028-S
Ø3.5	30mm	TI-1473-030-S	1473-030-S
Ø3.5	32mm	TI-1473-032-S	1473-032-S
Ø3.5	34mm	TI-1473-034-S	1473-034-S
Ø3.5	36mm	TI-1473-036-S	1473-036-S
Ø3.5	38mm	TI-1473-038-S	1473-038-S
Ø3.5	40mm	TI-1473-040-S	1473-040-S
Ø3.5	42mm	TI-1473-042-S	1473-042-S
Ø3.5	44mm	TI-1473-044-S	1473-044-S
Ø3.5	46mm	TI-1473-046-S	1473-046-S
Ø3.5	48mm	TI-1473-048-S	1473-048-S
Ø3.5	50mm	TI-1473-050-S	1473-050-S
Ø3.5	52mm	TI-1473-052-S	1473-052-S
Ø3.5	54mm	TI-1473-054-S	1473-054-S
Ø3.5	56mm	TI-1473-056-S	1473-056-S
Ø3.5	58mm	TI-1473-058-S	1473-058-S
Ø3.5	60mm	TI-1473-060-S	1473-060-S
Ø3.5	65mm	TI-1473-065-S	1473-065-S
Ø3.5	70mm	TI-1473-070-S	1473-070-S
Ø3.5	75mm	TI-1473-075-S	1473-075-S
Ø3.5	80mm	TI-1473-080-S	1473-080-S
Ø3.5	85mm	TI-1473-085-S	1473-085-S
Ø3.5	90mm	TI-1473-090-S	1473-090-S
Ø3.5	95mm	TI-1473-095-S	1473-095-S



3.5mm Cortical Screw, Self-Tapping, (Star Head)





Dia	Length	Titanium	Stainless Steel
Ø3.5	10mm	11-002-010TI	11-002-010SS
Ø3.5	12mm	11-002-012TI	11-002-012SS
Ø3.5	14mm	11-002-014TI	11-002-014SS
Ø3.5	16mm	11-002-016TI	11-002-016SS
Ø3.5	18mm	11-002-018TI	11-002-018SS
Ø3.5	20mm	11-002-020TI	11-002-020SS
Ø3.5	22mm	11-002-022TI	11-002-022SS
Ø3.5	24mm	11-002-024TI	11-002-024SS
Ø3.5	26mm	11-002-026TI	11-002-026SS
Ø3.5	28mm	11-002-028TI	11-002-028SS
Ø3.5	30mm	11-002-030TI	11-002-030SS
Ø3.5	32mm	11-002-032TI	11-002-032SS
Ø3.5	34mm	11-002-034TI	11-002-034SS
Ø3.5	36mm	11-002-036TI	11-002-036SS
Ø3.5	38mm	11-002-038TI	11-002-038SS
Ø3.5	40mm	11-002-040TI	11-002-040SS
Ø3.5	42mm	11-002-042TI	11-002-042SS
Ø3.5	44mm	11-002-044TI	11-002-044SS
Ø3.5	46mm	11-002-046TI	11-002-046SS
Ø3.5	48mm	11-002-048TI	11-002-048SS
Ø3.5	50mm	11-002-050TI	11-002-050SS
Ø3.5	55mm	11-002-055TI	11-002-055SS
Ø3.5	60mm	11-002-060TI	11-002-060SS
Ø3.5	65mm	11-002-065TI	11-002-065SS
Ø3.5	70mm	11-002-070TI	11-002-070SS
Ø3.5	75mm	11-002-075TI	11-002-075SS
Ø3.5	80mm	11-002-080TI	11-002-080SS
Ø3.5	85mm	11-002-085TI	11-002-085SS
Ø3.5	90mm	11-002-090TI	11-002-090SS
Ø3.5	95mm	11-002-095TI	11-002-095SS

STERILE



Dia	Length	Titanium	Stainless Steel
Ø3.5	10mm	11-002-010TI-S	11-002-010SS-S
Ø3.5	12mm	11-002-012TI-S	11-002-012SS-S
Ø3.5	14mm	11-002-014TI-S	11-002-014SS-S
Ø3.5	16mm	11-002-016TI-S	11-002-016SS-S
Ø3.5	18mm	11-002-018TI-S	11-002-018SS-S
Ø3.5	20mm	11-002-020TI-S	11-002-020SS-S
Ø3.5	22mm	11-002-022TI-S	11-002-022SS-S
Ø3.5	24mm	11-002-024TI-S	11-002-024SS-S
Ø3.5	26mm	11-002-026TI-S	11-002-026SS-S
Ø3.5	28mm	11-002-028TI-S	11-002-028SS-S
Ø3.5	30mm	11-002-030TI-S	11-002-030SS-S
Ø3.5	32mm	11-002-032TI-S	11-002-032SS-S
Ø3.5	34mm	11-002-034TI-S	11-002-034SS-S
Ø3.5	36mm	11-002-036TI-S	11-002-036SS-S
Ø3.5	38mm	11-002-038TI-S	11-002-038SS-S
Ø3.5	40mm	11-002-040TI-S	11-002-040SS-S
Ø3.5	42mm	11-002-042TI-S	11-002-042SS-S
Ø3.5	44mm	11-002-044TI-S	11-002-044SS-S
Ø3.5	46mm	11-002-046TI-S	11-002-046SS-S
Ø3.5	48mm	11-002-048TI-S	11-002-048SS-S
Ø3.5	50mm	11-002-050TI-S	11-002-050SS-S
Ø3.5	55mm	11-002-055TI-S	11-002-055SS-S
Ø3.5	60mm	11-002-060TI-S	11-002-060SS-S
Ø3.5	65mm	11-002-065TI-S	11-002-065SS-S
Ø3.5	70mm	11-002-070TI-S	11-002-070SS-S
Ø3.5	75mm	11-002-075TI-S	11-002-075SS-S
Ø3.5	80mm	11-002-080TI-S	11-002-080SS-S
Ø3.5	85mm	11-002-085TI-S	11-002-085SS-S
Ø3.5	90mm	11-002-090TI-S	11-002-090SS-S
Ø3.5	95mm	11-002-095TI-S	11-002-095SS-S



The background of the page features a light gray, wavy, topographic-like pattern that flows across the top and bottom edges, framing the central text.

INSTRUMENTS


7-012-01

Conical Extraction Screw, Ø3.5mm, For 2.7/3.5mm
Av-wiselock System



7-012-02

Drill Bit With Quick Coupling With Stopper, Ø2.8mmx Length 165mm,
For 2.7/3.5m Av-wiselock System



7-012-03

Star Screwdriver Shaft, T8, For 2.7/3.5mm Av-wiselock System




7-012-04

Depth Gauge Measuring Upto 60mm For 2.7/3.5mm
Av-wiselock System



7-012-05

Drill Bit With Quick Coupling, Ø2.0mm X Length 140mm ,
For 2.7/3.5mm Av-wiselock System




7-012-06

Bone Tap Quick Coupling For Cortical Screws, Ø3.5mm,
For 2.7/3.5mm Av-wiselock System



7-012-07

Screwdriver Shaft, T15, For 2.7/3.5mm Av-wiselock System



7-012-08

Threaded Drill Sleeve For Ø2.8mm Drill Bit - 2.7/3.5mm Av-wiselock System



7-012-09

Self-centering Double Drill Guide, Ø2.5/3.5mm, For 2.7/3.5mm Av-wiselock System



7-012-10

Torque Limiting Attachment, 1.2nm, For 2.7/3.5mm Av-wiselock System



7-012-11

Torque Limiting Attachment, 1.5nm, For 2.7/3.5mm Av-wiselock System



7-012-12

Drill Bit With Quick Coupling End, Ø2.5mm X Length 180mm - 2.7/3.5mm Av-wiselock System



7-012-13

Drill Bit With Quick Coupling End, Ø3.5mm X Length 110mm
- 2.7/3.5mm Av-wiselock System

**7-012-14**

2.0mm Variable Angle Guide For 2.7/3.5mm Av-wiselock System

**7-012-15**

Self-centering Double Drill Guide, Ø2.0/2.7mm, For 2.7/3.5mm
Av-wiselock System

**7-012-16**

2.8mm Variable Angle Double Drill Guide

**7-012-17**

Handle With Quick Coupling For 2.7/3.5mm Av-wiselock System

**7-012-18**

Handle For Torque Limiting Attachment, For 2.7/3.5mm
Av-wiselock System



7-012-19

2.5nm Torque Limiting Handle For 2.7/3.5mm Av-wiselock System



7-012-20

Depth Gauge Measuring Upto 100mm, For 2.7/3.5mm Av-wiselock System




7-012-21

Percutaneous Depth Gauge Adaptor For 3.5mm Screws - 2.7/3.5mm Av-wiselock System



7-012-22

Kirschner Wire With One Sided Trocar Tip Ø1.6mm X Length 150mm



7-012-23

Periarticular Reduction Forcep, Small, For 2.7/3.5mm Av-wiselock System



7-012-24

Periarticular Reduction Forcep, Medium, For 2.7/3.5mm Av-wiselock System



7-012-25

Large Compression Forcep With Speed Lock For 2.7/3.5mm
Av-wiselock System




7-012-26

Large Distraction Forcep With Speed Lock For 2.7/3.5mm
Av-wiselock System



7-012-27

3.5mm Post For Av-wiselock Locking Hole



7-012-28

2.8mm Compression Wire, 10mm Thread For 2.7/3.5mm Av-wiselock System



7-012-29

2.8mm Compression Wire, 15mm Thread For 2.7/3.5mm Av-wiselock System



7-012-30

2.8mm Compression Wire, 20mm Thread For 2.7/3.5mm Av-wiselock System



7-012-31

2.8mm Compression Wire, 25mm Thread For 2.7/3.5mm Av-wiselock System



7-012-32

2.8mm Compression Wire, 30mm Thread For 2.7/3.5mm Av-wiselock System



7-012-33

2.8mm Compression Wire, 35mm Thread For 2.7/3.5mm Av-wiselock System



7-012-34

2.8mm Compression Wire, 40mm Thread For 2.7/3.5mm Av-wiselock System



7-012-35

2.8mm Compression Wire, 45mm Thread For 2.7/3.5mm Av-wiselock System



7-012-36

2.8mm Compression Wire, 50mm Thread For 2.7/3.5mm Av-wiselock System



7-012-37

2.8mm Compression Wire, 55mm Thread For 2.7/3.5mm
Av-wiselock System




7-012-38

2.8mm Compression Wire, 60mm Thread For 2.7/3.5mm
Av-wiselock System




7-012-39

Guide Wire, Ø2.0mm X Length 150mm




7-012-40

Guide Wire With Threaded Tip, Ø2.0mm X Thread Length 10mm
X Length 150mm



7-012-41

Threaded Drill Sleeve For Drill Bit 2.0mm - 2.7/3.5mm
Av-wiselock System



7-012-42

Tap, Ha 2.7mm For 2.7/3.5mm Av-wiselock System



7-012-43

Screwdriver Star Head, T8, For 2.7/3.5mm Av-wiselock System



7-012-44

T-handle With Quick Coupling For 2.7/3.5mm Av-wiselock System



7-012-45

Self-centering Bone Holding Forcep, Speed Lock, 190mm
For 2.7/3.5mm Av-wiselock System

7-012-46

Reduction Forcep, Serrated Jaws, Speed Lock, 160mm
For 2.7/3.5mm Av-wiselock System

7-012-47

Reduction Forcep, Pointed, Ratchet Lock, 180mm
For 2.7/3.5mm Av-wiselock System

7-012-48

Bending Iron 1 For 2.7/3.5mm Av-wiselock System




7-012-49

Bending Iron 2 For 2.7/3.5mm Av-wiselock System



7-012-50

Ø1.2mm Kirschner Wire For 2.7/3.5mm Av-wiselock System



7-012-51

Hohmann Retractor, 6.5mm For 2.7/3.5mm Av-wiselock System



7-012-52

Hohmann Retractor, 8.5mm For 2.7/3.5mm Av-wiselock System



7-012-53

Hohmann Retractor, 15.5mm For 2.7/3.5mm Av-wiselock System



7-012-54

Periosteal Elevator With Silicon Handle, Curved, 6mm
For 2.7/3.5mm Av-wiselock System


7-012-55

Screw Holding Forcep For 2.7/3.5mm Av-wiselock System



7-012-56

Countersink For 2.7/3.5mm Av-wiselock System



7-012-57

Screw Holding Sleeve For 2.7mm Cortical Screw



7-012-60

Screw Holding Sleeve For 2.7mm Locking Screw



7-012-61

Screw Holding Sleeve For 3.5mm Cortical Screw



7-012-62

Screw Holding Sleeve For 3.5mm Locking Screw



7-050-05

T- Handle For Threaded Drill Sleeve




7-012-58

Instrument Trays For 2.7/3.5mm Av-wiselock System



7-012-59

Container For 2.7/3.5mm Av-wiselock System Instrument Set



7-012 : 2.7/3.5mm AV-Wiselock System Instrument Set

ORDER INFORMATION

7-012	2.7/3.5mm AV-Wiselock System Instrument Set
-------	---

Code	Set Consisting Of	Unit
7-012-01	Conical Extraction Screw, Ø3.5mm, for 2.7/3.5mm AV-Wiselock System	1
7-012-02	Drill Bit with Quick Coupling with Stopper, Ø2.8mm x Length 165mm, for 2.7/3.5mm AV-Wiselock System	2
7-012-03	Star Screwdriver Shaft, T8, for 2.7/3.5mm AV-Wiselock System	2
7-012-04	Depth Gauge measuring upto 60mm for 2.7/3.5mm AV-Wiselock System	1
7-012-05	Drill Bit with Quick Coupling, Ø2.0mm x Length 140mm, for 2.7/3.5mm AV-Wiselock System	2
7-012-06	Bone Tap Quick Coupling for Cortical Screws, Ø3.5mm, for 2.7/3.5mm AV-Wiselock System	1
7-012-07	Screwdriver Shaft, T15, for 2.7/3.5mm AV-Wiselock System	2
7-012-08	Threaded Drill Sleeve for Ø2.8mm Drill Bit - 2.7/3.5mm AV-Wiselock System	3
7-012-09	Self-Centering Double Drill Guide, Ø2.5/3.5mm, for 2.7/3.5mm AV-Wiselock System	1
7-012-10	Torque Limiting Attachment, 1.2Nm, for 2.7/3.5mm AV-Wiselock System	1
7-012-11	Torque Limiting Attachment, 1.5Nm, for 2.7/3.5mm AV-Wiselock System	1
7-012-12	Drill Bit with Quick Coupling End, Ø2.5mm x Length 180mm - 2.7/3.5mm AV-Wiselock System	2
7-012-13	Drill Bit with Quick Coupling End, Ø3.5mm x Length 110mm - 2.7/3.5mm AV-Wiselock System	1
7-012-14	2.0mm Variable Angle Guide for 2.7/3.5mm AV-Wiselock System	1
7-012-15	Self-Centering Double Drill Guide, Ø2.0/2.7mm, for 2.7/3.5mm AV-Wiselock System	1
7-012-16	2.8mm Variable Angle Double Drill Guide	1
7-012-17	Handle with Quick Coupling for 2.7/3.5mm AV-Wiselock System	1
7-012-18	Handle for Torque Limiting Attachment, for 2.7/3.5mm AV-Wiselock System	1
7-012-19	2.5Nm Torque Limiting Handle for 2.7/3.5mm AV-Wiselock System	1
7-012-20	Depth Gauge measuring upto 100mm, for 2.7/3.5mm AV-Wiselock System	1
7-012-21	Percutaneous Depth Gauge Adaptor for 3.5mm Screws - 2.7/3.5mm AV-Wiselock System	1
7-012-22	Kirschner Wire with One Sided Trocar Tip Ø1.6mm x Length 150mm	5
7-012-23	Periarticular Reduction Forcep, Small, for 2.7/3.5mm AV-Wiselock System	1
7-012-24	Periarticular Reduction Forcep, Medium, for 2.7/3.5mm AV-Wiselock System	1
7-012-25	Large Compression Forcep with Speed Lock for 2.7/3.5mm AV-Wiselock System	1
7-012-26	Large Distraction Forcep with Speed Lock for 2.7/3.5mm AV-Wiselock System	1
7-012-27	3.5mm Post for AV-Wiselock Locking Hole	2
7-012-28	2.8mm Compression Wire, 10mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-29	2.8mm Compression Wire, 15mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-30	2.8mm Compression Wire, 20mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-31	2.8mm Compression Wire, 25mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-32	2.8mm Compression Wire, 30mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-33	2.8mm Compression Wire, 35mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-34	2.8mm Compression Wire, 40mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-35	2.8mm Compression Wire, 45mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-36	2.8mm Compression Wire, 50mm Thread for 2.7/3.5mm AV-Wiselock System	2

7-012**2.7/3.5mm AV-Wiselock System Instrument Set**

Code	Set Consisting Of	Unit
7-012-37	2.8mm Compression Wire, 55mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-38	2.8mm Compression Wire, 60mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-39	Guide Wire, Ø2.0mm x Length 150mm	5
7-012-40	Guide Wire with Threaded Tip, Ø2.0mm x Thread Length 10mm x Length 150mm	2
7-012-41	Threaded Drill Sleeve for Drill Bit 2.0mm - 2.7/3.5mm AV-Wiselock System	3
7-012-42	Tap, HA 2.7mm for 2.7/3.5mm AV-Wiselock System	1
7-012-43	Screwdriver Star Head, T8, for 2.7/3.5mm AV-Wiselock System	1
7-012-44	T-Handle with Quick Coupling for 2.7/3.5mm AV-Wiselock System	1
7-012-45	Self-Centering Bone Holding Forcep, Speed Lock, 190mm for 2.7/3.5mm AV-Wiselock System	2
7-012-46	Reduction Forcep, Serrated Jaws, Speed Lock, 160mm for 2.7/3.5mm AV-Wiselock System	1
7-012-47	Reduction Forcep, Pointed, Ratchet Lock, 180mm for 2.7/3.5mm AV-Wiselock System	1
7-012-48	Bending Iron 1 for 2.7/3.5mm AV-Wiselock System	1
7-012-49	Bending Iron 2 for 2.7/3.5mm AV-Wiselock System	1
7-012-50	Ø1.2mm Kirschner Wire for 2.7/3.5mm AV-Wiselock System	3
7-012-51	Hohmann Retractor, 6.5mm for 2.7/3.5mm AV-Wiselock System	1
7-012-52	Hohmann Retractor, 8.5mm for 2.7/3.5mm AV-Wiselock System	1
7-012-53	Hohmann Retractor, 15.5mm for 2.7/3.5mm AV-Wiselock System	2
7-012-54	Periosteal Elevator with Silicon Handle, Curved, 6mm for 2.7/3.5mm AV-Wiselock System	1
7-012-55	Screw Holding Forcep for 2.7/3.5mm AV-Wiselock System	1
7-012-56	Countersink for 2.7/3.5mm AV-Wiselock System	1
7-012-57	Screw Holding Sleeve for 2.7mm Cortical Screw for 2.7/3.5mm AV-Wiselock System	1
7-012-58	Instrument Trays for 2.7/3.5mm AV-Wiselock System	1
7-012-59	Container for 2.7/3.5mm AV-Wiselock System Instrument Set	1
7-012-60	Screw Holding Sleeve for 2.7mm Locking Screw	1
7-012-61	Screw Holding Sleeve for 3.5mm Cortical Screw	1
7-012-62	Screw Holding Sleeve for 3.5mm Locking Screw	3
7-050-05	T- Handle for Threaded Drill Sleeve	1

GOING
GROWING



STRONG
STRONGER