

Auxein

EXTERNAL FIXATORS



NEXT-LEVEL CARE

ABOUT US

Auxein Medical is an integrated, research based, orthopaedic implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavour to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in India and serves customers in over 75 countries worldwide.



IMPACTING HEALTHCARE GLOBALLY

Our Reach & Influence

500+

Employees Worldwide

75+

Operations in Countries

150+

Distributors Worldwide

20_{M+}
Patients Cured

200+

FDA Approved Products

250+

Machines Worldwide

3000+

Innovative Orthopedic
Products & Solutions

Auxein is dedicated to advancing medical technology through innovative orthopedic solutions, including **Trauma, Arthroscopy, and Reconstructive Procedures**, with a commitment to improving patient outcomes worldwide. Our approach integrates advanced research, precision engineering, and high-quality manufacturing to develop reliable implants and surgical instruments that address the evolving needs of modern orthopedic practice.


OUR ACHIEVEMENTS



OVER THE YEARS

Disclaimer: Not all products in this catalogue are covered by European CE, US FDA, or other certifications. Please contact our Sales Team for product-specific certification details.






At Auxein, it is our goal to provide the best in class quality and value to every individual. Consistency is a virtue we pursue with a passion. From the procurement of raw materials to the manufacturing of our range of 1200+ orthopaedic implant products, we employ world-class systems and technology at par with the highest standards of globally recognized evaluation systems. The External Fixator System is part of our orthopedic trauma portfolio designed for temporary or definitive stabilization of fractures and deformities. Developed with clinical input, the system enables stable external fixation through percutaneous pin or wire placement, allowing fracture alignment, compression or distraction, and soft tissue management while preserving the biological environment.

The background of the page is a light cream color, decorated with thin, grey, wavy lines that resemble topographical map contours. These lines are concentrated in the top right and bottom left corners, creating a sense of depth and movement.

INDEX




Extra Caliber External Fixator	10-19
Dynamic Axial Fixator	20-33
Small Dynamic Axial Fixator	34-39
Fixation For HTO	40-47
Pertrochanteric Fixator	48-51
Hybrid Fixation System and Ring Fixator	52-67
Adult Limb Reconstruction System	68-76
Pediatric Limb Reconstruction System	78-85
Pelvic Fixator	86-93
Elbow Fixators	94-101
Dynamic Wrist Fixator	102-109
Radiolucent Wrist Fixator	110-111
Mini Fixators	112-117
Minirail System	118-121
Bone Screw	122-129
Bonofix 5 External Fixation System	130-133
Bonofix 8 External Fixation System	134-143
AO III External Fixation System	144-149
LineFix and Calles Fixator	150-153
AO Type External Fixators	154-159
Ilizarov Ring Fixators	160-177
Distal Radius Fixators	178-179



EXTRA CALIBER EXTERNAL FIXATOR




Extra - Caliber Fixator Kits Complete with Universal Allen Wrench and Clamp Cap



EF 160
Radiolucent Tibia Kit, Long




EF161
Radiolucent Tibia Kit, Short



EF162
Hybrid Kit, Long



EF163
Hybrid Kit, Short



EF164
Radiolucent Ankle Kit, Long




EF165
Radiolucent Ankle Kit, Short



EF949
Dynamic Compression-Distractor Unit (a)



COMPONENTS



EF 23
Supplementary Screw Holder
(Bar and Clamps), Sterile



EF24
Supplementary Screw Holder
(Bar, Clamps and Fixator Body
Attachment)



EF25
Supplementary Screw Holder
Bar Radiolucent Length 150mm

EF26
Supplementary Screw Holder
Bar Radiolucent Length 300mm



EF27
Supplementary Screw
Holder Clamp



EF28
Post, Length 50mm

EF29
Post, Length 100mm

EF30
Post, Length 150mm



TEMPLATE AND INSTRUMENTATION



EF33

RadioLucent Ankle Pin Guide, Sterile



Screw Guides


EF44 - Length 45mm

EF45 - Length 60mm

EF46 - Length 80mm

EF47 - Length 100mm

EF53 - Length 160mm



EF57

Ruland Pilot Wire Guide Ø 2.0mm
Length 75mm

EF58

Ruland Pilot Wire Guide Ø 2.0mm
Length 115mm



Drill Guides - Ø 3.2mm

EF62 - Length 20mm

EF63 - Length 40mm

EF64 - Length 80mm

Drill Guides - Ø 4.8mm

EF65 - Length 40mm

EF66 - Length 60mm

EF67 - Length 80mm

EF68 - Length 140mm



INSTRUMENTS



EF69

Hammer

EF70

Hammer Interchangeable Head

EF71

Hammer Locking Nut



EF72

Tapered Trocar



EF76

Extra-Caliber Manipulation Forceps
(1 pair)



EF77

Universal Allen Wrench



EF78

Allen Wrench 3mm



INSTRUMENTS



Drill Bit Ø 3.2mm

EF86 - Length 140mm

EF87 - Length 200mm

Drill Bit Ø 4.8mm

EF88 - Length 180mm

EF89 - Length 240mm

EF90 - Length 280mm



Cannulated Bits Ø3.2mm

EF94 -Length 150mm


EF84 -Length 300mm

Cannulated Bits Ø4.8mm

EF96 -Length 150mm

EF97 -Length 250mm

EF98 -Length 280mm



EF82

Drill Bit Stop Units Ø3.2mm

EF83

Drill Bit Stop Units Ø4.8mm

X-Wires without Olive

EF99 - Ø1.5mm, Length 250mm

EF100 - Ø2.0mm, Length 150mm




EF108- T-wrench



EF01

Extra-Caliber Radiolucent Tibia/Ankle Instruments Sterilization Box, Empty




Codes	Set consisting of:	Units
EF01-T	Tray for Extra-Caliber Radiolucent Tibia/Ankle Instruments Sterile	2
EF80	Allen Wrench 6mm	1
EF87	Drill Bit Ø3.2mm x Length 200mm	1
EF88	Drill Bit Ø4.8mm x Length 180mm	1
EF89	Drill Bit Ø4.8mm x Length 240mm	1
EF94	Cannulated Drill Bit Ø3.2mm x Length 150mm	1
EF63	Drill Guide Ø3.2mm Length 40mm	1
EF64	Drill Guide Ø3.2mm Length 80mm	1
EF65	Drill Guide Ø4.8mm Length 40mm	1
EF67	Drill Guide Ø4.8mm Length 80mm	1
EF45	Screw Guide, Length 60mm	4
EF47	Screw Guide, Length 100mm	4
EF78	Allen Wrench 3mm	2
EF108	T-Wrench	1
EF72	Tapered Trocar	1
EF76	Extra-caliber Manipulation Forceps (1 pair)	1
EF24	Supplementary Screw Holder (Bar, Clamp and Fixator Body Attachment)	1
EF33	RadioLucent Ankle Pin Guide, Sterile	1
EF99	X-Wires without Olive Ø1.5mm x Length 250mm	2
EF740	Screw Connector	1
EF83	Drill Bit Stop Unit Ø4.8mm	1
EF82	Drill Bit Stop Unit Ø3.2mm	1



EF02

Hybrid Instruments Steribox, Empty




Codes	Set consisting of:	Units
EF02-T	Tray for Hybrid Instruments Sterile	2
EF46	Screw Guide, Length 80mm	3
EF66	Drill Guide Ø4.8mm Length 60mm	1
EF89	Drill Bit Ø4.8mm x Length 240mm	1
EF80	Allen Wrench 6mm	1
EF72	Tapered Trocar	1
EF37	Wire Tensioner	1
EF38	Wire Cutter	1
EF26	Supplementary Screw Holder Bar Radiolucent Length 300mm	2
EF36	Wire Guide	1
EF28	Post, Length 50mm	2
EF440	Two-Hole X-Wire Securing Pin	2
EF445	Three-Hole Wire Clamp Slider Unit	4
EF104	X-Wire with lateral Olive Ø2mm x Length 350mm	4
EF101	X-Wire without Olive Ø2mm x Length 400mm	4
EF114	X-Wires with Central Olive Ø2mm x Length 400mm	2
EF118	Open End Wrench Ø10mm	2
EF441	Independent Wire Clamp	2
EF444	Single Screw Clamp	2
EF27	Supplementary Screw Holder Clamp	4
EF78	Allen Wrench 3mm	1
EF83	Drill Bit Stop Unit Ø4.8mm	1
EF207	Hybrid Coupling with Ball-joint	1



EF03

Extra-Caliber Hybrid System Ring Steribox, Empty



Codes	Set consisting of:	Units
EF03-T	Tray for Extra-Caliber Hybrid System Ring Sterile	1
EF429	2/3 Radiolucent Ring Ø150mm	2
EF435	1/3 Ring, Inner Ø150mm (Not illustrated)	2
EF430	2/3 Radiolucent Ring Ø175mm	2
EF436	1/3 Ring, Inner Ø175mm (Not illustrated)	2
EF431	2/3 Radiolucent Ring Ø190mm	2
EF437	1/3 Ring, Inner Ø190mm (Not illustrated)	2
EF433	Ring Locking Screw	2







DYNAMIC AXIAL FIXATOR

**EF135**

D.A.F. Fixator Kit, Short
Complete With appropriate
Compression-Distraction Unit
and Allen Wrench 6mm

**EF140**

D.A.F. Fixator Kit, Standard
Complete With appropriate
Compression-Distraction Unit
and Allen Wrench 6mm

**EF141**

D.A.F. Fixator Kit, Long
Complete With appropriate
Compression-Distraction Unit
and Allen Wrench 6mm



BALL-JOINTED MODULES




EF202

D.A.F. Straight Clamp



EF203

T-Clamp



EF204

Ball-Jointed Modules
For Independent Screw Placement
For iliac crest applications.

Consisting of:

*1 Bar with Ball-Joint and 2 independent Screw Clamp.
Enables screws to be placed freehand in chosen positions
on the iliac crest.*



EF205

Elbow Extended Range Clamp




EF206

Radiolucent Ankle Clamp



BALL-JOINTED MODULES



EF207
Hybrid Coupling with Ball-Joint



EF208
Metaphyseal Clamp



EF19
Torbay-Garches Clamp



EF209
Trauma-Garches Clamp




EF210
Angled Clamp



EF211
Double-Coupling Clamp



NON BALL-JOINTED MODULES



EF24

Supplementary Screw Holder
Bar, Clamps and Fixator Body Attachment



EF27

Supplementary Screw Holder Clamp



EF156

Self-aligning Articulated Body Assembly Kit
Consisting of: 1 Self-aligning Articulated Body,
1 Standard Body Male Component,
1 Central Body Locking Nut, 2 Bushes,
2 Cams, 1 ProCallus T-Clamp,
1 Micromovement Locking Nut,
1 Compression-Distraction Unit Standard.
See section "B. Dynamic Axial Fixator" for components.



COMPONENTS / REPLACEMENT




D.A.F. Central Bodies

Complete With 2 Cams, 2 Bushes,
1 D.A.F. Micromovement Locking Nut
and 1 D.A.F. Central Body Locking Nut

EF214 - Standard

EF221 - Long

EF222 - Short



D.A.F. Central Body Male Components

EF757 - Standard

EF759 - Long

BF761 - Short



D.A.F. Central Body Female Components

EF758 - Standard

EF760 - Long

EF762 - Short



EF215

D.A.F. Micromovement
Locking Nut



EF216

D.A.F. Central Body Locking Bolt



EF217

Bush



EF218

Cam



EF219

Clamp Screw for Micrometric
Swivelling Clamp and Multiplaner Clamp



EF220

RadioLucent Ankle Clamp
Articulation Locking Nut



COMPONENTS / REPLACEMENT



Compression-Distractor Units

EF227

Standard, Extends to 4cm

EF228

Long, Extends to 8cm

EF229

Short, Extends to 2.75cm

**EF278**

Actuators Standard, Extends to 4cm

EF279

Actuators Long, Extends to 8cm

EF280

Actuators Short, Extends to 2.75cm



TEMPLATES

**EF21**

RadioLucent Ankle Pin Guide

**EF33**RadioLucent Ankle
Pin Guide, Sterile



INSTRUMENTATION



Screw Guides

EF43 - Length 32mm


EF44 - Length 45mm

EF45 - Length 60mm

EF46 - Length 80mm

EF47 - Length 100mm

EF53 - Length 160mm



Ruland Pilot Wire Guide
Ø2.0mm

EF57 - Length 75mm

EF58 -Length 115mm



Drill Guides Ø2.9mm

EF61 - Length 20mm

Drill Guides Ø3.2mm

EF62 - Length 20mm

EF63 - Length 40mm

EF64 - Length 80mm

Drill Guides Ø4.8mm

EF65 - Length 40mm

EF66 - Length 60mm


EF67 - Length 80mm

EF68 - Length 140mm



EF72

Tapered Trocar



EF80

Allen Wrench 6 mm



INSTRUMENTATION

**EF69**

Hammer

EF70

Hammer Interchangeable Head

EF71

Hammer Locking Nut

**EF78**

Allen Wrench 3mm




Drill Bit Ø2.9mm

EF91 - Length 140mm

Drill Bit Ø3.2mm

EF86 - Length 140mm**EF87** - Length 200mm

Drill Bit Ø4.8mm

EF88 - Length 180mm**EF89** - Length 240mm**EF90** - Length 280mm

Drill Bit Stop Units

EF736 - Ø2.7mm**EF735** - Ø2.9mm**EF82** - Ø3.2mm**EF83** - Ø4.8mm



INSTRUMENTATION



Cannulated Drill Bits Ø3.2mm

EF94 - Length 150mm

EF84 - Length 300mm

Cannulated Drill Bits Ø4.8mm

EF96 - Length 150mm

EF97 - Length 250mm

EF98 - Length 280mm


X-wires without olive

EF99 - Ø1.5mm, Length 250mm

EF100 - Ø2.0mm, Length 150mm



INSTRUMENTATION



EF108- T-wrench



EF109
Torque wrench 6 mm



EF110
Fixator Self-Locking
Manipulation Forceps

Assembly Kit, Consisting of:

Screw Guides

EF45 -Length 60mm

EF47 -Length 100mm

Drill Guides Ø3.2mm

EF63 -Length 40mm


EF64 -Length 80mm

Drill Guides Ø4.8mm

EF65 -length 40mm

EF67 -Length 80mm

EF72 -Tapered Trocar

**EF04****Instruments Steribox, Empty**


Codes	Set consisting of:	Units
EF04-T	Tray for Instruments Steribox	1
EF80	Allen Wrench 6mm	2
EF87	Drill Bit Ø3.2mm x Length 200mm	2
EF88	Drill Bit Ø4.8mm x Length 180mm	2
EF89	Drill Bit Ø4.8mm x Length 240mm	2
EF78	Allen Wrench 3mm	2
EF108	T-Wrench	1
EF63	Drill Guide Ø3.2mm Length 40mm	1
EF64	Drill Guide Ø3.2mm Length 80mm	1
EF65	Drill Guide Ø4.8mm Length 40mm	1
EF67	Drill Guide Ø4.8mm Length 80mm	1
EF45	Screw Guide, Length 60mm	6
EF47	Screw Guide, Length 100mm	6
EF109	Torque Wrench 6mm	1
EF83	Drill Bit Stop Unit Ø4.8mm	1
EF82	Drill Bit Stop Unit Ø3.2mm	1
EF72	Tapered Trocar	1

EF05**D.A.F. System with Actuator Steribox, Empty**

Codes	Set consisting of:	Units
EF05-T	Tray for D.A.F. System with Actuator Steribox	1
EF135	D.A.F. Fixator Kit, Short Complete With appropriate Compression-Distraction Unit and Allen Wrench 6mm	1
EF140	D.A.F. Fixator Kit, Standard Complete With appropriate Compression-Distraction Unit and Allen Wrench 6mm	1
EF141	D.A.F. Fixator Kit, Long Complete With appropriate Compression-Distraction Unit and Allen Wrench 6mm	1
EF203	T-Clamp	1
EF108	T-Wrench	1
EF740	Screw Connector	1
EF24	Supplementary Screw Holder (Bar, Clamp and Fixator Body Attachment)	1
EF80	Allen Wrench 6mm	2



SMALL DYNAMIC AXIAL FIXATOR



Ef133

Small D.A.F. Kit
Complete with Compression-
Distraction unit and Allen Wrench 5mm

BALL -JOINTED MODULES



EF212

Small D.A.F. Straight Clamp



EF213

Small D.A.F. T-Clamp

NON BALL-JOINTED MODULES



EF27

Supplementary Screw Holder Clamp



COMPONENTS / REPLACEMENTS

**EF292**

Small D.A.F. Central Body
Locking Nut

**EF218**

Cam

**EF293**


Small D.A.F. Bush

**EF265**

Clamp Screw

**EF296**

Small D.A.F. Clamp Cover

**EF269**

Compression-Distractor Unit for
Small Fixator, Extends to 5cm

**EF27**

Supplementary Screw
Holder Clamp



INSTRUMENTATION



Cannulated Drill Bit Ø3.2mm
EF94 - Length 150mm

Cannulated Drill Bit Ø4.8mm.
EF96 - Length 150mm



Screw Guides

EF43 -Length 32mm

EF44 -Length 45mm

EF45 -Length 60mm



Drill Guides Ø2.9mm

EF61 -Length 20mm

Drill Guides Ø3.2mm

EF62 -Length 20mm

EF63 -Length 40mm

Drill Guides Ø4.8mm

EF65 -Length 40mm

EF66 -Length 60mm



EF72

Tapered Trocar



EF78

Allen Wrench 3mm




INSTRUMENTATION




Drill Bit Ø2.9mm
EF91 - Length 140mm

Drill Bit Ø3.2mm
EF86 - Length 140mm

Drill Bit Ø4.8mm
EF88 - Length 180mm



Drill Bit Stop Units
EF736 - Ø2.7mm
EF735 - Ø2.9mm
EF82 - Ø3.2mm
EF83 - Ø4.8mm



EF79
Allen wrench 5mm



EF69
Hammer

EF70
Hammer Interchangeable Head

EF71
Hammer Locking Nut



INSTRUMENTATION

X-Wires Without Olive

EF99 - Ø1.5mm, Length 250mm

EF100 - Ø2.0mm, Length 150mm



EF108-T-wrench



EF116
Torque Wrench 5mm

EF133 - Small D.A.F.

Complete with
Compression-
Distraction unit and
Allen Wrench 5mm

Screw Guides

EF43- Length 32mm

EF45- Length 60mm

Drill Guide Ø3.2mm

EF62 - Length 20mm

EF63 - Length 40mm

Drill Guide Ø4.8mm

EF65 - Length 40mm

EF72 -Tapered Trocar

EF110
Fixator Self-Locking
Manipulation Forceps





FIXATION FOR HTO

HIGH TIBIAL OSTEOTOMY

**EF156**

Self-aligning Articulated Body Assembly kit

*Consisting of: 1 Self-aligning Articulated Body,
1 Standard Body Male Component, 1 Central Body Locking Nut,
2 Bushes, 2 Cams, 1 ProCallus T-Clamp, 1 Micromovement Locking Nut,
1 Compression-Distractor Unit Standard.*

See section "B. Dynamic Axial Fixator" for components.

**EF151**

Fixator Assembly with LRS Rail
and L.R.S. OF-Garches T-Clamp

*Consisting of: 1 OF-Garches T-Clamp,
1 Of-Garches Standard Compression-Distractor Unit,
1 L.R.S. Rail short, 1 L.R.S. Straight Central Clamp.*



COMPONENTS / REPLACEMENTS



EF249
L.R.S. OF-Garches T-Clamp



EF299
OF-Garches Standard T-Clamp



EF306
OF-Garches Compression Distraction
Units, Standard,
Extends to 5.5cm

EF307
OF-Garches Compression Distraction
Units, Long,
Extends to 10cm



EF300
OF-Garches T-Clamp Axis
Locking Nut



EF301
D.A.F. Central Body Locking Nut



EF265
Clamp Screw



EF253
Clamp Cover



INSTRUMENTATION



Screw Guides

EF43 - Length 32mm


EF44 - Length 45mm

EF45 - Length 60mm

EF46 - Length 80mm

EF47 - Length 100mm

EF53 - Length 160mm



EF57

Ruland Pilot wire Guides Ø2.0mm, Length 75mm

EF58

Ruland Pilot wire Guides Ø2.0mm, Length 115mm



Drill Guides; Ø 3.2mm

EF62 - Length 20mm

EF63 -Length 40mm

EF64 - Length 80mm

Drill Guides; Ø 4.8mm

EF65 - Length 40mm

EF66 -Length 60mm

EF67 -Length 80mm

EF68 -Length 140mm



INSTRUMENTATION



EF69
Hammer

EF70
Hammer Interchangeable Head

EF71
Hammer Locking Nut



EF72
Tapered Trocar



EF78
Allen Wrench 3mm



EF80
Allen wrench 6mm




INSTRUMENTATION



Cannulated Drill Bit Ø 3.2mm
EF86 -Length 140mm
EF87 -Length 200mm

Cannulated Drill Bit Ø 4.8mm
EF88 - Length 180mm
EF89 - Length 240mm
EF90 - Length 280mm



EF82
Drill Bit Stop Unit Ø3.2m

EF83
Drill Bit Stop Unit Ø4.8mm



EF116
Torque Wrench 5mm



INSTRUMENTATION



Cannulated Drill Bit Ø3.2mm
EF94 -Length 150mm
EF84 -Length 300mm

Cannulated Drill Bit Ø4.8mm
EF96 -Length 150mm
EF97 -Length 250mm
EF98 - Length 280mm



X-Wires without Olive
EF99 -Ø1.5mm, Length 250mm
EF100 - Ø2.0mm, Length 150mm




EF108 T-wrench





PERTROCHANTERIC FIXATOR



EF142
Pertrochanteric Fixator Complete
with Allen Wrench 5mm

Components / Replacements



EF122
Pertrochanteric Fixator
Angulation Locking Screw



EF123
Pertrochanteric Fixator Distal
Rotation locking Screw



EF124
Pertrochanteric Fixator
Rotation Locking Screw



INSTRUMENTATION



EF731
Pertrochanteric Screw Guide,
Length 150mm.



EF53
Screw Guide, Length 160mm.

EF68
Drill Guide Ø4.8mm
Length 140mm



EF73
Trocar Ø 4.8mm.
Length 240mm.




EF74
Trocar 6mm, Length 240mm.



EF78
Allen wrench 3mm



INSTRUMENTATION



The background of the page features a series of thin, light gray wavy lines that create a sense of movement and depth. These lines are concentrated in the top right and bottom left corners, leaving the central area relatively clear for the text.

HYBRID FIXATION SYSTEM AND RING FIXATOR



EF166

Hybrid Fixator Assembly

Consisting of: 1 2/3 Ring, 2 Two-Hole X-Wire Securing Pins, 2 Three-Hole Wire Clamp Slider Unit, 1 ProCallus Central Body, 1 ProCallus Straight Clamp, 1 Hybrid Coupling with Ball-Joint, 2 Posts, 2 Supplementary Screw Holder Bars, 4 Supplementary Screw Holder Clamps.

EF167

Sheffield Ring Fixator Assembly for Trauma

Consisting of: 2 2/3 Rings, 2 Two-Hole Securing Pins, 2 Three-Hole Wire Clamp Slider Units, 1 Sheffield Clamp, 3 Reduction Units.



EF168

Sheffield Ring Fixator Assembly for Reconstruction

Consisting of: 2 2/3 Ring, 2 1/3 Rings, 4 Locking Screws, 1 Sheffield Clamp, 2 Hinge Clamps, 4 Threaded Bars and 8 Nuts and Washer, 1 Compression-Distraction Clicker Unit, 2 Extension Plates and 2 Single Screw Clamps.



COMPONENTS / REPLACEMENTS



2/3 Rings

EF428 - Inner Ø125mm

EF429 - Inner Ø150mm

EF430 - Inner Ø175mm

EF431 - Inner Ø190mm

EF432 - Inner Ø220mm

2/3 RadioLucent Rings

(Not illustrated)

EF451 - Inner Ø125mm

EF452 - Inner Ø150mm

EF453 - Inner Ø175mm

EF454 - Inner Ø190mm

EF455 - Inner Ø220mm



1/3 Rings

(Not illustrated)

EF434 - Inner Ø125mm

EF435 - Inner Ø150mm

EF436 - Inner Ø175mm

EF437 - Inner Ø190mm

EF438 - Inner Ø220mm

1/3 RadioLucent Rings(Not illustrated)

EF456 - Inner Ø125mm

EF457 - Inner Ø150mm

EF458 - Inner Ø175mm

EF459 - Inner Ø190mm

EF460 - Inner Ø220mm




EF433

Ring Locking Screw



COMPONENTS / REPLACEMENTS



EF439
Foot Ring Inner Ø150mm, Aluminium



EF445
Three-Hole Wire Clamp Slider Unit



EF440
Two-Hole X-Wire Securing Pin



EF441
Independent Wire Clamp



EF207
Hybrid Coupling with Ball-Joint




EF450
Ring Clamp



EF23
Supplementary Screw Holder
Bar and Clamp, sterile



COMPONENTS / REPLACEMENTS



EF27
Supplementary Screw
Holder Clamp



Posts
EF28 - Length 50mm
EF29 - Length 100mm
EF30 - Length 150mm




EF444
Single Screw Clamp



EF446
Hinge Clamp

COMPONENTS / REPLACEMENTS




EF447
Extension Plate



EF448
Translation Rotation Hinge



EF942
Counter Nut



Threaded Bars
with nuts and washers (set of 3)
EF31 - Length 80mm
EF32 - Length 120mm
EF745- Length 160mm
EF42 - Length 200mm



Compression-Distraction Clicker Units
EF288 - Standard, Extends to 4cm.
EF297 - Long, Extends to 8cm.



Bolts
with nuts and washers (set of 3)
EF443 - Length 60mm



COMPONENTS / REPLACEMENTS



EF449

Reduction Unit

3 units are required to connect 2 rings.



Hybrid Fixator Assembly Kits

EF304 - Standard

Consisting of:

EF214 - D.A.F. Standard Central Body

EF227 - Compression-Distractor Unit Standard

Extends to 4cm

EF202 - D.A.F. Straight Clamp

EF207 - Hybrid Coupling with Ball-Joint.

EF305 - Short

Consisting of:

EF222 - D.A.F. Short Central Body

EF229 - Compression-Distractor Unit Short

Extends to 2.75cm

EF202 - D.A.F. Straight Clamp

EF207 - Hybrid Coupling with Ball-Joint.



RING FIXATOR KITS

EF20	Standard Trauma Kit
------	---------------------

Codes	Set consisting of:	Units
EF37	Wire Tensioner	1
EF36	Wire Guide	1
EF153	Cortical Screw, Ø6-5mm Shaft Ø6mm, Thread 40mm Length 130mm	6
EF155	Cortical Screw, Ø6-5mm Shaft Ø6mm, Thread 40mm Length 150mm	6
EF89	Drill Bit Kit Ø4.8mm x Length 240mm	1
EF80	Allen Wrench 6mm	1
EF46	Screw Guide, Length 80mm	3
EF47	Screw Guide, Length 100mm	3
EF67	Drill Guide Ø4.8mm Length 80mm	1
EF38	Wire Cutter	1
EF440	Two-Hole X-Wire Securing Pin	2
EF445	Three-Hole Wire Clamp Slider Unit	3
EF450	Ring Clamp	1
EF447	Extension Plate	2
EF103	X-Wire with lateral Olive Ø2mm x Length 310mm	4
EF104	X-Wire with lateral Olive Ø2mm x Length 350mm	4
EF105	X-Wire with lateral Olive Ø2mm x Length 400mm	4
EF114	X-Wires with Central Olive Ø2mm x Length 400mm	4
EF115	X-Wires with Central Olive Ø2mm x Length 450mm	4
EF101	X-Wire without Olive Ø2mm x Length 400mm	4
EF102	X-Wire without Olive Ø2mm x Length 450mm	4
EF428	2/3 Ring Inner Ø125mm	2
EF429	2/3 Ring Inner Ø150mm	2
EF430	2/3 Ring Inner Ø175mm	2
EF431	2/3 Ring Inner Ø190mm	2
EF432	2/3 Ring Inner Ø220mm	2
EF434	1/3 Ring, Inner Ø125mm (Not illustrated)	2
EF435	1/3 Ring, Inner Ø150mm (Not illustrated)	2
EF436	1/3 Ring, Inner Ø175mm (Not illustrated)	2
EF437	1/3 Ring, Inner Ø190mm (Not illustrated)	2
EF438	1/3 Ring, Inner Ø220mm (Not illustrated)	2
EF433	Ring Locking Screw	4
EF446	Hinge Clamp	2
EF733	Open End Wrench, 10mm	2
EF444	Single Screw Clamp	2
EF449	Reduction Unit (3 Unit are required to connect 2 Ring)	3
EF77	Universal Allen Wrench	1



RING FIXATOR KITS

EF693	Lengthening and Correction of Deformities Kit
--------------	--

Codes	Set consisting of:	Units
EF288	Compression-Distractor Clicker Unit, Standard, Extends to 4cm.	3
EF297	Compression-Distractor Clicker Unit, Long, Extends to 8cm.	3
EF153	Cortical Screw, Ø6-5mm Shaft Ø6mm, Thread 40mm Length 130mm	3
EF155	Cortical Screw, Ø6-5mm Shaft Ø6mm, Thread 40mm Length 150mm	3
EF445	Three-Hole Wire Clamp Slider Unit	3
EF440	Two-Hole X-Wire Securing Pin	2
EF450	Ring Clamp	1
EF447	Extension Plate	4
EF433	Ring Locking Screw	4
EF446	Hinge Clamp	4
EF448	Translation Rotation Hinge	3
EF444	Single Screw Clamp	4

EF920	Cross Joint and Foot Fixation System
--------------	---

Codes	Set consisting of:	Units
EF288	Compression-Distractor Clicker Unit, Standard, Extends to 4cm.	3
EF297	Compression-Distractor Clicker Unit, Long, Extends to 8cm.	3
EF153	Cortical Screw, Ø6-5mm Shaft Ø6mm, Thread 40mm Length 130mm	3
EF155	Cortical Screw, Ø6-5mm Shaft Ø6mm, Thread 40mm Length 150mm	3
EF445	Three-Hole Wire Clamp Slider Unit	3
EF440	Two-Hole X-Wire Securing Pin	2
EF450	Ring Clamp	1
EF447	Extension Plate	2
EF433	Ring Locking Screw	2
EF446	Hinge Clamp	2
EF444	Single Screw Clamp	2

INSTRUMENTATION



Screw Guides
EF46 - Length 80mm
EF47 - Length 100mm




Drill Guide Ø4.8mm
EF67 - Length 80mm



EF78
Allen Wrench 3mm



Drill Bit Ø4.8mm
EF88 - Length 180mm
EF89 - Length 240mm



Drill Bit Stop Unit
EF83 - Ø4.8mm



EF79
Allen Wrench, 5mm

EF80
Allen Wrench, 6mm



EF733
Open End Wrench 10mm



INSTRUMENTATION



EF77

Universal Allen Wrench




EF36

Wire Guide



EF37

Wire Tensioner



EF38


Wire Cutter

EF39

Wire Cutter Re-use Kit
Consisting of 2 blades and
2 rubbers
(not illustrate)

EF07

Hybrid Fixation System Sterilization Box, Empty





EF07	Hybrid Fixation System Sterilization Box, Empty
-------------	--

Codes	Set consisting of:	Units
EF07-T	Tray for Hybrid Fixation System Sterilization Box	3
EF429	2/3 Ring Inner Ø150mm	2
EF435	1/3 Ring, Inner Ø150mm (Not illustrated)	2
EF431	2/3 Ring Inner Ø190mm	2
EF437	1/3 Ring, Inner Ø190mm (Not illustrated)	2
EF433	Ring Locking Screw	8
EF440	Two-Hole X-Wire Securing Pin	8
EF445	Three-Hole Wire Clamp Slider Unit	8

Lower Tray


Codes	Set consisting of:	Units
EF36	Wire Guide	1
EF37	Wire Tensioner	1
EF304	Hybrid Fixator Standard Assembly Kit Consisting of: 1 D.A.F. Standard Central Body, 1 Compression-Distraction Unit Standard, 1 D.A.F. Straight Clamp, 1 Hybrid Coupling with Ball-Joint.	1
EF38	Wire Cutter	1
EF27	Supplementary Screw Holder Clamp	8

Upper Tray

Codes	Set consisting of:	Units
EF80	Allen Wrench 6mm	2
EF79	Allen Wrench 5mm	2
EF78	Allen Wrench 3mm	2
EF441	Independent Wire Clamp	2
EF28	Post, Length 50mm	3
EF29	Post, Length 100mm	3
EF30	Post, Length 150mm	3
EF105	X-Wire with lateral Olive Ø2mm x Length 400mm	8
EF114	X-Wires with Central Olive Ø2mm x Length 400mm	4
EF115	X-Wires with Central Olive Ø2mm x Length 450mm	4
EF101	X-Wire without Olive Ø2mm x Length 400mm	4
EF102	X-Wire without Olive Ø2mm x Length 450mm	4

EF08

Hybrid Fixation System Sterilization Box, Empty





EF08	Hybrid Fixation System Sterilization Box, Empty
-------------	--

Codes	Set consisting of:	Units
EF08-1	Tray for Hybrid Fixation System Sterilization Box.	3
EF432	2/3 Ring Inners Ø220mm	2
EF438	1/3 Ring, Inner Ø220mm (Not illustrated)	2
EF430	2/3 Ring Inners Ø175mm	2
EF436	1/3 Ring, Inner Ø175mm (Not illustrated)	2

Lower Tray

Codes	Set consisting of:	Units
EF450	Ring Clamp	2
EF444	Single Screw Clamp	2
EF449	Reduction Unit- 3 Unit are required to connect 2 Ring	3
EF46	Screw Guide, Length 80mm	3
EF47	Screw Guide, Length 100mm	3
EF67	Drill Guide Ø4.8mm Length 80mm	1
EF66	Drill Guide Ø4.8mm Length 60mm	1
EF739	Bolts 60mm (set of 3)	1
EF31	Threaded Bars with Nuts and Washers (set of 3), Length 80mm	1
EF32	Threaded Bars with Nuts and Washers (set of 3), Length 120mm	1
EF745	Threaded Bars with nuts and washers (set of 3), Length 160mm	1
EF42	Threaded Bars with Nuts and Washers (set of 3), Length 200mm	1

Upper Tray

Codes	Set consisting of:	Units
EF37	Wire Tensioner	1
EF433	Ring Locking screw	8
EF440	Two-Hole X-Wire Securing Pin	8
EF305	Hybrid Fixator Short Assembly Kit Consisting of: 1 D.A.F. Straight Clamp, 1 D.A.F. Short Central Body, 1 Compression- Distraction Unit Short, 1 Hybrid Coupling with Ball-Joint	1
EF449	Reduction Unit : 3 Unit are required to connect 2 Ring	3
EF445	Three-Hole Wire Clamp Slider Unit	8




The background of the page features a series of thin, light gray wavy lines that create a sense of movement and depth. These lines are concentrated in the top right and bottom left corners, leaving the central area relatively clear for the text.


ADULT LIMB RECONSTRUCTION SYSTEM

ADULT LIMB RECONSTRUCTION SYSTEM (L.R.S.) KITS


Complete with Compression-Distraction Units and Allen Wrench 6mm.



EF143 - Adult Limb Reconstruction System (L.R.S.) Kit, Long
With 2 L.R.S. Straight Outer Clamp
EF243 and 1 L.R.S. Straight
Central Clamp EF244



EF 144 - Adult Limb Reconstruction System (L.R.S.) Kit, Standard
With 2 L.R.S. Straight Outer Clamps
EF243 and 1 L.R.S. Straight
Central Clamp EF244



EF145 - Adult Limb Reconstruction System (L.R.S.) Kit, Short
With 2 L.R.S. Straight Outer Clamps
EF243 and 1 L.R.S. Straight
Central Clamp EF244

EF 146 - Adult Limb Reconstruction System (L.R.S.) Kit, Extra Short

COMPONENTS / REPLACEMENTS



L.R.S. Rails

EF239 - Long, Length 40cm.

EF240 - Standard, Length 30cm.

EF241 - Short, Length 23cm.

EF242 - Extra-Short Rail, Length 12cm.



EF243

L.R.S. Straight Outer Clamp



EF244

L.R.S. Straight Central Clamp



EF245

L.R.S. T-Clamp



EF246

L.R.S. Ball-Joint Coupling



EF247


L.R.S. Micrometric Swivelling
Clamp



EF248

L.R.S. Multiplanar Clamp

COMPONENTS / REPLACEMENTS



EF249
L.R.S. OF-Garches T-Clamp



EF250
L.R.S. Hinged Rail Coupling



EF251
L.R.S. Sandwich Clamp with
Compression-Distraction Unit
Attachments



EF252
L.R.S. Dyna-Ring with
Compression-Distraction
Unit Attachment



EF253
Clamp Cover

COMPONENTS / REPLACEMENTS



EF219 - Clamp Screw
For Micrometric Swivelling
Clamp and Multiplanar Clamp.



EF300
OF-Garches T -Clamp Axis
Locking Nut



EF265
Clamp Screw



EF217
Bush



EF218
Cam



Compression-Distractor Unit
EF227 - Standard, Extends to 4cm.
EF228 - Long, Extends to 8cm.



Compression-Distractor Clicker Unit
EF288 - Standard, Extends to 4cm.
EF297 - Long, Extends to 8cm.



OF-Garches Compression-Distractor Unit
EF306 - Standard, Extends to 5.5cm.
EF307 - Long, Extends to 10cm.

INSTRUMENTATION



Ruland Pilot wire Guides
Ø2.0mm
EF57 - Length 75mm
EF58 - Length 115mm



Screw Guides
EF43 - Length 32mm
EF44 - Length 45mm
EF45 - Length 60mm
EF46 - Length 80mm
EF47 - Length 100mm
EF53 - Length 160mm



Drill Guides Ø2.9mm
EF61 - Length 20mm



Drill Guides Ø3.2mm
EF62 - Length 20mm
EF63 - Length 40mm
EF64 - Length 80mm



Drill Guides Ø4.8mm
EF65 - Length 40mm
EF66 - Length 60mm
EF67 - Length 80mm
EF68 - Length 140mm



EF78
Allen Wrench 3mm



EF80
Allen Wrench 6 mm



EF108 - T-wrench




INSTRUMENTATION



Drill Bit Ø2.9mm
EF91 - Length 140mm

Drill Bit Ø3.2 mm
EF86 - Length 140mm
EF87 - Length 200mm

Drill Bit Ø4.8mm
EF88 - Length 180mm
EF89 - Length 240mm
EF90 - Length 280mm



Drill Bit Stop Units

EF735 - Ø2.9mm

EF82 - Ø3.2mm

EF83 - Ø4.8mm



EF109
Torque wrench 6mm



EF72
Tapered Trocar




EF69
Hammer

EF70
Hammer Interchangeable Head

EF71
Hammer Locking Nut

EF10

L.R.S Instrument Sterilization Box, Empty





EF10	L.R.S Instrument Sterilization Box, Empty
-------------	--

Codes	Set consisting of:	Units
EF10-T	Tray for L.R.S. Instrument Sterilization Box	3
EF78	Allen Wrench 3mm	2
EF80	Allen Wrench 6mm	2
EF109	Torque Wrench 6mm	1
EF108	T-Wrench	2
EF72	Tapered Trocar	1
EF116	Torque Wrench 5mm	1
EF95	Cannulated Drill Bit Ø3.2mm x Length 200mm	2
EF98	Cannulated Drill Bit Ø4.8mm x Length 280mm	2
EF57	Ruland Pilot Wire Guide Ø2.0mm x Length 75mm	2
EF58	Ruland Pilot Wire Guide Ø2.0mm x Length 115mm	2
EF79	Allen Wrench 5mm	2
EF746	Allen Wrench 4mm	2
EF69	Hammer	1


Lower Tray

Codes	Set consisting of:	Units
EF65	Drill Guide Ø4.8mm Length 40mm	2
EF67	Drill Guide Ø4.8mm Length 80mm	2
EF63	Drill Guide Ø3.2mm Length 40mm	2
EF64	Drill Guide Ø3.2mm Length 80mm	2
EF68	Drill Guide Ø4.8mm Length 140mm	2
EF101	X-Wire without Olive Ø2mm x Length 400mm	5
EF89	Drill Bit Ø4.8mm x Length 240mm	2
EF87	Drill Bit Ø3.2mm x Length 200mm	2
EF90	Drill Bit Ø4.8mm x Length 280mm	2
EF83	Drill Bit Stop Unit Ø4.8mm	2
EF82	Drill Bit Stop Unit Ø3.2mm	2

Upper Tray

Codes	Set consisting of:	Units
EF45	Screw Guide, Length 60mm	10
EF47	Screw Guide, Length 100mm	10
EF53	Screw Guide, Length 160mm	10



The background of the page features a series of thin, light gray wavy lines that flow from the top right towards the bottom left, creating a sense of movement and depth.


PEDIATRIC LIMB RECONSTRUCTION SYSTEM

PEDIATRIC LIMB RECONSTRUCTION SYSTEM (L.R.S.) KITS


Complete with Compression-Distraction Units and Allen Wrench 5mm.



EF147- Pediatric Limb Reconstruction System (L.R.S.) Kit, Long With 2 Pediatric L.R.S. Straight Clamp (20010205)



EF148- Pediatric Limb Reconstruction System (L.R.S.) Kit, Standard With 2 Pediatric L.R.S. Straight Clamp (20010205)



EF149- Pediatric Limb Reconstruction System (L.R.S.) Kit, Short With 2 Pediatric L.R.S. Straight Clamp (20010205)



EF150 - Extra Short



COMPONENTS / REPLACEMENTS



Pediatric L.R.S. Rails

EF254 - Long, Length 25cm.

EF255 - Standard, Length 20cm.

EF256 - Short, Length 15cm.

EF734 - Extra-Short Length 10cm.



EF257

Pediatric L.R.S. Straight Clamp



EF258

Pediatric L.R.S. T-Clamp



EF259

Pediatric L.R.S. Ball-Joint Coupling



EF260

Pediatric L.R.S. Swivelling Clamp



EF261

Pediatric L.R.S. OF-Garches T-Clamp Kit



EF262

Pediatric L.R.S. Hinged Rail Coupling

COMPONENTS / REPLACEMENTS



EF263
Pediatric L.R.S. Dyna-Ring with
Compression-Distractor Unit Attachment



EF264
Pediatric L.R.S. locking Screw




EF265
Clamp Screw



EF266
Small Blue D.A.F. Bush



EF218
Cam



EF296
Small D.A.F. Clamp Cover



EF269 - Compression-Distractor Unit for Small Fixator, Extends to 5cm
EF270 - Pediatric L.R.S. Compression-Distractor Unit, Extends to 6.2cm



INSTRUMENTATION



Screw Guides

- EF43 - Length 32mm
- EF44 - Length 45mm
- EF45 - Length 60mm
- EF46 - Length 80mm
- EF47 - Length 100mm
- EF53 - Length 160mm



- Ruland Pilot Wire Guides
Ø2.0mm
- EF57 -Length 75mm
 - EF58 -Length 115mm




- Drill Guides Ø2.9mm
- EF61 - Length 20mm

- Drill Guides Ø3.2mm
- EF62 - Length 20mm
 - EF63 - Length 40mm
 - EF64 - Length 80mm

- Drill Guides Ø4.8mm
- EF65 - Length 40mm
 - EF66 - Length 60mm
 - EF67 - Length 80mm
 - EF68 - Length 140mm



- EF78
Allen Wrench 3mm



- EF79
Allen Wrench 5mm



INSTRUMENTATION



Drill Bit Ø2.9mm

EF91 - Length 140mm

Drill Bit Ø3.2mm

EF86 - Length 140mm


EF87 - Length 200mm

Drill Bit Ø4.8mm

EF88 - Length 180mm

EF89 - Length 240mm

EF90 - Length 280mm



Drill Bit Stop Units

EF735 - Ø2.9mm

EF82 - Ø3.2mm

EF83 - Ø4.8mm

INSTRUMENTATION



EF69
Hammer

EF70
Hammer Interchangeable Head

EF71
Hammer Locking Nut



EF116
Torque wrench 5mm



EF72
Tapered Trocar



EF108- T-wrench



EF12	L.R.S Instrument Sterilization Box, Empty
-------------	--

Codes	Set consisting of:	Units
EF12-T	Tray- L.R.S. Instrument Sterilization Box	3
EF78	Allen Wrench 3mm	2
EF80	Allen Wrench 6mm	2
EF109	Torque Wrench 6mm	1
EF108	T-Wrench	2
EF72	Tapered Trocar	1
EF116	Torque Wrench 5mm	1
EF95	Cannulated Drill Bit Ø3.2mm x Length 200mm	2
EF98	Cannulated Drill Bit, Ø4.8mm x Length 280mm	2
EF57	Ruland Pilot Wire Guide Ø2.0mm x Length 75mm	2
EF58	Ruland Pilot Wire Guide Ø2.0mm x Length 115mm	2
EF79	Allen Wrench 5mm	2
EF746	Allen Wrench 4mm	2
EF69	Hammer	1

Lower Tray


Codes	Set consisting of:	Units
EF65	Drill Guide Ø4.8mm Length 40mm	2
EF67	Drill Guide Ø4.8mm Length 80mm	2
EF63	Drill Guide Ø3.2mm Length 40mm	2
EF64	Drill Guide Ø3.2mm Length 80mm	2
EF68	Drill Guide Ø4.8mm Length 140mm	2
EF101	X-Wire without Olive Ø2mm x Length 400mm	5
EF89	Drill Bit Ø4.8mm x Length 240mm	2
EF87	Drill Bit Ø3.2mm x Length 200mm	2
EF90	Drill Bit Ø4.8mm x Length 280mm	2
EF83	Drill Bit Stop Unit Ø4.8mm	2
EF82	Drill Bit Stop Unit Ø3.2mm	2

Upper Tray

Codes	Set consisting of:	Units
EF45	Screw Guide, Length 60mm	10
EF47	Screw Guide, Length 100mm	10
EF53	Screw Guide, Length 160mm	10




PELVIC FIXATOR



EF157
Short Fixator with T-Clamp
*Consisting of: 2 Primary Links,
1 Pelvic Fixator Connecting Unit,
2 Cams, 2 Bushes and 2 T-Clamps.*



EF158
Long Fixator with Straight Clamp
*Consisting of: 2 Primary Links,
1 Pelvic Fixator Supplementary Link,
2 Pelvic Fixator Connecting Units,
2 Cams, 2 Bushes and 2 Straight Clamps.*



EF159
Long Fixator with Ball-Jointed Modules
for Independent Screw Placement
*Consisting of: 2 Primary Links,
1 Pelvic Fixator Supplementary Link
2 Pelvic Fixator Connecting Units,
2 Cams, 2 Bushes and 2 Ball-jointed
Modules for independent Screw Placement.*


BALL-JOINTED MODULES



EF202 - D.A.F. Straight Clamp



EF203 -T -Clamp



EF204

Ball-Jointed Modules for
Independent Screw Placement

For iliac crest applications.

*Consisting of 1 Bar with Ball-joint
and 2 independent Screw Clamps*

*Enables Screws to be placed freehand in
chosen positions on the iliac crest.*

COMPONENTS / REPLACEMENTS

EF237

Pelvic Fixator Long Body

Consisting of : 2 Primary Links,
1 Pelvic Fixator Supplementary Links,
2 Pelvic Fixator Connecting Units,
2 Cams and Bushes.

Maximum linear extension = 127mm.



EF238

Pelvic Fixator Short Body

Consisting of : 2 Primary Links,
1 Pelvic Fixator Connecting Unit,
2 Cams and Bushes.

Maximum linear extension = 77mm.



EF217
Bush



EF218
Cam



EF223
Pelvic Fixator
Connecting Unit



EF224
Primary Link



EF225
Pelvic Fixator Supplementary Link



EF226
Independent Screw Clamp

INSTRUMENTATION



EF41
Pelvic Screw Guide
Length 100mm



Screw Guides
EF45 - Length 60mm
EF47 - Length 100mm



Drill Guides Ø3.2mm
EF63 - Length 40mm
EF64 - Length 80mm



EF78
Allen Wrench 3mm



EF86
Drill Bit Ø3.2mm, Length 140mm



EF82
Drill Bit Stop Unit Ø3.2mm

INSTRUMENTATION



EF108
T-wrench



EF80
Allen Wrench 6 mm



EF72
Tapered Trocar



EF69
Hammer

EF70
Hammer Interchangeable Head


EF71
Hammer Locking Nut



EF109
Torque Wrench 6mm

EF13

Pelvic Fixator Sterilization Box, Empty



EF13	Pelvic Fixator Sterilization Box, Empty
-------------	--

Codes	Set consisting of:	Units
EF13-T	Tray for Pelvic Fixator Sterilization Box	2

Upper Tray


Codes	Set consisting of:	Units
EF755	Self-drilling Cortical Screw, Shaft Ø6mm, Thread Ø6-5mm, Total Length 220mm, Thread Length 50mm	6
EF754	Self-drilling Cortical Screw, Shaft Ø6mm, Thread Ø6-5mm, Total Length 150mm, Thread Length 50mm	6
EF45	Screw Guide, Length 60mm	4
EF47	Screw Guide, Length 100mm	4
EF41	Pelvic Screw Guide, Length 100mm	2
EF108	T-Wrench	1
EF72	Tapered Trocar	1
EF109	Torque Wrench 6mm	1
EF69	Hammer	1
EF80	Allen Wrench 6mm	2
EF86	Drill Bit Ø3.2mm x Length 140mm	2
EF63	Drill Guide Ø3.2mm x Length 40mm	2
	or	
EF64	Drill Guide Ø3.2mm x Length 80mm	2

Lower Tray

Codes	Set consisting of:	Units
EF237	Pelvic Fixator Long Body Consisting of: 2 Primary Links, 1 Pelvic Fixator Supplementary Links, 2 Pelvic Fixator Connecting Unit, 2 Cams and 2 Bushes. Maximum linear extension = 127mm	1
EF238	Pelvic Fixator Short Body Consisting of: 2 primary Links, 1 Pelvic Fixator Connecting Unit, 2 Cams and 2 Bushes Maximum linear extension = 77mm	1
EF204	Ball-Jointed Modules for Independent Screw Placement For iliac crest applications. Consisting of: 1 Bar with Ball-Joint and 2 Independent Screw Clamp, Enables Screws to be placed freehand in chosen positions on the iliac crest	2
EF226	Independent Screw Clamp	2
EF202	D.A.F. Straight Clamp	2
EF203	T-Clamp	2

The background of the page features a series of thin, light gray wavy lines that create a sense of movement and depth. These lines are concentrated in the top right and bottom left corners, leaving the central area relatively clear for the text.

ELBOW FIXATOR

**EF134**

Elbow Fixator with Straight Clamp

Consisting of: 2 Primary Links,
1 Elbow Fixator Central Connecting Unit,
2 Cams, 2 Bushes, 1 Elbow Fixator Small Distractor,
2 Straight Clamp, 1 Pair of Elbow Fixator
Compression-Distraction Unit Spacers,
1 Long Compression-Distraction Unit,
1 Allen Wrench 6mm and 1 Allen Wrench 3mm

**EF202**

D.A.F. Straight Clamp

**EF205**

Elbow Extended Range Clamp

COMPONENTS / REPLACEMENTS



EF230

Elbow Fixator Central Body

Consisting of : 2 Primary Links,
1 Elbow Fixator Central Connecting Unit,
1 Elbow Fixator Small Distraction,
2 Cams and 2 Bushes.



EF217

Bush



EF218

Cam



EF231

Elbow Fixator Central
Connecting Unit



EF232

Elbow Fixator Link
Locking Screw



EF233

Plate Locking Screw



Primary Link

EF224 - Length 130mm.



EF235

Elbow Fixator
Small Distractor



Short Primary Link

EF234- Length 100mm.



Compression-Distraction Unit

EF227 - Standard, Extends to 4cm

EF228 - Long, Extends to 8cm

EF229 - Short, Extends to 2.75cm



EF236

Elbow Fixator

Compression-Distraction Unit
Spacers (1 Pair)

INSTRUMENTATION



Screw Guide
EF45 -Length 60mm



Drill Guide Ø3.2mm
EF63 - Length 40mm

Drill Guide Ø4.8mm
EF65 - Length 40mm



Drill Bit Ø 3.2 mm
EF86 - Length 140 mm


Drill Bit Ø 4.8 mm
EF88 - Length 180 mm



EF78
Allen Wrench 3mm



EF108
T-Wrench



Drill Bit Stop Units
EF82 - Ø3.2mm
EF83 - Ø4.8mm



EF80
Allen Wrench 6mm

INSTRUMENTATION



EF116
Torque Wrench 5mm



X-Wire without Olive Ø2.0mm
EF100 - Length 150mm



EF72
Tapered Trocar

EF14**Elbow Fixator Sterilization Box, Empty**

Codes	Set consisting of:	Units
EF14-T	Tray for Elbow Fixator Sterilization Box, Empty	1
EF134	Elbow Fixator with Straight Clamp Consisting of: 2 Primary Links, 1 Elbow Fixator Central Connecting Unit, 2 Cams, 2 Bushes, 1 Elbow Fixator Small Distractor, 2 Straight Clamp, 1 Pair of Elbow Fixator Compression-Distractor Unit Spacers, 1 Long Compression-Distractor Unit, 1 Allen Wrench 6mm and 1 Allen Wrench 3mm	1
EF205	Elbow Extended Range Clamp	1
EF227	Compression-Distractor Unit, Standard, Extends to 4cm.	1
EF45	Screw Guide, Length 60mm	4
EF63	Drill Guide Ø3.2mm Length 40mm	2
EF65	Drill Guide Ø4.8mm Length 40mm	2
EF86	Drill Bit Ø3.2mm x Length 140mm	2
EF88	Drill Bit Ø4.8mm x Length 180mm	2
EF109	Torque Wrench 6mm	1
EF72	Tapered Trocar	1
EF108	T-Wrench	1
EF100	X-Wires without Olive Ø2.0mm x Length 150mm	4
EF749	Cortical Screw, Shaft Ø6mm, Thread Ø4.5-3.5mm, Total Length 120mm, Thread Length 20mm	4
EF747	Cortical Screw, Shaft Ø6mm, Thread Ø6-5mm, Total Length 110mm, Thread Length 30mm	4



2-022-000

Knee Fixator with Straight Clamp

Consisting of: 2 Primary Links,

1 Knee Fixator Central Connecting Unit,

2 Cams, 2 Bushes, 1 Knee Fixator Small Distractor,

2 Straight Clamp,

1 Pair of Knee Fixator Compression-Distraction Unit Spacers,

1 Long Compression-Distraction Unit,

1 Allen Wrench 6mm and 1 Allen Wrench 4mm





DYNAMIC WRIST FIXATOR



EF136

Dynamic Wrist Fixator

Consisting of: 1 Short Module with Straight Clamp, Cam and Bush,
1 Long Module with Straight Clamp, Cam and Bush and 1 Double Ball-joint.



EF137

Dynamic Wrist Fixator with

Compression-Distraction Module

Consisting of : 1 Long Module with Straight Clamp, Cam and Bush, 1 Compression-Distraction Module with Straight Clamp, Cam and Bush and 1 Double Ball-joint.



EF138

Dynamic Wrist Fixator for Extra-Articular Fractures and Corrective Osteotomies

Consisting of : 1 Short Module with T-Ciamp, Cam and Bush, 1 Compression-Distraction Module with Straight Clamp, Cam and Bush, and 1 Double Ball-joint.

COMPONENTS / REPLACEMENTS



EF271
Dynamic Wrist Fixator
Short Module
With Straight Clamp, Cam and Bush.



EF272
Dynamic Wrist Fixator
Short Module
With T-Clamp, Cam and Bush.



EF273
Dynamic Wrist Fixator
Long Module



EF274
Dynamic Wrist Fixator
Compression-Distraction Module




EF275
Dynamic Wrist Fixator
Straight Clamp




EF276
Dynamic Wrist Fixator
T-Clamp

COMPONENTS / REPLACEMENTS



EF277
Extra-Articular RadioLucent
Wrist Clamp (not illustrated)



EF282
Dynamic Wrist Fixator-Ulnar
Outrigger Clamp



EF283
Dynamic Wrist
Fixator - Double Ball-Joint



EF284
Dynamic
Wrist Fixator-Cam





EF285
Dynamic
Wrist Fixator - Bush



EF286
Dynamic Wrist Fixator -
Clamp Cover/ Anchoring Screw

INSTRUMENTATION



INSTRUMENTATION



EF72
Tapered Trocar




EF56
Pilot Wire Guide Ø 1.6mm



EF198
K-Wire Ø 1.5mm
Length 150mm.




EF111
T-wrench
For Clamp Screws.



EF112
T-wrench
For Bone Screws.

EF15	Penning Dynamic Wrist Fixator Sterilization Box, Empty
------	--




Codes	Set consisting of:	Units
EF15-T	Tray for Penning Dynamic Wrist Fixator Sterilization Box	1
EF136	Dynamic Wrist Fixator Consisting of: 1 Short Module with Straight Clamp, Cam and Bush, 1 Long Module with Straight Clamp, Cam and Bush and 1 Double Ball-Joint	2
EF138	Dynamic Wrist Fixator for Extra-Articular Fractures and Corrective Osteotomies Consisting of: 1 Short Module with T-Clamp, Cam and Bush, 1 Compression-Distracton Module with Straight Clamp, Cam and Bush, and 1 Double Ball-Joint	1
EF54	Template with Handle	2
EF283	Dynamic Wrist Fixator - Double Ball-Joint	2
EF285	Dynamic Wrist Fixator - Bush	4
EF284	Dynamic Wrist Fixator-Cam	4
EF92	Drill Bit Kit Ø2.7mm Length 127mm	2
EF59	Consisting of: 1 Drill Bit Stop Unit Ø2.9mm and 1 Allen Wrench 3mm Drill Guide Ø2.7mm	2
EF753	Penning Wrist Cortical Screw, Shaft Ø4mm, Thread Ø3.3-3mm, Total Length 80mm, Thread Length 35mm	5
EF738 EF112	Penning Wrist Cortical Screw, Shaft Ø4mm, Thread Ø3.3-3mm, Total Length 70mm, Thread Length 20mm	12
	T-Wrench, For Bone Screws	1
EF111	T-Wrench, For Clamp Screws	1
EF286	Penning Dynamic Wrist Fixator - Clamp Cover/ Anchoring Screw	4
EF72	Tapered Trocar	1
EF56	Pilot Wire Guide Ø1.6mm	2
EF276	Dynamic Wrist Fixator T-Clamp	1
EF282	Penning Dynamic Wrist Fixator-Ulnar Outrigger Clamp	1
EF277	Extra-Articular RadioLucent Wrist Clamp (not illustrated)	1





RADIOLUCENT WRIST FIXATOR

**EF139**

RadioLucent Wrist Fixator Kit, Sterile

Consisting of: 1 RadioLucent Wrist Fixator,

1 Compression-Distraction Unit,

1 Allen Wrench 4mm,

2 Drill Guide,

1 Drill Bit Ø2.7mm,

1 Screw Template,

1 T-Wrench,

2 Self-drilling Screws,

Shaft Ø4mm x Thread Ø3.3-3mm x Total Length 70mm x Thread Length 20mm,


2 Self-drilling Screws, Shaft Ø4mm x Thread Ø3.3-3mm x Total Length 80mm x

Thread Length 25mm and 4 Bone Screw Covers

The background of the page features a series of thin, light gray wavy lines that create a sense of movement and depth. These lines are concentrated in the top right and bottom left corners, while the center of the page is mostly clear, providing a space for the main text.

MINI FIXATOR


MINI FIXATOR



EF125

Mini Fixator Kits, Long


Consisting of : 1 Mini fixator Long Body,
2 Mini fixator Standard Clamps and
2 Mini fixator Compression-Distractive Nuts.



EF127

Mini Fixator Kits, Standard

Consisting of: 1 Mini fixator Standard Body,
2 Mini fixator Standard Clamps and
2 Mini fixator Compression-Distractive Nuts.




EF129

Mini Fixator Kits, Short

Consisting of: 1 Minifixator Short Body,
2 Minifixator Standard Clamps and
2 Minifixator Compression-Distractive Nuts.

COMPONENTS / REPLACEMENT



Mini Fixator Bodies

EF713- Long

1 Long Bar and 1 Standard Bar

Total Length 87mm

EF712 - Standard

2 Standard Bar Total Length 72mm

EF711 - Short

1 Standard Bar and 1 Short Bar

Total Length 62mm



EF714

Mini Fixator L-Clamp
(Left-Right)



EF715

Mini Fixator Standard Clamp
(1 pair)



EF716

Mini Fixator Compression-Distraction Nut



EF717

Mini Fixator Clamp Locking Screw



Pennig Mini Fixator Lengthening Bar Complete with Spacer

EF718 - Short, Length 80mm


EF719 - Standard, Length 100mm

EF720 - Long, Length 120mm



EF721

Minifixator Spacer




EF723

Threaded Wires Ø2mm

(Pack of 4) Total Length 100mm

Thread Length 15mm



EF722

Threaded Wires Ø1.6mm


(Pack of 4) Total Length 70mm

Thread Length 15mm

INSTRUMENTATION




EF78
Allen Wrench 3mm



EF40
Mini Fixator Reduction Forceps
(2 pair)




EF724
Threaded Wire Extractor



EF38
Wire Cutter

EF39
Wire Cutter Re-use Kit
Consisting of blades and 2 rubbers
(not illustrated).

EF16	Penning Mini Fixator Sterilization Box, Empty
------	---



Codes	Set consisting of:	Units
EF16-T	Tray for Pennig Mini Fixator sterilization Box	1
EF711	Mini Fixator Bodies Short (1 Standard Bar and 1 Short Bar, Total Length 62mm)	1
EF712	Mini Fixator Bodies Standard (2 Standard Bar, Total Length 72mm)	1
EF713	Mini Fixator Bodies Long (1 Long Bar and 1 Standard Bar; Total Length 87mm)	1
EF718	Mini Fixator Lengthening Bar Short Complete with Spacer, Length 80mm	1
EF719	Mini Fixator Lengthening Bar Standard Complete with Spacer, Length 100mm	1
EF720	Mini Fixator Lengthening Bar Long Complete with Spacer, Length 120mm	1
EF714	Mini Fixator L-Clamp (Left-Right)	1
EF715	Mini Fixator Standard Clamp (1 pair)	3
EF724	Threaded Wire Extractor	1
EF78	Allen Wrench 3mm	2
EF40	Mini Fixator Reduct on Forceps (2 pair)	1
EF723	Threaded Wires Ø2mm(Pack of 4) Total Length 100mm Thread Length 15mm	3
EF722	Threaded Wires Ø1.6mm(Pack of 4) Total Length 70mm Thread Length 15mm	2
EF716	Mini Fixator Compression-Distraction Nut	6
EF717	Mini Fixator Clamp Locking Screw	2



The background of the page features a series of thin, light gray wavy lines that create a sense of movement and depth. These lines are concentrated in the top right and bottom left corners, while the center of the page is mostly clear, providing a space for the main text.

MINIRAIL SYSTEM




EF126
MiniRail Lengtheners, Long
Allows Lengthening up to 135.5mm




EF128
MiniRail Lengtheners, Standard
Allows Lengthening up to 57.5mm



EF130
MiniRail Lengtheners, Short
Allows Lengthening up to 32.5mm



EF131
Articulated MiniRail Fixators, Vertical
Axis Allows articulation round an
axis vertical to the bone screws.



EF132
Articulated MiniRail Fixators, Horizontal
Axis Allows articulation round an axis
parallel to the bone screws.



EF154
Articulated MiniRail Fixators, Multiplanar
Axis Allows Articulation round an axis at
90 degrees to the bone screws.

COMPONENTS / REPLACEMENTS



EF708
MiniRail Lengthener
T-Clamp




EF710
MiniRail Lengthener/Fixator
Clamp Cover



EF709
MiniRail Lengthener/Fixator
Clamp Cover Screw

INSTRUMENTATION




EF78
Allen Wrench 3mm



EF111
T-Wrench
For Clamp screws




EF724
Threaded Wire Extractor



EF112
T-Wrench
For Bone screws.

EF17	MiniRail System Sterilization Box, Empty
------	--



Codes	Set consisting of:	Units
EF17-T	Tray for MiniRail System Sterilization Box	1
EF128	MiniRail Lengtheners, Standard, Allows Lengthening up to 57.5mm	1
EF126	MiniRail Lengtheners, Long, Allows Lengthening up to 135.5mm	1
EF130	MiniRail Lengtheners, Short, Allows Lengthening up to 32.5mm	1
EF132	Articulated MiniRail Fixators, Horizontal Axis Allows Articulation round an axis parallel to the bone Screws	1
EF131	Articulated MiniRail Fixators, Vertical Axis Allows Articulation round an axis vertical to the bone Screws	1
Ef112	T-Wrench, For Bone Screws	1
EF111	T-Wrench, For Clamp Screws	1
EF78	Allen Wrench 3mm	1
EF724	Threaded Wire Extractor	1
EF708	MiniRail Lengthener T-Clamp	2
EF751	Self-drilling Cortical Screws Shaft Ø3mm, Thread Ø2.5-2mm, Total Length 40mm, Thread Length 15mm	4
EF752	Self-drilling Cortical Screws Shaft Ø3mm, Thread Ø2.5-2mm, Total Length 45mm, Thread Length 20mm	4
EF722	Threaded Wires Ø1.6mm(Pack of 4) Total Length 70mm Thread Length 15mm	1
EF723	Threaded Wires Ø2mm(Pack of 4) Total Length 100mm Thread Length 15mm	1



BONE SCREWS

**Product Features:**

- **True 3D Fixation:** Provides multidirectional stability for enhanced construct strength
- **Simplified Surgical Technique:** Enables efficient and user-friendly implantation
- **Percutaneous Insertion Capability:** Allows pin insertion with minimal soft tissue disruption
- **U-Groove Design:** Increases bone-implant contact area to enhance fixation stability
- **Improved Fixation Security:** Reduces the risk of screw loosening under physiological loads

BONE SCREWS



Cortical Screw Shaft Ø6mm, Thread Ø6.0-5.0mm

Drill Bit Ø4.8mm

Thread Length 30mm

EF767 - Total Length 90mm
EF768 - Total Length 100mm
EF769 - Total Length 110mm
EF770 - Total Length 130mm
EF771 - Total Length 150mm

Thread Length 50mm

EF778 - Total Length 110mm
EF779 - Total Length 130mm
EF780 - Total Length 150mm
EF781 - Total Length 160mm
EF782 - Total Length 180mm
EF783 - Total Length 200mm

Thread Length 40mm

EF772 - Total Length 100mm
EF773 - Total Length 110mm
EF774 - Total Length 120mm
EF775 - Total Length 130mm
EF776 - Total Length 150mm
EF777 - Total Length 170mm

Thread Length 60mm

EF784 - Total Length 150mm
EF785 - Total Length 170mm
EF786 - Total Length 180mm
EF787 - Total Length 200mm

Thread Length 70mm

EF788 - Total Length 160mm

Thread Length 80mm

EF789 - Total Length 170mm

Thread Length 90mm

EF790 - Total Length 200mm



Cortical Screw Shaft Ø6mm, Thread Ø4.5-3.5mm

Drill Bit Ø3.2mm

Thread Length 20mm

EF791 - Total Length 70mm
EF792 - Total Length 80mm
EF793 - Total Length 100mm
EF794 - Total Length 120mm

Thread Length 30mm

EF795 - Total Length 80mm
EF796 - Total Length 100mm
EF797 - Total Length 120mm

Thread Length 40mm

EF798 - Total Length 100mm
EF799 - Total Length 120mm



Cortical Screw Shaft Ø6mm, Thread Ø3.5-3.2mm

Drill Bit Ø2.9mm

Thread Length 20mm

EF800 - Total Length 70mm

Thread Length 30mm

EF801 - Total Length 80mm



Wrist Cortical Screws Shaft Ø4mm, Thread Ø3.3-3.0mm

Drill Bit Ø2.7mm

Thread Length 20mm

EF802 - Total Length 70mm

Thread Length 35mm

EF803 - Total Length 80mm



Wrist Cortical Screws Shaft Ø4mm, Thread Ø3.0-2.5mm

Drill Bit Ø2.0mm

Thread Length 20mm

EF804 - Total Length 70mm

SELF-DRILLING BONE SCREWS



Self-Drilling Screws, Shaft Ø6mm, Thread Ø6.0-5.6mm

Thread Length 30mm

EF805 - Total Length 150mm
EF806 - Total Length 260mm

Thread Length 70mm

EF813 - Total Length 150mm
EF814 - Total Length 260mm

Thread Length 40mm

EF807 - Total Length 150mm
EF808 - Total Length 260mm

Thread Length 80mm

EF815 - Total Length 150mm
EF816 - Total Length 260mm

Thread Length 50mm

EF809 - Total Length 150mm
EF810 - Total Length 260mm

Thread Length 90mm

EF817 - Total Length 150mm
EF818 - Total Length 260mm

Thread Length 60mm

EF811 - Total Length 150mm
EF812 - Total Length 260mm

Self-Drilling Screws, Sterile Shaft Ø6mm, Thread Ø6.0-5.6mm

Thread Length 30mm

EF819 - Total Length 150mm
EF820 - Total Length 260mm

Thread Length 70mm

EF827 - Total Length 150mm
EF828 - Total Length 260mm

Thread Length 40mm

EF821 - Total Length 150mm
EF822 - Total Length 260mm

Thread Length 80mm

EF829 - Total Length 150mm
EF830 - Total Length 260mm

Thread Length 50mm

EF823 - Total Length 150mm
EF824 - Total Length 260mm

Thread Length 90mm

EF831 - Total Length 150mm
EF832 - Total Length 260mm

Thread Length 60mm

EF825 - Total Length 150mm
EF826 - Total Length 260mm



Self-Drilling Cortical Screws, Shaft Ø6mm, Thread Ø6.0-5.0mm

Thread Length 40mm

EF833 - Total Length 150mm
EF834 - Total Length 220mm

Thread Length 60mm

EF837 - Total Length 220mm

Thread Length 50mm

EF835 - Total Length 150mm
EF836 - Total Length 220mm

Thread Length 80mm

EF838 - Total Length 250mm

SELF-DRILLING BONE SCREWS



Self-Drilling Cortical Screws, Shaft Ø6mm, Thread Ø4.5-3.5mm

Thread Length 40mm EF839 - Total Length 120mm	Thread Length 60mm EF840 - Total Length 120mm
---	---



Wrist Self-Drilling Cortical Screws, Shaft Ø4mm, Thread Ø3.3-3.0mm

Drill Bit Ø2.0mm	
Thread Length 20mm EF841 - Total Length 60mm	Thread Length 30mm EF842 - Total Length 70mm EF843 - Total Length 100mm



Self-Drilling Cortical Screws, Shaft Ø3mm, Thread Ø3.0-2.5mm

For use with the M100 Range	
Thread Length 18mm EF844 - Total Length 50mm	Thread Length 25mm EF84 7 - Total Length 60mm EF848 - Total Length 70mm
Thread Length 20mm EF845 - Total Length 60mm EF846 - Total Length 70mm	Thread Length 30mm EF849 - Total Length 60mm EF850 - Total Length 70mm EF851 - Total Length 100mm



Self-Drilling Cortical Screws, Shaft Ø3mm, Thread Ø2.5-2.0mm

For use with the M100 Range	
Thread Length 15mm EF852 - Total Length 40mm	Thread Length 20mm EF853 - Total Length 45mm

BONE SCREWS (WITH HYDROXYAPATITE COATING)



Cortical Screws, Sterile Shaft Ø6mm, Thread Ø6.0-5.0mm

Drill Bit Ø4.8mm

Thread Length 30mm

EF854 - Total Length 90mm
 EF855 - Total Length 100mm
 EF856 - Total Length 110mm
 EF857 - Total Length 130mm
 EF858 - Total Length 150mm
 EF859 - Total Length 220mm

Thread Length 40mm

EF860 - Total Length 100mm
 EF861 - Total Length 110mm
 EF862 - Total Length 120mm
 EF863 - Total Length 130mm
 EF864 - Total Length 150mm
 EF865 - Total Length 170mm
 EF866 - Total Length 220mm

Thread Length 50mm

EF867 - Total Length 110mm
 EF868 - Total Length 130mm
 EF869 - Total Length 150mm
 EF870 - Total Length 160mm
 EF871 - Total Length 180mm
 EF872 - Total Length 200mm
 EF873 - Total Length 250mm

Thread Length 60mm

EF874 - Total Length 150mm
 EF875 - Total Length 170mm
 EF876 - Total Length 180mm
 EF877 - Total Length 200mm
 EF878 - Total Length 250mm

Thread Length 70mm

EF879 - Total Length 160mm

Thread Length 80mm

EF880 - Total Length 170mm

Thread Length 90mm

EF881 - Total Length 200mm

Cortical Screws, Sterile Shaft Ø4mm, Thread Ø3.3-3.0mm

Drill Bit Ø2.7mm

Thread Length 20mm

EF882 - Total Length 70mm

Thread Length 35mm

EF883 - Total Length 80mm

SELF-DRILLING BONE SCREWS (WITH HYDROXYAPATATITE COATING)



Self-Drilling Cortical Screws, Sterile Shaft Ø6mm, Thread Ø6.0-5.0mm

Thread Length 30mm
EF884 - Total Length 220mm
EF885 - Total Length 250mm

Thread Length 40mm
EF886 - Total Length 220mm

Thread Length 50mm
EF887 - Total Length 220mm

Thread Length 60mm
EF888 - Total Length 150mm

Thread Length 70mm
EF889 - Total Length 160mm

Thread Length 80mm
EF890 - Total Length 170mm

Thread Length 90mm
EF891 - Total Length 180mm



Bone Screws, Sterile Shaft Ø6mm, Thread Ø6.0-5.6mm

Thread Length 30mm
EF892 - Total Length 150mm
EF893 - Total Length 260mm

Thread Length 40mm
EF894 - Total Length 150mm
EF895 - Total Length 260mm

Thread Length 50mm
EF896 - Total Length 150mm
EF897 - Total Length 260mm


Thread Length 60mm
EF898 - Total Length 150mm
EF899 - Total Length 260mm

Thread Length 70mm
EF900 - Total Length 150mm
EF901 - Total Length 260mm

Thread Length 80mm
EF902 - Total Length 150mm
EF903 - Total Length 260mm

Thread Length 90mm
EF904 - Total Length 150mm
EF905 - Total Length 260mm





BONOFIX 5 EXTERNAL FIXATION SYSTEM

COMPONENTS



EF364 - Straight Post Ø5mm



EF366 - 30° Angled Post Ø5mm



Ø5mm Carbon Connecting Rod

EF528 - 65mm

EF529 - 100mm

EF530 - 150mm

EF531 - 200mm

EF532 - 250mm

EF533 - 300mm



Ø5mm Aluminium Connecting Rod

EF546 - 65mm

EF547 - 100mm

EF548 - 150mm

EF549 - 200mm

EF550 - 250mm

EF551 - 300mm



Ø3mm Self Drilling/ Self Tapping Pins

EF657 - Total Length 60mm

Thread Length 10mm

EF658 - Total Length 80mm

Thread Length 10mm

EF660 - Total Length 80mm

Thread Length 15mm

EF662 - Total Length 80mm

Thread Length 20mm

EF659 - Total Length 110mm

Thread Length 10mm

EF661 - Total Length 110mm

Thread Length 25mm



EF361

Compression-Distractor Tube Ø15mm,
Aluminium



COMPONENTS

Ø4mm Self Drilling/ Self Tapping Pins

EF663 - Total Length 90mm

Thread Length 10mm

EF667 - Total Length 90mm

Thread Length 20mm

EF671 - Total Length 90mm

Thread Length 30mm

EF664 - Total Length 120mm

Thread Length 30mm

EF668 - Total Length 120mm

Thread Length 35mm

EF665 - Total Length 150mm

Thread Length 40mm

EF669 - Total Length 150mm

Thread Length 50mm

EF666 - Total Length 180mm

Thread Length 40mm

EF670 - Total Length 180mm

Thread Length 50mm



INSTRUMENTATION



EF390
Stabilization/Reduction Wrench Ø5mm



EF391
Thumbwheel



EF392
5mm Wrench/Ø3-4mm Pin Driver



EF393
5mm Spanner Wrench



2-017-01
PEEK Rod to pin clamp, 5mm Rod
& 3-4mm pin



2-017-02
PEEK Rod to Rod clamp, 5mm



2-017-03
Elbow Joint device, 5mm




2-017-04
Wrist Joint device, 5mm



2-017-05
Pin to Rod clamp, mini, 3mm Rod
& 2mm pin



2-017-06
Rod to Rod clamp, mini, 3mm Rod



BONOFIX 8 EXTERNAL FIXATION SYSTEM

COMPONENTS



EF363
Straight Post Ø8mm



EF365
30° Angled Post Ø8mm



EF367
90° Angled Post Ø8mm



Ø8mm Carbon Connecting Rod

EF534 - 65mm
EF536 - 100mm
EF538 - 150mm
EF539 - 200mm
EF540 - 250mm
EF541 - 300mm
EF542 - 350mm
EF543 - 400mm
EF544 - 450mm
EF545 - 500mm



Ø8mm Aluminium Connecting Rod

EF552- 65mm
EF553 - 100mm
EF554 - 150mm
EF555 - 200mm
EF556 - 250mm
EF557 - 300mm
EF558 - 350mm
EF559 - 400mm
EF560 - 450mm
EF561 - 500mm




Ø8mm Stainless Steel Connecting Rod

EF501 - 65mm
EF502 - 100mm
EF503 - 150mm
EF504 - 200mm
EF505 - 250mm
EF506 - 300mm
EF507 - 350mm
EF508 - 400mm
EF509 - 450mm
EF510 - 500mm



COMPONENTS



Semi-Circular Curved Rod-Aluminium

EF595 - Small, 160mm (L)

EF596 - Medium, 180mm (L)

EF597 - Large, 200mm (L)



Ø5mm Self Drilling/Self Tapping Pin

EF672 - Total Length 120mm

Thread Length 30mm

EF677 - Total Length 120mm

Thread Length 35mm

EF673 - Total Length 150mm

Thread Length 40mm

EF678 - Total Length 150mm

Thread Length 50mm

EF674 - Total Length 180mm

Thread Length 35mm

EF679 - Total Length 180mm

Thread Length 50mm

EF682 - Total Length 180mm

Thread Length 60mm

EF675 - Total Length 200mm

Thread Length 50mm

EF680 - Total Length 200mm

Thread Length 60mm

EF676 - Total Length 250mm

Thread Length 50mm

EF681 - Total Length 250mm

Thread Length 70mm



Ø6mm Self Drilling/Self Tapping Pin

EF683 - Total Length 150mm

Thread Length 50mm

EF684 - Total Length 180mm

Thread Length 60mm

EF685 - Total Length 200mm

Thread Length 70mm

EF686 - Total Length 250mm

Thread Length 80mm

COMPONENTS



Cancellous Pins Ø6.0/5.0mm

EF638 - Total Length 120mm

Thread Length 35mm

EF639 - Total Length 150mm

Thread Length 40mm

EF642 - Total Length 250mm

Thread Length 50mm

EF640 - Total Length 180mm

Thread Length 60mm

EF641 - Total Length 250mm

Thread Length 70mm



Ø5mm Blunt Pins,

EF621 - Total Length 120mm

Thread Length 30mm

EF629 - Total Length 120mm

Thread Length 35mm

EF623 - Total Length 150mm

Thread Length 40mm

EF630 - Total Length 150mm

Thread Length 50mm

EF625 - Total Length 180mm

Thread Length 50mm

EF631 - Total Length 180mm

Thread Length 60mm

EF626 - Total Length 200mm

Thread Length 50mm

EF632 - Total Length 200mm

Thread Length 60mm

EF628 - Total Length 250mm

Thread Length 50mm

EF633 - Total Length 250mm

Thread Length 70mm



COMPONENTS



Ø6.0mm Blunt Pins

EF634 - Total Length 150mm

Thread Length 50mm

EF635 - Total Length 180mm

Thread Length 60mm

EF636 - Total Length 200mm

Thread Length 70mm

EF637 - Total Length 250mm

Thread Length 80mm



Trans

xing Pins Ø5.0/4.0mm

EF643 - Total Length 200mm

Thread Length 35mm

EF646 - Total Length 200mm

Thread Length 40mm

EF649 - Total Length 200mm

Thread Length 50mm

EF644 - Total Length 250mm

Thread Length 35mm

EF64 7 - Total Length 250mm

Thread Length 40mm

EF650 - Total Length 250mm

Thread Length 50mm

EF652 - Total Length 250mm

Thread Length 60mm

EF645 - Total Length 300mm

Thread Length 40mm

EF648 - Total Length 300mm

Thread Length 50mm


EF651 - Total Length 300mm

Thread Length 70mm

INSTRUMENTATION



EF912
Stabilization/Reduction Wrench Ø8mm



EF391
Thumbwheel



EF915
7mm Spanner Wrench



2-018-01
PEEK Rod to pin clamp, 8mm
Rod & 5-6mm pin



2-018-02
PEEK Rod to Rod clamp, 8mm



2-018-03
Ankle Joint device, 8mm




2-018-04
Knee Joint device, 8mm



EF359
Dynamination Tube Ø20mm

EF940	BonoFix 5 & 8 External Fixator Kit
-------	------------------------------------



EF940**BonoFix 5 & 8 External Fixator Kit**


Codes	Set consisting of:	Units
EF321	4-Hole Pin Clamp Ø3/4mm	2
EF327	Peri-Articular Pin Clamp Ø3mm	1
EF343	Rod-to-Rod Coupling,Rods or Posts of Ø5mm	4
EF335	Pin-to-Rod Coupling,Rods or Posts of Ø5mm, Pin of Ø3mm/Ø4mm	2
EF339	Tube-to-Rod Coupling,Tubes of Ø15mm, Rods or posts of Ø5mm	2
EF364	Straight Post Ø5mm	4
EF366	30° Angled Post Ø5mm	4
EF528	Ø5mm Carbon Connecting Rod 65mm	2
EF529	Ø5mm Carbon Connecting Rod 100mm	2
EF530	Ø5mm Carbon Connecting Rod 150mm	2
EF531	Ø5mm Carbon Connecting Rod 200mm	2
EF532	Ø5mm Carbon Connecting Rod 250mm	2
EF361	Compression-Distractor Tube Ø15mm	1
EF662	Ø3mm Self Drilling/ Self Tapping Pins, Total Length 80mm, Thread Length 15mm	8
EF664	Ø4mm Self Drilling/ Self Tapping Pins, Total Length 120mm, Thread Length 30mm	8
EF665	Ø4mm Self Drilling/ Self Tapping Pins, Total Length 150mm, Thread Length 40mm	8
EF666	Ø4mm Self Drilling/ Self Tapping Pins, Total Length 180mm, Thread Length 40mm	4
EF673	Ø5mm Self Drilling/Self Tapping Pin, Total Length 150mm, Thread Length 40mm	12
EF679	Ø5mm Self Drilling/Self Tapping Pin, Total Length 180mm, Thread Length 50mm	8
EF391	ThumbWheel	1
EF323	5-Hole Pin Clamp,Ø4mm,Ø5mm,Ø6mm	2
EF325	10-Hole Pin Clamp,Ø4mm,Ø5mm,Ø6mm	2
EF329	Pelvic Clamp, Ø4mm,Ø5mm,Ø6mm	2
EF341	Rod-to-Rod Coupling,Rods or Posts of Ø8mm	6
EF333	Pin-to-Rod Coupling,Rods or Posts of Ø8mm ;Pin of Ø4mm/Ø5mm	6
EF337	Tube-to-Rod Coupling,Tubes of Ø20mm;Rods or posts of Ø8mm	2
EF363	Straight Post Ø8mm	4
EF365	30° Angled Post Ø8mm	4
EF367	90° Angled Post Ø8mm	4
EF538	Ø8mm Carbon Connecting Rod 150mm	2
EF539	Ø8mm Carbon Connecting Rod 200mm	2
EF540	Ø8mm Carbon Connecting Rod 250mm	2
EF541	Ø8mm Carbon Connecting Rod 300mm	2
EF542	Ø8mm Carbon Connecting Rod 350mm	2
EF543	Ø8mm Carbon Connecting Rod 400mm	2
EF595	Small Semi-Circular Curved Rod-Aluminium Length 160mm	2
EF596	Medium Semi-Circular Curved Rod-Aluminium Length 180mm	2
EF597	Large Semi-Circular Curved Rod-Aluminium Length 200mm	2
EF359	Dynamination Tube Ø20mm	1
EF921	Spanner	1

**EF940****BonoFix 5 & 8 External Fixator Kit**


Codes	Set consisting of:	Units
EF922	Screw Guide(Small)	1
EF923	Screw Guide(Large)	1
EF924	Pressurizer	1
EF925	Hand Drill	1
EF926	Ø5mm Quick Coupling	1
EF927	Ø4mm Quick Coupling	1
EF928	Torque Wrench	1
EF912	Stabilization/Reduction Wrench Ø8mm	1
EF390	Stabilization/Reduction Wrench Ø5mm	1
EF931	Ø3.8mm Drill bit	1
EF932	Ø3.0mm Drill bit	1
EF934	Ø5mm T-Wrench	1
EF935	Ø7mm T-Wrench	1
EF936	Ø5mm T Quick Coupling	1
EF937	Ø4mm T Quick Coupling	1
EF938	Tray for BonoFix 5-8 External Fixator Kit	3
EF939	Container for BonoFix 5-8 External Fixator Kit	1

Sterilization Boxes

EF940	BonoFix 5 & 8 External Fixator Kit
-------	------------------------------------




Codes	Set consisting of:	Units
EF921	Spanner	1
EF922	Screw Guide(Small)	1
EF923	Screw Guide(Large)	1
EF924	Pressurizer	1
EF925	Hand Drill	1
EF926	Ø5mm Quick Coupling	1
EF927	Ø4mm Quick Coupling	1
EF928	Torque Wrench	1
EF929	Ø8mm Reduction Wrench	1
EF930	Ø5mm Reduction Wrench	1
EF931	Ø3.8mm Drill bit	1
EF932	Ø3.0mm Drill bit	1
EF391	Thumbwheel	1
EF934	Ø5mm T-Wrench	1
EF935	Ø7mm T-Wrench	1
EF936	Ø5mm T Quick Coupling	1
EF937	Ø4mm T Quick Coupling	1
EF938	Tray for BonoFix 5-8 External Fixator Kit	1
EF939	Container for BonoFix 5-8 External Fixator Kit	1

The background of the page features a series of thin, light gray wavy lines that create a sense of movement and depth. These lines are concentrated in the top right and bottom left corners, leaving the central area relatively clear for the text.

AO III EXTERNAL FIXATION SYSTEM

Distal Radius Fixator I



Distal Radius Fixator 1



COMPONENTS



Pin to Rod Coupling (Large)
 EF345 - Aluminium
 EF348 - Titanium
 Pin to Rod Coupling (Medium)
 EF346 - Aluminium
 EF349 - Titanium
 Pin to Rod Coupling (Small)
 EF347 - Aluminium
 EF350 - Titanium



Pin to Rod Coupling (Large)
 EF465 - Stainless Steel
 EF469- Titanium
 Pin to Rod Coupling (Medium)
 EF 471 - Stainless Steel
 EF 472 - Titanium
 Pin to Rod Coupling (Small)
 EF468- Stainless Steel
 EF470- Titanium



Tube To Tube Coupling (Large)
 EF464 - Stainless Steel
 EF473 - Titanium
 Tube To Tube Coupling {Medium)
 EF474 - Stainless Steel
 EF476 - Titanium
 Tube To Tube Coupling (Small)
 EF475 - Stainless Steel
 EF477 - Titanium



COMPONENTS



EF467
Pin To Tube Clamp



EF466
Cross Tube Clamp



EF351
Rod To Rod Coupling



EF352
Pin To Rod Coupling



EF535
Ø3mm Connecting Rod 75mm



EF537
Ø4mm Connecting Rod 120mm



COMPONENTS



Ø0.8mm Carbon Connecting Rods

EF534 - 65mm
EF536 - 100mm
EF538 - 150mm
EF539 - 200mm
EF540 - 250mm
EF541 - 300mm
EF542 - 350mm
EF543 - 400mm
EF544 - 450mm
EF545 - 500mm



Ø8mm Aluminium Connecting Rods

EF552 - 65mm
EF553 - 100mm
EF554 - 150mm
EF555 - 200mm
EF556 - 250mm
EF557 - 300mm
EF558 - 350mm
EF559 - 400mm
EF560 - 450mm
EF561 - 500mm



Ø8mm Stainless Steel Connecting Rods

EF501 - 65mm
EF502 - 100mm
EF503 - 150mm
EF504 - 200mm
EF505 - 250mm
EF506 - 300mm
EF507 - 350mm
EF508 - 400mm
EF509 - 450mm
EF510 - 500mm



COMPONENTS / REPLACEMENTS



Ø11mm Carbon Connecting Rods

EF562 - 65mm
 EF563 - 100mm
 EF564 - 150mm
 EF565 - 200mm
 EF566 - 250mm
 EF567 - 300mm
 EF568 - 350mm
 EF569 - 400mm
 EF570 - 450mm
 EF571 - 500mm



Ø11mm Aluminium Connecting Rods

EF578 - 65mm
 EF579 - 100mm
 EF580 - 150mm
 EF581 - 200mm
 EF582 - 250mm
 EF583 - 300mm
 EF584 - 350mm
 EF585 - 400mm
 EF586 - 450mm
 EF587 - 500mm



Ø11 mm Stainless Steel Connecting Rods

EF598 - 65mm
 EF599 - 100mm
 EF600 - 150mm
 EF601 - 200mm
 EF602 - 250mm
 EF603 - 300mm
 EF604 - 350mm
 EF605 - 400mm
 EF606 - 450mm
 EF607 - 500mm



Bone Pin

EF907 - 3.0mm
 EF908 - 4.0mm
 EF909 - 5.0mm
 EF910 - 6.0mm



K-Wire

EF911- Ø1.5mm




LINE FIX AND COLLES FIXATOR

LINE FIX


Fracture Region

Radius / Humerus / Tibia / Femur




LF01-001-130	Linefix set 13cm	
1609-085	Front Threaded Pins, Ø2.4mm x Thread 20mm x Length 80mm	8
1604-4.70-130	Connecting Rod, Aluminium, Ø4.70mm x Length 130mm	2
1603-3.10-110	Connecting Rod, Stainless Steel, Ø3.10mm x Length 110mm	2
BCM-001	Blue Clamp, Mini with Connecting Nut, Washer and Schanz Pin Fixation Bolt	12
1608-007	Wrench 7mm	1
LF01-001-150	Linefix set 15cm	
1609-030	Front Threaded Pins, Ø3.0mm, Thread Length 25mm, Length 100mm	8
1604-6.30-150	Connecting Rod, Aluminium, Ø6.30mm x Length 150mm	2
1603-3.10-150	Connecting Rod, Stainless Steel, Ø3.10mm x Length 150mm	2
BCS-001	Blue Clamp, Small with Connecting Nut, Washer and Schanz Pin Fixation Bolt	12
1608-010	Wrench 10mm	1
LF01-001-170	Linefix set 17cm	
1609-030	Front Threaded Pins, Ø3.0mm, Thread Length 25mm, Length 100mm	8
1604-6.30-170	Connecting Rod, Aluminium, Ø6.30mm x Length 170mm	2
1603-3.10-150	Connecting Rod, Stainless Steel, Ø3.10mm x Length 150mm	2
BCS-001	Blue Clamp, Small with Connecting Nut, Washer and Schanz Pin Fixation Bolt	12
1608-010	Wrench 10mm	1
LF01-001-200	Linefix set 20cm	
1609-086	Front Threaded Pins, Ø4.5mm, Thread Length 42mm, Length 180mm	8
1604-7.90-200	Connecting Rod, Aluminium, Ø7.90mm x Length 200mm	2
1603-4.75-200	Connecting Rod, Stainless Steel, Ø4.75mm x Length 200mm	2
BCS-001	Blue Clamp, Small with Connecting Nut, Washer and Schanz Pin Fixation Bolt	12
1608-010	Wrench 10mm	1
LF01-001-250	Linefix set 25cm	
1609-076	Front Threaded Pins, Ø5.0mm, Thread Length 42mm, Length 180mm	8
1604-7.90-250	Connecting Rod, Aluminium, Ø7.90mm x Length 250mm	2
1603-4.75-200	Connecting Rod, Stainless Steel, Ø4.75mm x Length 200mm	2
BCS-001	Blue Clamp, Small with Connecting Nut, Washer and Schanz Pin Fixation Bolt	12
1608-010	Wrench 10mm	1
LF01-001-300	Linefix set 30cm	
1609-076	Front Threaded Pins, Ø5.0mm, Thread Length 42mm, Length 180mm	8
1604-7.90-300	Connecting Rod, Aluminium, Ø7.90mm x Length 300mm	2
1603-4.75-200	Connecting Rod, Stainless Steel, Ø4.75mm x Length 200mm	2
BCS-001	Blue Clamp, Small with Connecting Nut, Washer and Schanz Pin Fixation Bolt	12
1608-010	Wrench 10mm	1
LF01-001-350	Linefix set 35cm	
1609-076	Front Threaded Pins, Ø5.0mm, Thread Length 42mm, Length 180mm	8
1604-7.90-350	Connecting Rod, Aluminium, Ø7.90mm x Length 350mm	2
1603-4.75-200	Connecting Rod, Stainless Steel, Ø4.75mm x Length 200mm	2
BCS-001	Blue Clamp, Small with Connecting Nut, Washer and Schanz Pin Fixation Bolt	12
1608-010	Wrench 10mm	1
LF01-001-400	Linefix set 40cm	
1609-076	Front Threaded Pins, Ø5.0mm, Thread Length 42mm, Length 180mm	8
1604-7.90-400	Connecting Rod, Aluminium, Ø7.90mm x Length 400mm	2
1603-4.75-200	Connecting Rod, Stainless Steel, Ø4.75mm x Length 200mm	2
BCS-001	Blue Clamp, Small with Connecting Nut, Washer and Schanz Pin Fixation Bolt	12
1608-010	Wrench 10mm	1

COLLES FIXATOR




CF01-001-145	Colles Fixatar, 145mm	
CF01.001	Colles Fixatar 145mm	1
1609-030-24	Front Threaded Pins, Ø3.0mm x Thread 24mm x Length 100mm	3
1609-029	Front Threaded Pins, Ø 2.5mm, Thread Length 20mm, Length 80mm	3
1608-014	Wrench 14mm	2
3766-004	Allen Key 4mm	1
CF0 1-001-165	Colles Fixatar 165mm	
CF01-002	Colles Fixatar 165mm	1
1609-030-24	Front Threaded Pins, Ø3.0mm x Thread 24mm x Length 100mm	3
1609-029	Front Threaded Pins, Ø 2.5mm, Thread Length 20mm, Length 80mm	3
1608-014	Wrench 14mm	2
3766-004	Allen Key 4mm	1



AO TYPE EXTERNAL FIXATION

FOREARM AND METACARPALS

1601-001M
Small Single Pin Clamp 4.0 x 2.5mm Modified



1602-000M
Small Connecting Clamp 4.0 x 4.0mm Modified



1603-000
Small Compression Distraction Device



Spanner



2163-03 7mm
2163-02 8mm

Front Threaded Pins




Codes	Dia(mm)	Length
1607-1.5-075	1.5	75mm
1607-2.0-075	2.0	75mm
1607-2.5-100	2.5	100mm

Front Threaded Pins



Codes	Description
1607-001 1607-002	Dia 3.5mm, Shaft 4.0mm, Length 125mm Thread Length- 14mm Thread Length- 20mm
1607-003 1607-004	Dia 3.5mm, Shaft 4.5mm, Length 125mm Thread Length- 18mm Thread Length- 22mm
1607-005	Cancellous Front Threaded Pins Dia 3.5mm, Shaft 4.0mm, Length 125mm

T-Wrench



2-023-000 7mm
1608-008 8mm

Carbon Fibre Rod - 4.0mm Dia.



Connecting Rod - 4.0mm Dia.



Codes	Length(mm)
1605-0100	100
1605-0150	150
1605-0200	200
1605-0250	250
1605-0300	300

Codes	Length(mm)
C1605-0100	100
C1605-0150	150
C1605-0200	200
C1605-0250	250
C1605-0300	300

COMPONENTS

1620-0005 Single Pin Clamp (Deluxe),
Stainless Steel



1620-0006 Double Pin Clamp,
Stainless Steel



1620-0008 Connecting Clamp,
Stainless Steel



1620-0009 Twin Adjustable Clamp
- Straight, Stainless Steel



1620-0011 Transverse Clamp, Stainless Steel



1620-0021 Self Holding Clamp,
Stainless Steel



1620-0007 Open Single Pin Clamp,
Stainless Steel



1620-0010 Twin Adjustable Clamp Curved,
Stainless Steel



1620-0012 Tube to Tube Clamp, Stainless Steel



1620-0018 Compression Distraction
Device, Stainless Steel



1620-005S Single Pin Clamp (Open Wise),
Stainless Steel

1620-0008a Open Connecting Clamp,
, Stainless Steel

INSTRUMENTATION

1620-0019 Schanz Pin Introducer,
Stainless Steel



2163-01 Spanner 11 mm, Stainless Steel



1620-0004 Triple Trocar Set with Handle
Extra Long, Stainless Steel



Tubular Rod 11mm Dia, Stainless Steel



Codes	Size
1620-0100	100mm(4")
1620-0150	150mm(6")
1620-0200	200mm(8")
1620-0250	250mm(10")
1620-0300	300mm(12")
1620-0350	350mm(14")
1620-0400	400mm(16")
1620-0450	450mm(18")

1620-0023 Box Spanner - 11 mm, Stainless Steel

1620-0020 Connecting Rod 5.0mm,
Stainless Steel



1608-011 T-Wrench 11mm, Stainless Steel



1620-0003 Triple Trocar Set with Handle,
Stainless Steel



1620-0002 Sleeve with Handle for 6.5mm
Cancellous Schanz Pin, Stainless Steel



Carbon Connecting Rod 11mm Dia,
Stainless Steel



Codes	Size
EF563	100mm
EF564	150mm
EF565	200mm
EF566	250mm
EF567	300mm
EF568	350mm
EF569	400mm
EF570	450mm
EF571	500mm

BONE PINS

Cancellous Front Threaded Pins 6.5mm



Codes	Thread Length	Length
1610-001	16mm	225mm
1610-002	32mm	225mm
1610-003	50mm	225mm

Front Threaded Pins (Shanz Screws)



Codes	Dia(mm)	Length(mm)	Thread(mm)
1609-001	4.0	125	18
1609-002	4.0	175	25
1609-003	4.5	150	18
1609-004	4.5	150	25
1609-005	4.5	175	18
1609-006	4.5	175	25
1609-007	4.5	175	32
1609-008	4.5	175	40
1609-009	4.5	200	25
1609-010	4.5	200	32
1609-011	4.5	200	40
1609-012	4.5	225	25
1609-013	4.5	225	32
1609-014	4.5	225	40
1609-015	5.0	150	18
1609-016	5.0	150	25
1609-017	5.0	175	18
1609-018	5.0	175	25
1609-019	5.0	175	32
1609-020	5.0	175	40
1609-021	5.0	200	25
1609-022	5.0	200	32
1609-023	5.0	200	40
1609-024	5.0	225	25
1609-025	5.0	225	32
1609-026	5.0	225	40
1609-027	6.0	225	32
1609-028	6.0	225	40

EF348
Pin to Rod Coupling (Large), Titanium



3/4 Ring

Codes	Length
EF961	115mm
EF962	140mm
EF963	165mm
EF964	205mm



EF965
Ring to Rod Clamp, Stainless Steel



EF966
Wire Ring Clamp, Stainless Steel




The background of the page features a series of thin, light gray wavy lines that create a sense of movement and depth. These lines are concentrated in the top right and bottom left corners, leaving the central area relatively clear for the text.

ILIZAROV RING FIXATORS

COMPONENTS


ILIZAROV- Full Ring

Codes	Size(mm)
1611-080	80
1611-100	100
1611-120	120
1611-140	140
1611-160	160
1611-180	180
1611-200	200
1611-220	220
1611-240	240




ILIZAROV - 1/2 Ring

Codes	Size(mm)
1612-080	80
1612-100	100
1612-120	120
1612-140	140
1612-160	160
1612-180	180
1612-200	200
1612-220	220
1612-240	240



ILIZAROV - 5/8 Ring

Codes	Size(mm)
1613-100	100
1613-120	120
1613-140	140
1613-160	160
1613-180	180
1613-200	200
1613-220	220
1613-240	240
1613-260	260



COMPONENTS

Omega Ring

Codes	Size(mm)
1617-120	120
1617-140	140
1617-160	160
1617-180	180
1617-200	200



Carbon Fiber Ring



ID(mm)	OD(mm)	Full Ring	3/4 Ring	5/8 Ring	Half Ring	3/8 Ring	1/4 Ring
80	120	FR1618-01	3/4- 1618-01	5/8- 1618-01	HF- 1618-01	3/8- 1618-01	1/4- 1618-01
100	140	FR1618-02	3/4- 1618-02	5/8- 1618-02	HF- 1618-02	3/8- 1618-02	1/4- 1618-02
110	150	FR1618-03	3/4- 1618-03	5/8- 1618-03	HF- 1618-03	3/8- 1618-03	1/4- 1618-03
120	160	FR1618-04	3/4- 1618-04	5/8- 1618-04	HF- 1618-04	3/8- 1618-04	1/4- 1618-04
130	170	FR1618-05	3/4- 1618-05	5/8- 1618-05	HF- 1618-05	3/8- 1618-05	1/4- 1618-05
140	180	FR1618-06	3/4- 1618-06	5/8- 1618-06	HF- 1618-06	3/8- 1618-06	1/4- 1618-06
150	190	FR1618-07	3/4- 1618-07	5/8- 1618-07	HF- 1618-07	3/8- 1618-07	1/4- 1618-07
160	200	FR1618-08	3/4- 1618-08	5/8- 1618-08	HF- 1618-08	3/8- 1618-08	1/4- 1618-08
180	220	FR1618-09	3/4- 1618-09	5/8- 1618-09	HF- 1618-09	3/8- 1618-09	1/4- 1618-09
200	240	FR1618-10	3/4- 1618-10	5/8- 1618-10	HF- 1618-10	3/8- 1618-10	1/4- 1618-10
220	260	FR1618-11	3/4- 1618-11	5/8- 1618-11	HF- 1618-11	3/8- 1618-11	1/4- 1618-11
240	280	x	x	5/8- 1618-12	HF- 1618-12	x	x

Italian Femoral Arche 90° Italian Femoral Arche 120°



Stainless Steel

90°		120°	
Codes	Size (mm)	Codes	Size (mm)
1621-100	100	1622-100	100
1621-120	120	1622-120	120
1621-140	140	1622-140	140
1621-160	160	1622-160	160
1621-180	180	1622-180	180
1621-200	200	1622-200	200
1621-220	220	1622-220	220
1621-240	240	1622-240	240
1621-260	260	1622-260	260



Carbon Fibre

90°		120°	
Codes	Size (mm)	Codes	Size (mm)
1621-100-C	100	1622-100-C	100
1621-120-C	120	1622-120-C	120
1621-140-C	140	1622-140-C	140
1621-160-C	160	1622-160-C	160
1621-180-C	180	1622-180-C	180
1621-200-C	200	1622-200-C	200
1621-220-C	220	1622-220-C	220
1621-240-C	240	1622-240-C	240
1621-260-C	260	1622-260-C	260



Threaded Rods



Codes	Size(mm)
1614-075	75
1614-100	100
1614-125	125
1614-150	150
1614-200	200
1614-250	250
1614-300	300
1614-350	350
1614-400	400
1614-450	450

Short Connection Plates



Codes	Holes
1625-02	2
1625-03	3
1625-04	4
1625-05	5
1625-06	6
1625-07	7
1625-08	8
1625-09	9
1625-10	10

Threaded Rods Slotted



Codes	Size(mm)
1615-075	75
1615-100	100
1615-125	125
1615-150	150

Graduated Telescopic Rod



Codes	Size(mm)
1623-100	100
1623-150	150

Carbon Fibre Connecting Plate

Codes	Holes
1616-02	2
1616-03	3
1616-04	4
1616-05	5
1616-06	6

Curved Plate



Codes	Holes
1633-03	3
1633-04	4
1633-05	5
1633-06	6

Twisted Plate



Codes	Holes
1632-02	2
1632-03	3
1632-04	4

Connection Plates with Threaded Ends



Codes	Holes
1626-05	5
1626-07	7
1626-09	9
1626-11	11

Post Male



Codes	Holes
1630-02	2
1630-03	3
1630-04	4

Post Female



Codes	Holes
1631-02	2
1631-03	3
1631-04	4

Telescopic Rod



Codes	Holes
1624-100	100
1624-150	150
1624-200	200
1624-250	250
1624-300	300

Threaded Sockets



Codes	Holes
1634-020	20
1634-030	30
1634-040	40
1634-060	60

1627-000 Hinge- Male



1628-000 Hinge - 90°



1629-000 Hinge - Female



Multi Pin Fixation Clamp



Codes	Slots
1640-000-2	2
1640-000-3	3
1640-000-4	4

1641-000 Single Pin Fixation Bolt



1642-000
Single Pin Fixation Bolt with Washer
- Deluxe



1644-000 Wire Fixation Bolt Cannulated



Ilizarov Wire Bayonet Point



Codes	Ø Dia	Length
1636-01	1.5mm	400
1636-02	1.8mm	400

Ilizarov Wire Trocar Point



Codes	Ø Dia	Length
1637-01	1.5mm	400
1637-02	1.8mm	400

Olive Wire Bayonet Point



Codes	Ø Dia	Length
1638-01	1.5mm	400
1638-02	1.8mm	400

Olive Wire Trocar Point




Codes	Ø Dia	Length
1639-01	1.5mm	400
1639-02	1.8mm	400

1645-000 Wire Fixation Bolt Slotted



1648-000 Wire Buckle



1647-000 Oblique Support with Thread



1646-000 Oblique Support without Thread



1649-000 Connical Washer (Pair)



Connection Bolt



Codes	Length(mm)
1650-000-10	10
1650-000-15	15
1650-000-20	20
1650-000-025	25
1650-000-030	30

1651-000 Nut- 10mm



1652-000 Square Nut



1653-000 Nylon Insert Nut



1654-000 Flat Sided Washer



1655-000 Slotted Washer



1656-000 Spacing Washer




1643-01 Wrenchocube One Hole



1643-02 Wrenchocube Two Hole



1643-03 Wrenchocube Three Hole**1643-04 Wrenchocube Four Hole****1656-000A Repox Device****1657-000 Washer Hanger****1658-000 Angulated Wrench 10mm Cannulated****1660-000 Wire Tensioner Dynamometer****Wire Tensioner**

Codes	Machanical
1660-001	1
7-045-08	2

Corticotome

Codes	Size(mm)
1659-005	5
1659-007	7
1659-010	10
1659-015	15

Illizarov Foot Ring, Large

Codes	Length(mm)
1698-220	120
1698-240	140
1698-260	160
1698-280	180
1698-300	200
1698-320	220

Codes	Discription
1663-000	Box Spanner 10mm
1664-000	Cannulated Bolt with Female Thread in Head
2163-04	Spanner 10mm
EF952	Universal Joint for Illizarov

Illizarov Foot Ring, Small

Codes	Length(mm)
1698-220S	120
1698-240S	140
1698-300S	200
1698-320S	220

Illizarov Foot Ring

Codes	Size(mm)	
2-019-01SS	140	Short
2-019-02SS	140	Long
2-020-01SS	160	Short
2-020-02SS	160	Long

SET COMPOSITION

2-016-000	Ilizarov Set- Small- 80mm
------------------	----------------------------------

Codes	Set consisting of:	Units
1612-08	ILIZAROV - 1/2 Ring, Length 80mm, Stainless Steel	8
1644-000	Wire Fixation Bolt, Cannulated	20
1650-000-10	Connection Bolt- 10mm	10
1624-100	Telescopic Rod, Length 100mm	3
1624-150	Telescopic Rod, Length 150mm	3
1651-000	Nut- 10mm	100
1614-100	Threaded Rod, Length 100mm	6
1614-150	Threaded Rod, Length 150mm	6
1614-200	Threaded Rod, Length 200mm	6
1625-03	Short Connection Plate, 3 Holes	4
1625-04	Short Connection Plate, 4 Holes	4
1630-02	Post Male, 2 Holes	2
1631-02	Post Female, 2 Holes	2
1656-000	Spacing Washer	20
1627-000	Hinge- Male	3
1629-000	Hinge - Female	3
1634-020	Threaded Socket, Length 20mm	3
1634-030	Threaded Socket, Length 30mm	3
1636-02	Ilizarov Wire Bayonet Point, Ø1.8mm, Length 400mm, Stainless Steel	10
1638-02	Olive Wire Bayonet Point, Ø1.8mm, Length 400mm, Stainless Steel	2
A10-036	Aluminium box for Ilizarov Set-Small	1



SET COMPOSITION

5005-000	Ilizarov Set- Small- 100mm
-----------------	-----------------------------------

Codes	Set consisting of:	Units
1612-100	ILIZAROV - 1/2 Ring, Length 100mm, Stainless Steel	8
1644-000	Wire Fixation Bolt, Cannulated	20
1650-000-10	Connection Bolt- 10mm	10
1624-100	Telescopic Rod, Length 100mm	3
1624-150	Telescopic Rod, Length 150mm	3
1651-000	Nut- 10mm	100
1614-100	Threaded Rod, Length 100mm	6
1614-150	Threaded Rod, Length 150mm	6
1614-200	Threaded Rod, Length 200mm	6
1625-03	Short Connection Plate, 3 Holes	4
1625-04	Short Connection Plate, 4 Holes	4
1630-02	Post Male, 2 Holes	2
1631-02	Post Female, 2 Holes	2
1656-000	Spacing Washer	20
1627-000	Hinge- Male	3
1629-000	Hinge - Female	3
1634-020	Threaded Socket, Length 20mm	3
1634-030	Threaded Socket, Length 30mm	3
1636-02	Ilizarov Wire Bayonet Point, Ø1.8mm, Length 400mm, Stainless Steel	10
1638-02	Olive Wire Bayonet Point, Ø1.8mm, Length 400mm, Stainless Steel	2
A10-036	Aluminium box for Ilizarov Set-Small	1

SET COMPOSITION

2-006-000	Ilizarov Set- Small- 120mm
------------------	-----------------------------------

Codes	Set consisting of:	Units
1612-12	ILIZAROV - 1/2 Ring, Length 120mm, Stainless Steel	8
1644-000	Wire Fixation Bolt, Cannulated	20
1650-000-10	Connection Bolt- 10mm	10
1624-100	Telescopic Rod, Length 100mm	3
1624-150	Telescopic Rod, Length 150mm	3
1651-000	Nut- 10mm	100
1614-100	Threaded Rod, Length 100mm	6
1614-150	Threaded Rod, Length 150mm	6
1614-200	Threaded Rod, Length 200mm	6
1625-03	Short Connection Plate, 3 Holes	4
1625-04	Short Connection Plate, 4 Holes	4
1630-02	Post Male, 2 Holes	2
1631-02	Post Female, 2 Holes	2
1656-000	Spacing Washer	20
1627-000	Hinge- Male	3
1629-000	Hinge - Female	3
1634-020	Threaded Socket, Length 20mm	3
1634-030	Threaded Socket, Length 30mm	3
1636-02	Ilizarov Wire Bayonet Point, Ø1.8mm, Length 400mm, Stainless Steel	10
1638-02	Olive Wire Bayonet Point, Ø1.8mm, Length 400mm, Stainless Steel	2
A10-036	Aluminium box for Ilizarov Set-Small	1

SET COMPOSITION

5006-000	Ilizarov Set- Medium- 140mm(Carbon Fiber Items)
-----------------	--

Codes	Set consisting of:	Units
1612-140	ILIZAROV - 1/2 Ring, Length 140mm, Stainless Steel	8
1644-000	Wire Fixation Bolt, Cannulated	20
1645-000	Wire Fixation Bolt Slotted	20
1650-000-10	Connection Bolt- 10mm	10
1624-100	Telescopic Rod, Length 100mm	3
1624-150	Telescopic Rod, Length 150mm	3
1651-000	Nut- 10mm	100
1614-100	Threaded Rod, Length 100mm	6
1614-150	Threaded Rod, Length 150mm	6
1614-200	Threaded Rod, Length 200mm	6
1625-03	Short Connection Plate, 3 Holes	4
1625-04	Short Connection Plate, 4 Holes	4
1625-05	Short Connection Plate, 5 Holes	2
1625-06	Short Connection Plate, 6 Holes	2
1625-07	Short Connection Plate, 7 Holes	2
1630-03	Post Male, 3 Holes	2
1631-03	Post Female, 3 Holes	2
1656-000	Spacing Washer	20
1627-000	Hinge- Male	3
1629-000	Hinge - Female	3
1634-030	Threaded Socket, Length 30mm	6
1636-02	Ilizarov Wire Bayonet Point, Ø1.8mm, Length 400mm, Stainless Steel	10
1638-02	Olive Wire Bayonet Point, Ø1.8mm, Length 400mm, Stainless Steel	2
A10-037	Aluminium box for Ilizarov Set-Medium	1

SET COMPOSITION

2-007-000	Ilizarov Set- Medium- 160mm
------------------	------------------------------------

Codes	Set consisting of:	Units
1612-160	ILIZAROV - 1/2 Ring, Length 160mm, Stainless Steel	8
1644-000	Wire Fixation Bolt, Cannulated	20
1645-000	Wire Fixation Bolt Slotted	20
1650-000-10	Connection Bolt- 10mm	10
1624-100	Telescopic Rod, Length 100mm	3
1624-150	Telescopic Rod, Length 150mm	3
1651-000	Nut- 10mm	100
1614-100	Threaded Rod, Length 100mm	6
1614-150	Threaded Rod, Length 150mm	6
1614-200	Threaded Rod, Length 200mm	6
1625-03	Short Connection Plate, 3 Holes	4
1625-04	Short Connection Plate, 4 Holes	4
1625-05	Short Connection Plate, 5 Holes	2
1625-06	Short Connection Plate, 6 Holes	2
1625-07	Short Connection Plate, 7 Holes	2
1630-03	Post Male, 3 Holes	2
1631-03	Post Female, 3 Holes	2
1656-000	Spacing Washer	20
1627-000	Hinge- Male	3
1629-000	Hinge - Female	3
1634-030	Threaded Socket, Length 30mm	6
1636-02	Ilizarov Wire Bayonet Point, Ø1.8mm, Length 400mm, Stainless Steel	10
1638-02	Olive Wire Bayonet Point, Ø1.8mm, Length 400mm, Stainless Steel	2
A10-037	Aluminium box for Ilizarov Set-Medium	1

SET COMPOSITION

2-008-000	Ilizarov Set- Medium- 180mm
------------------	------------------------------------

Codes	Set consisting of:	Units
1612-180	ILIZAROV - 1/2 Ring, Length 180mm, Stainless Steel	8
1644-000	Wire Fixation Bolt, Cannulated	20
1645-000	Wire Fixation Bolt Slotted	20
1650-000-10	Connection Bolt- 10mm	10
1624-100	Telescopic Rod, Length 100mm	3
1624-150	Telescopic Rod, Length 150mm	3
1651-000	Nut- 10mm	100
1614-100	Threaded Rod, Length 100mm	6
1614-150	Threaded Rod, Length 150mm	6
1614-200	Threaded Rod, Length 200mm	6
1625-03	Short Connection Plate, 3 Holes	4
1625-04	Short Connection Plate, 4 Holes	4
1625-05	Short Connection Plate, 5 Holes	2
1625-06	Short Connection Plate, 6 Holes	2
1625-07	Short Connection Plate, 7 Holes	2
1630-03	Post Male, 3 Holes	2
1631-03	Post Female, 3 Holes	2
1656-000	Spacing Washer	20
1627-000	Hinge- Male	3
1629-000	Hinge - Female	3
1634-030	Threaded Socket, Length 30mm	6
1636-02	Ilizarov Wire Bayonet Point, Ø1.8mm, Length 400mm, Stainless Steel	10
1638-02	Olive Wire Bayonet Point, Ø1.8mm, Length 400mm, Stainless Steel	2
A10-037	Aluminium box for Ilizarov Set-Medium	1

SET COMPOSITION

2-021-000	Ilizarov Set- Medium- 200mm
------------------	------------------------------------

Codes	Set consisting of:	Units
1612-200	ILIZAROV - 1/2 Ring, Length 200mm, Stainless Steel	8
1613-200	ILIZAROV - 5/8 Ring, Length 200mm, Stainless Steel	2
1622-120	Italian Femoral Arche 120°, 120mm, Stainless Steel	2
1644-000	Wire Fixation Bolt, Cannulated	20
1645-000	Wire Fixation Bolt Slotted	20
1650-000-10	Connection Bolt- 10mm	10
1624-100	Telescopic Rod, Length 100mm	3
1624-150	Telescopic Rod, Length 150mm	3
1624-200	Telescopic Rod, Length 200mm	3
1651-000	Nut- 10mm	100
1614-150	Threaded Rod, Length 150mm	6
1614-200	Threaded Rod, Length 200mm	6
1614-250	Threaded Rod, Length 250mm	6
1625-03	Short Connection Plate, 3 Holes	4
1625-04	Short Connection Plate, 4 Holes	4
1625-05	Short Connection Plate, 5 Holes	2
1625-06	Short Connection Plate, 6 Holes	2
1625-07	Short Connection Plate, 7 Holes	2
1630-04	Post Male, 4 Holes	2
1631-04	Post Female, 4 Holes	2
1656-000	Spacing Washer	20
1627-000	Hinge- Male	3
1629-000	Hinge - Female	3
1634-040	Threaded Socket, Length 40mm	6
1609-010	Front Threaded Pin (Shanz Screw), Ø4.5mm, Thread Length 32mm, Total Length 200mm, Stainless steel	4
1641-000	Single Pin Fixation Bolt	4
1647-000	Oblique Support with Thread	2
1636-02	Ilizarov Wire Bayonet Point, Ø1.8mm, Length 400mm, Stainless Steel	10
1638-02	Olive Wire Bayonet Point, Ø1.8mm, Length 400mm, Stainless Steel	10
A10-038	Aluminium box for Ilizarov Set-Large	1

SET COMPOSITION

2-009-000	Ilizarov Set- Large- 220mm
------------------	-----------------------------------

Codes	Set consisting of:	Units
1612-220	ILIZAROV - 1/2 Ring, Length 220mm, Stainless Steel	8
1613-220	ILIZAROV - 5/8 Ring, Length 220mm, Stainless Steel	2
1622-120	Italian Femoral Arche 120°, 120mm, Stainless Steel	2
1644-000	Wire Fixation Bolt, Cannulated	20
1645-000	Wire Fixation Bolt Slotted	20
1650-000-10	Connection Bolt- 10mm	10
1624-100	Telescopic Rod, Length 100mm	3
1624-150	Telescopic Rod, Length 150mm	3
1624-200	Telescopic Rod, Length 200mm	3
1651-000	Nut- 10mm	100
1614-150	Threaded Rod, Length 150mm	6
1614-200	Threaded Rod, Length 200mm	6
1614-250	Threaded Rod, Length 250mm	6
1625-03	Short Connection Plate, 3 Holes	4
1625-04	Short Connection Plate, 4 Holes	4
1625-05	Short Connection Plate, 5 Holes	2
1625-06	Short Connection Plate, 6 Holes	2
1625-07	Short Connection Plate, 7 Holes	2
1630-04	Post Male, 4 Holes	2
1631-04	Post Female, 4 Holes	2
1656-000	Spacing Washer	20
1627-000	Hinge- Male	3
1629-000	Hinge - Female	3
1634-040	Threaded Socket, Length 40mm	6
1609-010	Front Threaded Pin (Shanz Screw), Ø4.5mm, Thread Length 32mm, Total Length 200mm, Stainless steel	4
1641-000	Single Pin Fixation Bolt	4
1647-000	Oblique Support with Thread	2
1636-02	Ilizarov Wire Bayonet Point, Ø1.8mm, Length 400mm, Stainless Steel	10
1638-02	Olive Wire Bayonet Point, Ø1.8mm, Length 400mm, Stainless Steel	10
A10-038	Aluminium box for Ilizarov Set-Large	1

SET COMPOSITION

2-010-000	Ilizarov Set- Large- 240mm(Carbon Fiber Items)
------------------	---

Codes	Set consisting of:	Units
1612-240	ILIZAROV - 1/2 Ring, Length 240mm, Stainless Steel	8
1613-240	ILIZAROV - 5/8 Ring, Length 240mm, Stainless Steel	6
1622-120	Italian Femoral Arche 120°, 120mm, Stainless Steel	2
1644-000	Wire Fixation Bolt, Cannulated	20
1645-000	Wire Fixation Bolt Slotted	20
1650-000-10	Connection Bolt- 10mm	10
1624-100	Telescopic Rod, Length 100mm	3
1624-150	Telescopic Rod, Length 150mm	3
1624-200	Telescopic Rod, Length 200mm	3
1651-000	Nut- 10mm	100
1614-150	Threaded Rod, Length 150mm	6
1614-200	Threaded Rod, Length 200mm	6
1614-250	Threaded Rod, Length 250mm	6
1625-03	Short Connection Plate, 3 Holes	4
1625-04	Short Connection Plate, 4 Holes	4
1625-05	Short Connection Plate, 5 Holes	2
1625-06	Short Connection Plate, 6 Holes	2
1625-07	Short Connection Plate, 7 Holes	2
1630-04	Post Male, 4 Holes	2
1631-04	Post Female, 4 Holes	2
1656-000	Spacing Washer	20
1627-000	Hinge- Male	3
1629-000	Hinge - Female	3
1634-040	Threaded Socket, Length 40mm	6
1609-010	Front Threaded Pin (Shanz Screw), Ø4.5mm, Thread Length 32mm, Total Length 200mm, Stainless steel	4
1641-000	Single Pin Fixation Bolt	4
1647-000	Oblique Support with Thread	2
1636-02	Ilizarov Wire Bayonet Point, Ø1.8mm, Length 400mm, Stainless Steel	10
1638-02	Olive Wire Bayonet Point, Ø1.8mm, Length 400mm, Stainless Steel	10
A10-038	Aluminium box for Ilizarov Set-Large	1





DISTAL RADIUS FIXATOR


3723-000 Distal Radius Fixator



3724-000 Jig for Distal Radius Fixator



3727-000 Half Pin Introducer 4mm



3726-000 Trocar



Tapered Pins



Codes	Size(mm)
3725-01	4x50Lx18TL
3725-02	4x60Lx20TL
3725-03	4x70Lx20TL
3725-04	4x80Lx20TL
3725-05	4x100Lx20TL

3728-000 Drill Sleeve 4 x 2.5mm



**3730-000 Drill Bit Plain Shank/Jacob
Chuck End Ø2.5mm x Length 200mm**



3729-000 Screw Guide 6 x 4mm



3731-000 Allen Key Universal A/F 3mm



EXTERNAL FIXATOR SYSTEM

The content of this catalogue is intended for healthcare professionals only. All product information, illustrations, descriptions, and specifications provided here are for general informational purposes and are subject to change without notice. Auxein Medical Pvt. Ltd. reserves the right to improve or discontinue products at any time without prior notification.

This document does not constitute medical advice. Surgeons and medical practitioners must exercise their own professional judgment when using any medical devices or instruments, including those provided by Auxein.

CLINICAL USE

The implants and instruments in this catalogue are intended for use by qualified orthopedic trauma surgeons trained in fracture fixation and intramedullary nailing procedures. Improper selection of the implant, incorrect surgical technique, misuse, or insufficient surgical experience may result in implant failure or other complications.

REGULATORY COMPLIANCE CERTIFIED

ISO 13485:2016

CDSCO


MDSAP

DISCLAIMER: Not all products in this catalogue are covered by European CE, US FDA, or other certifications. Please contact our Sales Team for product-specific certification details.

INTELLECTUAL PROPERTY

All trademarks, trade names, brand names, product names, and logos mentioned in this catalogue are the property of Auxein Medical Pvt. Ltd. Unauthorized use or reproduction is strictly prohibited.

V1-2026



SCAN TO KNOW MORE

INDIA
AUXEIN MEDICAL PVT.LTD.

PLOT NO. 168-169-170, PHASE-4,
KUNDLI INDUSTRIAL AREA,
HSIIDC, SECTOR-57,
SONEPAT-131028, HARYANA
TEL: +91 99106 43638
FAX: +91 86077 70197
EMAIL: INFO@AUXEIN.COM

MÉXICO
AUXEIN MÉXICO, S.A DE CV

TEPIC 39 INT 801, COLONIA
ROMA SUR,
ALCADÍA CUAUHTÉMOC, CDMX
MÉXICO, C.P 06760
TEL: +521 55 7261 0318
EMAIL: INFO@AUXEIN.MX

USA
AUXEIN INC.

1500 NW 89TH COURT,
SUITE 107-108
DORAL, FLORIDA 33172
TEL. +1 305 395 6062
E FAX: +1 305 395 62 62
EMAIL: USOFFICE@AUXEIN.COM

DUBAI UAE
AUXEIN MEDICAL

TRADING FZCO OUD METHA
OFFICES BUILDING B BLOCK,
1ST FLOOR, OFFICE NO.26,
NEAR HEALTHCARE CITY,
PO BOX-25572
TEL. +971 438 71333
EMAIL: INFO@AUXEIN.COM

WWW.AUXEIN.COM

AMPL-CTG-06-01C Rev.01 Dated : 08-06-2026