

Auxein

CANNULATED SCREW SYSTEM



NEXT-LEVEL CARE

ABOUT US

Auxein Medical is an integrated, research based, orthopaedic implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavour to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in India and serves customers in over 75 countries worldwide.



IMPACTING HEALTHCARE GLOBALLY

Our Reach & Influence

500+

Employees Worldwide

75+

Operations in Countries

150+

Distributors Worldwide

20^{M+}

Patients Cured

200+

FDA Approved Products

250+

Machines Worldwide

3000+

Products Orthopedic Solutions

Auxein is dedicated to advancing medical technology through innovative orthopedic solutions, including **Trauma, Arthroscopy, and Reconstructive Procedures**, with a commitment to improving patient outcomes worldwide. Our approach integrates advanced research, precision engineering, and high-quality manufacturing to develop reliable implants and surgical instruments that address the evolving needs of modern orthopedic practice.


OUR ACHIEVEMENTS



OVER THE YEARS

Disclaimer: Not all products in this catalogue are covered by European CE, US FDA, or other certifications. Please contact our Sales Team for product-specific certification details.







At Auxein, it is our goal to provide the best in class quality and value to every individual. Consistency is a virtue we pursue with a passion. From the procurement of raw materials to the manufacturing of our range of 1200+ orthopaedic implant products, we employee world-class systems and technology at par with the highest standards of globally recognized evaluation systems. At Auxein, it is our goal to provide the best in class quality and value to every individual. Consistency is a virtue we pursue with a passion. From the procurement of raw materials to the manufacturing of our range of 1200+ orthopaedic implant products, we employee world-class systems and technology at par with the highest standards of globally recognized evaluation systems. The Cannulated Screw System is designed for precise fixation in orthopedic trauma procedures. Developed with input from global experts, the system allows guided screw placement over a guide wire, facilitating accurate positioning and stable fixation while supporting compression across fracture sites.

The background of the page features a series of thin, light gray wavy lines that create a sense of movement and depth. These lines are concentrated in the top right and bottom left corners, while the center of the page is mostly clear, providing a space for the main text.

INDEX



3.0mm Cannulated Cancellous Screw	10-19
3.5mm Cannulated Cortical Screw	20-33
4.0mm Cannulated Cancellous Screw	34-45
4.5mm Cannulated Cortical Screw	46-55
6.5mm Cannulated Screw	56-67
7.0mm Cannulated Cancellous Screw	68-81
7.3mm Cannulated Cancellous Screw	82-96



CANNULATED SCREW SYSTEM

3.0mm | 3.5mm | 4.0mm | 4.5mm | 6.5mm | 7.0mm | 7.3mm

Cannulated Screw System

3.0mm | 3.5mm | 4.0mm | 4.5mm | 6.5mm | 7.0mm | 7.3mm

Sharp Self-drilling and self-tapping Speed Tip thread design:

Sharp Self-Drilling and Self-Tapping Speed Tip Thread Design:

Functionally engineered for efficient cutting with immediate bone engagement. The triangular tip design allows simultaneous drilling and tapping during insertion, enabling controlled compression of bone tissue and enhancing pull-out strength and fixation stability.



The background of the page features a light gray, abstract pattern of thin, wavy lines that flow across the top and bottom edges, creating a sense of movement and depth. The central area is a solid, light gray color.

3.0MM CANNULATED CANCELLOUS SCREW



3.0mm Cannulated Cancellous Screw, Short Thread



Length	Titanium	Stainless Steel
10mm	TI-123.010	123.010
12mm	TI-123.012	123.012
14mm	TI-123.014	123.014
16mm	TI-123.016	123.016
18mm	TI-123.018	123.018
20mm	TI-123.020	123.020
22mm	TI-123.022	123.022
24mm	TI-123.024	123.024
26mm	TI-123.026	123.026
28mm	TI-123.028	123.028
30mm	TI-123.030	123.030
32mm	TI-123.032	123.032
34mm	TI-123.034	123.034
36mm	TI-123.036	123.036
38mm	TI-123.038	123.038
40mm	TI-123.040	123.040
42mm	TI-123.042	123.042
44mm	TI-123.044	123.044
46mm	TI-123.046	123.046
48mm	TI-123.048	123.048
50mm	TI-123.050	123.050




STERILE

Length	Titanium	Stainless Steel
10mm	TI-123.010-S	123.010-S
12mm	TI-123.012-S	123.012-S
14mm	TI-123.014-S	123.014-S
16mm	TI-123.016-S	123.016-S
18mm	TI-123.018-S	123.018-S
20mm	TI-123.020-S	123.020-S
22mm	TI-123.022-S	123.022-S
24mm	TI-123.024-S	123.024-S
26mm	TI-123.026-S	123.026-S
28mm	TI-123.028-S	123.028-S
30mm	TI-123.030-S	123.030-S
32mm	TI-123.032-S	123.032-S
34mm	TI-123.034-S	123.034-S
36mm	TI-123.036-S	123.036-S
38mm	TI-123.038-S	123.038-S
40mm	TI-123.040-S	123.040-S
42mm	TI-123.042-S	123.042-S
44mm	TI-123.044-S	123.044-S
46mm	TI-123.046-S	123.046-S
48mm	TI-123.048-S	123.048-S
50mm	TI-123.050-S	123.050-S

3.0mm Cannulated Cancellous Screw, Long Thread



Length	Titanium	Stainless Steel
14mm	TI-123.214	123.214
16mm	TI-123.216	123.216
18mm	TI-123.218	123.218
20mm	TI-123.220	123.220
22mm	TI-123.222	123.222
24mm	TI-123.224	123.224
26mm	TI-123.226	123.226
28mm	TI-123.228	123.228
30mm	TI-123.230	123.230
32mm	TI-123.232	123.232
34mm	TI-123.234	123.234
36mm	TI-123.236	123.236
38mm	TI-123.238	123.238
40mm	TI-123.240	123.240
42mm	TI-123.242	123.242
44mm	TI-123.244	123.244
46mm	TI-123.246	123.246
48mm	TI-123.248	123.248
50mm	TI-123.250	123.250



STERILE

Length	Titanium	Stainless Steel
14mm	TI-123.214-S	123.214-S
16mm	TI-123.216-S	123.216-S
18mm	TI-123.218-S	123.218-S
20mm	TI-123.220-S	123.220-S
22mm	TI-123.222-S	123.222-S
24mm	TI-123.224-S	123.224-S
26mm	TI-123.226-S	123.226-S
28mm	TI-123.228-S	123.228-S
30mm	TI-123.230-S	123.230-S
32mm	TI-123.232-S	123.232-S
34mm	TI-123.234-S	123.234-S
36mm	TI-123.236-S	123.236-S
38mm	TI-123.238-S	123.238-S
40mm	TI-123.240-S	123.240-S
42mm	TI-123.242-S	123.242-S
44mm	TI-123.244-S	123.244-S
46mm	TI-123.246-S	123.246-S
48mm	TI-123.248-S	123.248-S
50mm	TI-123.250-S	123.250-S




Washer for Small Fragment Screw




Titanium	Stainless Steel
TI-100.110	100.110

STERILE

Titanium	Stainless Steel
TI-100.110-S	100.110-S



The background of the page features a series of thin, light gray, wavy lines that flow from the top and bottom edges towards the center, creating a sense of movement and depth. These lines are more densely packed in some areas, particularly near the top and bottom, and more spread out in the middle where the text is located.

INSTRUMENTS



2190-2.0-150

Cannulated Drill Bit with Quick Coupling End, Ø2.0mm x Length 150mm



3767-1.1

Double Drill Guide, Ø2.0/1.1mm



3433-01

Cannulated Screwdriver, Hex 2.5mm



3446-000

Fiber Handle with Quick Coupling, Cannulated



3430-03

Cannulated Countersink with Quick Coupling for Ø3.0mm Screw



3436-09

Stylet Ø1.1mm






3426-03

Guide Wire, Ø1.1mm x Length 150mm




3438-17

Direct Measuring Device for Cannulated Screw Ø3.0mm



3425-000

Holding Clip for Washer



3415-000

Screw Holding Forcep



7-120-03

Tap for 3.0mm Cannulated Cancellous Screw



7-120-01

Instrument Tray for 3.0mm Cannulated Cancellous Screw Instrument Set





7-120-02

Container for 3.0mm Cannulated Cancellous Screw Instrument Set



ORDER INFORMATION

2335-000

3.0mm Cannulated Cancellous Screw Instrument Set



Code	Set Consisting Of	Units
2190-2.0-150	Cannulated Drill Bit with Quick Coupling End, Ø2.0mm x Length 150mm	2
3767-1.1	Double Drill Guide, Ø2.0/1.1mm	1
3433-01	Cannulated Screwdriver, Hex 2.5mm	1
3446-000	Fiber Handle with Quick Coupling, Cannulated	1
3430-03	Cannulated Countersink with Quick Coupling for Ø3.0mm Screw	1
3436-09	Stylet Ø1.1mm	1
3426-03	Guide Wire, Ø1.1mm x Length 150mm	10
3438-17	Direct Measuring Device for Cannulated Screw Ø3.0mm	1
3425-000	Holding Clip for Washer	1
3415-000	Screw Holding Forcep	1
7-120-03	Tap for 3.0mm Cannulated Cancellous Screw	1
7-120-01	Instrument Tray for 3.0mm Cannulated Cancellous Screw Instrument Set	1
7-120-02	Container for 3.0mm Cannulated Cancellous Screw Instrument Set	1

The background of the page features a light gray, abstract pattern of thin, wavy lines that flow across the top and bottom edges, creating a sense of movement and depth. The central area is a solid, light gray color.

3.5MM CANNULATED CORTICAL SCREW

3.5mm Cannulated Cortical Screw, Short Thread



Length	Titanium	Stainless Steel
10mm	TI-121.010	121.010
12mm	TI-121.012	121.012
14mm	TI-121.014	121.014
16mm	TI-121.016	121.016
18mm	TI-121.018	121.018
20mm	TI-121.020	121.020
22mm	TI-121.022	121.022
24mm	TI-121.024	121.024
26mm	TI-121.026	121.026
28mm	TI-121.028	121.028
30mm	TI-121.030	121.030
32mm	TI-121.032	121.032
34mm	TI-121.034	121.034
36mm	TI-121.036	121.036
38mm	TI-121.038	121.038
40mm	TI-121.040	121.040
42mm	TI-121.042	121.042
44mm	TI-121.044	121.044
46mm	TI-121.046	121.046
48mm	TI-121.048	121.048
50mm	TI-121.050	121.050




STERILE

Length	Titanium	Stainless Steel
10mm	TI-121.010-S	121.010-S
12mm	TI-121.012-S	121.012-S
14mm	TI-121.014-S	121.014-S
16mm	TI-121.016-S	121.016-S
18mm	TI-121.018-S	121.018-S
20mm	TI-121.020-S	121.020-S
22mm	TI-121.022-S	121.022-S
24mm	TI-121.024-S	121.024-S
26mm	TI-121.026-S	121.026-S
28mm	TI-121.028-S	121.028-S
30mm	TI-121.030-S	121.030-S
32mm	TI-121.032-S	121.032-S
34mm	TI-121.034-S	121.034-S
36mm	TI-121.036-S	121.036-S
38mm	TI-121.038-S	121.038-S
40mm	TI-121.040-S	121.040-S
42mm	TI-121.042-S	121.042-S
44mm	TI-121.044-S	121.044-S
46mm	TI-121.046-S	121.046-S
48mm	TI-121.048-S	121.048-S
50mm	TI-121.050-S	121.050-S

3.5mm Cannulated Cortical Screw, Full Thread



Length	Titanium	Stainless Steel
10mm	TI-121.210	121.210
12mm	TI-121.212	121.212
14mm	TI-121.214	121.214
16mm	TI-121.216	121.216
18mm	TI-121.218	121.218
20mm	TI-121.220	121.220
22mm	TI-121.222	121.222
24mm	TI-121.224	121.224
26mm	TI-121.226	121.226
28mm	TI-121.228	121.228
30mm	TI-121.230	121.230
32mm	TI-121.232	121.232
34mm	TI-121.234	121.234
36mm	TI-121.236	121.236
38mm	TI-121.238	121.238
40mm	TI-121.240	121.240
42mm	TI-121.242	121.242
44mm	TI-121.244	121.244
46mm	TI-121.246	121.246
48mm	TI-121.248	121.248
50mm	TI-121.250	121.250



STERILE

Length	Titanium	Stainless Steel
10mm	TI-121.210-S	121.210-S
12mm	TI-121.212-S	121.212-S
14mm	TI-121.214-S	121.214-S
16mm	TI-121.216-S	121.216-S
18mm	TI-121.218-S	121.218-S
20mm	TI-121.220-S	121.220-S
22mm	TI-121.222-S	121.222-S
24mm	TI-121.224-S	121.224-S
26mm	TI-121.226-S	121.226-S
28mm	TI-121.228-S	121.228-S
30mm	TI-121.230-S	121.230-S
32mm	TI-121.232-S	121.232-S
34mm	TI-121.234-S	121.234-S
36mm	TI-121.236-S	121.236-S
38mm	TI-121.238-S	121.238-S
40mm	TI-121.240-S	121.240-S
42mm	TI-121.242-S	121.242-S
44mm	TI-121.244-S	121.244-S
46mm	TI-121.246-S	121.246-S
48mm	TI-121.248-S	121.248-S
50mm	TI-121.250-S	121.250-S

3.5mm Cannulated Cortical Screw, Short Thread, Self-Tapping



Length	Titanium	Stainless Steel
10mm	TI-121.310	121.310
12mm	TI-121.312	121.312
14mm	TI-121.314	121.314
16mm	TI-121.316	121.316
18mm	TI-121.318	121.318
20mm	TI-121.320	121.320
22mm	TI-121.322	121.322
24mm	TI-121.324	121.324
26mm	TI-121.326	121.326
28mm	TI-121.328	121.328
30mm	TI-121.330	121.330
32mm	TI-121.332	121.332
34mm	TI-121.334	121.334
36mm	TI-121.336	121.336
38mm	TI-121.338	121.338
40mm	TI-121.340	121.340
42mm	TI-121.342	121.342
44mm	TI-121.344	121.344
46mm	TI-121.346	121.346
48mm	TI-121.348	121.348
50mm	TI-121.350	121.350




STERILE

Length	Titanium	Stainless Steel
10mm	TI-121.310-S	121.310-S
12mm	TI-121.312-S	121.312-S
14mm	TI-121.314-S	121.314-S
16mm	TI-121.316-S	121.316-S
18mm	TI-121.318-S	121.318-S
20mm	TI-121.320-S	121.320-S
22mm	TI-121.322-S	121.322-S
24mm	TI-121.324-S	121.324-S
26mm	TI-121.326-S	121.326-S
28mm	TI-121.328-S	121.328-S
30mm	TI-121.330-S	121.330-S
32mm	TI-121.332-S	121.332-S
34mm	TI-121.334-S	121.334-S
36mm	TI-121.336-S	121.336-S
38mm	TI-121.338-S	121.338-S
40mm	TI-121.340-S	121.340-S
42mm	TI-121.342-S	121.342-S
44mm	TI-121.344-S	121.344-S
46mm	TI-121.346-S	121.346-S
48mm	TI-121.348-S	121.348-S
50mm	TI-121.350-S	121.350-S

3.5mm Cannulated Cortical Screw, Full Thread, Self-Tapping

Length	Titanium	Stainless Steel
10mm	TI-121.410	121.410
12mm	TI-121.412	121.412
14mm	TI-121.414	121.414
16mm	TI-121.416	121.416
18mm	TI-121.418	121.418
20mm	TI-121.420	121.420
22mm	TI-121.422	121.422
24mm	TI-121.424	121.424
26mm	TI-121.426	121.426
28mm	TI-121.428	121.428
30mm	TI-121.430	121.430
32mm	TI-121.432	121.432
34mm	TI-121.434	121.434
36mm	TI-121.436	121.436
38mm	TI-121.438	121.438
40mm	TI-121.440	121.440
42mm	TI-121.442	121.442
44mm	TI-121.444	121.444
46mm	TI-121.446	121.446
48mm	TI-121.448	121.448
50mm	TI-121.450	121.450



Length	Titanium	Stainless Steel
10mm	TI-121.410-S	121.410-S
12mm	TI-121.412-S	121.412-S
14mm	TI-121.414-S	121.414-S
16mm	TI-121.416-S	121.416-S
18mm	TI-121.418-S	121.418-S
20mm	TI-121.420-S	121.420-S
22mm	TI-121.422-S	121.422-S
24mm	TI-121.424-S	121.424-S
26mm	TI-121.426-S	121.426-S
28mm	TI-121.428-S	121.428-S
30mm	TI-121.430-S	121.430-S
32mm	TI-121.432-S	121.432-S
34mm	TI-121.434-S	121.434-S
36mm	TI-121.436-S	121.436-S
38mm	TI-121.438-S	121.438-S
40mm	TI-121.440-S	121.440-S
42mm	TI-121.442-S	121.442-S
44mm	TI-121.444-S	121.444-S
46mm	TI-121.446-S	121.446-S
48mm	TI-121.448-S	121.448-S
50mm	TI-121.450-S	121.450-S


Washer for Small Fragment Screw




Titanium	Stainless Steel
TI-100.110	100.110

STERILE

Titanium	Stainless Steel
TI-100.110-S	100.110-S



The background of the page features a series of thin, light gray, wavy lines that flow across the top and bottom edges, creating a sense of movement and depth. The central area is a solid, light gray color.

INSTRUMENTS

7-121-03

Cannulated Drill Bit Quick Coupling, Ø2.7mm x Length 160mm



3435-01

Stopper for Drill Bit Ø2.7mm



3435-02

Cannulated Drill Bit with Quick Coupling, Ø3.5mm x Length 160mm



3767-2.7

Double Drill Guide, Ø2.7/1.2 mm



3767-3.5

Double Drill Guide, Ø2.5/3.5mm



3434-000

Double Drill Guide, Ø3.5/2.7mm




3433-01

Cannulated Screwdriver, Hex 2.5mm



3446-000

Fiber Handle With Quick Coupling, Cannulated



3430-04

Cannulated Countersink With Quick Coupling For Ø3.5mm Screw



1485-1.2

Stylet, Ø1.2mm




3426-01

Guide Wire With Threaded Trocar Tip, Ø1.2mm X Thread Length 5mm X Length 150mm



3438-19

Direct Measuring Device For Ø3.5mm Screw




3415-000

Screw Holding Forcep



3425-000

Holding Clip For Washer



7-056-02

Wire Sleeve for Guide Wire Ø2.7/1.2mm




1489-01

Holding Sleeve for Cannulated Screwdriver, Hex 2.5mm




7-121-02

Bone Tap for 3.5mm Cannulated Cortical Screw



7-121-01

Instrument Tray for 3.0mm Cannulated Cancellous Screw Instrument Set





A10-045

Container for 3.5mm Cannulated Cortical Instruments Set





ORDER INFORMATION

2336-000**3.5mm Cannulated Cortical Screw Instrument Set**

Code	Set Consisting Of	Units
7-121-03	Cannulated Drill Bit Quick Coupling, Ø2.7mm x Length 160mm	2
3435-01	Stopper for Drill Bit Ø2.7mm	1
3435-02	Cannulated Drill Bit with Quick Coupling, Ø3.5mm x Length 160mm	2
3767-2.7	Double Drill Guide, Ø2.7/1.2 mm	1
3767-3.5	Double Drill Guide, Ø2.5/3.5mm	1
3434-000	Double Drill Guide, Ø3.5/2.7mm	1
3433-01	Cannulated Screwdriver, Hex 2.5mm	1
3446-000	Fiber Handle with Quick coupling, Cannulated	1
3430-04	Cannulated Countersink with Quick Coupling for Ø3.5mm Screw	1
1485-1.2	Stylet, Ø1.2mm	1
3426-01	Guide Wire with Threaded Trocar Tip, Ø1.2mm x Thread Length 5mm x Length 150mm	6
3438-19	Direct Measuring Device for Ø3.5mm Screw	1
3415-000	Screw Holding Forcep	1
3425-000	Holding Clip for Washer	1
7-056-02	Wire Sleeve for Guide Wire Ø2.7/1.2mm	1
1489-01	Holding Sleeve for Cannulated Screwdriver, Hex 2.5mm	1
7-121-02	Bone Tap for 3.5mm Cannulated Cortical Screw	1
7-121-01	Instrument Tray for 3.5mm Cannulated Cortical Screw Instrument Set	1
A10-045	Container for 3.5mm Cannulated Cortical Instruments Set	1

The background of the page features a series of thin, light gray wavy lines that flow from the top left towards the bottom right, creating a sense of movement and depth. These lines are more densely packed in some areas, particularly near the top and bottom edges, and more spread out in the center where the text is located.

4.0MM CANNULATED CANCELLOUS SCREW



4.0mm Cannulated Cancellous Screw, Short Thread, Self-Tapping



STERILE

	Length	Titanium	Stainless Steel
*	10mm	TI-124.210	124.210
*	12mm	TI-124.212	124.212
*	14mm	TI-124.214	124.214
	16mm	TI-124.216	124.216
	18mm	TI-124.218	124.218
	20mm	TI-124.220	124.220
	22mm	TI-124.222	124.222
	24mm	TI-124.224	124.224
	26mm	TI-124.226	124.226
	28mm	TI-124.228	124.228
	30mm	TI-124.230	124.230
	32mm	TI-124.232	124.232
	34mm	TI-124.234	124.234
	36mm	TI-124.236	124.236
	38mm	TI-124.238	124.238
	40mm	TI-124.240	124.240
	42mm	TI-124.242	124.242
	44mm	TI-124.244	124.244
	46mm	TI-124.246	124.246
	48mm	TI-124.248	124.248
	50mm	TI-124.250	124.250
	55mm	TI-124.255	124.255
	60mm	TI-124.260	124.260
*	65mm	TI-124.265	124.265
*	70mm	TI-124.270	124.270



	Length	Titanium	Stainless Steel
*	10mm	TI-124.210-S	124.210-S
*	12mm	TI-124.212-S	124.212-S
*	14mm	TI-124.214-S	124.214-S
	16mm	TI-124.216-S	124.216-S
	18mm	TI-124.218-S	124.218-S
	20mm	TI-124.220-S	124.220-S
	22mm	TI-124.222-S	124.222-S
	24mm	TI-124.224-S	124.224-S
	26mm	TI-124.226-S	124.226-S
	28mm	TI-124.228-S	124.228-S
	30mm	TI-124.230-S	124.230-S
	32mm	TI-124.232-S	124.232-S
	34mm	TI-124.234-S	124.234-S
	36mm	TI-124.236-S	124.236-S
	38mm	TI-124.238-S	124.238-S
	40mm	TI-124.240-S	124.240-S
	42mm	TI-124.242-S	124.242-S
	44mm	TI-124.244-S	124.244-S
	46mm	TI-124.246-S	124.246-S
	48mm	TI-124.248-S	124.248-S
	50mm	TI-124.250-S	124.250-S
	55mm	TI-124.255-S	124.255-S
	60mm	TI-124.260-S	124.260-S
*	65mm	TI-124.265-S	124.265-S
*	70mm	TI-124.270-S	124.270-S

* Sizes not available in Screw Caddy



4.0mm Cannulated Cancellous Screw, Full Thread, Self-Tapping



	Length	Titanium	Stainless Steel
*	10mm	TI-125.210	125.210
*	12mm	TI-125.212	125.212
*	14mm	TI-125.214	125.214
	16mm	TI-125.216	125.216
	18mm	TI-125.218	125.218
	20mm	TI-125.220	125.220
	22mm	TI-125.222	125.222
	24mm	TI-125.224	125.224
	26mm	TI-125.226	125.226
	28mm	TI-125.228	125.228
	30mm	TI-125.230	125.230
	32mm	TI-125.232	125.232
	34mm	TI-125.234	125.234
	36mm	TI-125.236	125.236
	38mm	TI-125.238	125.238
	40mm	TI-125.240	125.240
	42mm	TI-125.242	125.242
	44mm	TI-125.244	125.244
	46mm	TI-125.246	125.246
	48mm	TI-125.248	125.248
	50mm	TI-125.250	125.250
	55mm	TI-125.255	125.255
	60mm	TI-125.260	125.260
*	65mm	TI-125.265	125.265
*	70mm	TI-125.270	125.270



STERILE


	Length	Titanium	Stainless Steel
*	10mm	TI-125.210-S	125.210-S
*	12mm	TI-125.212-S	125.212-S
*	14mm	TI-125.214-S	125.214-S
	16mm	TI-125.216-S	125.216-S
	18mm	TI-125.218-S	125.218-S
	20mm	TI-125.220-S	125.220-S
	22mm	TI-125.222-S	125.222-S
	24mm	TI-125.224-S	125.224-S
	26mm	TI-125.226-S	125.226-S
	28mm	TI-125.228-S	125.228-S
	30mm	TI-125.230-S	125.230-S
	32mm	TI-125.232-S	125.232-S
	34mm	TI-125.234-S	125.234-S
	36mm	TI-125.236-S	125.236-S
	38mm	TI-125.238-S	125.238-S
	40mm	TI-125.240-S	125.240-S
	42mm	TI-125.242-S	125.242-S
	44mm	TI-125.244-S	125.244-S
	46mm	TI-125.246-S	125.246-S
	48mm	TI-125.248-S	125.248-S
	50mm	TI-125.250-S	125.250-S
	55mm	TI-125.255-S	125.255-S
	60mm	TI-125.260-S	125.260-S
*	65mm	TI-125.265-S	125.265-S
*	70mm	TI-125.270-S	125.270-S

* Sizes not available in Screw Caddy

Washer for Small Fragment Screw



Titanium	Stainless Steel
TI-100.110	100.110




STERILE

Titanium	Stainless Steel
TI-100.110-S	100.110-S

**ST-005-4.0**

Screw Caddy for 4.0mm Cannulated Cancellous Screw(Tray Without Screw)



4.0mm Cannulated
Cancellous Screw, Short
Thread, Self-Tapping

4.0mm Cannulated
Cancellous Screw, Full
Thread, Self-Tapping

Unit	Length	Titanium	Stainless Steel
3	16mm	TI-124.216	124.216
3	18mm	TI-124.218	124.218
3	20mm	TI-124.220	124.220
3	22mm	TI-124.222	124.222
3	24mm	TI-124.224	124.224
3	26mm	TI-124.226	124.226
3	28mm	TI-124.228	124.228
3	30mm	TI-124.230	124.230
3	32mm	TI-124.232	124.232
3	34mm	TI-124.234	124.234
3	36mm	TI-124.236	124.236
3	38mm	TI-124.238	124.238
3	40mm	TI-124.240	124.240
3	42mm	TI-124.242	124.242
3	44mm	TI-124.244	124.244
3	46mm	TI-124.246	124.246
3	48mm	TI-124.248	124.248
3	50mm	TI-124.250	124.250
3	55mm	TI-124.255	124.255
3	60mm	TI-124.260	124.260

Unit	Length	Titanium	Stainless Steel
3	16mm	TI-125.216	125.216
3	18mm	TI-125.218	125.218
3	20mm	TI-125.220	125.220
3	22mm	TI-125.222	125.222
3	24mm	TI-125.224	125.224
3	26mm	TI-125.226	125.226
3	28mm	TI-125.228	125.228
3	30mm	TI-125.230	125.230
3	32mm	TI-125.232	125.232
3	34mm	TI-125.234	125.234
3	36mm	TI-125.236	125.236
3	38mm	TI-125.238	125.238
3	40mm	TI-125.240	125.240
3	42mm	TI-125.242	125.242
3	44mm	TI-125.244	125.244
3	46mm	TI-125.246	125.246
3	48mm	TI-125.248	125.248
3	50mm	TI-125.250	125.250
3	55mm	TI-125.255	125.255
3	60mm	TI-125.260	125.260

The background of the page is a light cream color. It features abstract, wavy, topographic-style lines in a light gray color. These lines are concentrated in the top-left and bottom-right corners, creating a sense of depth and movement. The lines flow and curve, resembling contour lines on a map or perhaps the ripples of water.

INST **INSTRUMENTS**



451-1.2-180E

Guide Wire, Ø1.2mm x Length 180mm



3426-01-180

Guide Wire with Threaded Trocar Tip, Ø1.2mm x Thread Length 5mm x Length 180mm



3767-1.20

Double Drill Guide, Ø2.7/1.2mm



1484-4.0

Cannulated Countersink with T-Handle for Ø4.0mm Screw



1486-2.7

Protection Sleeve for Drill Bit Ø2.8mm



7-056-02


Wire Sleeve for Guide Wire Ø2.7/1.2mm





7-056-01

Screw Holding Forcep for Ø4.0mm Screw



2103-2.7-160

Cannulated Drill Bit Plain Shank/Jacob Chuck,
Ø2.7mm x Length 160mm

3435-01

Stopper for Drill Bit Ø2.7mm



3433-01

Cannulated Screwdriver, Hex 2.5mm




1489-01

Holding Sleeve for Cannulated Screwdriver, Hex 2.5mm



1487-2.5

Allen Key, Hex 2.5mm, for Ø4.0mm Screw





2103-11

Cannulated T-Handle Tap for Ø4.0mm Screw



3431-01

Direct Measuring Device for Ø4.0mm Screw



1488-01

Extraction Screw with T- Handle for Ø4.0mm Screw




1485-1.2

Stylet, Ø1.2mm




ST-005-4.0

Screw Caddy for 4.0mm Cannulated Cancellous Screw



A10-009

Container with Lid for 4.0mm Cannulated Cancellous Screw Instrument set






ORDER INFORMATION

2337-000

4.0mm Cannulated Cancellous Screw Instrument Set



Code	Set Consisting Of	Unit
451-1.2-180E	Guide Wire, Ø1.2mm x Length 180mm	3
3426-01-180	Guide Wire with Threaded Trocar Tip, Ø1.2mm x Thread Length 5mm x Length 180mm	3
3767-1.20	Double Drill Guide, Ø2.7/1.2mm	1
1484-4.0	Cannulated Countersink with T-Handle for Ø4.0mm Screw	1
1486-2.7	Protection Sleeve for Drill Bit Ø2.8mm	1
7-056-02	Wire Sleeve for Guide Wire Ø2.7/1.2mm	1
7-056-01	Screw Holding Forcep for Ø4.0mm Screw	1
2103-2.7-160	Cannulated Drill Bit Plain Shank/Jacob Chuck, Ø2.7mm x Length 160mm	2
3435-01	Stopper for Drill Bit Ø2.7mm	1
3433-01	Cannulated Screwdriver, Hex 2.5mm	2
1489-01	Holding Sleeve for Cannulated Screwdriver, Hex 2.5mm	1
1487-2.5	Allen Key, Hex 2.5mm, for Ø4.0mm Screw	1
2103-11	Cannulated T-Handle Tap for Ø4.0mm Screw	1
3431-01	Direct Measuring Device for Ø4.0mm Screw	1
1488-01	Extraction Screw with T- Handle for Ø4.0mm Screw	1
1485-1.2	Stylet, Ø1.2mm	1
ST-005-4.0	Screw Caddy for 4.0mm Cannulated Cancellous Screw	1
A10-009	Container with Lid for 4.0mm Cannulated Cancellous Screw Instrument set	1



The background of the page features a light gray, abstract pattern of thin, wavy lines that flow across the top and bottom edges, creating a sense of movement and depth. The central area is a solid, light gray color.

4.5MM CANNULATED CORTICAL SCREW



4.5mm Cannulated Cortical Screw, Short Thread



Length	Titanium	Stainless Steel
20mm	TI-122.020	122.020
22mm	TI-122.022	122.022
24mm	TI-122.024	122.024
26mm	TI-122.026	122.026
28mm	TI-122.028	122.028
30mm	TI-122.030	122.030
32mm	TI-122.032	122.032
34mm	TI-122.034	122.034
36mm	TI-122.036	122.036
38mm	TI-122.038	122.038
40mm	TI-122.040	122.040
42mm	TI-122.042	122.042
44mm	TI-122.044	122.044
46mm	TI-122.046	122.046
48mm	TI-122.048	122.048
50mm	TI-122.050	122.050
52mm	TI-122.052	122.052
54mm	TI-122.054	122.054
56mm	TI-122.056	122.056
58mm	TI-122.058	122.058
60mm	TI-122.060	122.060
62mm	TI-122.062	122.062
64mm	TI-122.064	122.064
66mm	TI-122.066	122.066
68mm	TI-122.068	122.068
70mm	TI-122.070	122.070
72mm	TI-122.072	122.072
74mm	TI-122.074	122.074
76mm	TI-122.076	122.076
78mm	TI-122.078	122.078
80mm	TI-122.080	122.080


STERILE

Length	Titanium	Stainless Steel
20mm	TI-122.020-S	122.020-S
22mm	TI-122.022-S	122.022-S
24mm	TI-122.024-S	122.024-S
26mm	TI-122.026-S	122.026-S
28mm	TI-122.028-S	122.028-S
30mm	TI-122.030-S	122.030-S
32mm	TI-122.032-S	122.032-S
34mm	TI-122.034-S	122.034-S
36mm	TI-122.036-S	122.036-S
38mm	TI-122.038-S	122.038-S
40mm	TI-122.040-S	122.040-S
42mm	TI-122.042-S	122.042-S
44mm	TI-122.044-S	122.044-S
46mm	TI-122.046-S	122.046-S
48mm	TI-122.048-S	122.048-S
50mm	TI-122.050-S	122.050-S
52mm	TI-122.052-S	122.052-S
54mm	TI-122.054-S	122.054-S
56mm	TI-122.056-S	122.056-S
58mm	TI-122.058-S	122.058-S
60mm	TI-122.060-S	122.060-S
62mm	TI-122.062-S	122.062-S
64mm	TI-122.064-S	122.064-S
66mm	TI-122.066-S	122.066-S
68mm	TI-122.068-S	122.068-S
70mm	TI-122.070-S	122.070-S
72mm	TI-122.072-S	122.072-S
74mm	TI-122.074-S	122.074-S
76mm	TI-122.076-S	122.076-S
78mm	TI-122.078-S	122.078-S
80mm	TI-122.080-S	122.080-S



4.5mm Cannulated Cortical Screw, Full Thread



Length	Titanium	Stainless Steel
20mm	TI-122.220	122.220
22mm	TI-122.222	122.222
24mm	TI-122.224	122.224
26mm	TI-122.226	122.226
28mm	TI-122.228	122.228
30mm	TI-122.230	122.230
32mm	TI-122.232	122.232
34mm	TI-122.234	122.234
36mm	TI-122.236	122.236
38mm	TI-122.238	122.238
40mm	TI-122.240	122.240
42mm	TI-122.242	122.242
44mm	TI-122.244	122.244
46mm	TI-122.246	122.246
48mm	TI-122.248	122.248
50mm	TI-122.250	122.250
52mm	TI-122.252	122.252
54mm	TI-122.254	122.254
56mm	TI-122.256	122.256
58mm	TI-122.258	122.258
60mm	TI-122.260	122.260
62mm	TI-122.262	122.262
64mm	TI-122.264	122.264
66mm	TI-122.266	122.266
68mm	TI-122.268	122.268
70mm	TI-122.270	122.270
72mm	TI-122.272	122.272
74mm	TI-122.274	122.274
76mm	TI-122.276	122.276
78mm	TI-122.278	122.278
80mm	TI-122.280	122.280



STERILE

Length	Titanium	Stainless Steel
20mm	TI-122.220-S	122.220-S
22mm	TI-122.222-S	122.222-S
24mm	TI-122.224-S	122.224-S
26mm	TI-122.226-S	122.226-S
28mm	TI-122.228-S	122.228-S
30mm	TI-122.230-S	122.230-S
32mm	TI-122.232-S	122.232-S
34mm	TI-122.234-S	122.234-S
36mm	TI-122.236-S	122.236-S
38mm	TI-122.238-S	122.238-S
40mm	TI-122.240-S	122.240-S
42mm	TI-122.242-S	122.242-S
44mm	TI-122.244-S	122.244-S
46mm	TI-122.246-S	122.246-S
48mm	TI-122.248-S	122.248-S
50mm	TI-122.250-S	122.250-S
52mm	TI-122.252-S	122.252-S
54mm	TI-122.254-S	122.254-S
56mm	TI-122.256-S	122.256-S
58mm	TI-122.258-S	122.258-S
60mm	TI-122.260-S	122.260-S
62mm	TI-122.262-S	122.262-S
64mm	TI-122.264-S	122.264-S
66mm	TI-122.266-S	122.266-S
68mm	TI-122.268-S	122.268-S
70mm	TI-122.270-S	122.270-S
72mm	TI-122.272-S	122.272-S
74mm	TI-122.274-S	122.274-S
76mm	TI-122.276-S	122.276-S
78mm	TI-122.278-S	122.278-S
80mm	TI-122.280-S	122.280-S



4.5mm Cannulated Cortical Screw, Short Thread, Self-Tapping



Length	Titanium	Stainless Steel
20mm	TI-122.120	122.120
22mm	TI-122.122	122.122
24mm	TI-122.124	122.124
26mm	TI-122.126	122.126
28mm	TI-122.128	122.128
30mm	TI-122.130	122.130
32mm	TI-122.132	122.132
34mm	TI-122.134	122.134
36mm	TI-122.136	122.136
38mm	TI-122.138	122.138
40mm	TI-122.140	122.140
42mm	TI-122.142	122.142
44mm	TI-122.144	122.144
46mm	TI-122.146	122.146
48mm	TI-122.148	122.148
50mm	TI-122.150	122.150
52mm	TI-122.152	122.152
54mm	TI-122.154	122.154
56mm	TI-122.156	122.156
58mm	TI-122.158	122.158
60mm	TI-122.160	122.160
62mm	TI-122.162	122.162
64mm	TI-122.164	122.164
66mm	TI-122.166	122.166
68mm	TI-122.168	122.168
70mm	TI-122.170	122.170
72mm	TI-122.172	122.172
74mm	TI-122.174	122.174
76mm	TI-122.176	122.176
78mm	TI-122.178	122.178
80mm	TI-122.180	122.180


STERILE

Length	Titanium	Stainless Steel
20mm	TI-122.120-S	122.120-S
22mm	TI-122.122-S	122.122-S
24mm	TI-122.124-S	122.124-S
26mm	TI-122.126-S	122.126-S
28mm	TI-122.128-S	122.128-S
30mm	TI-122.130-S	122.130-S
32mm	TI-122.132-S	122.132-S
34mm	TI-122.134-S	122.134-S
36mm	TI-122.136-S	122.136-S
38mm	TI-122.138-S	122.138-S
40mm	TI-122.140-S	122.140-S
42mm	TI-122.142-S	122.142-S
44mm	TI-122.144-S	122.144-S
46mm	TI-122.146-S	122.146-S
48mm	TI-122.148-S	122.148-S
50mm	TI-122.150-S	122.150-S
52mm	TI-122.152-S	122.152-S
54mm	TI-122.154-S	122.154-S
56mm	TI-122.156-S	122.156-S
58mm	TI-122.158-S	122.158-S
60mm	TI-122.160-S	122.160-S
62mm	TI-122.162-S	122.162-S
64mm	TI-122.164-S	122.164-S
66mm	TI-122.166-S	122.166-S
68mm	TI-122.168-S	122.168-S
70mm	TI-122.170-S	122.170-S
72mm	TI-122.172-S	122.172-S
74mm	TI-122.174-S	122.174-S
76mm	TI-122.176-S	122.176-S
78mm	TI-122.178-S	122.178-S
80mm	TI-122.180-S	122.180-S



4.5mm Cannulated Cortical Screw, Full Thread, Self-Tapping



Length	Titanium	Stainless Steel
20mm	TI-122.320	122.320
22mm	TI-122.322	122.322
24mm	TI-122.324	122.324
26mm	TI-122.326	122.326
28mm	TI-122.328	122.328
30mm	TI-122.330	122.330
32mm	TI-122.332	122.332
34mm	TI-122.334	122.334
36mm	TI-122.336	122.336
38mm	TI-122.338	122.338
40mm	TI-122.340	122.340
42mm	TI-122.342	122.342
44mm	TI-122.344	122.344
46mm	TI-122.346	122.346
48mm	TI-122.348	122.348
50mm	TI-122.350	122.350
52mm	TI-122.352	122.352
54mm	TI-122.354	122.354
56mm	TI-122.356	122.356
58mm	TI-122.358	122.358
60mm	TI-122.360	122.360
62mm	TI-122.362	122.362
64mm	TI-122.364	122.364
66mm	TI-122.366	122.366
68mm	TI-122.368	122.368
70mm	TI-122.370	122.370
72mm	TI-122.372	122.372
74mm	TI-122.374	122.374
76mm	TI-122.376	122.376
78mm	TI-122.378	122.378
80mm	TI-122.380	122.380




STERILE

Length	Titanium	Stainless Steel
20mm	TI-122.320-S	122.320-S
22mm	TI-122.322-S	122.322-S
24mm	TI-122.324-S	122.324-S
26mm	TI-122.326-S	122.326-S
28mm	TI-122.328-S	122.328-S
30mm	TI-122.330-S	122.330-S
32mm	TI-122.332-S	122.332-S
34mm	TI-122.334-S	122.334-S
36mm	TI-122.336-S	122.336-S
38mm	TI-122.338-S	122.338-S
40mm	TI-122.340-S	122.340-S
42mm	TI-122.342-S	122.342-S
44mm	TI-122.344-S	122.344-S
46mm	TI-122.346-S	122.346-S
48mm	TI-122.348-S	122.348-S
50mm	TI-122.350-S	122.350-S
52mm	TI-122.352-S	122.352-S
54mm	TI-122.354-S	122.354-S
56mm	TI-122.356-S	122.356-S
58mm	TI-122.358-S	122.358-S
60mm	TI-122.360-S	122.360-S
62mm	TI-122.362-S	122.362-S
64mm	TI-122.364-S	122.364-S
66mm	TI-122.366-S	122.366-S
68mm	TI-122.368-S	122.368-S
70mm	TI-122.370-S	122.370-S
72mm	TI-122.372-S	122.372-S
74mm	TI-122.374-S	122.374-S
76mm	TI-122.376-S	122.376-S
78mm	TI-122.378-S	122.378-S
80mm	TI-122.380-S	122.380-S

Washer for Large Fragment Screw




Titanium	Stainless Steel
TI-100.120	100.120



STERILE

Titanium	Stainless Steel
TI-100.120-S	100.120-S

The background of the page is a light gray color. It features abstract, wavy, topographic-style lines in a slightly darker gray. These lines are concentrated in the top-left and bottom-right corners, creating a sense of depth and movement. The central area of the page is mostly blank, providing a clear space for the text.

INSTRUMENTS



3435-04

Cannulated Drill Bit with Quick Coupling, Ø4.5mm x Length 170mm



3435-03

Cannulated Drill Bit with Quick Coupling, Ø3.2mm x Length 170mm



2100-1.5-160

Drill Bit Quick Coupling End, Ø1.5mm x Length 160mm



3767-4.5

Double Drill Guide, Ø4.5/3.2mm



3767-3.2

Double Drill Guide, Ø3.2/1.6mm



3432-01

Bone Tap with Quick Coupling for Ø4.5mm Cannulated Screw, Quick Coupling





3433-02

Cannulated Screwdriver, Hex 3.5mm



3446-000

Fiber Handle With Quick Coupling, Cannulated



3430-02

Cannulated Countersink With Quick Coupling For Ø4.5mm Screw



3427-02

Stylet, Ø1.6mm




3426-02

Guide Wire With Threaded Trocar Tip Ø1.6mmxthread Length 14mmx Length 240mm



3431-02


Direct Measuring Device For Ø4.5mm Screw





3425-000

Holding Clip For Washer




3415-000

Screw Holding Forcep



7-122-01

Instrument Tray for 4.5mm Cannulated Cortical Screw Instruments Set



A10-010

Container for 4.5mm Cannulated Cortical Screw Instruments Set






ORDER INFORMATION

2338-000

4.5mm Cannulated Cortical Screw Instrument Set



Code	Set Consisting Of	Units
3435-04	Cannulated Drill Bit with Quick Coupling, Ø4.5mm x Length 170mm	1
3435-03	Cannulated Drill Bit with Quick Coupling, Ø3.2mm x Length 170mm	2
2100-1.5-160	Drill Bit Quick Coupling End, Ø1.5mm x Length 160mm	1
3767-4.5	Double Drill Guide, Ø4.5/3.2mm	1
3767-3.2	Double Drill Guide, Ø3.2/1.6mm	1
3432-01	Bone Tap with Quick Coupling for Ø4.5mm Cannulated Screw, Quick Coupling	1
3433-02	Cannulated Screwdriver, Hex 3.5mm	1
3446-000	Fiber Handle with Quick Coupling, Cannulated	1
3430-02	Cannulated Countersink with Quick Coupling for Ø4.5mm Screw	1
3427-02	Stylet, Ø1.6mm	1
3426-02	Guide Wire with Threaded Trocar Tip Ø1.6mm x Thread Length 14mm x Length 240mm	10
3431-02	Direct Measuring Device for Ø4.5mm Screw	1
3425-000	Holding Clip For Washer	1
3415-000	Screw Holding Forcep	1
7-122-01	Instrument Tray for 4.5mm Cannulated Cortical Screw Instruments Set	1
A10-010	Container for 4.5mm Cannulated Cortical Screw Instruments Set	1

The background of the page features a series of thin, light gray wavy lines that flow from the top left towards the bottom right, creating a sense of movement and depth. These lines are more densely packed in some areas, particularly near the top and bottom edges, and more spread out in the center where the text is located.

6.5MM CANNULATED SCREW



6.5mm Cannulated Screw, 16mm Thread, Self-Drilling



Length	Titanium	Stainless Steel
30mm	TI-132.030	132.030
35mm	TI-132.035	132.035
40mm	TI-132.040	132.040
45mm	TI-132.045	132.045
50mm	TI-132.050	132.050
55mm	TI-132.055	132.055
60mm	TI-132.060	132.060
65mm	TI-132.065	132.065
70mm	TI-132.070	132.070
75mm	TI-132.075	132.075
80mm	TI-132.080	132.080
85mm	TI-132.085	132.085
90mm	TI-132.090	132.090
95mm	TI-132.095	132.095
100mm	TI-132.100	132.100
105mm	TI-132.105	132.105
110mm	TI-132.110	132.110
115mm	TI-132.115	132.115
120mm	TI-132.120	132.120
125mm	TI-132.125	132.125
130mm	TI-132.130	132.130
135mm	TI-132.135	132.135
140mm	TI-132.140	132.140
145mm	TI-132.145	132.145
150mm	TI-132.150	132 150


STERILE

Length	Titanium	Stainless Steel
30mm	TI-132.030-S	132.030-S
35mm	TI-132.035-S	132.035-S
40mm	TI-132.040-S	132.040-S
45mm	TI-132.045-S	132.045-S
50mm	TI-132.050-S	132.050-S
55mm	TI-132.055-S	132.055-S
60mm	TI-132.060-S	132.060-S
65mm	TI-132.065-S	132.065-S
70mm	TI-132.070-S	132.070-S
75mm	TI-132.075-S	132.075-S
80mm	TI-132.080-S	132.080-S
85mm	TI-132.085-S	132.085-S
90mm	TI-132.090-S	132.090-S
95mm	TI-132.095-S	132.095-S
100mm	TI-132.100-S	132.100-S
105mm	TI-132.105-S	132.105-S
110mm	TI-132.110-S	132.110-S
115mm	TI-132.115-S	132.115-S
120mm	TI-132.120-S	132.120-S
125mm	TI-132.125-S	132.125-S
130mm	TI-132.130-S	132.130-S
135mm	TI-132.135-S	132.135-S
140mm	TI-132.140-S	132.140-S
145mm	TI-132.145-S	132.145-S
150mm	TI-132.150-S	132 150-S



6.5mm Cannulated Screw, 32mm Thread, Self-Drilling



Length	Titanium	Stainless Steel
45mm	TI-132.245	132.245
50mm	TI-132.250	132.250
55mm	TI-132.255	132.255
60mm	TI-132.260	132.260
65mm	TI-132.265	132.265
70mm	TI-132.270	132.270
75mm	TI-132.275	132.275
80mm	TI-132.280	132.280
85mm	TI-132.285	132.285
90mm	TI-132.290	132.290
95mm	TI-132.295	132.295
100mm	TI-132.300	132.300
105mm	TI-132.305	132.305
110mm	TI-132.310	132.310
115mm	TI-132.315	132.315
120mm	TI-132.320	132.320
125mm	TI-132.325	132.325
130mm	TI-132.330	132.330
135mm	TI-132.335	132.335
140mm	TI-132.340	132.340
145mm	TI-132.345	132.345
150mm	TI-132.350	132.350



STERILE

Length	Titanium	Stainless Steel
45mm	TI-132.245-S	132.245-S
50mm	TI-132.250-S	132.250-S
55mm	TI-132.255-S	132.255-S
60mm	TI-132.260-S	132.260-S
65mm	TI-132.265-S	132.265-S
70mm	TI-132.270-S	132.270-S
75mm	TI-132.275-S	132.275-S
80mm	TI-132.280-S	132.280-S
85mm	TI-132.285-S	132.285-S
90mm	TI-132.290-S	132.290-S
95mm	TI-132.295-S	132.295-S
100mm	TI-132.300-S	132.300-S
105mm	TI-132.305-S	132.305-S
110mm	TI-132.310-S	132.310-S
115mm	TI-132.315-S	132.315-S
120mm	TI-132.320-S	132.320-S
125mm	TI-132.325-S	132.325-S
130mm	TI-132.330-S	132.330-S
135mm	TI-132.335-S	132.335-S
140mm	TI-132.340-S	132.340-S
145mm	TI-132.345-S	132.345-S
150mm	TI-132.350-S	132.350-S



6.5mm Cannulated Screw, Full Threaded, Self-Drilling



Length	Titanium	Stainless Steel
20mm	TI-133.020	133.020
25mm	TI-133.025	133.025
30mm	TI-133.030	133.030
35mm	TI-133.035	133.035
40mm	TI-133.040	133.040
45mm	TI-133.045	133.045
50mm	TI-133.050	133.050
55mm	TI-133.055	133.055
60mm	TI-133.060	133.060
65mm	TI-133.065	133.065
70mm	TI-133.070	133.070
75mm	TI-133.075	133.075
80mm	TI-133.080	133.080
85mm	TI-133.085	133.085
90mm	TI-133.090	133.090
95mm	TI-133.095	133.095
100mm	TI-133.100	133.100
105mm	TI-133.105	133.105
110mm	TI-133.110	133.110
115mm	TI-133.115	133.115
120mm	TI-133.120	133.120
125mm	TI-133.125	133.125
130mm	TI-133.130	133.130




STERILE

Length	Titanium	Stainless Steel
20mm	TI-133.020-S	133.020-S
25mm	TI-133.025-S	133.025-S
30mm	TI-133.030-S	133.030-S
35mm	TI-133.035-S	133.035-S
40mm	TI-133.040-S	133.040-S
45mm	TI-133.045-S	133.045-S
50mm	TI-133.050-S	133.050-S
55mm	TI-133.055-S	133.055-S
60mm	TI-133.060-S	133.060-S
65mm	TI-133.065-S	133.065-S
70mm	TI-133.070-S	133.070-S
75mm	TI-133.075-S	133.075-S
80mm	TI-133.080-S	133.080-S
85mm	TI-133.085-S	133.085-S
90mm	TI-133.090-S	133.090-S
95mm	TI-133.095-S	133.095-S
100mm	TI-133.100-S	133.100-S
105mm	TI-133.105-S	133.105-S
110mm	TI-133.110-S	133.110-S
115mm	TI-133.115-S	133.115-S
120mm	TI-133.120-S	133.120-S
125mm	TI-133.125-S	133.125-S
130mm	TI-133.130-S	133.130-S

Washer for Large Fragment Screw




Titanium	Stainless Steel
TI-100.120	100.120



STERILE

Titanium	Stainless Steel
TI-100.120-S	100.120-S

The background of the page features a series of thin, light gray wavy lines that flow from the top left towards the bottom right, creating a sense of movement and depth. These lines are more densely packed in some areas, particularly near the top and bottom edges, and more spread out in the center.

INSTRUMENTS



7-133-01

Guide Wire, Ø2.0mm x Length 250mm for 6.5mm Cannulated Cancellous Screw



7-133-02

Guide Wire with Threaded Tip, Ø2.0mm x Thread Length 8mm x Length 250mm for 6.5mm Cannulated Cancellous Screw



7-133-03

Cannulated Screwdriver, Hex 3.5mm for 6.5mm Cannulated Cancellous Screw



7-133-04

Holding Sleeve for Cannulated Screwdriver, Hex 3.5mm for 6.5mm Cannulated Cancellous Screw



7-133-05

Wire Sleeve Bolt Key for 6.5mm Cannulated Cancellous Screw



7-133-06

Wire Sleeve for Guide Wire, Ø2.0mm for 6.5mm Cannulated Cancellous Screw





7-133-07

Double Drill Guide, Ø4.5/2.0mm for 6.5mm Cannulated Cancellous Screw



7-133-08

Parallel Guide for Guide Wire, Ø2.0mm for 6.5mm Cannulated Cancellous Screw



7-133-09

Direct Measuring Device for Cannulated Cancellous Screw for 6.5mm
Cannulated Cancellous Screw

7-133-10

Screw Holding Forcep for Ø6.5mm Screw



7-133-11

Cannulated Drill Bit Plain Shank/Jacob Chuck End, Ø4.5mm x Length 250mm
for 6.5mm Cannulated Cancellous Screw

7-133-12

Stopper for Drill Bit Ø4.5mm for 6.5mm Cannulated Cancellous Screw





7-133-13

Allen Key, Hex 3.0mm, for Ø6.5mm Screw



7-133-14

Trocár Ø2.0mm, for Ø6.5mm Screw



7-133-15

Guide Sleeve Ø4.5/2.0mm for 6.5mm Cannulated Cancellous Screw



7-133-16

Drill Sleeve Ø8.0/4.5mm for 6.5mm Cannulated Cancellous Screw



7-133-17

Cannulated Countersink with T-Handle for Ø6.5mm Screw



7-133-18

Extraction Screw, Ø5.0mm





7-133-19

Protection Sleeve for Drill Bit Ø4.5mm for 6.5mm
Cannulated Cancellous Screw



7-133-20

Cannulated Tap with T-Handle for Ø6.5mm Screw




7-133-21

Stylet Ø2.0mm for 6.5mm Cannulated Cancellous Screw




7-133-22

Screw Caddy for 6.5mm Cannulated Cancellous Screw




7-133-23

Instrument Tray for 6.5mm Cannulated Cancellous Screw Instrument Set



7-133-24

Container for 6.5mm Cannulated Cancellous Screw Instrument Set





7-133

6.5mm Cannulated Cancellous Instrument Set





ORDER INFORMATION

7-133

6.5mm Cannulated Cancellous Screw Instrument Set

Code	Set Consisting Of	Unit
451-2.0-250E	Guide Wire, Ø2.0mm x Length 250mm	3
GW-2.0-250	Guide Wire with Threaded Tip, Ø2.0mm x Thread Length 8mm x Length 250mm	3
3433-02	Cannulated Screwdriver, Hex 3.5mm	2
1490-01	Holding Sleeve for Cannulated Screwdriver, Hex 3.5mm	1
1482-01	Wire Sleeve Bolt Key	1
1481-01	Wire Sleeve for Guide Wire, Ø2.0mm	1
1480-4.5-2.0	Double Drill Guide, Ø4.5/2.0mm	1
3436-013	Parallel Guide for Guide Wire, Ø2.0mm	1
3436-014	Direct Measuring Device for Cannulated Cancellous Screw	1
7-057-01	Screw Holding Forcep for Ø7.0mm Screw	1
2191-4.5-250	Cannulated Drill Bit Plain Shank/Jacob Chuck End, Ø4.5mm x Length 250mm	1
3435-000	Stopper for Drill Bit Ø4.5mm	1
7-057-02	Allen Key, Hex 3.0mm, for Ø7.0mm Screw	1
1491-01	Trocar Ø2.0mm, for Ø7.0mm Screw	1
1492-01	Guide Sleeve Ø4.5/2.0mm	1
1493-01	Drill Sleeve Ø8.0/4.5mm	1
1483-01	Cannulated Countersink with T-Handle for Ø7.0mm Screw	1
3443-38-5.0	Extraction Screw, Ø5.0mm for Large Fragment	1
3436-012	Protection Sleeve for Drill Bit Ø4.5mm	1
3436-011	Cannulated Tap with T-Handle for Ø7.0mm Screw	1
3436-016	Stylet Ø2.0mm	1
ST-005-7.0	Screw Caddy for 7.0mm Cannulated Cancellous Screw	1
7-057-03	Instrument Tray for 7.0mm Cannulated Cancellous Screw Instrument Set	2
A10-012	Container For 7.0mm Cannulated Cancellous Screw Instrument Set	1

4.10 Injection process filling phase



5.00

0.48 sec

P_{max}

5.4 bar



8

7

6

5

4

3

2

1

0



s

5.0

10.0 mm

v

12.0

14.0 %

p

55

55 bar

v [%]

110.0

88.0

66.0

44.0

22.0

0.0

mm

mm

160 P [bar]

128

96

64

32

0

80.0 s [mm]

0.0

16.0

32.0

48.0

64.0

80.0

mm

mm



1.00





7.0MM CANNULATED CANCELLOUS SCREW



7.0mm Cannulated Cancellous Screw, 16mm Thread, Self-Drilling



Length	Titanium	Stainless Steel
30mm	TI-129.230	129.230
35mm	TI-129.235	129.235
40mm	TI-129.240	129.240
45mm	TI-129.245	129.245
50mm	TI-129.250	129.250
55mm	TI-129.255	129.255
60mm	TI-129.260	129.260
65mm	TI-129.265	129.265
70mm	TI-129.270	129.270
75mm	TI-129.275	129.275
80mm	TI-129.280	129.280
85mm	TI-129.285	129.285
90mm	TI-129.290	129.290
95mm	TI-129.295	129.295
100mm	TI-129.300	129.300
105mm	TI-129.305	129.305
110mm	TI-129.310	129.310
* 115mm	TI-129.315	129.315
* 120mm	TI-129.320	129.320



STERILE

Length	Titanium	Stainless Steel
30mm	TI-129.230-S	129.230-S
35mm	TI-129.235-S	129.235-S
40mm	TI-129.240-S	129.240-S
45mm	TI-129.245-S	129.245-S
50mm	TI-129.250-S	129.250-S
55mm	TI-129.255-S	129.255-S
60mm	TI-129.260-S	129.260-S
65mm	TI-129.265-S	129.265-S
70mm	TI-129.270-S	129.270-S
75mm	TI-129.275-S	129.275-S
80mm	TI-129.280-S	129.280-S
85mm	TI-129.285-S	129.285-S
90mm	TI-129.290-S	129.290-S
95mm	TI-129.295-S	129.295-S
100mm	TI-129.300-S	129.300-S
105mm	TI-129.305-S	129.305-S
110mm	TI-129.310-S	129.310-S
* 115mm	TI-129.315-S	129.315-S
* 120mm	TI-129.320-S	129.320-S

7.0mm Cannulated Cancellous Screw, 32mm Thread, Self-Drilling



Length	Titanium	Stainless Steel
35mm	TI-130.235	130.235
40mm	TI-130.240	130.240
45mm	TI-130.245	130.245
50mm	TI-130.250	130.250
55mm	TI-130.255	130.255
60mm	TI-130.260	130.260
65mm	TI-130.265	130.265
70mm	TI-130.270	130.270
75mm	TI-130.275	130.275
80mm	TI-130.280	130.280
85mm	TI-130.285	130.285
90mm	TI-130.290	130.290
95mm	TI-130.295	130.295
100mm	TI-130.300	130.300
105mm	TI-130.305	130.305
110mm	TI-130.310	130.310
* 115mm	TI-130.315	130.315
* 120mm	TI-130.320	130.320



STERILE


Length	Titanium	Stainless Steel
35mm	TI-130.235-S	130.235-S
40mm	TI-130.240-S	130.240-S
45mm	TI-130.245-S	130.245-S
50mm	TI-130.250-S	130.250-S
55mm	TI-130.255-S	130.255-S
60mm	TI-130.260-S	130.260-S
65mm	TI-130.265-S	130.265-S
70mm	TI-130.270-S	130.270-S
75mm	TI-130.275-S	130.275-S
80mm	TI-130.280-S	130.280-S
85mm	TI-130.285-S	130.285-S
90mm	TI-130.290-S	130.290-S
95mm	TI-130.295-S	130.295-S
100mm	TI-130.300-S	130.300-S
105mm	TI-130.305-S	130.305-S
110mm	TI-130.310-S	130.310-S
* 115mm	TI-130.315-S	130.315-S
* 120mm	TI-130.320-S	130.320-S



7.0mm Cannulated Cancellous Screw, Full Thread, Self-Drilling



Length	Titanium	Stainless Steel
30mm	TI-131.230	131.230
35mm	TI-131.235	131.235
40mm	TI-131.240	131.240
45mm	TI-131.245	131.245
50mm	TI-131.250	131.250
55mm	TI-131.255	131.255
60mm	TI-131.260	131.260
65mm	TI-131.265	131.265
70mm	TI-131.270	131.270
75mm	TI-131.275	131.275
80mm	TI-131.280	131.280
85mm	TI-131.285	131.285
90mm	TI-131.290	131.290
95mm	TI-131.295	131.295
100mm	TI-131.300	131.300
105mm	TI-131.305	131.305
110mm	TI-131.310	131.310
* 115mm	TI-131.315	131.315
* 120mm	TI-131.320	131.320




STERILE

Length	Titanium	Stainless Steel
30mm	TI-131.230-S	131.230-S
35mm	TI-131.235-S	131.235-S
40mm	TI-131.240-S	131.240-S
45mm	TI-131.245-S	131.245-S
50mm	TI-131.250-S	131.250-S
55mm	TI-131.255-S	131.255-S
60mm	TI-131.260-S	131.260-S
65mm	TI-131.265-S	131.265-S
70mm	TI-131.270-S	131.270-S
75mm	TI-131.275-S	131.275-S
80mm	TI-131.280-S	131.280-S
85mm	TI-131.285-S	131.285-S
90mm	TI-131.290-S	131.290-S
95mm	TI-131.295-S	131.295-S
100mm	TI-131.300-S	131.300-S
105mm	TI-131.305-S	131.305-S
110mm	TI-131.310-S	131.310-S
* 115mm	TI-131.315-S	131.315-S
* 120mm	TI-131.320-S	131.320-S

Washer for Large Fragment Screw



Titanium	Stainless Steel
TI-100.120	100.120





STERILE

Titanium	Stainless Steel
TI-100.120-S	100.120-S



ST-005-7.0 : Screw Caddy for 7.0mm Cannulated Cancellous Screw (Tray without Screw)



The background of the page features a series of thin, light gray wavy lines that flow from the top left towards the bottom right, creating a sense of movement and depth. These lines are more densely packed in some areas, particularly near the top and bottom edges, and more spread out in the center.

INSTRUMENTS




451-2.0-250E

Guide Wire, Ø2.0mm x Length 250mm



GW-2.0-250

Guide Wire with Threaded Tip, Ø2.0mm x Thread Length 8mm x Length 250mm



3433-02

Cannulated Screwdriver, Hex 3.5mm



1490-01

Holding Sleeve for Cannulated Screwdriver, Hex 3.5mm



1482-01

Wire Sleeve Bolt Key



1481-01

Wire Sleeve for Guide Wire, Ø2.0mm





1480-4.5-2.0

Double Drill Guide, Ø4.5/2.0mm



3436-013

Parallel Guide for Guide Wire, Ø2.0mm



3436-014

Direct Measuring Device for Cannulated Cancellous Screw



7-057-01

Screw Holding Forcep for Ø7.0mm Screw



2191-4.5-250

Cannulated Drill Bit Plain Shank/Jacob Chuck End, Ø4.5mm x Length 250mm



3435-000

Stopper for Drill Bit Ø4.5mm





7-057-02

Allen Key, Hex 3.0mm, for Ø7.0mm Screw



1491-01

Trocár Ø2.0mm, for Ø7.0mm Screw



1492-01

Guide Sleeve Ø4.5/2.0mm



1493-01

Drill Sleeve Ø8.0/4.5mm



1483-01

Cannulated Countersink with T-Handle for Ø7.0mm Screw



3443.38-5.0

Extraction Screw, Ø5.0mm





3436-012

Protection Sleeve for Drill Bit Ø4.5mm



3436-011

Cannulated Tap with T-Handle for Ø7.0mm Screw




3436-016

Stylet Ø2.0mm




ST-005-7.0

Screw Caddy for 7.0mm Cannulated Cancellous Screw




7-057-03

Instrument Tray for 7.0mm Cannulated Cancellous Screw Instrument Set



A10-012

Container For 7.0mm Cannulated Cancellous Screw Instrument Set





IS006 : 7.0mm Cannulated Cancellous Instrument Set





ORDER INFORMATION

IS006	7.0mm Cannulated Cancellous Screw Instrument Set
-------	--

Code	Set Consisting Of	Unit
451-2.0-250E	Guide Wire, Ø2.0mm x Length 250mm	3
GW-2.0-250	Guide Wire with Threaded Tip, Ø2.0mm x Thread Length 8mm x Length 250mm	3
3433-02	Cannulated Screwdriver, Hex 3.5mm	2
1490-01	Holding Sleeve for Cannulated Screwdriver, Hex 3.5mm	1
1482-01	Wire Sleeve Bolt Key	1
1481-01	Wire Sleeve for Guide Wire, Ø2.0mm	1
1480-4.5-2.0	Double Drill Guide, Ø4.5/2.0mm	1
3436-013	Parallel Guide for Guide Wire, Ø2.0mm	1
3436-014	Direct Measuring Device for Cannulated Cancellous Screw	1
7-057-01	Screw Holding Forcep for Ø7.0mm Screw	1
2191-4.5-250	Cannulated Drill Bit Plain Shank/Jacob Chuck End, Ø4.5mm x Length 250mm	1
3435-000	Stopper for Drill Bit Ø4.5mm	1
7-057-02	Allen Key, Hex 3.0mm, for Ø7.0mm Screw	1
1491-01	Trocar Ø2.0mm, for Ø7.0mm Screw	1
1492-01	Guide Sleeve Ø4.5/2.0mm	1
1493-01	Drill Sleeve Ø8.0/4.5mm	1
1483-01	Cannulated Countersink with T-Handle for Ø7.0mm Screw	1
3443-38-5.0	Extraction Screw, Ø5.0mm	1
3436-012	Protection Sleeve for Drill Bit Ø4.5mm	1
3436-011	Cannulated Tap with T-Handle for Ø7.0mm Screw	1
3436-016	Stylet Ø2.0mm	1
ST-005-7.0	Screw Caddy for 7.0mm Cannulated Cancellous Screw	1
7-057-03	Instrument Tray for 7.0mm Cannulated Cancellous Screw Instrument Set	2
A10-012	Container For 7.0mm Cannulated Cancellous Screw Instrument Set	1



The background of the page features a light gray, abstract pattern of wavy, concentric lines that resemble topographical contours or liquid ripples. These lines are more densely packed in the upper and lower corners, creating a sense of depth and movement, while the center of the page is relatively clear to accommodate the text.

7.3MM CANNULATED CANCELLOUS SCREW



7.3mm Cannulated Cancellous Screw, 16mm Thread, Self-Drilling



Length	Titanium	Stainless Steel
30	TI-134.030	134.030
35	TI-134.035	134.035
40	TI-134.040	134.040
45	TI-134.045	134.045
50	TI-134.050	134.050
55	TI-134.055	134.055
60	TI-134.060	134.060
65	TI-134.065	134.065
70	TI-134.070	134.070
75	TI-134.075	134.075
80	TI-134.080	134.080
85	TI-134.085	134.085
90	TI-134.090	134.090
95	TI-134.095	134.095
100	TI-134.100	134.100
105	TI-134.105	134.105
110	TI-134.110	134.110
* 115	TI-134.115	134.115
* 120	TI-134.120	134.120
* 125	TI-134.125	134.125
* 130	TI-134.130	134.130
* 135	TI-134.135	134.135
* 140	TI-134.140	134.140
* 145	TI-134.145	134.145
* 150	TI-134.150	134.150



STERILE

Length	Titanium	Stainless Steel
30	TI-134.030-S	134.030-S
35	TI-134.035-S	134.035-S
40	TI-134.040-S	134.040-S
45	TI-134.045-S	134.045-S
50	TI-134.050-S	134.050-S
55	TI-134.055-S	134.055-S
60	TI-134.060-S	134.060-S
65	TI-134.065-S	134.065-S
70	TI-134.070-S	134.070-S
75	TI-134.075-S	134.075-S
80	TI-134.080-S	134.080-S
85	TI-134.085-S	134.085-S
90	TI-134.090-S	134.090-S
95	TI-134.095-S	134.095-S
100	TI-134.100-S	134.100-S
105	TI-134.105-S	134.105-S
110	TI-134.110-S	134.110-S
* 115	TI-134.115-S	134.115-S
* 120	TI-134.120-S	134.120-S
* 125	TI-134.125-S	134.125-S
* 130	TI-134.130-S	134.130-S
* 135	TI-134.135-S	134.135-S
* 140	TI-134.140-S	134.140-S
* 145	TI-134.145-S	134.145-S
* 150	TI-134.150-S	134.150-S


* Sizes not available in Screw Caddy



7.3mm Cannulated Cancellous Screw, 32mm Thread, Self-Drilling


STERILE

Length	Titanium	Stainless Steel
35	TI-134.235	134.235
40	TI-134.240	134.240
45	TI-134.245	134.245
50	TI-134.250	134.250
55	TI-134.255	134.255
60	TI-134.260	134.260
65	TI-134.265	134.265
70	TI-134.270	134.270
75	TI-134.275	134.275
80	TI-134.280	134.280
85	TI-134.285	134.285
90	TI-134.290	134.290
95	TI-134.295	134.295
100	TI-134.300	134.300
105	TI-134.305	134.305
110	TI-134.310	134.310
* 115	TI-134.315	134.315
* 120	TI-134.320	134.320
* 125	TI-134.325	134.325
* 130	TI-134.330	134.330
* 135	TI-134.335	134.335
* 140	TI-134.340	134.340
* 145	TI-134.345	134.345
* 150	TI-134.350	134.350



Length	Titanium	Stainless Steel
35	TI-134.235-S	134.235-S
40	TI-134.240-S	134.240-S
45	TI-134.245-S	134.245-S
50	TI-134.250-S	134.250-S
55	TI-134.255-S	134.255-S
60	TI-134.260-S	134.260-S
65	TI-134.265-S	134.265-S
70	TI-134.270-S	134.270-S
75	TI-134.275-S	134.275-S
80	TI-134.280-S	134.280-S
85	TI-134.285-S	134.285-S
90	TI-134.290-S	134.290-S
95	TI-134.295-S	134.295-S
100	TI-134.300-S	134.300-S
105	TI-134.305-S	134.305-S
110	TI-134.310-S	134.310-S
* 115	TI-134.315-S	134.315-S
* 120	TI-134.320-S	134.320-S
* 125	TI-134.325-S	134.325-S
* 130	TI-134.330-S	134.330-S
* 135	TI-134.335-S	134.335-S
* 140	TI-134.340-S	134.340-S
* 145	TI-134.345-S	134.345-S
* 150	TI-134.350-S	134.350-S


* Sizes not available in Screw Caddy



7.3mm Cannulated Cancellous Screw, Full Thread, Self-Drilling



Length	Titanium	Stainless Steel
20	TI-148.020	148.020
25	TI-148.025	148.025
30	TI-148.030	148.030
35	TI-148.035	148.035
40	TI-148.040	148.040
45	TI-148.045	148.045
50	TI-148.050	148.050
55	TI-148.055	148.055
60	TI-148.060	148.060
65	TI-148.065	148.065
70	TI-148.070	148.070
75	TI-148.075	148.075
80	TI-148.080	148.080
85	TI-148.085	148.085
90	TI-148.090	148.090
95	TI-148.095	148.095
100	TI-148.100	148.100
105	TI-148.105	148.105
110	TI-148.110	148.110
* 115	TI-148.115	148.115
* 120	TI-148.120	148.120
* 125	TI-148.125	148.125
* 130	TI-148.130	148.130




STERILE

Length	Titanium	Stainless Steel
20	TI-148.020	148.020
25	TI-148.025	148.025
30	TI-148.030	148.030
35	TI-148.035	148.035
40	TI-148.040	148.040
45	TI-148.045	148.045
50	TI-148.050	148.050
55	TI-148.055	148.055
60	TI-148.060	148.060
65	TI-148.065	148.065
70	TI-148.070	148.070
75	TI-148.075	148.075
80	TI-148.080	148.080
85	TI-148.085	148.085
90	TI-148.090	148.090
95	TI-148.095	148.095
100	TI-148.100	148.100
105	TI-148.105	148.105
110	TI-148.110	148.110
* 115	TI-148.115	148.115
* 120	TI-148.120	148.120
* 125	TI-148.125	148.125
* 130	TI-148.130	148.130

Washer for Large Fragment Screw



Titanium	Stainless Steel
TI-100.120	100.120




STERILE


Titanium	Stainless Steel
TI-100.120-S	100.120-S

* Sizes not available in Screw Caddy



ST-005-7.3 : Screw Caddy for 7.3mm Cannulated Cancellous Screw (Tray without Screw)



The background of the page features a light gray, wavy, topographic-like pattern that flows from the top and bottom edges towards the center, framing the central text.

INSTRUMENTS




GW-2.8-250T

Guide Wire with Threaded Tip, Ø2.8mm x Thread Length
10mm x Length 250mm



GW-2.8-250

Guide Wire, Ø2.8mm x Length 250mm



1532-01

Drill Bit Cannulated Plain Shank/Jacob Chuck End, Ø5.0mm x Length 250mm



1532-0S

Stopper for Drill Bit Ø5.0mm



1539-01

Allen Key, Hex 3mm, for Ø7.3mm Screw



3436-014

Direct Measuring Device for Cannulated Cancellous Screw





1538-01

Cannulated Countersink with T-Handle for Ø7.3mm Screw



1543-01

Wire Sleeve Ø10.0/2.9mm



1542-01

Drill Guide for Drill Bit Ø5.0mm



1541-01

Guide Sleeve Ø5.0/2.9mm



1537-01

Protection Sleeve for Ø10.0/2.9mm



1482-01

Wire Sleeve Bolt Key






1544-01

Wire Sleeve, Adjustable



1534-01

Cannulated Tap with T-Handle, for Ø7.3mm Screw




1533-01

Parallel Guide for Guide Wire Ø2.8mm



7-058-01

Screw Holding Forcep for Ø7.3mm Screw



1478-2.5

Stylet, Ø2.5mm



3433-03

Hexagonal Screwdriver, Cannulation Ø2.9mm, Hex 4.0mm






1540-01

Hexagonal Screwdriver, Hex 4.0mm for Ø7.3mm Screw



1540-0S

Holding Sleeve for Cannulated Hexagonal Screwdriver, Hex 4.0mm




1545-01

Extraction Screw with T-Handle, for Ø7.3mm Screw




ST-005-7.3

Screw Caddy for 7.3mm Cannulated Cancellous Screw



IS006-7.3T

Instrument Tray for 7.3mm Cannulated Screw Instrument Set




IS006-7.3B

Container For 7.3mm Cannulated Cancellous Screw Instrument Set





IS006-7.3	7.3mm Cannulated Cancellous Screw Instrument Set
-----------	--





ORDER INFORMATION

IS006-7.3

7.3mm Cannulated Cancellous Screw Instrument Set

Code	Set Consisting Of	Unit
GW-2.8-250T	Guide Wire with Threaded Tip, Ø2.8mm x Thread Length 10mm x Length 250mm	3
GW-2.8-250	Guide Wire, Ø2.8mm x Length 250mm	3
1532-01	Drill Bit Cannulated Plain Shank/Jacob Chuck End, Ø5.0mm x Length 250mm	2
1532-0S	Stopper for Drill Bit Ø5.0mm	1
1539-01	Allen Key, Hex 3mm, for Ø7.3mm Screw	1
3436-014	Direct Measuring Device for Cannulated Cancellous Screw	1
1538-01	Cannulated Countersink with T-Handle for Ø7.3mm Screw	1
1543-01	Wire Sleeve Ø10.0/2.9mm	1
1542-01	Drill Guide for Drill Bit Ø5.0mm	1
1541-01	Guide Sleeve Ø5.0/2.9mm	1
1537-01	Protection Sleeve for Ø10.0/2.9mm	1
1482-01	Wire Sleeve Bolt Key	1
1544-01	Wire Sleeve, Adjustable	1
1534-01	Cannulated Tap with T-Handle, for Ø7.3mm Screw	1
1533-01	Parallel Guide for Guide Wire Ø2.8mm	1
7-058-01	Screw Holding Forcep for Ø7.3mm Screw	1
1478-2.5	Stylet, Ø2.5mm	1
3433-03	Hexagonal Screwdriver, Cannulation Ø2.9mm, Hex 4.0mm	1
1540-01	Hexagonal Screwdriver, Hex 4.0mm for Ø7.3mm Screw	1
1540-0S	Holding Sleeve for Cannulated Hexagonal Screwdriver, Hex 4.0mm	1
1545-01	Extraction Screw with T-Handle, for Ø7.3mm Screw	1
ST-005-7.3	Screw Caddy for 7.3mm Cannulated Cancellous Screw	1
IS006-7.3T	Instrument Tray for 7.3mm Cannulated Screw Instrument Set	2
IS006-7.3B	Container For 7.3mm Cannulated Cancellous Screw Instrument Set	1

CANNULATED SCREW SYSTEM

The content of this catalogue is intended for healthcare professionals only. All product information, illustrations, descriptions, and specifications provided here are for general informational purposes and are subject to change without notice. Auxein Medical Pvt. Ltd. reserves the right to improve or discontinue products at any time without prior notification.

This document does not constitute medical advice. Surgeons and medical practitioners must exercise their own professional judgment when using any medical devices or instruments, including those provided by Auxein.

CLINICAL USE

The implants and instruments in this catalogue are intended for use by qualified orthopedic surgeons trained in fracture fixation procedures. Proper implant selection and surgical technique are essential for safe and effective use. Improper selection, misuse, or insufficient surgical experience may result in complications.

REGULATORY COMPLIANCE

US FDA 510(k)

CE Mark (Class IIb)

MDSAP (Medical Device Single Audit Program)

ISO 13485:2016 Quality Management System


CDSCO – Indian Medical Device Rules, 2017

EU MDR- European Union Medical Device Regulation, 2017

DISCLAIMER: Not all products in this catalogue are covered by European CE, US FDA, or other certifications. Please contact our Sales Team for product-specific certification details.

INTELLECTUAL PROPERTY

All trademarks, trade names, brand names, product names, and logos mentioned in this catalogue are the property of Auxein Medical Pvt. Ltd. Unauthorized use or reproduction is strictly prohibited.



SCAN TO KNOW MORE

INDIA
AUXEIN MEDICAL PVT.LTD.

PLOT NO. 168-169-170, PHASE-4,
KUNDLI INDUSTRIAL AREA,
HSIIDC, SECTOR-57,
SONEPAT-131028, HARYANA
TEL: +91 99106 43638
FAX: +91 86077 70197
EMAIL: INFO@AUXEIN.COM

MÉXICO
AUXEIN MÉXICO, S.A DE CV

TEPIC 39 INT 801, COLONIA
ROMA SUR,
ALCADÍA CUAUHTÉMOC, CDMX
MÉXICO, C.P 06760
TEL: +521 55 7261 0318
EMAIL: INFO@AUXEIN.MX

USA
AUXEIN INC.

1500 NW 89TH COURT,
SUITE 107-108
DORAL, FLORIDA 33172
TEL. +1 305 395 6062
E FAX: +1 305 395 62 62
EMAIL: USOFFICE@AUXEIN.COM

DUBAI UAE
AUXEIN MEDICAL

TRADING FZCO OUD METHA
OFFICES BUILDING B BLOCK,
1ST FLOOR, OFFICE NO.26,
NEAR HEALTHCARE CITY,
PO BOX-25572
TEL. +971 438 71333
EMAIL: INFO@AUXEIN.COM

WWW.AUXEIN.COM