


Auxein

TRAUMATOLOGY



NEXT-LEVEL CARE

ABOUT US

Auxein Medical is an integrated, research based, orthopaedic implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavour to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in India and serves customers in over 75 countries worldwide.



IMPACTING HEALTHCARE GLOBALLY

Our Reach & Influence

500+

Employees Worldwide

75+

Operations in Countries

150+

Distributors Worldwide

20^{M+}

Patients Cured

200+

FDA Approved Products

250+

Machines Worldwide

3000+

Orthopedic Solutions

Auxein is dedicated to advancing medical technology through innovative orthopedic solutions, including **Trauma, Arthroscopy, and Reconstructive Procedures**, with a commitment to improving patient outcomes worldwide. Our approach integrates advanced research, precision engineering, and high-quality manufacturing to develop reliable implants and surgical instruments that address the evolving needs of modern orthopedic practice.


OUR ACHIEVEMENTS



OVER THE YEARS

Disclaimer: Not all products in this catalogue are covered by European CE, US FDA, or other certifications. Please contact our Sales Team for product-specific certification details.






At Auxein, it is our goal to provide the best in class quality and value to every individual. Consistency is a virtue we pursue with a passion. From the procurement of raw materials to the manufacturing of our range of 1200+ orthopaedic implant products, we employ world-class systems and technology at par with the highest standards of globally recognized evaluation systems. The Traumatology System provides comprehensive solutions for fracture fixation in orthopedic trauma. It includes Mini Fragment, Small Fragment, and Large Fragment Plating Systems, along with Staples, Pins & Wires, offering versatile fixation options for the management of a wide range of fractures. Developed with input from global experts, the system supports effective and reliable fracture management.



INDEX



Mini Fragment	08-27
Small Fragment	28-57
Large Fragment	58-90
Angle Blade Plates	92-103
Staples, Pins & Wires	104-115

The background of the page features a light gray, wavy, topographic-like pattern that flows from the top left towards the bottom right, creating a sense of movement and depth. The lines are thin and closely spaced in some areas, while more spread out in others, giving it a fluid, organic appearance.

MINI FRAGMENT



2.0mm Mini "L" Plate



Left



Right



Left Direction

Right Direction

Holes	Length(mm)	Titanium	Stainless Steel	Titanium	Stainless Steel
2	18.60mm	TI-906.100L	906.100L	TI-906.100R	906.100R
4	30.20mm	TI-906.100-04L	906.100-04L	TI-906.100-04R	906.100-04R
6	41.80mm	TI-906.100-06L	906.100-06L	TI-906.100-06R	906.100-06R
8	53.40mm	TI-906.100-08L	906.100-08L	TI-906.100-08R	906.100-08R

STERILE

Left Direction

Right Direction

Holes	Length(mm)	Titanium	Stainless Steel	Titanium	Stainless Steel
2	18.60mm	TI-906.100L-S	906.100L-S	TI-906.100R-S	906.100R-S
4	30.20mm	TI-906.100-04L-S	906.100-04L-S	TI-906.100-04R-S	906.100-04R-S
6	41.80mm	TI-906.100-06L-S	906.100-06L-S	TI-906.100-06R-S	906.100-06R-S
8	53.40mm	TI-906.100-08L-S	906.100-08L-S	TI-906.100-08R-S	906.100-08R-S

2.0mm Mini Straight Plate

**STERILE**

Holes	Length	Titanium	Stainless Steel
2	12.40mm	TI-909.102	909.102
3	18.40mm	TI-909.103	909.103
4	24.40mm	TI-909.104	909.104
5	30.40mm	TI-909.105	909.105
6	36.40mm	TI-909.106	909.106
7	42.40mm	TI-909.107	909.107
8	48.40mm	TI-909.108	909.108

Holes	Length	Titanium	Stainless Steel
2	12.40mm	TI-909.102-S	909.102-S
3	18.40mm	TI-909.103-S	909.103-S
4	24.40mm	TI-909.104-S	909.104-S
5	30.40mm	TI-909.105-S	909.105-S
6	36.40mm	TI-909.106-S	909.106-S
7	42.40mm	TI-909.107-S	909.107-S
8	48.40mm	TI-909.108-S	909.108-S



2.0mm DCP Plate



Holes	Length	Titanium	Stainless Steel
4	22.00mm	TI-908.104	908.104
5	27.00mm	TI-908.105	908.105
6	32.00mm	TI-908.106	908.106
7	37.00mm	TI-908.107	908.107
8	42.00mm	TI-908.108	908.108

STERILE

Holes	Length	Titanium	Stainless Steel
4	22.00mm	TI-908.104-S	908.104-S
5	27.00mm	TI-908.105-S	908.105-S
6	32.00mm	TI-908.106-S	908.106-S
7	37.00mm	TI-908.107-S	908.107-S
8	42.00mm	TI-908.108-S	908.108-S

2.0mm Mini "T" Plate

Shaft with two holes & Head with two holes.



Holes	Length	Titanium	Stainless Steel
4	18.60mm	TI-905.100	905.100

STERILE

Holes	Length	Titanium	Stainless Steel
4	18.60mm	TI-905.100-S	905.100-S

2.7mm Mini "L" plate



Type	Length	Titanium	Stainless Steel
Left	32.80mm	TI-705.100L	705.100L
Right	32.80mm	TI-705.100R	705.100R

STERILE

Type	Length	Titanium	Stainless Steel
Left	32.80mm	TI-705.100L-S	705.100L-S
Right	32.80mm	TI-705.100R-S	705.100R-S

2.7mm Reconstruction Plate, Straight



Holes	Length	Titanium	Stainless Steel
2	16.20mm	TI-726.102	726.102
3	24.20mm	TI-726.103	726.103
4	32.20mm	TI-726.104	726.104
5	40.20mm	TI-726.105	726.105
6	48.20mm	TI-726.106	726.106
7	56.20mm	TI-726.107	726.107
8	64.20mm	TI-726.108	726.108
9	72.20mm	TI-726.109	726.109
10	80.20mm	TI-726.110	726.110
11	88.20mm	TI-726.111	726.111
12	96.20mm	TI-726.112	726.112
14	112.2mm	TI-726.114	726.114
16	128.2mm	TI-726.116	726.116
18	144.2mm	TI-726.118	726.118
20	160.2mm	TI-726.120	726.120
22	176.2mm	TI-726.122	726.122
24	192.2mm	TI-726.124	726.124

STERILE

Holes	Length	Titanium	Stainless Steel
2	16.20mm	TI-726.102-S	726.102-S
3	24.20mm	TI-726.103-S	726.103-S
4	32.20mm	TI-726.104-S	726.104-S
5	40.20mm	TI-726.105-S	726.105-S
6	48.20mm	TI-726.106-S	726.106-S
7	56.20mm	TI-726.107-S	726.107-S
8	64.20mm	TI-726.108-S	726.108-S
9	72.20mm	TI-726.109-S	726.109-S
10	80.20mm	TI-726.110-S	726.110-S
11	88.20mm	TI-726.111-S	726.111-S
12	96.20mm	TI-726.112-S	726.112-S
14	112.2mm	TI-726.114-S	726.114-S
16	128.2mm	TI-726.116-S	726.116-S
18	144.2mm	TI-726.118-S	726.118-S
20	160.2mm	TI-726.120-S	726.120-S
22	176.2mm	TI-726.122-S	726.122-S
24	192.2mm	TI-726.124-S	726.124-S

2.7mm "T" Plate



Length	Titanium	Stainless Steel
34.50mm	TI-707.100	707.100



STERILE

Length	Titanium	Stainless Steel
34.50mm	TI-707.100-S	707.100-S

2.7mm Dynamic Compression Plate



Holes	Length	Titanium	Stainless Steel
2	20.00mm	TI-719.102	719.102
3	28.00mm	TI-719.103	719.103
4	36.00mm	TI-719.104	719.104
5	44.00mm	TI-719.105	719.105
6	52.00mm	TI-719.106	719.106
7	60.00mm	TI-719.107	719.107
8	68.00mm	TI-719.108	719.108
9	76.00mm	TI-719.109	719.109
10	84.00mm	TI-719.110	719.110
12	100.0mm	TI-719.112	719.112

STERILE

Holes	Length	Titanium	Stainless Steel
2	20.00mm	TI-719.102-S	719.102-S
3	28.00mm	TI-719.103-S	719.103-S
4	36.00mm	TI-719.104-S	719.104-S
5	44.00mm	TI-719.105-S	719.105-S
6	52.00mm	TI-719.106-S	719.106-S
7	60.00mm	TI-719.107-S	719.107-S
8	68.00mm	TI-719.108-S	719.108-S
9	76.00mm	TI-719.109-S	719.109-S
10	84.00mm	TI-719.110-S	719.110-S
12	100.0mm	TI-719.112-S	719.112-S


2.7mm Quarter Tubular Plate



Holes	Length	Titanium	Stainless Steel
2	18.00mm	TI-701.102	701.102
3	26.00mm	TI-701.103	701.103
4	34.00mm	TI-701.104	701.104
5	42.00mm	TI-701.105	701.105
6	50.00mm	TI-701.106	701.106
7	58.00mm	TI-701.107	701.107
8	66.00mm	TI-701.108	701.108
9	74.00mm	TI-701.109	701.109
10	82.00mm	TI-701.110	701.110

STERILE

Holes	Length	Titanium	Stainless Steel
2	18.00mm	TI-701.102-S	701.102-S
3	26.00mm	TI-701.103-S	701.103-S
4	34.00mm	TI-701.104-S	701.104-S
5	42.00mm	TI-701.105-S	701.105-S
6	50.00mm	TI-701.106-S	701.106-S
7	58.00mm	TI-701.107-S	701.107-S
8	66.00mm	TI-701.108-S	701.108-S
9	74.00mm	TI-701.109-S	701.109-S
10	82.00mm	TI-701.110-S	701.110-S

The background of the page features a light gray, wavy, topographical-like pattern that flows from the top and bottom edges towards the center, creating a sense of movement and depth. The lines are thin and closely spaced in some areas, while more spread out in others.

SCREWS



1.5mm Cortical Screw (Hex Head)

Drill Bit: 1.1mm



Length	Titanium	Stainless Steel
6mm	TI-101.006	101.006
7mm	TI-101.007	101.007
8mm	TI-101.008	101.008
9mm	TI-101.009	101.009
10mm	TI-101.010	101.010
11mm	TI-101.011	101.011
12mm	TI-101.012	101.012
14mm	TI-101.014	101.014
16mm	TI-101.016	101.016
18mm	TI-101.018	101.018
20mm	TI-101.020	101.020



STERILE

Length	Titanium	Stainless Steel
6mm	TI-101.006-S	101.006-S
7mm	TI-101.007-S	101.007-S
8mm	TI-101.008-S	101.008-S
9mm	TI-101.009-S	101.009-S
10mm	TI-101.010-S	101.010-S
11mm	TI-101.011-S	101.011-S
12mm	TI-101.012-S	101.012-S
14mm	TI-101.014-S	101.014-S
16mm	TI-101.016-S	101.016-S
18mm	TI-101.018-S	101.018-S
20mm	TI-101.020-S	101.020-S

1.5mm Cortical Screw, Self -Tapping, (Hex Head)

Drill Bit: 1.1mm



Length	Titanium	Stainless Steel
6mm	TI-101.006	101.006
7mm	TI-101.007	101.007
8mm	TI-101.008	101.008
9mm	TI-101.009	101.009
10mm	TI-101.010	101.010
11mm	TI-101.011	101.011
12mm	TI-101.012	101.012
14mm	TI-101.014	101.014
16mm	TI-101.016	101.016
18mm	TI-101.018	101.018
20mm	TI-101.020	101.020



STERILE

Length	Titanium	Stainless Steel
6mm	TI-101.006-S	101.006-S
7mm	TI-101.007-S	101.007-S
8mm	TI-101.008-S	101.008-S
9mm	TI-101.009-S	101.009-S
10mm	TI-101.010-S	101.010-S
11mm	TI-101.011-S	101.011-S
12mm	TI-101.012-S	101.012-S
14mm	TI-101.014-S	101.014-S
16mm	TI-101.016-S	101.016-S
18mm	TI-101.018-S	101.018-S
20mm	TI-101.020-S	101.020-S



2.0mm Cortical Screw (Hex Head)

Drill Bit: 1.5mm



Length	Titanium	Stainless Steel
6mm	TI-102.006	102.006
7mm	TI-102.007	102.007
8mm	TI-102.008	102.008
9mm	TI-102.009	102.009
10mm	TI-102.010	102.010
12mm	TI-102.012	102.012
14mm	TI-102.014	102.014
16mm	TI-102.016	102.016
18mm	TI-102.018	102.018
20mm	TI-102.020	102.020



STERILE

Length	Titanium	Stainless Steel
6mm	TI-102.006-S	102.006-S
7mm	TI-102.007-S	102.007-S
8mm	TI-102.008-S	102.008-S
9mm	TI-102.009-S	102.009-S
10mm	TI-102.010-S	102.010-S
12mm	TI-102.012-S	102.012-S
14mm	TI-102.014-S	102.014-S
16mm	TI-102.016-S	102.016-S
18mm	TI-102.018-S	102.018-S
20mm	TI-102.020-S	102.020-S

2.0mm Cortical Screw, Self -Tapping, (Hex Head)

Drill Bit: 1.5mm



Length	Titanium	Stainless Steel
6mm	TI-102.306	102.306
7mm	TI-102.307	102.307
8mm	TI-102.308	102.308
9mm	TI-102.309	102.309
10mm	TI-102.310	102.310
12mm	TI-102.312	102.312
14mm	TI-102.314	102.314
16mm	TI-102.316	102.316
18mm	TI-102.318	102.318
20mm	TI-102.320	102.320



STERILE

Length	Titanium	Stainless Steel
6mm	TI-102.306-S	102.306-S
7mm	TI-102.307-S	102.307-S
8mm	TI-102.308-S	102.308-S
9mm	TI-102.309-S	102.309-S
10mm	TI-102.310-S	102.310-S
12mm	TI-102.312-S	102.312-S
14mm	TI-102.314-S	102.314-S
16mm	TI-102.316-S	102.316-S
18mm	TI-102.318-S	102.318-S
20mm	TI-102.320-S	102.320-S



2.7mm Cortical Screw(Hex Head)

Drill Bit: 2.0mm



Length	Titanium	Stainless Steel
6mm	TI-103.006	103.006
8mm	TI-103.008	103.008
10mm	TI-103.010	103.010
12mm	TI-103.012	103.012
14mm	TI-103.014	103.014
16mm	TI-103.016	103.016
18mm	TI-103.018	103.018
20mm	TI-103.020	103.020
22mm	TI-103.022	103.022
24mm	TI-103.024	103.024
26mm	TI-103.026	103.026
28mm	TI-103.028	103.028
30mm	TI-103.030	103.030



STERILE

Length	Titanium	Stainless Steel
6mm	TI-103.006-S	103.006-S
8mm	TI-103.008-S	103.008-S
10mm	TI-103.010-S	103.010-S
12mm	TI-103.012-S	103.012-S
14mm	TI-103.014-S	103.014-S
16mm	TI-103.016-S	103.016-S
18mm	TI-103.018-S	103.018-S
20mm	TI-103.020-S	103.020-S
22mm	TI-103.022-S	103.022-S
24mm	TI-103.024-S	103.024-S
26mm	TI-103.026-S	103.026-S
28mm	TI-103.028-S	103.028-S
30mm	TI-103.030-S	103.030-S

2.7mm Cortical Screw, Self -Tapping(Hex Head)

Drill Bit: 2.0mm



Length	Titanium	Stainless Steel
6mm	TI-103.206	103.206
8mm	TI-103.208	103.208
10mm	TI-103.210	103.210
12mm	TI-103.212	103.212
14mm	TI-103.214	103.214
16mm	TI-103.216	103.216
18mm	TI-103.218	103.218
20mm	TI-103.220	103.220
22mm	TI-103.222	103.222
24mm	TI-103.224	103.224
26mm	TI-103.226	103.226
28mm	TI-103.228	103.228
30mm	TI-103.230	103.230



STERILE

Length	Titanium	Stainless Steel
6mm	TI-103.206-S	103.206-S
8mm	TI-103.208-S	103.208-S
10mm	TI-103.210-S	103.210-S
12mm	TI-103.212-S	103.212-S
14mm	TI-103.214-S	103.214-S
16mm	TI-103.216-S	103.216-S
18mm	TI-103.218-S	103.218-S
20mm	TI-103.220-S	103.220-S
22mm	TI-103.222-S	103.222-S
24mm	TI-103.224-S	103.224-S
26mm	TI-103.226-S	103.226-S
28mm	TI-103.228-S	103.228-S
30mm	TI-103.230-S	103.230-S

Mini Fragment- Implant Box

7-044-32

Implant Box for 1.5mm Mini Fragment System



7-044-34


Implant Box for 2.0mm Mini Fragment System



7-044-35

Implant Box for 2.7mm Mini Fragment System



The background of the page features a light gray, wavy, topographic-like pattern that flows across the top and bottom edges, framing the central text.

INSTRUMENTS

Instruments for 1.5mm System

7-044-01

DRILL BIT, Ø1.1MM X LENGTH 60MM

**7-044-02**

DOUBLE DRILL GUIDE, Ø1.5/1.1MM

**7-044-03**

TAP Ø1.5MM, QUICK COUPLING



Instruments for 2.0mm System

7-044-05

DRILL BIT, Ø1.5MM X LENGTH 80MM

**7-044-06**

DOUBLE DRILL GUIDE, Ø2.0/1.5MM

**7-044-07**

TAP Ø2.0MM, QUICK COUPLING





Instruments for 2.7mm System

7-044-08

DRILL BIT, Ø2.0MM X LENGTH 80MM

**7-044-09**

DOUBLE DRILL GUIDE, Ø2.7/2.0MM FOR MINI FRAGMENT

**7-044-10**

TAP Ø2.7MM, QUICK COUPLING

**7-044-12**

COUNTERSINK, FOR SCREWS 2.7MM



General Instruments


7-044-13

T-HANDLE WITH QUICK COUPLING FOR MINI FRAGMENT



7-044-14


Mini Plate and Bone Holding Forcep

**7-044-15**

Mini Screw Holding Forcep

**7-044-16**

Mini Plate Bender (Pair)

**7-044-17**

Mini Periosteal Elevator, Curved, 3mm

**7-044-18**

Mini Periosteal Elevator, Curved, 5mm

**7-044-19**

Mini Reduction Forcep - Pointed



7-044-20

Mini Reduction Forcep - Serrated

**7-044-21**

Mini Hohmann Retractor, 6.0mm

**7-044-22**

Mini Hohmann Retractor, 8.5mm

**7-044-23**

Mini Hohmann Retractor, 15mm

**7-044-24**

Sharp Hook

**7-044-25**

Quick Coupling Handle, Straight



7-044-26

Depth Gauge measuring upto 50mm for Mini Fragment

**7-044-27**

Skin Hook, Single Pronged

**7-044-28**

Skin Hook, Double Pronged

**7-044-29**

Plate Cutter & Bending Forcep

**7-044-30**

Plate Holding Forcep

**7-044-31**

Plate Cutter



7-044-04

Screwdriver Shaft with Holding Sleeve, Hex 1.5mm, for 1.5mm Screw

**7-044-11**

Screwdriver Shaft with Holding Sleeve, Hex 1.5mm, for 2.0mm Screw

**7-044-36**

Screwdriver Shaft with Holding Sleeve, Hex 2.0mm, for 2.7mm Screw

**7-044-32**

Implant Box for 1.5mm Mini Fragment System

**7-044-34**

Implant Box for 2.0mm Mini Fragment System


**7-044-35**

Implant Box for 2.7mm Mini Fragment System



7-044-33

Container for Mini Fragment Instrument Set



ORDER INFORMATION

7-044**Mini Fragment Instrument Set**



ORDER INFORMATION

7-044**Mini Fragment Instrument Set**

Code	Set Consisting Of	Unit
	Instruments for 1.5mm System	
7-044-01	Drill Bit, Ø1.1mm x Length 60mm	3
7-044-02	Double Drill Guide, Ø1.5/1.1mm	1
7-044-03	Tap Ø1.5mm, Quick Coupling	1
	Instruments for 2.0mm System	
7-044-05	Drill Bit, Ø1.5mm x Length 80mm	3
7-044-06	Double Drill Guide, Ø2.0/1.5mm	1
7-044-07	Tap Ø2.0mm, Quick Coupling	1
	Instruments for 2.7mm System	
7-044-08	Drill Bit, Ø2.0mm x Length 80mm	3
7-044-09	Double Drill Guide, Ø2.7/2.0mm for Mini Fragment	1
7-044-10	Tap Ø2.7mm, Quick Coupling	1
7-044-12	Countersink, for Screws 2.7mm	1
	General Instrument	
7-044-13	T-Handle with Quick Coupling for Mini Fragment	1
7-044-14	Mini Plate and Bone Holding Forcep	1
7-044-15	Mini Screw Holding Forcep	1
7-044-16	Mini Plate Bender (Pair)	1
7-044-17	Mini Periosteal Elevator, Curved, 3mm	1
7-044-18	Mini Periosteal Elevator, Curved, 5mm	1
7-044-19	Mini Reduction Forcep - Pointed	1
7-044-20	Mini Reduction Forcep - Serrated	1
7-044-21	Mini Hohmann Retractor, 6.0mm	1
7-044-22	Mini Hohmann Retractor, 8.5mm	1
7-044-23	Mini Hohmann Retractor, 15mm	1
7-044-24	Sharp Hook	1
7-044-25	Quick Coupling Handle, Straight	1
7-044-26	Depth Gauge measuring upto 50mm for Mini Fragment	1
7-044-27	Skin Hook, Single Pronged	1
7-044-28	Skin Hook, Double Pronged	1
7-044-29	Plate Cutter & Bending Forcep	1
7-044-30	Plate Holding Forcep	1
7-044-31	Plate Cutter	1
7-044-04	Screwdriver Shaft with Holding Sleeve, Hex 1.5mm, for 1.5mm Screw	1
7-044-11	Screwdriver Shaft with Holding Sleeve, Hex 1.5mm, for 2.0mm Screw	1
7-044-36	Screwdriver Shaft with Holding Sleeve, Hex 2.0mm, for 2.7mm Screw	1
7-044-32	Implant Box for 1.5mm Mini Fragment System	1
7-044-34	Implant Box for 2.0mm Mini Fragment System	1
7-044-35	Implant Box for 2.7mm Mini Fragment System	1
7-044-33	Container for Mini Fragment Instrument Set	1

The background of the page features a light gray, wavy, topographical-like pattern that flows from the top left towards the bottom right, creating a sense of movement and depth. The lines are thin and closely spaced in some areas, while more spread out in others.

SMALL FRAGMENT



3.5mm One Third Tubular Plate

Used as a buttress or neutralization plate for fixation of fractures of the fibula, metatarsals, and metacarpals. Compatible with 3.5 mm and 4.0 mm screws, allowing stable fixation and compression when required.



Holes	Titanium	Stainless Steel
2	TI-702.102	702.102
3	TI-702.103	702.103
4	TI-702.104	702.104
5	TI-702.105	702.105
6	TI-702.106	702.106
7	TI-702.107	702.107
8	TI-702.108	702.108
9	TI-702.109	702.109
10	TI-702.110	702.110
11	TI-702.111	702.111
12	TI-702.112	702.112

STERILE

Holes	Titanium	Stainless Steel
2	TI-702.102-S	702.102-S
3	TI-702.103-S	702.103-S
4	TI-702.104-S	702.104-S
5	TI-702.105-S	702.105-S
6	TI-702.106-S	702.106-S
7	TI-702.107-S	702.107-S
8	TI-702.108-S	702.108-S
9	TI-702.109-S	702.109-S
10	TI-702.110-S	702.110-S
11	TI-702.111-S	702.111-S
12	TI-702.112-S	702.112-S



3.5mm Small Dynamic Compression Plate

For Ulna and Radius, fixation with 3.5mm & 4.0mm Screws



Holes	Titanium	Stainless Steel
2	TI-720.102	720.102
3	TI-720.103	720.103
4	TI-720.104	720.104
5	TI-720.105	720.105
6	TI-720.106	720.106
7	TI-720.107	720.107
8	TI-720.108	720.108
9	TI-720.109	720.109
10	TI-720.110	720.110
11	TI-720.111	720.111
12	TI-720.112	720.112
13	TI-720.113	720.113
14	TI-720.114	720.114

STERILE

Holes	Titanium	Stainless Steel
2	TI-720.102-S	720.102-S
3	TI-720.103-S	720.103-S
4	TI-720.104-S	720.104-S
5	TI-720.105-S	720.105-S
6	TI-720.106-S	720.106-S
7	TI-720.107-S	720.107-S
8	TI-720.108-S	720.108-S
9	TI-720.109-S	720.109-S
10	TI-720.110-S	720.110-S
11	TI-720.111-S	720.111-S
12	TI-720.112-S	720.112-S
13	TI-720.113-S	720.113-S
14	TI-720.114-S	720.114-S





3.5mm One-Third Tubular Plate with Collar

3.5mm One-Third Tubular Plate with Collar Used as a buttress or neutralization plate for fixation of fractures of the fibula, metatarsals, and metacarpals. Compatible with 3.5 mm and 4.0 mm screws, allowing stable fixation and compression when required.



Holes	Titanium	Stainless Steel
2	TI-703.102	703.102
3	TI-703.103	703.103
4	TI-703.104	703.104
5	TI-703.105	703.105
6	TI-703.106	703.106
7	TI-703.107	703.107
8	TI-703.108	703.108
9	TI-703.109	703.109
10	TI-703.110	703.110
11	TI-703.111	703.111
12	TI-703.112	703.112

STERILE

Holes	Titanium	Stainless Steel
2	TI-703.102-S	703.102-S
3	TI-703.103-S	703.103-S
4	TI-703.104-S	703.104-S
5	TI-703.105-S	703.105-S
6	TI-703.106-S	703.106-S
7	TI-703.107-S	703.107-S
8	TI-703.108-S	703.108-S
9	TI-703.109-S	703.109-S
10	TI-703.110-S	703.110-S
11	TI-703.111-S	703.111-S
12	TI-703.112-S	703.112-S



3.5mm Small Limited Contact Dynamic Compression Plate

For Ulna and Radius, fixation with 3.5mm & 4.0mm Screws



Holes	Titanium	Stainless Steel
3	TI-721.103	721.103
4	TI-721.104	721.104
5	TI-721.105	721.105
6	TI-721.106	721.106
7	TI-721.107	721.107
8	TI-721.108	721.108
9	TI-721.109	721.109
10	TI-721.110	721.110
11	TI-721.111	721.111
12	TI-721.112	721.112
13	TI-721.113	721.113
14	TI-721.114	721.114

STERILE

Holes	Titanium	Stainless Steel
3	TI-721.103-S	721.103-S
4	TI-721.104-S	721.104-S
5	TI-721.105-S	721.105-S
6	TI-721.106-S	721.106-S
7	TI-721.107-S	721.107-S
8	TI-721.108-S	721.108-S
9	TI-721.109-S	721.109-S
10	TI-721.110-S	721.110-S
11	TI-721.111-S	721.111-S
12	TI-721.112-S	721.112-S
13	TI-721.113-S	721.113-S
14	TI-721.114-S	721.114-S





3.5mm Small "T" Plate, Right Angled, (3 Head Holes)

For Distal Radius, Can be used for Left & Right Radius. Fixtation with 3.5mm & 4.0mm Screws.



Holes	Titanium	Stainless Steel
3	TI-708.103	708.103
4	TI-708.104	708.104
5	TI-708.105	708.105
6	TI-708.106	708.106
7	TI-708.107	708.107
8	TI-708.108	708.108

STERILE

Holes	Titanium	Stainless Steel
3	TI-708.103-S	708.103-S
4	TI-708.104-S	708.104-S
5	TI-708.105-S	708.105-S
6	TI-708.106-S	708.106-S
7	TI-708.107-S	708.107-S
8	TI-708.108-S	708.108-S



3.5mm Small "T" Plate, Right Angled, (4 Head Holes)

For the Left & Right Radius.
Fixtation with 3.5mm & 4.0mm Screws.



Holes	Titanium	Stainless Steel
3	TI-710.103	710.103
4	TI-710.104	710.104
5	TI-710.105	710.105
6	TI-710.106	710.106
7	TI-710.107	710.107
8	TI-710.108	710.108

STERILE

Holes	Titanium	Stainless Steel
3	TI-710.103-S	710.103-S
4	TI-710.104-S	710.104-S
5	TI-710.105-S	710.105-S
6	TI-710.106-S	710.106-S
7	TI-710.107-S	710.107-S
8	TI-710.108-S	710.108-S



3.5mm Hook Plate



Holes	Titanium	Stainless Steel
4	TI-729.104	729.104
5	TI-729.105	729.105
6	TI-729.106	729.106
7	TI-729.107	729.107
8	TI-729.108	729.108
10	TI-729.110	729.110

STERILE

Holes	Titanium	Stainless Steel
4	TI-729.104-S	729.104-S
5	TI-729.105-S	729.105-S
6	TI-729.106-S	729.106-S
7	TI-729.107-S	729.107-S
8	TI-729.108-S	729.108-S
10	TI-729.110-S	729.110-S





3.5mm Small "T" Plate, Oblique Angled


For Distal Radius, Can be used for Left & Right Radius.
Fixtation with 3.5mm & 4.0mm Screws



Left Direction			Right Direction	
Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
3	TI-711.103L	711.103L	TI-711.103R	711.103R
4	TI-711.104L	711.104L	TI-711.104R	711.104R
5	TI-711.105L	711.105L	TI-711.105R	711.105R
6	TI-711.106L	711.106L	TI-711.106R	711.106R
7	TI-711.107L	711.107L	TI-711.107R	711.107R
8	TI-711.108L	711.108L	TI-711.108R	711.108R

STERILE

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
3	TI-711.103L-S	711.103L-S	TI-711.103R-S	711.103R-S
4	TI-711.104L-S	711.104L-S	TI-711.104R-S	711.104R-S
5	TI-711.105L-S	711.105L-S	TI-711.105R-S	711.105R-S
6	TI-711.106L-S	711.106L-S	TI-711.106R-S	711.106R-S
7	TI-711.107L-S	711.107L-S	TI-711.107R-S	711.107R-S
8	TI-711.108L-S	711.108L-S	TI-711.108R-S	711.108R-S



3.5mm Pelvic Reconstruction Plate, Straight



STERILE

Holes	Titanium	Stainless Steel	Holes	Titanium	Stainless Steel
5	TI-727.105S	727.105S	5	TI-727.105S-S	727.105S-S
6	TI-727.106S	727.106S	6	TI-727.106S-S	727.106S-S
7	TI-727.107S	727.107S	7	TI-727.107S-S	727.107S-S
8	TI-727.108S	727.108S	8	TI-727.108S-S	727.108S-S
9	TI-727.109S	727.109S	9	TI-727.109S-S	727.109S-S
10	TI-727.110S	727.110S	10	TI-727.110S-S	727.110S-S
11	TI-727.111S	727.111S	11	TI-727.111S-S	727.111S-S
12	TI-727.112S	727.112S	12	TI-727.112S-S	727.112S-S
13	TI-727.113S	727.113S	13	TI-727.113S-S	727.113S-S
14	TI-727.114S	727.114S	14	TI-727.114S-S	727.114S-S
15	TI-727.115S	727.115S	15	TI-727.115S-S	727.115S-S
16	TI-727.116S	727.116S	16	TI-727.116S-S	727.116S-S
17	TI-727.117S	727.117S	17	TI-727.117S-S	727.117S-S
18	TI-727.118S	727.118S	18	TI-727.118S-S	727.118S-S
19	TI-727.119S	727.119S	19	TI-727.119S-S	727.119S-S
20	TI-727.120S	727.120S	20	TI-727.120S-S	727.120S-S
21	TI-727.121S	727.121S	21	TI-727.121S-S	727.121S-S
22	TI-727.122S	727.122S	22	TI-727.122S-S	727.122S-S



3.5mm Pelvic Reconstruction Plate, Curved



Holes	Titanium	Stainless Steel
5	TI-727.105C	727.105C
6	TI-727.106C	727.106C
7	TI-727.107C	727.107C
8	TI-727.108C	727.108C
9	TI-727.109C	727.109C
10	TI-727.110C	727.110C
11	TI-727.111C	727.111C
12	TI-727.112C	727.112C
13	TI-727.113C	727.113C
14	TI-727.114C	727.114C
16	TI-727.116C	727.116C
18	TI-727.118C	727.118C

STERILE

Holes	Titanium	Stainless Steel
5	TI-727.105C-S	727.105C-S
6	TI-727.106C-S	727.106C-S
7	TI-727.107C-S	727.107C-S
8	TI-727.108C-S	727.108C-S
9	TI-727.109C-S	727.109C-S
10	TI-727.110C-S	727.110C-S
11	TI-727.111C-S	727.111C-S
12	TI-727.112C-S	727.112C-S
13	TI-727.113C-S	727.113C-S
14	TI-727.114C-S	727.114C-S
16	TI-727.116C-S	727.116C-S
18	TI-727.118C-S	727.118C-S

3.5mm Semi Tubular Plate



Holes	Titanium	Stainless Steel
2	TI-704.102	704.102
3	TI-704.103	704.103
4	TI-704.104	704.104
5	TI-704.105	704.105
6	TI-704.106	704.106
7	TI-704.107	704.107
8	TI-704.108	704.108
9	TI-704.109	704.109
10	TI-704.110	704.110
12	TI-704.112	704.112

STERILE

Holes	Titanium	Stainless Steel
2	TI-704.102-S	704.102-S
3	TI-704.103-S	704.103-S
4	TI-704.104-S	704.104-S
5	TI-704.105-S	704.105-S
6	TI-704.106-S	704.106-S
7	TI-704.107-S	704.107-S
8	TI-704.108-S	704.108-S
9	TI-704.109-S	704.109-S
10	TI-704.110-S	704.110-S
12	TI-704.112-S	704.112-S






3.5mm Calcaneal Plate



Size	Titanium	Stainless Steel
60mm	TI-717.16	717.16
70mm	TI-717.17	717.17



STERILE

Size	Titanium	Stainless Steel
60mm	TI-717.16-S	717.16-S
70mm	TI-717.17-S	717.17-S

Cervical Spine H-Plate



Holes	Holes Spacing	Titanium	Stainless Steel
5	16mm	TI-912-01	912-01
5	18mm	TI-912-02	912-02
5	21mm	TI-912-03	912-03
8	16+18mm	TI-912-04	912-04
8	18+21mm	TI-912-05	912-05
11	16+18+18mm	TI-912-06	912-06



STERILE

Holes	Holes Spacing	Titanium	Stainless Steel
5	16mm	TI-912-01	912-01
5	18mm	TI-912-02	912-02
5	21mm	TI-912-03	912-03
8	16+18mm	TI-912-04	912-04
8	18+21mm	TI-912-05	912-05
11	16+18+18mm	TI-912-06	912-06

3.5mm Cloverleaf Plate



Holes	Titanium	Stainless Steel
3	TI-715.103	715.103
4	TI-715.104	715.104
5	TI-715.105	715.105
6	TI-715.106	715.106
7	TI-715.107	715.107
8	TI-715.108	715.108
9	TI-715.109	715.109
10	TI-715.110	715.110
11	TI-715.111	715.111
12	TI-715.112	715.112
13	TI-715.113	715.113
14	TI-715.114	715.114



STERILE

Holes	Titanium	Stainless Steel
3	TI-715.103-S	715.103-S
4	TI-715.104-S	715.104-S
5	TI-715.105-S	715.105-S
6	TI-715.106-S	715.106-S
7	TI-715.107-S	715.107-S
8	TI-715.108-S	715.108-S
9	TI-715.109-S	715.109-S
10	TI-715.110-S	715.110-S
11	TI-715.111-S	715.111-S
12	TI-715.112-S	715.112-S
13	TI-715.113-S	715.113-S
14	TI-715.114-S	715.114-S



3.5mm Clavicle Hook Plate



Left Direction				Right Direction	
Holes	Hook Depth	Titanium	Stainless Steel	Titanium	Stainless Steel
4	12mm	TI-751.4.12L	751.4.12L	TI-751.4.12R	751.4.12R
4	15mm	TI-751.4.15L	751.4.15L	TI-751.4.15R	751.4.15R
4	18mm	TI-751.4.18L	751.4.18L	TI-751.4.18R	751.4.18R
5	12mm	TI-751.5.12L	751.5.12L	TI-751.5.12R	751.5.12R
5	15mm	TI-751.5.15L	751.5.15L	TI-751.5.15R	751.5.15R
5	18mm	TI-751.5.18L	751.5.18L	TI-751.5.18R	751.5.18R
6	12mm	TI-751.6.12L	751.6.12L	TI-751.6.12R	751.6.12R
6	15mm	TI-751.6.15L	751.6.15L	TI-751.6.15R	751.6.15R
6	18mm	TI-751.6.18L	751.6.18L	TI-751.6.18R	751.6.18R
7	12mm	TI-751.7.12L	751.7.12L	TI-751.7.12R	751.7.12R
7	15mm	TI-751.7.15L	751.7.15L	TI-751.7.15R	751.7.15R
7	18mm	TI-751.7.18L	751.7.18L	TI-751.7.18R	751.7.18R

STERILE

Left Direction				Right Direction	
Holes	Hook Depth	Titanium	Stainless Steel	Titanium	Stainless Steel
4	12mm	TI-751.4.12L-S	751.4.12L-S	TI-751.4.12R-S	751.4.12R-S
4	15mm	TI-751.4.15L-S	751.4.15L-S	TI-751.4.15R-S	751.4.15R-S
4	18mm	TI-751.4.18L-S	751.4.18L-S	TI-751.4.18R-S	751.4.18R-S
5	12mm	TI-751.5.12L-S	751.5.12L-S	TI-751.5.12R-S	751.5.12R-S
5	15mm	TI-751.5.15L-S	751.5.15L-S	TI-751.5.15R-S	751.5.15R-S
5	18mm	TI-751.5.18L-S	751.5.18L-S	TI-751.5.18R-S	751.5.18R-S
6	12mm	TI-751.6.12L-S	751.6.12L-S	TI-751.6.12R-S	751.6.12R-S
6	15mm	TI-751.6.15L-S	751.6.15L-S	TI-751.6.15R-S	751.6.15R-S
6	18mm	TI-751.6.18L-S	751.6.18L-S	TI-751.6.18R-S	751.6.18R-S
7	12mm	TI-751.7.12L-S	751.7.12L-S	TI-751.7.12R-S	751.7.12R-S
7	15mm	TI-751.7.15L-S	751.7.15L-S	TI-751.7.15R-S	751.7.15R-S
7	18mm	TI-751.7.18L-S	751.7.18L-S	TI-751.7.18R-S	751.7.18R-S

3.5mm Y-Reconstruction Plate



STERILE

Holes	Size	Titanium	Stainless Steel
6	Small	TI-716.106S	716.106S
9	Medium	TI-716.109M	716.109M
12	Large	TI-716.112L	716.112L

Holes	Size	Titanium	Stainless Steel
6	Small	TI-716.106S-S	716.106S-S
9	Medium	TI-716.109M-S	716.109M-S
12	Large	TI-716.112L-S	716.112L-S

Lecestre Plate - Modified, 2 Head Holes



Left Direction

Right Direction

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
5	TI-800.005L	800.005L	TI-800.005R	800.005R
7	TI-800.007L	800.007L	TI-800.007R	800.007R
9	TI-800.009L	800.009L	TI-800.009R	800.009R
10	TI-800.010L	800.010L	TI-800.010R	800.010R
11	TI-800.011L	800.011L	TI-800.011R	800.011R
13	TI-800.013L	800.013L	TI-800.013R	800.013R

STERILE

Left Direction

Right Direction

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
5	TI-800.005L-S	800.005L-S	TI-800.005R-S	800.005R-S
7	TI-800.007L-S	800.007L-S	TI-800.007R-S	800.007R-S
9	TI-800.009L-S	800.009L-S	TI-800.009R-S	800.009R-S
10	TI-800.010L-S	800.010L-S	TI-800.010R-S	800.010R-S
11	TI-800.011L-S	800.011L-S	TI-800.011R-S	800.011R-S
13	TI-800.013L-S	800.013L-S	TI-800.013R-S	800.013R-S

Proximal Ulna Plate



STERILE

Holes	Size	Titanium	Stainless Steel
3	Left	TI-CP05.03L	CP05.03L
3	Right	TI-CP05.03R	CP05.03R
4	Left	TI-CP05.04L	CP05.04L
4	Right	TI-CP05.04R	CP05.04R
5	Left	TI-CP05.05L	CP05.05L
5	Right	TI-CP05.05R	CP05.05R
6	Left	TI-CP05.06L	CP05.06L
6	Right	TI-CP05.06R	CP05.06R

Holes	Size	Titanium	Stainless Steel
3	Left	TI-CP05.03L-S	CP05.03L-S
3	Right	TI-CP05.03R-S	CP05.03R-S
4	Left	TI-CP05.04L-S	CP05.04L-S
4	Right	TI-CP05.04R-S	CP05.04R-S
5	Left	TI-CP05.05L-S	CP05.05L-S
5	Right	TI-CP05.05R-S	CP05.05R-S
6	Left	TI-CP05.06L-S	CP05.06L-S
6	Right	TI-CP05.06R-S	CP05.06R-S

Proximal Humeral Plate



Holes	Titanium	Stainless Steel
4	TI-CP03.04	CP03.04
5	TI-CP03.05	CP03.05
6	TI-CP03.06	CP03.06
7	TI-CP03.07	CP03.07
8	TI-CP03.08	CP03.08
9	TI-CP03.09	CP03.09
10	TI-CP03.10	CP03.10

STERILE

Holes	Titanium	Stainless Steel
4	TI-CP03.04-S	CP03.04-S
5	TI-CP03.05-S	CP03.05-S
6	TI-CP03.06-S	CP03.06-S
7	TI-CP03.07-S	CP03.07-S
8	TI-CP03.08-S	CP03.08-S
9	TI-CP03.09-S	CP03.09-S
10	TI-CP03.10-S	CP03.10-S

Distal Radius Plate



Holes	Titanium	Stainless Steel
2	TI-CP06.02	CP06.02
3	TI-CP06.03	CP06.03
4	TI-CP06.04	CP06.04
5	TI-CP06.05	CP06.05
6	TI-CP06.06	CP06.06

STERILE

Holes	Titanium	Stainless Steel
2	TI-CP06.02-S	CP06.02-S
3	TI-CP06.03-S	CP06.03-S
4	TI-CP06.04-S	CP06.04-S
5	TI-CP06.05-S	CP06.05-S
6	TI-CP06.06-S	CP06.06-S

Y-Fork Plate



Holes	Size	Titanium	Stainless Steel
4L	Left	TI-CP04.04L	CP04.04L
4R	Right	TI-CP04.04R	CP04.04R
5L	Left	TI-CP04.05L	CP04.05L
5R	Right	TI-CP04.05R	CP04.05R
6L	Left	TI-CP04.06L	CP04.06L
6R	Right	TI-CP04.06R	CP04.06R

STERILE

Holes	Size	Titanium	Stainless Steel
4L	Left	TI-CP04.04L-S	CP04.04L-S
4R	Right	TI-CP04.04R-S	CP04.04R-S
5L	Left	TI-CP04.05L-S	CP04.05L-S
5R	Right	TI-CP04.05R-S	CP04.05R-S
6L	Left	TI-CP04.06L-S	CP04.06L-S
6R	Right	TI-CP04.06R-S	CP04.06R-S



Distal Fibular Plate



Left Direction

Right Direction

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
3	TI-DF01.03L	DF01.03L	TI-DF01.03R	DF01.03R
4	TI-DF01.04L	DF01.04L	TI-DF01.04R	DF01.04R
5	TI-DF01.05L	DF01.05L	TI-DF01.05R	DF01.05R
6	TI-DF01.06L	DF01.06L	TI-DF01.06R	DF01.06R
7	TI-DF01.07L	DF01.07L	TI-DF01.07R	DF01.07R
8	TI-DF01.08L	DF01.08L	TI-DF01.08R	DF01.08R
9	TI-DF01.09L	DF01.09L	TI-DF01.09R	DF01.09R

STERILE

Left Direction

Right Direction

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
5	TI-DF01.03L-S	DF01.03L-S	TI-DF01.03R-S	DF01.03R-S
7	TI-DF01.04L-S	DF01.04L-S	TI-DF01.04R-S	DF01.04R-S
9	TI-DF01.05L-S	DF01.05L-S	TI-DF01.05R-S	DF01.05R-S
10	TI-DF01.06L-S	DF01.06L-S	TI-DF01.06R-S	DF01.06R-S
11	TI-DF01.07L-S	DF01.07L-S	TI-DF01.07R-S	DF01.07R-S
13	TI-DF01.08L-S	DF01.08L-S	TI-DF01.08R-S	DF01.08R-S
15	TI-DF01.09L-S	DF01.09L-S	TI-DF01.09R-S	DF01.09R-S



Distal Lateral Humeral Plate



Left Direction

Right Direction

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
3	TI-CP02.03R	CP02.03R	TI-CP02.03L	CP02.03L
4	TI-CP02.04R	CP02.04R	TI-CP02.04L	CP02.04L
5	TI-CP02.05R	CP02.05R	TI-CP02.05L	CP02.05L
6	TI-CP02.06R	CP02.06R	TI-CP02.06L	CP02.06L

STERILE

Left Direction


Right Direction

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
3	TI-CP02.03R-S	CP02.03R-S	TI-CP02.03L-S	CP02.03L-S
4	TI-CP02.04R-S	CP02.04R-S	TI-CP02.04L-S	CP02.04L-S
5	TI-CP02.05R-S	CP02.05R-S	TI-CP02.05L-S	CP02.05L-S
6	TI-CP02.06R-S	CP02.06R-S	TI-CP02.06L-S	CP02.06L-S



ORDER INFORMATION

ST-001NW**Screw Caddy For Small fragment Instrument Set**

The background of the page features a light gray, wavy, topographical-like pattern that flows from the top and bottom edges towards the center, framing the central text.

SCREWS



3.5mm Cortical Screw (Hex Head)

Drill Bit : 2.5mm



Length	Titanium	Stainless Steel
10mm	TI-104.010	104.010
12mm	TI-104.012	104.012
14mm	TI-104.014	104.014
16mm	TI-104.016	104.016
18mm	TI-104.018	104.018
20mm	TI-104.020	104.020
22mm	TI-104.022	104.022
24mm	TI-104.024	104.024
26mm	TI-104.026	104.026
28mm	TI-104.028	104.028
30mm	TI-104.030	104.030
32mm	TI-104.032	104.032
34mm	TI-104.034	104.034
36mm	TI-104.036	104.036
38mm	TI-104.038	104.038
40mm	TI-104.040	104.040
42mm	TI-104.042	104.042
44mm	TI-104.044	104.044
46mm	TI-104.046	104.046
48mm	TI-104.048	104.048
50mm	TI-104.050	104.050
55mm	TI-104.055	104.055
60mm	TI-104.060	104.060
65mm	TI-104.065	104.065
70mm	TI-104.070	104.070
75mm	TI-104.075	104.075
80mm	TI-104.080	104.080
85mm	TI-104.085	104.085
90mm	TI-104.090	104.090



STERILE

Length	Titanium	Stainless Steel
10mm	TI-104.010-S	104.010-S
12mm	TI-104.012-S	104.012-S
14mm	TI-104.014-S	104.014-S
16mm	TI-104.016-S	104.016-S
18mm	TI-104.018-S	104.018-S
20mm	TI-104.020-S	104.020-S
22mm	TI-104.022-S	104.022-S
24mm	TI-104.024-S	104.024-S
26mm	TI-104.026-S	104.026-S
28mm	TI-104.028-S	104.028-S
30mm	TI-104.030-S	104.030-S
32mm	TI-104.032-S	104.032-S
34mm	TI-104.034-S	104.034-S
36mm	TI-104.036-S	104.036-S
38mm	TI-104.038-S	104.038-S
40mm	TI-104.040-S	104.040-S
42mm	TI-104.042-S	104.042-S
44mm	TI-104.044-S	104.044-S
46mm	TI-104.046-S	104.046-S
48mm	TI-104.048-S	104.048-S
50mm	TI-104.050-S	104.050-S
55mm	TI-104.055-S	104.055-S
60mm	TI-104.060-S	104.060-S
65mm	TI-104.065-S	104.065-S
70mm	TI-104.070-S	104.070-S
75mm	TI-104.075-S	104.075-S
80mm	TI-104.080-S	104.080-S
85mm	TI-104.085-S	104.085-S
90mm	TI-104.090-S	104.090-S



3.5mm Cortical Screw, Self-Tapping (Hex Head)

Drill Bit : 2.5mm



Length	Titanium	Stainless Steel
10mm	TI-104.210	104.210
12mm	TI-104.212	104.212
14mm	TI-104.214	104.214
16mm	TI-104.216	104.216
18mm	TI-104.218	104.218
20mm	TI-104.220	104.220
22mm	TI-104.222	104.222
24mm	TI-104.224	104.224
26mm	TI-104.226	104.226
28mm	TI-104.228	104.228
30mm	TI-104.230	104.230
32mm	TI-104.232	104.232
34mm	TI-104.234	104.234
36mm	TI-104.236	104.236
38mm	TI-104.238	104.238
40mm	TI-104.240	104.240
42mm	TI-104.242	104.242
44mm	TI-104.244	104.244
46mm	TI-104.246	104.246
48mm	TI-104.248	104.248
50mm	TI-104.250	104.250
55mm	TI-104.255	104.255
60mm	TI-104.260	104.260
65mm	TI-104.265	104.265
70mm	TI-104.270	104.270
75mm	TI-104.275	104.275
80mm	TI-104.280	104.280
85mm	TI-104.285	104.285
90mm	TI-104.290	104.290



STERILE

Length	Titanium	Stainless Steel
10mm	TI-104.210-S	104.210-S
12mm	TI-104.212-S	104.212-S
14mm	TI-104.214-S	104.214-S
16mm	TI-104.216-S	104.216-S
18mm	TI-104.218-S	104.218-S
20mm	TI-104.220-S	104.220-S
22mm	TI-104.222-S	104.222-S
24mm	TI-104.224-S	104.224-S
26mm	TI-104.226-S	104.226-S
28mm	TI-104.228-S	104.228-S
30mm	TI-104.230-S	104.230-S
32mm	TI-104.232-S	104.232-S
34mm	TI-104.234-S	104.234-S
36mm	TI-104.236-S	104.236-S
38mm	TI-104.238-S	104.238-S
40mm	TI-104.240-S	104.240-S
42mm	TI-104.242-S	104.242-S
44mm	TI-104.244-S	104.244-S
46mm	TI-104.246-S	104.246-S
48mm	TI-104.248-S	104.248-S
50mm	TI-104.250-S	104.250-S
55mm	TI-104.255-S	104.255-S
60mm	TI-104.260-S	104.260-S
65mm	TI-104.265-S	104.265-S
70mm	TI-104.270-S	104.270-S
75mm	TI-104.275-S	104.275-S
80mm	TI-104.280-S	104.280-S
85mm	TI-104.285-S	104.285-S
90mm	TI-104.290-S	104.290-S



3.5mm Cancellous Screw, Short Thread

Drill Bit : 2.5mm



Length	Titanium	Stainless Steel
10mm	TI-108.010S	108.010S
12mm	TI-108.012S	108.012S
14mm	TI-108.014S	108.014S
16mm	TI-108.016S	108.016S
18mm	TI-108.018S	108.018S
20mm	TI-108.020S	108.020S
22mm	TI-108.022S	108.022S
24mm	TI-108.024S	108.024S
26mm	TI-108.026S	108.026S
28mm	TI-108.028S	108.028S
30mm	TI-108.030S	108.030S
32mm	TI-108.032S	108.032S
34mm	TI-108.034S	108.034S
36mm	TI-108.036S	108.036S
38mm	TI-108.038S	108.038S
40mm	TI-108.040S	108.040S
42mm	TI-108.042S	108.042S
44mm	TI-108.044S	108.044S
46mm	TI-108.046S	108.046S
48mm	TI-108.048S	108.048S
50mm	TI-108.050S	108.050S
52mm	TI-108.052S	108.052S
54mm	TI-108.054S	108.054S
56mm	TI-108.056S	108.056S
58mm	TI-108.058S	108.058S
60mm	TI-108.060S	108.060S



STERILE

Length	Titanium	Stainless Steel
10mm	TI-108.010S-S	108.010S-S
12mm	TI-108.012S-S	108.012S-S
14mm	TI-108.014S-S	108.014S-S
16mm	TI-108.016S-S	108.016S-S
18mm	TI-108.018S-S	108.018S-S
20mm	TI-108.020S-S	108.020S-S
22mm	TI-108.022S-S	108.022S-S
24mm	TI-108.024S-S	108.024S-S
26mm	TI-108.026S-S	108.026S-S
28mm	TI-108.028S-S	108.028S-S
30mm	TI-108.030S-S	108.030S-S
32mm	TI-108.032S-S	108.032S-S
34mm	TI-108.034S-S	108.034S-S
36mm	TI-108.036S-S	108.036S-S
38mm	TI-108.038S-S	108.038S-S
40mm	TI-108.040S-S	108.040S-S
42mm	TI-108.042S-S	108.042S-S
44mm	TI-108.044S-S	108.044S-S
46mm	TI-108.046S-S	108.046S-S
48mm	TI-108.048S-S	108.048S-S
50mm	TI-108.050S-S	108.050S-S
52mm	TI-108.052S-S	108.052S-S
54mm	TI-108.054S-S	108.054S-S
56mm	TI-108.056S-S	108.056S-S
58mm	TI-108.058S-S	108.058S-S
60mm	TI-108.060S-S	108.060S-S



3.5mm Cancellous Screw, Full Thread

Drill Bit : 2.5mm



Length	Titanium	Stainless Steel
10mm	TI-108.010	108.010
12mm	TI-108.012	108.012
14mm	TI-108.014	108.014
16mm	TI-108.016	108.016
18mm	TI-108.018	108.018
20mm	TI-108.020	108.020
22mm	TI-108.022	108.022
24mm	TI-108.024	108.024
26mm	TI-108.026	108.026
28mm	TI-108.028	108.028
30mm	TI-108.030	108.030
32mm	TI-108.032	108.032
34mm	TI-108.034	108.034
36mm	TI-108.036	108.036
38mm	TI-108.038	108.038
40mm	TI-108.040	108.040
42mm	TI-108.042	108.042
44mm	TI-108.044	108.044
46mm	TI-108.046	108.046
48mm	TI-108.048	108.048
50mm	TI-108.050	108.050
52mm	TI-108.052	108.052
54mm	TI-108.054	108.054
56mm	TI-108.056	108.056
58mm	TI-108.058	108.058
60mm	TI-108.060	108.060



STERILE

Length	Titanium	Stainless Steel
10mm	TI-108.010-S	108.010-S
12mm	TI-108.012-S	108.012-S
14mm	TI-108.014-S	108.014-S
16mm	TI-108.016-S	108.016-S
18mm	TI-108.018-S	108.018-S
20mm	TI-108.020-S	108.020-S
22mm	TI-108.022-S	108.022-S
24mm	TI-108.024-S	108.024-S
26mm	TI-108.026-S	108.026-S
28mm	TI-108.028-S	108.028-S
30mm	TI-108.030-S	108.030-S
32mm	TI-108.032-S	108.032-S
34mm	TI-108.034-S	108.034-S
36mm	TI-108.036-S	108.036-S
38mm	TI-108.038-S	108.038-S
40mm	TI-108.040-S	108.040-S
42mm	TI-108.042-S	108.042-S
44mm	TI-108.044-S	108.044-S
46mm	TI-108.046-S	108.046-S
48mm	TI-108.048-S	108.048-S
50mm	TI-108.050-S	108.050-S
52mm	TI-108.052-S	108.052-S
54mm	TI-108.054-S	108.054-S
56mm	TI-108.056-S	108.056-S
58mm	TI-108.058-S	108.058-S
60mm	TI-108.060-S	108.060-S



4.0mm Cancellous Screw, Short Thread

Drill Bit : 2.5mm



Length	Titanium	Stainless Steel
12mm	TI-109.012	109.012
14mm	TI-109.014	109.014
16mm	TI-109.016	109.016
18mm	TI-109.018	109.018
20mm	TI-109.020	109.020
22mm	TI-109.022	109.022
24mm	TI-109.024	109.024
26mm	TI-109.026	109.026
28mm	TI-109.028	109.028
30mm	TI-109.030	109.030
32mm	TI-109.032	109.032
34mm	TI-109.034	109.034
36mm	TI-109.036	109.036
38mm	TI-109.038	109.038
40mm	TI-109.040	109.040
45mm	TI-109.045	109.045
50mm	TI-109.050	109.050
55mm	TI-109.055	109.055
60mm	TI-109.060	109.060



STERILE

Length	Titanium	Stainless Steel
12mm	TI-109.012-S	109.012-S
14mm	TI-109.014-S	109.014-S
16mm	TI-109.016-S	109.016-S
18mm	TI-109.018-S	109.018-S
20mm	TI-109.020-S	109.020-S
22mm	TI-109.022-S	109.022-S
24mm	TI-109.024-S	109.024-S
26mm	TI-109.026-S	109.026-S
28mm	TI-109.028-S	109.028-S
30mm	TI-109.030-S	109.030-S
32mm	TI-109.032-S	109.032-S
34mm	TI-109.034-S	109.034-S
36mm	TI-109.036-S	109.036-S
38mm	TI-109.038-S	109.038-S
40mm	TI-109.040-S	109.040-S
45mm	TI-109.045-S	109.045-S
50mm	TI-109.050-S	109.050-S
55mm	TI-109.055-S	109.055-S
60mm	TI-109.060-S	109.060-S



4.0mm Cancellous Screw, Full Thread

Drill Bit : 2.5mm



Length	Titanium	Stainless Steel
12mm	TI-110.012	110.012
14mm	TI-110.014	110.014
16mm	TI-110.016	110.016
18mm	TI-110.018	110.018
20mm	TI-110.020	110.020
22mm	TI-110.022	110.022
24mm	TI-110.024	110.024
26mm	TI-110.026	110.026
28mm	TI-110.028	110.028
30mm	TI-110.030	110.030
32mm	TI-110.032	110.032
34mm	TI-110.034	110.034
36mm	TI-110.036	110.036
38mm	TI-110.038	110.038
40mm	TI-110.040	110.040
45mm	TI-110.045	110.045
50mm	TI-110.050	110.050
55mm	TI-110.055	110.055
60mm	TI-110.060	110.060



STERILE

Length	Titanium	Stainless Steel
12mm	TI-110.012-S	110.012-S
14mm	TI-110.014-S	110.014-S
16mm	TI-110.016-S	110.016-S
18mm	TI-110.018-S	110.018-S
20mm	TI-110.020-S	110.020-S
22mm	TI-110.022-S	110.022-S
24mm	TI-110.024-S	110.024-S
26mm	TI-110.026-S	110.026-S
28mm	TI-110.028-S	110.028-S
30mm	TI-110.030-S	110.030-S
32mm	TI-110.032-S	110.032-S
34mm	TI-110.034-S	110.034-S
36mm	TI-110.036-S	110.036-S
38mm	TI-110.038-S	110.038-S
40mm	TI-110.040-S	110.040-S
45mm	TI-110.045-S	110.045-S
50mm	TI-110.050-S	110.050-S
55mm	TI-110.055-S	110.055-S
60mm	TI-110.060-S	110.060-S



3.5mm Malleolar Screw

Drill Bit : 2.5mm



Length	Titanium	Stainless Steel
25mm	TI-106.225	106.225
30mm	TI-106.230	106.230
35mm	TI-106.235	106.235
40mm	TI-106.240	106.240
45mm	TI-106.245	106.245
50mm	TI-106.250	106.250
55mm	TI-106.255	106.255
60mm	TI-106.260	106.260
65mm	TI-106.265	106.265
70mm	TI-106.270	106.270




STERILE

Length	Titanium	Stainless Steel
25mm	TI-106.225-S	106.225-S
30mm	TI-106.230-S	106.230-S
35mm	TI-106.235-S	106.235-S
40mm	TI-106.240-S	106.240-S
45mm	TI-106.245-S	106.245-S
50mm	TI-106.250-S	106.250-S
55mm	TI-106.255-S	106.255-S
60mm	TI-106.260-S	106.260-S
65mm	TI-106.265-S	106.265-S
70mm	TI-106.270-S	106.270-S

Washer for Small Fragment Screws




Titanium	Stainless Steel
TI-100.110	100.110



STERILE

Titanium	Stainless Steel
TI-100.110-S	100.110-S

The background of the page features a series of thin, light gray, wavy lines that flow from the top left towards the bottom right, creating a sense of movement and depth. These lines are more densely packed in some areas, particularly near the top and bottom edges, and more spread out in the center.

INSTRUMENTS

2100-2.5-112

Drill Bit with Quick Coupling End, Ø2.5mm x Length 112mm, for Small Fragment

**2100-3.5-112**

Drill Bit with Quick Coupling End, Ø3.5mm x Length 112mm, for Small Fragment

**2104-02**

Bone Tap Quick Coupling for Cortical Screws, Ø3.5mm, for Small Fragment

**2104-04**

Bone Tap, Quick Coupling for Cancellous Screws, Ø4.0mm, for Small Fragment

**3420-01**

Drill Sleeve Insert Ø3.5/2.5mm, for Small Fragment

**3400-01**

Countersink for Ø3.5/4.0mm Screws, Quick Coupling, for Small Fragment



3402-000

T-Handle with Quick Coupling, for Small Fragment

**1472-054**

Quick Coupling Shaft

**3767-3.5**

Double Drill Guide, Ø2.5/3.5mm

**1472-046**

Self-Centering Double Drill Guide, Ø2.5/3.5mm, for Small Fragment

**1472-038**

Drill Guide for Neutral and Loaded Position, DCP Ø2.5mm, for Small Fragment

**1472-042**

Drill Guide for Neutral and Loaded Position, LC-DCP Ø2.5mm, for Small Fragment



3408-01

Hexagonal Screwdriver Shaft - 2.5mm Tip, Quick Coupling, for Small Fragment

**3443-37**


Depth Gauge measuring upto 80mm, for Small Fragment

**3443-38-3.5**

Extraction Screw Ø3.5mm

**PWG-2.0**

2.0mm Wire Sleeve for Small Fragment

**451-1.2-150E**

Kirschner Wire with One Side Trocar Tip, Ø1.2mm x Length 150mm


**451-1.5-150E**

Kirschner Wire with One Side Trocar Tip, Ø1.5mm x Length 150mm



**451-2.0-150E**

Kirschner Wire with One Side Trocar Tip, Ø2.0mm x Length 150mm

**3426-01**

Guide Wire with Threaded Trocar Tip, Ø1.2mm x Thread Length 5mm x Length 150mm

**3426-05**

Guide Wire with Threaded Trocar Tip, Ø1.5mm x Thread Length 5mm x Length 150mm

**GW-2.0-150**

Guide Wire with Threaded Trocar Tip, Ø2.0mm x Thread Length 8mm x Length 150mm

**3406-03**

Hexagonal Screwdriver - 2.5mm Tip, for Small Fragment

**3406-03S**

Screw Holding Sleeve for 2.5mm Tip Screwdriver, for Small Fragment



3409-04L

Bending Iron, Left for Small Fragment (Non Locking)

**3409-04R**

Bending Iron, Right for Small Fragment (Non Locking)

**7-046-01**

Reduction Forcep, Pointed, Ratchet Lock, for Small Fragment

**7-046-02**

Reduction Forcep, Serrated Jaws, Speed Lock, for Small Fragment

**2106-190**

Self-Centering Bone Holding Forcep, Speed Lock, 190mm, for Small Fragment

**2150-1006**

Periosteal Elevator with Silicon Handle, Curved, 6mm, for Small Fragment



**1472-066**

Hohmann Retractor, 6.5mm, for Small Fragment

**1472-068**

Hohmann Retractor, 8.5mm, for Small Fragment

**2146-018**

Hohmann Retractor, 15.5mm, for Small Fragment

**BT-3.5-87**

3.5mm Bending Template, Small

**BT-3.5-114**

3.5mm Bending Template, Large

**7-046-03**

Screw Holding Forcep, for Small Fragment



1472-032

Articulated Tension Device

**1472-034**

Articulated Tension Device Key

**7-046-04**

Instrument trays for Small Fragment Instrument Set

**ST-001NW**


Screw Caddy for Small fragment Instrument System

**IC-2310**

Container for Small Fragment Instrument Set



ORDER INFORMATION

2310-000**3.5mm Small Fragment Instrument Set**



ORDER INFORMATION

7-133	3.5mm Small Fragment Instrument Set
--------------	--

CODE	SET CONSISTING OF	UNIT
2100-2.5-112	Drill Bit with Quick Coupling End, Ø2.5mm x Length 112mm, for Small Fragment	3
2100-3.5-112	Drill Bit with Quick Coupling End, Ø3.5mm x Length 112mm, for Small Fragment	2
2104-02	Bone Tap Quick Coupling for Cortical Screws, Ø3.5mm, for Small Fragment	2
2104-04	Bone Tap, Quick Coupling for Cancellous Screws, Ø4.0mm, for Small Fragment	1
3420-01	Drill Sleeve Insert Ø3.5/2.5mm, for Small Fragment	1
3400-01	Countersink for Ø3.5/4.0mm Screws, Quick Coupling, for Small Fragment	1
3402-000	T-Handle with Quick Coupling, for Small Fragment	1
1472-054	Quick Coupling Shaft	1
3767-3.5	Double Drill Guide, Ø2.5/3.5mm	1
1472-046	Self-Centering Double Drill Guide, Ø2.5/3.5mm, for Small Fragment	1
1472-038	Drill Guide for Neutral and Loaded Position, DCP Ø2.5mm, for Small Fragment	1
1472-042	Drill Guide for Neutral and Loaded Position, LC-DCP Ø2.5mm, for Small Fragment	1
3408-01	Hexagonal Screwdriver Shaft - 2.5mm Tip, Quick Coupling, for Small Fragment	1
3443-37	Depth Gauge measuring upto 80mm, for Small Fragment	1
3443-38-3.5	Extraction Screw Ø 3.5mm	1
PWG-2.0	2.0mm Wire Sleeve, for Small Fragment	1
451-1.2-150E	Kirschner Wire with One Side Trocar Tip, Ø1.2mm x Length 150mm	3
451-1.5-150E	Kirschner Wire with One Side Trocar Tip, Ø1.5mm x Length 150mm	3
451-2.0-150E	Kirschner Wire with One Side Trocar Tip, Ø2.0mm x Length 150mm	3
3426-01	Guide Wire with Threaded Trocar Tip, Ø1.2mm x Thread Length 5mm x Length 150mm	3
3426-05	Guide Wire with Threaded Trocar Tip, Ø1.5mm x Thread Length 5mm x Length 150mm	3
GW-2.0-150	Guide Wire with Threaded Trocar Tip, Ø2.0mm x Thread Length 8mm x Length 150mm	3
3406-03	Hexagonal Screwdriver - 2.5mm Tip, for Small Fragment	1
3406-03S	Screw Holding Sleeve for 2.5mm Tip Screwdriver, for Small Fragment	1
3409-04L	Bending Iron, Left for Small Fragment (Non Locking)	1
3409-04R	Bending Iron, Right for Small Fragment (Non Locking)	1
7-046-01	Reduction Forcep, Pointed, Ratchet Lock, for Small Fragment	1
7-046-02	Reduction Forcep, Serrated Jaws, Speed Lock, for Small Fragment	1
2106-190	Self-Centering Bone Holding Forcep, Speed Lock, 190mm, for Small Fragment	2
2150-1006	Periosteal Elevator with Silicon Handle, Curved, 6mm, for Small Fragment	1
1472-066	Hohmann Retractor, 6.5mm, for Small Fragment	1
1472-068	Hohmann Retractor, 8.5mm, for Small Fragment	1
2146-018	Hohmann Retractor, 15.5mm, for Small Fragment	2
BT-3.5-87	3.5mm Bending Template, Small	1
BT-3.5-114	3.5mm Bending Template, Large	1
7-046-03	Screw Holding Forcep, for Small Fragment	1
1472-032	Articulated Tension Device	1
1472-034	Articulated Tension Device Key	1
7-046-04	Instrument trays for Small Fragment Instrument Set	2
ST-001NW	Screw Caddy for Small fragment Instrument System	1
IC-2310	Container for Small Fragment Instrument Set	1

The background of the page features a series of thin, light gray wavy lines that flow from the top left towards the bottom right, creating a sense of movement and depth. These lines are more densely packed in some areas, particularly near the top and bottom edges, and more spread out in the center.

LARGE FRAGMENT



4.5mm Narrow Dynamic Compression Plate



Holes	Titanium	Stainless Steel
4	TI-722.104	722.104
5	TI-722.105	722.105
6	TI-722.106	722.106
7	TI-722.107	722.107
8	TI-722.108	722.108
9	TI-722.109	722.109
10	TI-722.110	722.110
11	TI-722.111	722.111
12	TI-722.112	722.112
13	TI-722.113	722.113
14	TI-722.114	722.114
15	TI-722.115	722.115
16	TI-722.116	722.116
17	TI-722.117	722.117
18	TI-722.118	722.118

STERILE

Holes	Titanium	Stainless Steel
4	TI-722.104-S	722.104-S
5	TI-722.105-S	722.105-S
6	TI-722.106-S	722.106-S
7	TI-722.107-S	722.107-S
8	TI-722.108-S	722.108-S
9	TI-722.109-S	722.109-S
10	TI-722.110-S	722.110-S
11	TI-722.111-S	722.111-S
12	TI-722.112-S	722.112-S
13	TI-722.113-S	722.113-S
14	TI-722.114-S	722.114-S
15	TI-722.115-S	722.115-S
16	TI-722.116-S	722.116-S
17	TI-722.117-S	722.117-S
18	TI-722.118-S	722.118-S

4.5mm Narrow Limited Contact Dynamic Compression Plate



Holes	Titanium	Stainless Steel
4	TI-723.104	723.104
5	TI-723.105	723.105
6	TI-723.106	723.106
7	TI-723.107	723.107
8	TI-723.108	723.108
9	TI-723.109	723.109
10	TI-723.110	723.110
11	TI-723.111	723.111
12	TI-723.112	723.112
13	TI-723.113	723.113
14	TI-723.114	723.114
15	TI-723.115	723.115
16	TI-723.116	723.116
17	TI-723.117	723.117
18	TI-723.118	723.118

STERILE

Holes	Titanium	Stainless Steel
4	TI-723.104-S	723.104-S
5	TI-723.105-S	723.105-S
6	TI-723.106-S	723.106-S
7	TI-723.107-S	723.107-S
8	TI-723.108-S	723.108-S
9	TI-723.109-S	723.109-S
10	TI-723.110-S	723.110-S
11	TI-723.111-S	723.111-S
12	TI-723.112-S	723.112-S
13	TI-723.113-S	723.113-S
14	TI-723.114-S	723.114-S
15	TI-723.115-S	723.115-S
16	TI-723.116-S	723.116-S
17	TI-723.117-S	723.117-S
18	TI-723.118-S	723.118-S



4.5mm Broad Dynamic Compression Plate



Holes	Titanium	Stainless Steel
4	TI-724.104	724.104
5	TI-724.105	724.105
6	TI-724.106	724.106
7	TI-724.107	724.107
8	TI-724.108	724.108
9	TI-724.109	724.109
10	TI-724.110	724.110
11	TI-724.111	724.111
12	TI-724.112	724.112
13	TI-724.113	724.113
14	TI-724.114	724.114
15	TI-724.115	724.115
16	TI-724.116	724.116
17	TI-724.117	724.117
18	TI-724.118	724.118
19	TI-724.119	724.119
20	TI-724.120	724.120



STERILE

Holes	Titanium	Stainless Steel
4	TI-724.104-S	724.104-S
5	TI-724.105-S	724.105-S
6	TI-724.106-S	724.106-S
7	TI-724.107-S	724.107-S
8	TI-724.108-S	724.108-S
9	TI-724.109-S	724.109-S
10	TI-724.110-S	724.110-S
11	TI-724.111-S	724.111-S
12	TI-724.112-S	724.112-S
13	TI-724.113-S	724.113-S
14	TI-724.114-S	724.114-S
15	TI-724.115-S	724.115-S
16	TI-724.116-S	724.116-S
17	TI-724.117-S	724.117-S
18	TI-724.118-S	724.118-S
19	TI-724.119-S	724.119-S
20	TI-724.120-S	724.120-S

4.5mm Broad Limited Contact Dynamic Compression Plate



Holes	Titanium	Stainless Steel
4	TI-725.104	725.104
5	TI-725.105	725.105
6	TI-725.106	725.106
7	TI-725.107	725.107
8	TI-725.108	725.108
9	TI-725.109	725.109
10	TI-725.110	725.110
11	TI-725.111	725.111
12	TI-725.112	725.112
13	TI-725.113	725.113
14	TI-725.114	725.114
15	TI-725.115	725.115
16	TI-725.116	725.116
17	TI-725.117	725.117
18	TI-725.118	725.118
19	TI-725.119	725.119
20	TI-725.120	725.120



STERILE

Holes	Titanium	Stainless Steel
4	TI-725.104-S	725.104-S
5	TI-725.105-S	725.105-S
6	TI-725.106-S	725.106-S
7	TI-725.107-S	725.107-S
8	TI-725.108-S	725.108-S
9	TI-725.109-S	725.109-S
10	TI-725.110-S	725.110-S
11	TI-725.111-S	725.111-S
12	TI-725.112-S	725.112-S
13	TI-725.113-S	725.113-S
14	TI-725.114-S	725.114-S
15	TI-725.115-S	725.115-S
16	TI-725.116-S	725.116-S
17	TI-725.117-S	725.117-S
18	TI-725.118-S	725.118-S
19	TI-725.119-S	725.119-S
20	TI-725.120-S	725.120-S

4.5mm Narrow Lengthening Plate, 8 Holes



Lengthening	Length	Titanium	Stainless Steel
30mm	135mm	TI-755.030	755.030
40mm	145mm	TI-755.040	755.040
50mm	155mm	TI-755.050	755.050
60mm	165mm	TI-755.060	755.060
70mm	175mm	TI-755.070	755.070
80mm	185mm	TI-755.080	755.080

STERILE

Lengthening	Length	Titanium	Stainless Steel
30mm	135mm	TI-755.030-S	755.030-S
40mm	145mm	TI-755.040-S	755.040-S
50mm	155mm	TI-755.050-S	755.050-S
60mm	165mm	TI-755.060-S	755.060-S
70mm	175mm	TI-755.070-S	755.070-S
80mm	185mm	TI-755.080-S	755.080-S

4.5mm Broad Lengthening Plate, 8 Holes



Lengthening	Length	Titanium	Stainless Steel
30mm	135mm	TI-754.030	754.030
40mm	145mm	TI-754.040	754.040
50mm	155mm	TI-754.050	754.050
60mm	165mm	TI-754.060	754.060

STERILE

Lengthening	Length	Titanium	Stainless Steel
30mm	135mm	TI-754.030-S	754.030-S
40mm	145mm	TI-754.040-S	754.040-S
50mm	155mm	TI-754.050-S	754.050-S
60mm	165mm	TI-754.060-S	754.060-S

4.5mm Broad Lengthening Plate, 10 Holes



Lengthening	Length	Titanium	Stainless Steel
50mm	179mm	TI-753.050	753.050
60mm	189mm	TI-753.060	753.060
70mm	199mm	TI-753.070	753.070
80mm	209mm	TI-753.080	753.080
90mm	219mm	TI-753.090	753.090
100mm	229mm	TI-753.100	753.100
110mm	239mm	TI-753.110	753.110
120mm	249mm	TI-753.120	753.120

STERILE

Lengthening	Length	Titanium	Stainless Steel
50mm	179mm	TI-753.050-S	753.050-S
60mm	189mm	TI-753.060-S	753.060-S
70mm	199mm	TI-753.070-S	753.070-S
80mm	209mm	TI-753.080-S	753.080-S
90mm	219mm	TI-753.090-S	753.090-S
100mm	229mm	TI-753.100-S	753.100-S
110mm	239mm	TI-753.110-S	753.110-S
120mm	249mm	TI-753.120-S	753.120-S



4.5mm Semi Tubular Plate



Holes	Titanium	Stainless Steel
4	TI-704.204	704.204
5	TI-704.205	704.205
6	TI-704.206	704.206
7	TI-704.207	704.207
8	TI-704.208	704.208
9	TI-704.209	704.209
10	TI-704.210	704.210
11	TI-704.211	704.211
12	TI-704.212	704.212



STERILE

Holes	Titanium	Stainless Steel
4	TI-704.204-S	704.204-S
5	TI-704.205-S	704.205-S
6	TI-704.206-S	704.206-S
7	TI-704.207-S	704.207-S
8	TI-704.208-S	704.208-S
9	TI-704.209-S	704.209-S
10	TI-704.210-S	704.210-S
11	TI-704.211-S	704.211-S
12	TI-704.212-S	704.212-S

4.5mm Lateral Tibial Head Buttress Plate, Right

Fixation with 4.5 & 6.5mm Screws



Holes	Titanium	Stainless Steel
5	TI-733.105R	733.105R
7	TI-733.107R	733.107R
9	TI-733.109R	733.109R
11	TI-733.111R	733.111R
13	TI-733.113R	733.113R



STERILE

Holes	Titanium	Stainless Steel
5	TI-733.105R-S	733.105R-S
7	TI-733.107R-S	733.107R-S
9	TI-733.109R-S	733.109R-S
11	TI-733.111R-S	733.111R-S
13	TI-733.113R-S	733.113R-S

Lateral Tibial Head Buttress Plate, Left



Holes	Titanium	Stainless Steel
5	TI-733.105L	733.105L
7	TI-733.107L	733.107L
9	TI-733.109L	733.109L
11	TI-733.111L	733.111L
13	TI-733.113L	733.113L



STERILE

Holes	Titanium	Stainless Steel
5	TI-733.105L-S	733.105L-S
7	TI-733.107L-S	733.107L-S
9	TI-733.109L-S	733.109L-S
11	TI-733.111L-S	733.111L-S
13	TI-733.113L-S	733.113L-S

4.5mm Condylar Buttress Plate, Left



Holes	Titanium	Stainless Steel
7	TI-731.107L	731.107L
9	TI-731.109L	731.109L
11	TI-731.111L	731.111L
12	TI-731.112L	731.112L
13	TI-731.113L	731.113L
15	TI-731.115L	731.115L



STERILE

Holes	Titanium	Stainless Steel
7	TI-731.107L-S	731.107L-S
9	TI-731.109L-S	731.109L-S
11	TI-731.111L-S	731.111L-S
12	TI-731.112L-S	731.112L-S
13	TI-731.113L-S	731.113L-S
15	TI-731.115L-S	731.115L-S

4.5mm Condylar Buttress Plate, Right



Holes	Titanium	Stainless Steel
7	TI-731.107R	731.107R
9	TI-731.109R	731.109R
11	TI-731.111R	731.111R
12	TI-731.112R	731.112R
13	TI-731.113R	731.113R
15	TI-731.115R	731.115R



STERILE

Holes	Titanium	Stainless Steel
7	TI-731.107R-S	731.107R-S
9	TI-731.109R-S	731.109R-S
11	TI-731.111R-S	731.111R-S
12	TI-731.112R-S	731.112R-S
13	TI-731.113R-S	731.113R-S
15	TI-731.115R-S	731.115R-S

4.5mm "T" Plate



Holes	Titanium	Stainless Steel
3	TI-713.103	713.103
4	TI-713.104	713.104
5	TI-713.105	713.105
6	TI-713.106	713.106
7	TI-713.107	713.107
8	TI-713.108	713.108
9	TI-713.109	713.109
10	TI-713.110	713.110
11	TI-713.111	713.111
12	TI-713.112	713.112
13	TI-713.113	713.113
14	TI-713.114	713.114



STERILE

Holes	Titanium	Stainless Steel
3	TI-713.103-S	713.103-S
4	TI-713.104-S	713.104-S
5	TI-713.105-S	713.105-S
6	TI-713.106-S	713.106-S
7	TI-713.107-S	713.107-S
8	TI-713.108-S	713.108-S
9	TI-713.109-S	713.109-S
10	TI-713.110-S	713.110-S
11	TI-713.111-S	713.111-S
12	TI-713.112-S	713.112-S
13	TI-713.113-S	713.113-S
14	TI-713.114-S	713.114-S



4.5mm "T" Buttress Plate, Titanium



Holes	Titanium	Stainless Steel
4	TI-714.104	714.104
5	TI-714.105	714.105
6	TI-714.106	714.106
7	TI-714.107	714.107
8	TI-714.108	714.108
9	TI-714.109	714.109
10	TI-714.110	714.110
11	TI-714.111	714.111
12	TI-714.112	714.112
13	TI-714.113	714.113
14	TI-714.114	714.114

STERILE

Holes	Titanium	Stainless Steel
4	TI-714.104-S	714.104-S
5	TI-714.105-S	714.105-S
6	TI-714.106-S	714.106-S
7	TI-714.107-S	714.107-S
8	TI-714.108-S	714.108-S
9	TI-714.109-S	714.109-S
10	TI-714.110-S	714.110-S
11	TI-714.111-S	714.111-S
12	TI-714.112-S	714.112-S
13	TI-714.113-S	714.113-S
14	TI-714.114-S	714.114-S

4.5mm "L" Buttress Plate, Left



Holes	Titanium	Stainless Steel
3	TI-706.103L	706.103L
4	TI-706.104L	706.104L
5	TI-706.105L	706.105L
6	TI-706.106L	706.106L
7	TI-706.107L	706.107L
8	TI-706.108L	706.108L
9	TI-706.109L	706.109L
10	TI-706.110L	706.110L
11	TI-706.111L	706.111L
12	TI-706.112L	706.112L
13	TI-706.113L	706.113L
14	TI-706.114L	706.114L

STERILE

Holes	Titanium	Stainless Steel
3	TI-706.103L-S	706.103L-S
4	TI-706.104L-S	706.104L-S
5	TI-706.105L-S	706.105L-S
6	TI-706.106L-S	706.106L-S
7	TI-706.107L-S	706.107L-S
8	TI-706.108L-S	706.108L-S
9	TI-706.109L-S	706.109L-S
10	TI-706.110L-S	706.110L-S
11	TI-706.111L-S	706.111L-S
12	TI-706.112L-S	706.112L-S
13	TI-706.113L-S	706.113L-S
14	TI-706.114L-S	706.114L-S

4.5mm "L" Buttress Plate, Right



Holes	Titanium	Stainless Steel
3	TI-706.103R	706.103R
4	TI-706.104R	706.104R
5	TI-706.105R	706.105R
6	TI-706.106R	706.106R
7	TI-706.107R	706.107R
8	TI-706.108R	706.108R
9	TI-706.109R	706.109R
10	TI-706.110R	706.110R
11	TI-706.111R	706.111R
12	TI-706.112R	706.112R
13	TI-706.113R	706.113R
14	TI-706.114R	706.114R

STERILE

Holes	Titanium	Stainless Steel
3	TI-706.103R-S	706.103R-S
4	TI-706.104R-S	706.104R-S
5	TI-706.105R-S	706.105R-S
6	TI-706.106R-S	706.106R-S
7	TI-706.107R-S	706.107R-S
8	TI-706.108R-S	706.108R-S
9	TI-706.109R-S	706.109R-S
10	TI-706.110R-S	706.110R-S
11	TI-706.111R-S	706.111R-S
12	TI-706.112R-S	706.112R-S
13	TI-706.113R-S	706.113R-S
14	TI-706.114R-S	706.114R-S

4.5mm Spoon Plate



Holes	Titanium	Stainless Steel
5	TI-718.105	718.105
6	TI-718.10	718.106

STERILE

Holes	Titanium	Stainless Steel
5	TI-718.105-S	718.105-S
6	TI-718.106-S	718.106-S

4.5mm Hook Plate, 6 Holes

Compatible with 4.5 mm and 6.5 mm screws. The shaft features Dynamic Compression Plate (DCP) holes and a slot for a tension device. Used in cases where the blade of an angled hip plate has broken off in the proximal femur.



Holes	Titanium	Stainless Steel
6	TI-730.106	730.106

STERILE

Holes	Titanium	Stainless Steel
6	TI-730.106-S	730.106-S



4.5mm Reconstruction Plate, Straight

Do not bend more than 15°



Holes	Titanium	Stainless Steel
3	TI-728.103	728.103
4	TI-728.104	728.104
5	TI-728.105	728.105
6	TI-728.106	728.106
7	TI-728.107	728.107
8	TI-728.108	728.108
9	TI-728.109	728.109
10	TI-728.110	728.110
11	TI-728.111	728.111
12	TI-728.112	728.112
13	TI-728.113	728.113
14	TI-728.114	728.114
15	TI-728.115	728.115
16	TI-728.116	728.116
17	TI-728.117	728.117
18	TI-728.118	728.118
20	TI-728.120	728.120
22	TI-728.122	728.122



STERILE

Holes	Titanium	Stainless Steel
3	TI-728.103-S	728.103-S
4	TI-728.104-S	728.104-S
5	TI-728.105-S	728.105-S
6	TI-728.106-S	728.106-S
7	TI-728.107-S	728.107-S
8	TI-728.108-S	728.108-S
9	TI-728.109-S	728.109-S
10	TI-728.110-S	728.110-S
11	TI-728.111-S	728.111-S
12	TI-728.112-S	728.112-S
13	TI-728.113-S	728.113-S
14	TI-728.114-S	728.114-S
15	TI-728.115-S	728.115-S
16	TI-728.116-S	728.116-S
17	TI-728.117-S	728.117-S
18	TI-728.118-S	728.118-S
20	TI-728.120-S	728.120-S
22	TI-728.122-S	728.122-S

4.5mm Cobra Head Plate



Holes	Titanium	Stainless Steel
8	TI-735.108	735.108
9	TI-735.109	735.109
10	TI-735.110	735.110
11	TI-735.111	735.111
12	TI-735.112	735.112
13	TI-735.113	735.113
14	TI-735.114	735.114



STERILE

Holes	Titanium	Stainless Steel
8	TI-735.108-S	735.108-S
9	TI-735.109-S	735.109-S
10	TI-735.110-S	735.110-S
11	TI-735.111-S	735.111-S
12	TI-735.112-S	735.112-S
13	TI-735.113-S	735.113-S
14	TI-735.114-S	735.114-S



4.5mm Roy Camille Plate



Holes	Titanium	Stainless Steel
5	TI-804.005	804.005
6	TI-804.006	804.006
7	TI-804.007	804.007
8	TI-804.008	804.008
9	TI-804.009	804.009
10	TI-804.010	804.010
11	TI-804.011	804.011
12	TI-804.012	804.012
13	TI-804.013	804.013
14	TI-804.014	804.014

STERILE

Holes	Titanium	Stainless Steel
5	TI-804.005-S	804.005-S
6	TI-804.006-S	804.006-S
7	TI-804.007-S	804.007-S
8	TI-804.008-S	804.008-S
9	TI-804.009-S	804.009-S
10	TI-804.010-S	804.010-S
11	TI-804.011-S	804.011-S
12	TI-804.012-S	804.012-S
13	TI-804.013-S	804.013-S
14	TI-804.014-S	804.014-S

Sherman Plate



Holes	Titanium	Stainless Steel
3	TI-802.003	802.003
4	TI-802.004	802.004
5	TI-802.005	802.005
6	TI-802.006	802.006
7	TI-802.007	802.007
8	TI-802.008	802.008
9	TI-802.009	802.009
10	TI-802.010	802.010
11	TI-802.011	802.011
12	TI-802.012	802.012
13	TI-802.013	802.013
14	TI-802.014	802.014

STERILE

Holes	Titanium	Stainless Steel
3	TI-802.003-S	802.003-S
4	TI-802.004-S	802.004-S
5	TI-802.005-S	802.005-S
6	TI-802.006-S	802.006-S
7	TI-802.007-S	802.007-S
8	TI-802.008-S	802.008-S
9	TI-802.009-S	802.009-S
10	TI-802.010-S	802.010-S
11	TI-802.011-S	802.011-S
12	TI-802.012-S	802.012-S
13	TI-802.013-S	802.013-S
14	TI-802.014-S	802.014-S



4.5mm Proximal Tibia Medial Plate



Left Direction

Right Direction

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
3	TI-763.103L	763.103L	TI-763.103R	763.103R
5	TI-763.105L	763.105L	TI-763.105R	763.105R
7	TI-763.107L	763.107L	TI-763.107R	763.107R
9	TI-763.109L	763.109L	TI-763.109R	763.109R
11	TI-763.111L	763.111L	TI-763.111R	763.111R
13	TI-763.113L	763.113L	TI-763.113R	763.113R

STERILE

Left Direction

Right Direction

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
3	TI-763.103L-S	763.103L-S	TI-763.103R-S	763.103R-S
5	TI-763.105L-S	763.105L-S	TI-763.105R-S	763.105R-S
7	TI-763.107L-S	763.107L-S	TI-763.107R-S	763.107R-S
9	TI-763.109L-S	763.109L-S	TI-763.109R-S	763.109R-S
11	TI-763.111L-S	763.111L-S	TI-763.111R-S	763.111R-S
13	TI-763.113L-S	763.113L-S	TI-763.113R-S	763.113R-S



4.5mm Proximal Tibia Plate with Round Holes



Left Direction

Right Direction

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
4	TI-MP03.04L	MP03.04L	TI-MP03.04R	MP03.04R
5	TI-MP03.05L	MP03.05L	TI-MP03.05R	MP03.05R
6	TI-MP03.06L	MP03.06L	TI-MP03.06R	MP03.06R
7	TI-MP03.07L	MP03.07L	TI-MP03.07R	MP03.07R
8	TI-MP03.08L	MP03.08L	TI-MP03.08R	MP03.08R
9	TI-MP03.09L	MP03.09L	TI-MP03.09R	MP03.09R
10	TI-MP03.10L	MP03.10L	TI-MP03.10R	MP03.10R

STERILE

Left Direction

Right Direction

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
4	TI-MP03.04L-S	MP03.04L-S	TI-MP03.04R-S	MP03.04R-S
5	TI-MP03.05L-S	MP03.05L-S	TI-MP03.05R-S	MP03.05R-S
6	TI-MP03.06L-S	MP03.06L-S	TI-MP03.06R-S	MP03.06R-S
7	TI-MP03.07L-S	MP03.07L-S	TI-MP03.07R-S	MP03.07R-S
8	TI-MP03.08L-S	MP03.08L-S	TI-MP03.08R-S	MP03.08R-S
9	TI-MP03.09L-S	MP03.09L-S	TI-MP03.09R-S	MP03.09R-S
10	TI-MP03.10L-S	MP03.10L-S	TI-MP03.10R-S	MP03.10R-S





4.5mm Proximal Lateral Tibia Plate



Left Direction

Right Direction

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
3	TI-736.103L	736.103L	TI-736.103R	736.103R
5	TI-736.105L	736.105L	TI-736.105R	736.105R
7	TI-736.107L	736.107L	TI-736.107R	736.107R
9	TI-736.109L	736.109L	TI-736.109R	736.109R
11	TI-736.111L	736.111L	TI-736.111R	736.111R
13	TI-736.113L	736.113L	TI-736.113R	736.113R

STERILE

Left Direction

Right Direction

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
3	TI-736.103L-S	736.103L-S	TI-736.103R-S	736.103R-S
5	TI-736.105L-S	736.105L-S	TI-736.105R-S	736.105R-S
7	TI-736.107L-S	736.107L-S	TI-736.107R-S	736.107R-S
9	TI-736.109L-S	736.109L-S	TI-736.109R-S	736.109R-S
11	TI-736.111L-S	736.111L-S	TI-736.111R-S	736.111R-S
13	TI-736.113L-S	736.113L-S	TI-736.113R-S	736.113R-S



4.5mm Proximal Femur Plate



Left Direction

Right Direction

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
3	TI-764.103L	764.103L	TI-764.103R	764.103R
5	TI-764.105L	764.105L	TI-764.105R	764.105R
7	TI-764.107L	764.107L	TI-764.107R	764.107R
9	TI-764.109L	764.109L	TI-764.109R	764.109R
11	TI-764.111L	764.111L	TI-764.111R	764.111R
13	TI-764.113L	764.113L	TI-764.113R	764.113R

STERILE

Left Direction

Right Direction

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
3	TI-764.103L-S	764.103L-S	TI-764.103R-S	764.103R-S
5	TI-764.105L-S	764.105L-S	TI-764.105R-S	764.105R-S
7	TI-764.107L-S	764.107L-S	TI-764.107R-S	764.107R-S
9	TI-764.109L-S	764.109L-S	TI-764.109R-S	764.109R-S
11	TI-764.111L-S	764.111L-S	TI-764.111R-S	764.111R-S
13	TI-764.113L-S	764.113L-S	TI-764.113R-S	764.113R-S





4.5mm Distal Tibia Plate



Left Direction

Right Direction

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
3	TI-737.103L	737.103L	TI-737.103R	737.103R
5	TI-737.105L	737.105L	TI-737.105R	737.105R
7	TI-737.107L	737.107L	TI-737.107R	737.107R
9	TI-737.109L	737.109L	TI-737.109R	737.109R
11	TI-737.111L	737.111L	TI-737.111R	737.111R
13	TI-737.113L	737.113L	TI-737.113R	737.113R

STERILE

Left Direction

Right Direction

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
5	TI-737.103L-S	737.103L-S	TI-737.103R-S	737.103R-S
7	TI-737.105L-S	737.105L-S	TI-737.105R-S	737.105R-S
9	TI-737.107L-S	737.107L-S	TI-737.107R-S	737.107R-S
10	TI-737.109L-S	737.109L-S	TI-737.109R-S	737.109R-S
11	TI-737.111L-S	737.111L-S	TI-737.111R-S	737.111R-S
13	TI-737.113L-S	737.113L-S	TI-737.113R-S	737.113R-S



4.5mm Distal Lateral Tibia Plate with Round Holes



STERILE

Holes	Type	Titanium	Stainless Steel
5	Left	TI-CP01.05L	CP01.05L
5	Right	TI-CP01.05R	CP01.05R
7	Left	TI-CP01.07L	CP01.07L
7	Right	TI-CP01.07R	CP01.07R
9	Left	TI-CP01.09L	CP01.09L
9	Right	TI-CP01.09R	CP01.09R
11	Left	TI-CP01.11L	CP01.11L
11	Right	TI-CP01.11R	CP01.11R
13	Left	TI-CP01.13L	CP01.13L
13	Right	TI-CP01.13R	CP01.13R

Holes	Type	Titanium	Stainless Steel
5	Left	TI-CP01.05L-S	CP01.05L-S
5	Right	TI-CP01.05R-S	CP01.05R-S
7	Left	TI-CP01.07L-S	CP01.07L-S
7	Right	TI-CP01.07R-S	CP01.07R-S
9	Left	TI-CP01.09L-S	CP01.09L-S
9	Right	TI-CP01.09R-S	CP01.09R-S
11	Left	TI-CP01.11L-S	CP01.11L-S
11	Right	TI-CP01.11R-S	CP01.11R-S
13	Left	TI-CP01.13L-S	CP01.13L-S
13	Right	TI-CP01.13R-S	CP01.13R-S



4.5mm Distal Tibia Lateral Plate



Left Direction

Right Direction

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
3	TI-737.203L	737.203L	TI-737.203R	737.203R
5	TI-737.205L	737.205L	TI-737.205R	737.205R
7	TI-737.207L	737.207L	TI-737.207R	737.207R
9	TI-737.209L	737.209L	TI-737.209R	737.209R
11	TI-737.211L	737.211L	TI-737.211R	737.211R
13	TI-737.213L	737.213L	TI-737.213R	737.213R

STERILE

Left Direction

Right Direction

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
5	TI-737.203L-S	737.203L-S	TI-737.203R-S	737.203R-S
7	TI-737.205L-S	737.205L-S	TI-737.205R-S	737.205R-S
9	TI-737.207L-S	737.207L-S	TI-737.207R-S	737.207R-S
10	TI-737.209L-S	737.209L-S	TI-737.209R-S	737.209R-S
11	TI-737.211L-S	737.211L-S	TI-737.211R-S	737.211R-S
13	TI-737.213L-S	737.213L-S	TI-737.213R-S	737.213R-S

4.5mm Distal Femur Plate



Left Direction

Right Direction

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
3	TI-737.303L	737.303L	TI-737.303R	737.303R
5	TI-737.305L	737.305L	TI-737.305R	737.305R
7	TI-737.307L	737.307L	TI-737.307R	737.307R
9	TI-737.309L	737.309L	TI-737.309R	737.309R
11	TI-737.311L	737.311L	TI-737.311R	737.311R
13	TI-737.313L	737.313L	TI-737.313R	737.313R

STERILE

Left Direction

Right Direction

Holes	Titanium	Stainless Steel	Titanium	Stainless Steel
5	TI-737.303L-S	737.303L-S	TI-737.303R-S	737.303R-S
7	TI-737.305L-S	737.305L-S	TI-737.305R-S	737.305R-S
9	TI-737.307L-S	737.307L-S	TI-737.307R-S	737.307R-S
10	TI-737.309L-S	737.309L-S	TI-737.309R-S	737.309R-S
11	TI-737.311L-S	737.311L-S	TI-737.311R-S	737.311R-S
13	TI-737.313L-S	737.313L-S	TI-737.313R-S	737.313R-S






4.5mm Fibular Distal Tibia Plate




Holes	Type	Titanium	Stainless Steel
7	Left	TI-CP12.07L	CP12.07L
9	Left	TI-CP12.09L	CP12.09L
11	Left	TI-CP12.11L	CP12.11L
13	Left	TI-CP12.13L	CP12.13L
7	Right	TI-CP12.07R	CP12.07R
9	Right	TI-CP12.09R	CP12.09R
11	Right	TI-CP12.11R	CP12.11R
13	Right	TI-CP12.13R	CP12.13R

STERILE

Holes	Type	Titanium	Stainless Steel
7	Left	CP12.07L-S	TI-CP12.07L-S
9	Left	CP12.09L-S	TI-CP12.09L-S
11	Left	CP12.11L-S	TI-CP12.11L-S
13	Left	CP12.13L-S	TI-CP12.13L-S
7	Right	CP12.07R-S	TI-CP12.07R-S
9	Right	CP12.09R-S	TI-CP12.09R-S
11	Right	CP12.11R-S	TI-CP12.11R-S
13	Right	CP12.13R-S	TI-CP12.13R-S



The background of the page features a series of thin, light gray wavy lines that flow from the top left towards the bottom right, creating a sense of movement and depth. These lines are more densely packed in some areas, particularly near the top and bottom edges, and more spread out in the center.

SCREWS



4.5mm Cortical Screw, (Hex Head)

Drill Bit : 3.2mm



STERILE

Length	Titanium	Stainless Steel
12	TI-105.012	105.012
14	TI-105.014	105.014
16	TI-105.016	105.016
18	TI-105.018	105.018
20	TI-105.020	105.020
22	TI-105.022	105.022
24	TI-105.024	105.024
26	TI-105.026	105.026
28	TI-105.028	105.028
30	TI-105.030	105.030
32	TI-105.032	105.032
34	TI-105.034	105.034
36	TI-105.036	105.036
38	TI-105.038	105.038
40	TI-105.040	105.040
42	TI-105.042	105.042
44	TI-105.044	105.044
46	TI-105.046	105.046
48	TI-105.048	105.048
50	TI-105.050	105.050
52	TI-105.052	105.052
54	TI-105.054	105.054
56	TI-105.056	105.056
58	TI-105.058	105.058
60	TI-105.060	105.060
62	TI-105.062	105.062
64	TI-105.064	105.064
66	TI-105.066	105.066
68	TI-105.068	105.068
70	TI-105.070	105.070
72	TI-105.072	105.072
74	TI-105.074	105.074
76	TI-105.076	105.076
78	TI-105.078	105.078
80	TI-105.080	105.080



Length	Titanium	Stainless Steel
12	TI-105.012-S	105.012-S
14	TI-105.014-S	105.014-S
16	TI-105.016-S	105.016-S
18	TI-105.018-S	105.018-S
20	TI-105.020-S	105.020-S
22	TI-105.022-S	105.022-S
24	TI-105.024-S	105.024-S
26	TI-105.026-S	105.026-S
28	TI-105.028-S	105.028-S
30	TI-105.030-S	105.030-S
32	TI-105.032-S	105.032-S
34	TI-105.034-S	105.034-S
36	TI-105.036-S	105.036-S
38	TI-105.038-S	105.038-S
40	TI-105.040-S	105.040-S
42	TI-105.042-S	105.042-S
44	TI-105.044-S	105.044-S
46	TI-105.046-S	105.046-S
48	TI-105.048-S	105.048-S
50	TI-105.050-S	105.050-S
52	TI-105.052-S	105.052-S
54	TI-105.054-S	105.054-S
56	TI-105.056-S	105.056-S
58	TI-105.058-S	105.058-S
60	TI-105.060-S	105.060-S
62	TI-105.062-S	105.062-S
64	TI-105.064-S	105.064-S
66	TI-105.066-S	105.066-S
68	TI-105.068-S	105.068-S
70	TI-105.070-S	105.070-S
72	TI-105.072-S	105.072-S
74	TI-105.074-S	105.074-S
76	TI-105.076-S	105.076-S
78	TI-105.078-S	105.078-S
80	TI-105.080-S	105.080-S



4.5mm Cortical Screw, Self-Tapping, (Hex Head)

Drill Bit : 3.2mm



STERILE

Length	Titanium	Stainless Steel
12	TI-105.212	105.212
14	TI-105.214	105.214
16	TI-105.216	105.216
18	TI-105.218	105.218
20	TI-105.220	105.220
22	TI-105.222	105.222
24	TI-105.224	105.224
26	TI-105.226	105.226
28	TI-105.228	105.228
30	TI-105.230	105.230
32	TI-105.232	105.232
34	TI-105.234	105.234
36	TI-105.236	105.236
38	TI-105.238	105.238
40	TI-105.240	105.240
42	TI-105.242	105.242
44	TI-105.244	105.244
46	TI-105.246	105.246
48	TI-105.248	105.248
50	TI-105.250	105.250
52	TI-105.252	105.252
54	TI-105.254	105.254
56	TI-105.256	105.256
58	TI-105.258	105.258
60	TI-105.260	105.260
62	TI-105.262	105.262
64	TI-105.264	105.264
66	TI-105.266	105.266
68	TI-105.268	105.268
70	TI-105.270	105.270
72	TI-105.272	105.272
74	TI-105.274	105.274
76	TI-105.276	105.276
78	TI-105.278	105.278
80	TI-105.280	105.280
85	TI-105.285	105.285
90	TI-105.290	105.290



Length	Titanium	Stainless Steel
12	TI-105.212-S	105.212-S
14	TI-105.214-S	105.214-S
16	TI-105.216-S	105.216-S
18	TI-105.218-S	105.218-S
20	TI-105.220-S	105.220-S
22	TI-105.222-S	105.222-S
24	TI-105.224-S	105.224-S
26	TI-105.226-S	105.226-S
28	TI-105.228-S	105.228-S
30	TI-105.230-S	105.230-S
32	TI-105.232-S	105.232-S
34	TI-105.234-S	105.234-S
36	TI-105.236-S	105.236-S
38	TI-105.238-S	105.238-S
40	TI-105.240-S	105.240-S
42	TI-105.242-S	105.242-S
44	TI-105.244-S	105.244-S
46	TI-105.246-S	105.246-S
48	TI-105.248-S	105.248-S
50	TI-105.250-S	105.250-S
52	TI-105.252-S	105.252-S
54	TI-105.254-S	105.254-S
56	TI-105.256-S	105.256-S
58	TI-105.258-S	105.258-S
60	TI-105.260-S	105.260-S
62	TI-105.262-S	105.262-S
64	TI-105.264-S	105.264-S
66	TI-105.266-S	105.266-S
68	TI-105.268-S	105.268-S
70	TI-105.270-S	105.270-S
72	TI-105.272-S	105.272-S
74	TI-105.274-S	105.274-S
76	TI-105.276-S	105.276-S
78	TI-105.278-S	105.278-S
80	TI-105.280-S	105.280-S
85	TI-105.285-S	105.285-S
90	TI-105.290-S	105.290-S



4.5mm Malleolar Screw

Drill Bit : 3.2mm



Length	Titanium	Stainless Steel
25	TI-106.025	106.025
30	TI-106.030	106.030
35	TI-106.035	106.035
40	TI-106.040	106.040
45	TI-106.045	106.045
50	TI-106.050	106.050
55	TI-106.055	106.055
60	TI-106.060	106.060
65	TI-105.065	105.065
70	TI-106.070	106.070
75	TI-106.075	106.075
80	TI-106.080	106.080
85	TI-106.085	106.085
90	TI-106.090	106.090
95	TI-106.095	106.095
100	TI-106.100	106.100



STERILE

Length	Titanium	Stainless Steel
25	TI-106.025-S	106.025-S
30	TI-106.030-S	106.030-S
35	TI-106.035-S	106.035-S
40	TI-106.040-S	106.040-S
45	TI-106.045-S	106.045-S
50	TI-106.050-S	106.050-S
55	TI-106.055-S	106.055-S
60	TI-106.060-S	106.060-S
65	TI-105.065-S	105.065-S
70	TI-106.070-S	106.070-S
75	TI-106.075-S	106.075-S
80	TI-106.080-S	106.080-S
85	TI-106.085-S	106.085-S
90	TI-106.090-S	106.090-S
95	TI-106.095-S	106.095-S
100	TI-106.100-S	106.100-S

Washer for 4.5mm to 7mm Screws



Titanium	Stainless Steel
TI-100.120	100.120

STERILE

Titanium	Stainless Steel
TI-100.120-S	100.120-S



6.5mm Cancellous Screw, 16mm Thread

Drill Bit : 3.2mm



Length	Titanium	Stainless Steel
25	TI-111.025	111.025
30	TI-111.030	111.030
35	TI-111.035	111.035
40	TI-111.040	111.040
45	TI-111.045	111.045
50	TI-111.050	111.050
55	TI-111.055	111.055
60	TI-111.060	111.060
65	TI-111.065	111.065
70	TI-111.070	111.070
75	TI-111.075	111.075
80	TI-111.080	111.080
85	TI-111.085	111.085
90	TI-111.090	111.090
95	TI-111.095	111.095
100	TI-111.100	111.100
105	TI-111.105	111.105
110	TI-111.110	111.110
115	TI-111.115	111.115
120	TI-111.120	111.120

**STERILE**

Length	Titanium	Stainless Steel
25	TI-111.025-S	111.025-S
30	TI-111.030-S	111.030-S
35	TI-111.035-S	111.035-S
40	TI-111.040-S	111.040-S
45	TI-111.045-S	111.045-S
50	TI-111.050-S	111.050-S
55	TI-111.055-S	111.055-S
60	TI-111.060-S	111.060-S
65	TI-111.065-S	111.065-S
70	TI-111.070-S	111.070-S
75	TI-111.075-S	111.075-S
80	TI-111.080-S	111.080-S
85	TI-111.085-S	111.085-S
90	TI-111.090-S	111.090-S
95	TI-111.095-S	111.095-S
100	TI-111.100-S	111.100-S
105	TI-111.105-S	111.105-S
110	TI-111.110-S	111.110-S
115	TI-111.115-S	111.115-S
120	TI-111.120-S	111.120-S

6.5mm Cancellous Screw, 32mm Thread

Drill Bit : 3.2mm



Length	Titanium	Stainless Steel
35	TI-112.035	112.035
40	TI-112.040	112.040
45	TI-112.045	112.045
50	TI-112.050	112.050
55	TI-112.055	112.055
60	TI-112.060	112.060
65	TI-112.065	112.065
70	TI-112.070	112.070
75	TI-112.075	112.075
80	TI-112.080	112.080
85	TI-112.085	112.085
90	TI-112.090	112.090
95	TI-112.095	112.095
100	TI-112.100	112.100
105	TI-112.105	112.105
110	TI-112.110	112.110
115	TI-112.115	112.115
120	TI-112.120	112.120

**STERILE**

Length	Titanium	Stainless Steel
35	TI-112.035-S	112.035-S
40	TI-112.040-S	112.040-S
45	TI-112.045-S	112.045-S
50	TI-112.050-S	112.050-S
55	TI-112.055-S	112.055-S
60	TI-112.060-S	112.060-S
65	TI-112.065-S	112.065-S
70	TI-112.070-S	112.070-S
75	TI-112.075-S	112.075-S
80	TI-112.080-S	112.080-S
85	TI-112.085-S	112.085-S
90	TI-112.090-S	112.090-S
95	TI-112.095-S	112.095-S
100	TI-112.100-S	112.100-S
105	TI-112.105-S	112.105-S
110	TI-112.110-S	112.110-S
115	TI-112.115-S	112.115-S
120	TI-112.120-S	112.120-S



6.5mm Cancellous Screw, Full Thread

Drill Bit : 3.2mm




Length	Titanium	Stainless Steel
25	TI-113.025	113.025
30	TI-113.030	113.030
35	TI-113.035	113.035
40	TI-113.040	113.040
45	TI-113.045	113.045
50	TI-113.050	113.050
55	TI-113.055	113.055
60	TI-113.060	113.060
65	TI-113.065	113.065
70	TI-113.070	113.070
75	TI-113.075	113.075
80	TI-113.080	113.080
85	TI-113.085	113.085
90	TI-113.090	113.090
95	TI-113.095	113.095
100	TI-113.100	113.100
105	TI-113.105	113.105
110	TI-113.110	113.110
115	TI-113.115	113.115
120	TI-113.120	113.120




STERILE

Length	Titanium	Stainless Steel
25	TI-113.025-S	113.025-S
30	TI-113.030-S	113.030-S
35	TI-113.035-S	113.035-S
40	TI-113.040-S	113.040-S
45	TI-113.045-S	113.045-S
50	TI-113.050-S	113.050-S
55	TI-113.055-S	113.055-S
60	TI-113.060-S	113.060-S
65	TI-113.065-S	113.065-S
70	TI-113.070-S	113.070-S
75	TI-113.075-S	113.075-S
80	TI-113.080-S	113.080-S
85	TI-113.085-S	113.085-S
90	TI-113.090-S	113.090-S
95	TI-113.095-S	113.095-S
100	TI-113.100-S	113.100-S
105	TI-113.105-S	113.105-S
110	TI-113.110-S	113.110-S
115	TI-113.115-S	113.115-S
120	TI-113.120-S	113.120-S

ORDER INFORMATION

ST-002NW**Screw Caddy for Large Fragment Instrument Set**

The background of the page features a light gray, wavy, topographical-like pattern that flows across the top and bottom edges, framing the central text.


INSTRUMENTS

2100-2.0-100

Drill Bit with Quick Coupling End, Ø2.0mm x Length 100mm, for Large Fragment

**2100-3.2-145**

Drill Bit with Quick Coupling End, Ø3.2mm x Length 145mm, for Large Fragment

**2100-4.5-145**

Drill Bit with Quick Coupling End, Ø4.5mm x Length 145mm, for Large Fragment

**2103-07**

Bone Tap Quick Coupling for Ø4.5mm Cortical Screws, for Large Fragment

**2103-09**

Bone Tap Quick Coupling for Ø6.5mm Cancellous Screws, for Large Fragment

**3767-4.5**

Double Drill Guide, Ø4.5/3.2mm



1472-058

Self-Centering Double Drill Guide, Ø4.5/3.2mm, for Large Fragment

**3767-6.5**

Double Drill Guide, Ø6.5/3.2mm

**3767-01**

Drill Sleeve Ø6.5/4.5mm for Large Fragment

**3412-030**

Drill Guide For Neutral and Loaded Position, DCP Ø4.5/3.2mm, for Large Fragment

**1472-040**

Drill Guide For Neutral and Loaded Position, LC-DCP Ø3.2/3.2mm, for Large Fragment

**1472-050**

Pointed Drill Guide, Ø6.5mm, for Large Fragment



3420-020

Insert Drill Sleeve Ø4.5/3.2mm, for Large Fragment

**3443-38-5.0**

Extraction Screw, Ø5.0mm

**3402-000**

T-Handle With Quick Coupling

**1472-054**

Quick Coupling Shaft

**451-2.0-150E**

Kirschner Wire with One Side Trocar Tip, Ø2.0mm x Length 150mm

**GW-2.0-150**

Guide Wire with Threaded Trocar Tip, Ø2.0mm x Thread Length 8mm x Length 150mm



3409-010L

Large Bending Iron, Left, for Large Fragment

**3409-010R**

Large Bending Iron, Right, for Large Fragment

**BT-4.5-120**

Bending Template, Small, for Large Fragment

**BT-4.5-155**

Bending Template, Medium, for Large Fragment

**BT-4.5-212**

Bending Template, Large, for Large Fragment

**3406-04**

Screwdriver, Hex 3.5mm



1472-056

Screwdriver Holding Sleeve for Cortical and Cancellous Screws

**3408-04**

Hexagonal Screwdriver Shaft - 3.5mm Tip, Quick Coupling, Short

**3400-02**

T-Handle Countersink for Ø4.5/6.5mm Screws, for Large Fragment

**3443-400**

Depth Gauge measuring upto 110mm

**2139-000**

Sharp Hook for Large Fragment

**1472-032**

Articulated Tension Device



1472-034

Articulated Tension Device Key

**2150-1008**

Periosteal Elevator, Round, 8mm

**2150-1016**

Periosteal Elevator, Flat, 16mm

**1472-070**

Hohmann Retractor, 44mm, for Large Fragment

**2146-016**

Hohmann Retractor, 16mm, for Large Fragment

**2106-1210**

Serrated Reduction Forcep, Speed Lock, for Large Fragment



2106-2210

Pointed Reduction Forcep, Ratchet Lock, for Large Fragment

**2105-000**

Self-Centering Bone Holding Forcep, Speed Lock

**7-047-01**

Screw Holding forcep for Large Fragment

**ST-002NW**

Screw Caddy for Large Fragment Instrument System

**7-047-02**


Instrument trays for Large Fragment Instrument Set

**7-047-03**

Container for Large Fragment Instrument Set



ORDER INFORMATION

2309-000**Instrument Set For Large Fragment**



ORDER INFORMATION

2309-000	Instrument Set For Large Fragment
-----------------	--


CODE	SET CONSISTING OF	UNIT
2100-2.0-100	Drill Bit with Quick Coupling End, Ø2.0mm x Length 100mm, for Large Fragment	2
2100-3.2-145	Drill Bit with Quick Coupling End, Ø3.2mm x Length 145mm, for Large Fragment	3
2100-4.5-145	Drill Bit with Quick Coupling End, Ø4.5mm x Length 145mm, for Large Fragment	2
2103-07	Bone Tap Quick Coupling for Ø4.5mm Cortical Screws, for Large Fragment	3
2103-09	Bone Tap Quick Coupling for Ø6.5mm Cancellous Screws, for Large Fragment	1
3767-4.5	Double Drill Guide, Ø4.5/3.2mm	1
1472-058	Self-Centering Double Drill Guide, Ø4.5/3.2mm, for Large Fragment	1
3767-6.5	Double Drill Guide, Ø6.5/3.2mm	1
3767-01	Drill Sleeve Ø6.5/4.5mm for Large Fragment	1
3412-030	Drill Guide For Neutral and Loaded Position, DCP Ø4.5/3.2mm, for Large Fragment	1
1472-040	Drill Guide For Neutral and Loaded Position, LC-DCP Ø3.2/3.2mm, for Large Fragment	1
1472-050	Pointed Drill Guide, Ø6.5mm, for Large Fragment	1
3420-020	Insert Drill Sleeve Ø4.5/3.2mm, for Large Fragment	1
3443-38-5.0	Extraction Screw, Ø5.0mm	1
3402-000	T-Handle With Quick Coupling	1
1472-054	Quick Coupling Shaft	1
451-2.0-150E	Kirschner Wire with One Side Trocar Tip, Ø2.0mm x Length 150mm	3
GW-2.0-150	Guide Wire with Threaded Trocar Tip, Ø2.0mm x Thread Length 8mm x Length 150mm	3
3409-010L	Large Bending Iron, Left, for Large Fragment	1
3409-010R	Large Bending Iron, Right, for Large Fragment	1
BT-4.5-120	Bending Template, Small, for Large Fragment	1
BT-4.5-155	Bending Template, Medium, for Large Fragment	1
BT-4.5-212	Bending Template, Large, for Large Fragment	1
3406-04	Screwdriver, Hex 3.5mm	1
1472-056	Screwdriver Holding Sleeve for Cortical and Cancellous Screws	1
3408-04	Hexagonal Screwdriver Shaft - 3.5mm Tip, Quick Coupling, Short	1
3400-02	T-Handle Countersink for Ø4.5/6.5mm Screws, for Large Fragment	1
3443-400	Depth Gauge measuring upto 110mm	1
2139-000	Sharp Hook for Large Fragment	1
1472-032	Articulated Tension Device	1
1472-034	Articulated Tension Device Key	1
2150-1008	Periosteal Elevator, Round, 8mm	1
2150-1016	Periosteal Elevator, Flat, 16mm	1
1472-070	Hohmann Retractor, 44mm, for Large Fragment	1
2146-016	Hohmann Retractor, 16mm, for Large Fragment	1
2106-1210	Serrated Reduction Forcep, Speed Lock, for Large Fragment	1
2106-2210	Pointed Reduction Forcep, Ratchet Lock, for Large Fragment	1
2105-000	Self-Centering Bone Holding Forcep, Speed Lock	2
7-047-01	Screw Holding forcep for Large Fragment	1
ST-002NW	Screw Caddy for Large Fragment Instrument System	1
7-047-02	Instrument trays for Large Fragment Instrument Set	3
7-047-03	Container for Large Fragment Instrument Set	1



The background of the page features a light gray, wavy, topographical-style line pattern. These lines flow from the top left and bottom left corners towards the right side, creating a sense of movement and depth. The lines are thin and closely spaced in some areas, while more spread out in others.

Angle Blade Plates

95° Condylar Blade Plate



Blade Length

	Titanium				Stainless Steel			
Holes	50mm	60mm	70mm	80mm	50mm	60mm	70mm	80mm
5	TI-765.5005	TI-765.6005	TI-765.7005	TI-765.8005	765.5005	765.6005	765.7005	765.8005
6	TI-765.5006	TI-765.6006	TI-765.7006	TI-765.8006	765.5006	765.6006	765.7006	765.8006
7	TI-765.5007	TI-765.6007	TI-765.7007	TI-765.8007	765.5007	765.6007	765.7007	765.8007
8	TI-765.5008	TI-765.6008	TI-765.7008	TI-765.8008	765.5008	765.6008	765.7008	765.8008
9	TI-765.5009	TI-765.6009	TI-765.7009	TI-765.8009	765.5009	765.6009	765.7009	765.8009
10	TI-765.5010	TI-765.6010	TI-765.7010	TI-765.8010	765.5010	765.6010	765.7010	765.8010
11	TI-765.5011	TI-765.6011	TI-765.7011	TI-765.8011	765.5011	765.6011	765.7011	765.8011
12	TI-765.5012	TI-765.6012	TI-765.7012	TI-765.8012	765.5012	765.6012	765.7012	765.8012
14	TI-765.5014	TI-765.6014	TI-765.7014	TI-765.8014	765.5014	765.6014	765.7014	765.8014
16	TI-765.5016	TI-765.6016	TI-765.7016	TI-765.8016	765.5016	765.6016	765.7016	765.8016
18	TI-765.5018	TI-765.6018	TI-765.7018	TI-765.8018	765.5018	765.6018	765.7018	765.8018

STERILE

Blade Length

	Titanium				Stainless Steel			
Holes	50mm	60mm	70mm	80mm	50mm	60mm	70mm	80mm
5	TI-765.5005-S	TI-765.6005-S	TI-765.7005-S	TI-765.8005-S	765.5005-S	765.6005-S	765.7005-S	765.8005-S
6	TI-765.5006-S	TI-765.6006-S	TI-765.7006-S	TI-765.8006-S	765.5006-S	765.6006-S	765.7006-S	765.8006-S
7	TI-765.5007-S	TI-765.6007-S	TI-765.7007-S	TI-765.8007-S	765.5007-S	765.6007-S	765.7007-S	765.8007-S
8	TI-765.5008-S	TI-765.6008-S	TI-765.7008-S	TI-765.8008-S	765.5008-S	765.6008-S	765.7008-S	765.8008-S
9	TI-765.5009-S	TI-765.6009-S	TI-765.7009-S	TI-765.8009-S	765.5009-S	765.6009-S	765.7009-S	765.8009-S
10	TI-765.5010-S	TI-765.6010-S	TI-765.7010-S	TI-765.8010-S	765.5010-S	765.6010-S	765.7010-S	765.8010-S
11	TI-765.5011-S	TI-765.6011-S	TI-765.7011-S	TI-765.8011-S	765.5011-S	765.6011-S	765.7011-S	765.8011-S
12	TI-765.5012-S	TI-765.6012-S	TI-765.7012-S	TI-765.8012-S	765.5012-S	765.6012-S	765.7012-S	765.8012-S
14	TI-765.5014-S	TI-765.6014-S	TI-765.7014-S	TI-765.8014-S	765.5014-S	765.6014-S	765.7014-S	765.8014-S
16	TI-765.5016-S	TI-765.6016-S	TI-765.7016-S	TI-765.8016-S	765.5016-S	765.6016-S	765.7016-S	765.8016-S
18	TI-765.5018-S	TI-765.6018-S	TI-765.7018-S	TI-765.8018-S	765.5018-S	765.6018-S	765.7018-S	765.8018-S

130° Angle Blade Plate



Blade Length

	Stainless Steel								
Holes	50mm	60mm	70mm	75mm	80mm	85mm	90mm	95mm	100mm
4	766.5004	766.6004-S	766.7004-S	766.7504-S	766.8004-S	766.8504-S	766.9004-S	766.9504-S	766.10004-S
5	766.5005	766.6005-S	766.7005-S	766.7505-S	766.8005-S	766.8505-S	766.9005-S	766.9505-S	766.10005-S
6	766.5006	766.6006-S	766.7006-S	766.7506-S	766.8006-S	766.8506-S	766.9006-S	766.9506-S	766.10006-S
7	766.5007	766.6007-S	766.7007-S	766.7507-S	766.8007-S	766.8507-S	766.9007-S	766.9507-S	766.10007-S
8	766.5008	766.6008-S	766.7008-S	766.7508-S	766.8008-S	766.8508-S	766.9008-S	766.9508-S	766.10008-S
9	766.5009	766.6009-S	766.7009-S	766.7509-S	766.8009-S	766.8509-S	766.9009-S	766.9509-S	766.10009-S
10	766.5010	766.6010-S	766.7010-S	766.7510-S	766.8010-S	766.8510-S	766.9010-S	766.9510-S	766.10010-S
11	766.5011	766.6011-S	766.7011-S	766.7511-S	766.8011-S	766.8511-S	766.9011-S	766.9511-S	766.10011-S
12	766.5012	766.6012-S	766.7012-S	766.7512-S	766.8012-S	766.8512-S	766.9012-S	766.9512-S	766.10012-S

STERILE

Blade Length

	Stainless Steel								
Holes	50mm	60mm	70mm	75mm	80mm	85mm	90mm	95mm	100mm
4	766.5004-S	766.6004-S	766.7004-S	766.7504-S	766.8004-S	766.8504-S	766.9004-S	766.9504-S	766.10004-S
5	766.5005-S	766.6005-S	766.7005-S	766.7505-S	766.8005-S	766.8505-S	766.9005-S	766.9505-S	766.10005-S
6	766.5006-S	766.6006-S	766.7006-S	766.7506-S	766.8006-S	766.8506-S	766.9006-S	766.9506-S	766.10006-S
7	766.5007-S	766.6007-S	766.7007-S	766.7507-S	766.8007-S	766.8507-S	766.9007-S	766.9507-S	766.10007-S
8	766.5008-S	766.6008-S	766.7008-S	766.7508-S	766.8008-S	766.8508-S	766.9008-S	766.9508-S	766.10008-S
9	766.5009-S	766.6009-S	766.7009-S	766.7509-S	766.8009-S	766.8509-S	766.9009-S	766.9509-S	766.10009-S
10	766.5010-S	766.6010-S	766.7010-S	766.7510-S	766.8010-S	766.8510-S	766.9010-S	766.9510-S	766.10010-S
11	766.5011-S	766.6011-S	766.7011-S	766.7511-S	766.8011-S	766.8511-S	766.9011-S	766.9511-S	766.10011-S
12	766.5012-S	766.6012-S	766.7012-S	766.7512-S	766.8012-S	766.8512-S	766.9012-S	766.9512-S	766.10012-S

130° Angle Blade Plate



Blade Length


	Titanium								
Holes	50mm	60mm	70mm	75mm	80mm	85mm	90mm	95mm	100mm
4	TI-766.5004	TI-766.6004	TI-766.7004	TI-766.7504	TI-766.8004	TI-766.8504	TI-766.9504	TI-766.9004	TI-766.10004
5	TI-766.5005	TI-766.6005	TI-766.7005	TI-766.7505	TI-766.8005	TI-766.8505	TI-766.9505	TI-766.9005	TI-766.10005
6	TI-766.5006	TI-766.6006	TI-766.7006	TI-766.7506	TI-766.8006	TI-766.8506	TI-766.9506	TI-766.9006	TI-766.10006
7	TI-766.5007	TI-766.6007	TI-766.7007	TI-766.7507	TI-766.8007	TI-766.8507	TI-766.9507	TI-766.9007	TI-766.10007
8	TI-766.5008	TI-766.6008	TI-766.7008	TI-766.7508	TI-766.8008	TI-766.8508	TI-766.9508	TI-766.9008	TI-766.10008
9	TI-766.5009	TI-766.6009	TI-766.7009	TI-766.7509	TI-766.8009	TI-766.8509	TI-766.9509	TI-766.9009	TI-766.10009
10	TI-766.5010	TI-766.6010	TI-766.7010	TI-766.7510	TI-766.8010	TI-766.8510	TI-766.9510	TI-766.9010	TI-766.10010
11	TI-766.5011	TI-766.6011	TI-766.7011	TI-766.7511	TI-766.8011	TI-766.8511	TI-766.9511	TI-766.9011	TI-766.10011
12	TI-766.5012	TI-766.6012	TI-766.7012	TI-766.7512	TI-766.8012	TI-766.8512	TI-766.9512	TI-766.9012	TI-766.10012

STERILE

Blade Length

	Titanium								
Holes	50mm	60mm	70mm	75mm	80mm	85mm	90mm	95mm	100mm
4	TI-766.5004-S	TI-766.6004-S	TI-766.7004-S	TI-766.7504-S	TI-766.8004-S	TI-766.8504-S	TI-766.9504-S	TI-766.9004-S	TI-766.10004-S
5	TI-766.5005-S	TI-766.6005-S	TI-766.7005-S	TI-766.7505-S	TI-766.8005-S	TI-766.8505-S	TI-766.9505-S	TI-766.9005-S	TI-766.10005-S
6	TI-766.5006-S	TI-766.6006-S	TI-766.7006-S	TI-766.7506-S	TI-766.8006-S	TI-766.8506-S	TI-766.9506-S	TI-766.9006-S	TI-766.10006-S
7	TI-766.5007-S	TI-766.6007-S	TI-766.7007-S	TI-766.7507-S	TI-766.8007-S	TI-766.8507-S	TI-766.9507-S	TI-766.9007-S	TI-766.10007-S
8	TI-766.5008-S	TI-766.6008-S	TI-766.7008-S	TI-766.7508-S	TI-766.8008-S	TI-766.8508-S	TI-766.9508-S	TI-766.9008-S	TI-766.10008-S
9	TI-766.5009-S	TI-766.6009-S	TI-766.7009-S	TI-766.7509-S	TI-766.8009-S	TI-766.8509-S	TI-766.9509-S	TI-766.9009-S	TI-766.10009-S
10	TI-766.5010-S	TI-766.6010-S	TI-766.7010-S	TI-766.7510-S	TI-766.8010-S	TI-766.8510-S	TI-766.9510-S	TI-766.9010-S	TI-766.10010-S
11	TI-766.5011-S	TI-766.6011-S	TI-766.7011-S	TI-766.7511-S	TI-766.8011-S	TI-766.8511-S	TI-766.9511-S	TI-766.9011-S	TI-766.10011-S
12	TI-766.5012-S	TI-766.6012-S	TI-766.7012-S	TI-766.7512-S	TI-766.8012-S	TI-766.8512-S	TI-766.9512-S	TI-766.9012-S	TI-766.10012-S

90°/100° Angle Blade Plate for Adults



Holes	Blade Length	Displacement	Angle	Titanium	Stainless Steel
4	40mm	10mm	90°	TI-767.94010	767.94010
4	50mm	10mm	90°	TI-767.95010	767.95010
4	60mm	10mm	90°	TI-767.96010	767.96010
4	70mm	10mm	90°	TI-767.97010	767.97010
4	40mm	15mm	90°	TI-767.94015	767.94015
4	50mm	15mm	90°	TI-767.95015	767.95015
4	60mm	15mm	90°	TI-767.96015	767.96015
4	70mm	15mm	90°	TI-767.97015	767.97015
4	50mm	20mm	90°	TI-767.95020	767.95020
4	60mm	20mm	90°	TI-767.96020	767.96020
4	70mm	20mm	90°	TI-767.97020	767.97020
4	50mm	10mm	100°	TI-767.05010	767.05010
4	60mm	10mm	100°	TI-767.06010	767.06010
4	70mm	10mm	100°	TI-767.07010	767.07010

STERILE

Holes	Blade Length	Displacement	Angle	Titanium	Stainless Steel
4	40mm	10mm	90°	TI-767.94010-S	767.94010-S
4	50mm	10mm	90°	TI-767.95010-S	767.95010-S
4	60mm	10mm	90°	TI-767.96010-S	767.96010-S
4	70mm	10mm	90°	TI-767.97010-S	767.97010-S
4	40mm	15mm	90°	TI-767.94015-S	767.94015-S
4	50mm	15mm	90°	TI-767.95015-S	767.95015-S
4	60mm	15mm	90°	TI-767.96015-S	767.96015-S
4	70mm	15mm	90°	TI-767.97015-S	767.97015-S
4	50mm	20mm	90°	TI-767.95020-S	767.95020-S
4	60mm	20mm	90°	TI-767.96020-S	767.96020-S
4	70mm	20mm	90°	TI-767.97020-S	767.97020-S
4	50mm	10mm	100°	TI-767.05010-S	767.05010-S
4	60mm	10mm	100°	TI-767.06010-S	767.06010-S
4	70mm	10mm	100°	TI-767.07010-S	767.07010-S



110°/120°/130° Angle Blade Plate for Intertrochanteric Femoral Osteotomies in Adults




Blade Length	Angle	Titanium	Stainless Steel
60mm	110°	TI-768.1160	768.1160
75mm	110°	TI-768.1175	768.1175
85mm	110°	TI-768.1185	768.1185
60mm	120°	TI-768.1260	768.1260
65mm	120°	TI-768.1265	768.1265
70mm	120°	TI-768.1270	768.1270
75mm	120°	TI-768.1275	768.1275
80mm	120°	TI-768.1280	768.1280
85mm	120°	TI-768.1285	768.1285
65mm	130°	TI-768.1365	768.1365
70mm	130°	TI-768.1370	768.1370
75mm	130°	TI-768.1375	768.1375
80mm	130°	TI-768.1380	768.1380
85mm	130°	TI-768.1385	768.1385
90mm	130°	TI-768.1390	768.1390
95mm	130°	TI-768.1395	768.1395
100mm	130°	TI-768.1400	768.1400



STERILE

Blade Length	Angle	Titanium	Stainless Steel
60mm	110°	TI-768.1160-S	768.1160-S
75mm	110°	TI-768.1175-S	768.1175-S
85mm	110°	TI-768.1185-S	768.1185-S
60mm	120°	TI-768.1260-S	768.1260-S
65mm	120°	TI-768.1265-S	768.1265-S
70mm	120°	TI-768.1270-S	768.1270-S
75mm	120°	TI-768.1275-S	768.1275-S
80mm	120°	TI-768.1280-S	768.1280-S
85mm	120°	TI-768.1285-S	768.1285-S
65mm	130°	TI-768.1365-S	768.1365-S
70mm	130°	TI-768.1370-S	768.1370-S
75mm	130°	TI-768.1375-S	768.1375-S
80mm	130°	TI-768.1380-S	768.1380-S
85mm	130°	TI-768.1385-S	768.1385-S
90mm	130°	TI-768.1390-S	768.1390-S
95mm	130°	TI-768.1395-S	768.1395-S
100mm	130°	TI-768.1400-S	768.1400-S

The background of the page features a series of thin, light gray wavy lines that flow from the top left towards the bottom right, creating a sense of movement and depth. These lines are more densely packed in some areas, particularly near the top and bottom edges, and more spread out in the center.

INSTRUMENTS

2100-2.0-100

Drill Bit with Quick Coupling End, Ø2.0mm x Length 100mm, for Large Fragment



2100-3.2-145

Drill Bit with Quick Coupling End, Ø3.2mm x Length 145mm, for Large Fragment



2100-4.5-145

Drill Bit with Quick Coupling End, Ø4.5mm x Length 145mm, for Large Fragment



3767-4.5

Double Drill Guide, Ø4.5/3.2mm



3441-20

Drill Guide 4.5mm for Neutral and Loaded Position



3437-09

Triple Drill Guide 130°



3437-04

Slotted Hammer



3437-01

Chisel Guide with adjustable Angle



3437-02

Seating Chisel



3437-05

Impactor for Angle Blade Plate



3437-03

Inserter-Extractor



3437-07

Triangular Positioning Plate



3437-08

Quadrangular Positioning Plate



3437-06

Condylar Plate Guide



451-2.0-150E

Kirschner Wire with One Side Trocar Tip, Ø2.0mm x Length 150mm



A10-016

Aluminium Box For Angle Blade Plate Instrument Set

ORDER INFORMATION

2313-000	Angled Blade Plate Instrument Set
----------	-----------------------------------


Code	Set Consisting Of	Unit
2100-2.0-100	Drill Bit with Quick Coupling End, Ø2.0mm x Length 100mm, for Large Fragment	1
2100-3.2-145	Drill Bit with Quick Coupling End, Ø3.2mm x Length 145mm, for Large Fragment	1
2100-4.5-145	Drill Bit with Quick Coupling End, Ø4.5mm x Length 145mm, for Large Fragment	1
3767-4.5	Double Drill Guide, Ø4.5/3.2mm	1
3441-20	Drill Guide 4.5mm for Neutral and Loaded Position	1
3437-09	Triple Drill Guide 130°	1
3437-04	Slotted Hammer	1
3437-01	Chisel Guide with adjustable Angle	1
3437-02	Seating Chisel	1
3437-05	Impactor for Angle Blade Plate	1
3437-03	Insertor-Extractor	1
3437-07	Triangular Positioning Plate	1
3437-08	Quadrangular Positioning Plate	1
3437-06	Condylar Plate Guide	1
451-2.0-150E	Kirschner Wire with One Side Trocar Tip, Ø2.0mm x Length 150mm	10
A10-016	Aluminium Box For Angle Blade Plate Instrument Set	1



The background of the page features a series of thin, light gray wavy lines that flow from the top left towards the bottom right, creating a sense of movement and depth. These lines are more densely packed in some areas, particularly near the top and bottom edges, and more sparse in the center where the text is located.

STAPLES, PINS & WIRES

Staple with Offset




A	B	C	Titanium
24mm	36mm	5mm	462-24-36-5
24mm	40mm	10mm	462-24-40-10
24mm	44mm	14mm	462-24-44-14

STERILE

A	B	C	Titanium
24mm	36mm	5mm	462-24-36-5-S
24mm	40mm	10mm	462-24-40-10-S
24mm	44mm	14mm	462-24-44-14-S

Three Legged Staple




A	B	C	Titanium
26mm	34mm	5mm	463-26-34-5

STERILE

A	B	C	Titanium
26mm	34mm	5mm	463-26-34-5-S

Ligament Staple



A	X	Y	Z	Titanium
6mm	11mm	25mm	3.7mm	464-6-11-25
8mm	13mm	25mm	3.7mm	464-8-13-25
11mm	16mm	25mm	3.7mm	464-11-16-25
16mm	21mm	25mm	3.7mm	464-16-21-25
8mm	13mm	20mm	3.7mm	464-8-13-20

STERILE

A	X	Y	Z	Titanium
6mm	11mm	25mm	3.7mm	464-6-11-25-S
8mm	13mm	25mm	3.7mm	464-8-13-25-S
11mm	16mm	25mm	3.7mm	464-11-16-25-S
16mm	21mm	25mm	3.7mm	464-16-21-25-S
8mm	13mm	20mm	3.7mm	464-8-13-20-S

Staple-Coventry



Description	Codes
5mm Step	460.005
10mm Step	460.010
15mm Step	460.015

STERILE

Description	Codes
5mm Step	460.005-S
10mm Step	460.010-S
15mm Step	460.015-S

Staple U - Type, Stainless Steel



Dia	15mm	20mm	25mm
1.5mm	9-028-01SS	9-028-05SS	9-028-09SS
2.0mm	9-028-02SS	9-028-06SS	9-028-10SS
2.5mm	9-028-03SS	9-028-07SS	9-028-11SS
3.0mm	9-028-04SS	9-028-08SS	9-028-12SS

STERILE

Dia	15mm	20mm	25mm
1.5mm	9-028-01SS-S	9-028-05SS-S	9-028-09SS-S
2.0mm	9-028-02SS-S	9-028-06SS-S	9-028-10SS-S
2.5mm	9-028-03SS-S	9-028-07SS-S	9-028-11SS-S
3.0mm	9-028-04SS-S	9-028-08SS-S	9-028-12SS-S

Staple U-Type, Titanium




Dia	15mm	20mm	25mm
1.5mm	459-1.5-015	459-1.5-020	459-1.5-025
2.0mm	459-2.0-015	459-2.0-020	459-2.0-025
2.5mm	459-2.5-015	459-2.5-020	459-2.5-025
3.0mm	459-3.0-015	459-3.0-020	459-3.0-025

STERILE

Dia	15mm	20mm	25mm
1.5mm	459-1.5-015-S	459-1.5-020-S	459-1.5-025-S
2.0mm	459-2.0-015-S	459-2.0-020-S	459-2.0-025-S
2.5mm	459-2.5-015-S	459-2.5-020-S	459-2.5-025-S
3.0mm	459-3.0-015-S	459-3.0-020-S	459-3.0-025-S

Straight Staple



A	B	Titanium
25mm	10mm	461-25-10
25mm	16mm	461-25-16
25mm	20mm	461-25-20
25mm	22mm	461-25-22
25mm	26mm	461-25-26
25mm	36mm	461-25-36

STERILE

A	B	Titanium
25mm	10mm	461-25-10-S
25mm	16mm	461-25-16-S
25mm	20mm	461-25-20-S
25mm	22mm	461-25-22-S
25mm	26mm	461-25-26-S
25mm	36mm	461-25-36-S

Steinmann Pins



Size (mm)	Codes
Ø2.0mm x 230mm	453-2.0-230
Ø2.3mm x 230mm	453-2.3-230
Ø2.0mm x 300mm	453-2.0-300
Ø2.5mm x 300mm	453-2.5-300
Ø3.0mm x 150mm	453-3.0-150
Ø3.0mm x 175mm	453-3.0-175
Ø3.0mm x 200mm	453-3.0-200
Ø3.0mm x 225mm	453-3.0-225
Ø3.5mm x 150mm	453-3.5-150
Ø3.5mm x 175mm	453-3.5-175
Ø3.5mm x 200mm	453-3.5-200
Ø3.5mm x 225mm	453-3.5-225
Ø4.0mm x 150mm	453-4.0-150
Ø4.0mm x 175mm	453-4.0-175
Ø4.0mm x 200mm	453-4.0-200
Ø4.0mm x 225mm	453-4.0-225
Ø4.5mm x 150mm	453-4.5-150
Ø4.5mm x 175mm	453-4.5-175
Ø4.5mm x 200mm	453-4.5-200
Ø4.5mm x 225mm	453-4.5-225
Ø5.0mm x 150mm	453-5.0-150
Ø5.0mm x 175mm	453-5.0-175
Ø5.0mm x 200mm	453-5.0-200
Ø5.0mm x 225mm	453-5.0-225


STERILE

Size (mm)	Codes
Ø2.0mm x 230mm	453-2.0-230-S
Ø2.3mm x 230mm	453-2.3-230-S
Ø2.0mm x 300mm	453-2.0-300-S
Ø2.5mm x 300mm	453-2.5-300-S
Ø3.0mm x 150mm	453-3.0-150-S
Ø3.0mm x 175mm	453-3.0-175-S
Ø3.0mm x 200mm	453-3.0-200-S
Ø3.0mm x 225mm	453-3.0-225-S
Ø3.5mm x 150mm	453-3.5-150-S
Ø3.5mm x 175mm	453-3.5-175-S
Ø3.5mm x 200mm	453-3.5-200-S
Ø3.5mm x 225mm	453-3.5-225-S
Ø4.0mm x 150mm	453-4.0-150-S
Ø4.0mm x 175mm	453-4.0-175-S
Ø4.0mm x 200mm	453-4.0-200-S
Ø4.0mm x 225mm	453-4.0-225-S
Ø4.5mm x 150mm	453-4.5-150-S
Ø4.5mm x 175mm	453-4.5-175-S
Ø4.5mm x 200mm	453-4.5-200-S
Ø4.5mm x 225mm	453-4.5-225-S
Ø5.0mm x 150mm	453-5.0-150-S
Ø5.0mm x 175mm	453-5.0-175-S
Ø5.0mm x 200mm	453-5.0-200-S
Ø5.0mm x 225mm	453-5.0-225-S

Centrally Threaded Steinman Pin




Size (mm)	Codes
Ø2.0mm x 100mm	454-2.0-100
Ø3.5mm x 250mm	454-3.5-250
Ø3.5mm x 300mm	454-3.5-300
Ø4.0mm x 225mm	454-4.0-225
Ø4.0mm x 300mm	454-4.0-300
Ø4.5mm x 225mm	454-4.5-225
Ø4.5mm x 300mm	454-4.5-300
Ø5.0mm x 200mm	454-5.0-200
Ø5.0mm x 225mm	454-5.0-225
Ø5.0mm x 226mm	454-5.0-226
Ø5.0mm x 300mm	454-5.0-300
Ø6.0mm x 225mm	454-6.0-225



Size (mm)	Codes
Ø2.0mm x 100mm	454-2.0-100-S
Ø3.5mm x 250mm	454-3.5-250-S
Ø3.5mm x 300mm	454-3.5-300-S
Ø4.0mm x 225mm	454-4.0-225-S
Ø4.0mm x 300mm	454-4.0-300-S
Ø4.5mm x 225mm	454-4.5-225-S
Ø4.5mm x 300mm	454-4.5-300-S
Ø5.0mm x 200mm	454-5.0-200-S
Ø5.0mm x 225mm	454-5.0-225-S
Ø5.0mm x 226mm	454-5.0-226-S
Ø5.0mm x 300mm	454-5.0-300-S
Ø6.0mm x 225mm	454-6.0-225-S

Suture Wire



Size (mm)	Codes
1.60mm, 10 Meter	450.016-10
1.50mm, 10 Meter	450.017-10
1.25mm, 10 Meter	450.018-10
1.00mm, 10 Meter	450.019-10
0.90mm, 10 Meter	450.020-10
0.80mm, 10 Meter	450.021-10
0.70mm, 10 Meter	450.022-10
0.60mm, 10 Meter	450.023-10
0.55mm, 10 Meter	450.024-10
0.45mm, 10 Meter	450.026-10
0.35mm, 10 Meter	450.035-10
0.28mm, 10 Meter	450.028-10
0.30mm, 10 Meter	450.030-10
1.60mm, 5 Meter	450.016-5
1.50mm, 5 Meter	450.017-5
1.25mm, 5 Meter	450.018-5
1.00mm, 5 Meter	450.019-5
0.90mm, 5 Meter	450.020-5
0.80mm, 5 Meter	450.021-5
0.70mm, 5 Meter	450.022-5
0.60mm, 5 Meter	450.023-5
0.55mm, 5 Meter	450.024-5
0.45mm, 5 Meter	450.026-5
0.35mm, 5 Meter	450.035-5
0.28mm, 5 Meter	450.028-5
0.30mm, 5 Meter	450.030-5



Size (mm)	Codes
1.60mm, 10 Meter	450.016-10-S
1.50mm, 10 Meter	450.017-10-S
1.25mm, 10 Meter	450.018-10-S
1.00mm, 10 Meter	450.019-10-S
0.90mm, 10 Meter	450.020-10-S
0.80mm, 10 Meter	450.021-10-S
0.70mm, 10 Meter	450.022-10-S
0.60mm, 10 Meter	450.023-10-S
0.55mm, 10 Meter	450.024-10-S
0.45mm, 10 Meter	450.026-10-S
0.35mm, 10 Meter	450.035-10-S
0.28mm, 10 Meter	450.028-10-S
0.30mm, 10 Meter	450.030-10-S
1.60mm, 5 Meter	450.016-5-S
1.50mm, 5 Meter	450.017-5-S
1.25mm, 5 Meter	450.018-5-S
1.00mm, 5 Meter	450.019-5-S
0.90mm, 5 Meter	450.020-5-S
0.80mm, 5 Meter	450.021-5-S
0.70mm, 5 Meter	450.022-5-S
0.60mm, 5 Meter	450.023-5-S
0.55mm, 5 Meter	450.024-5-S
0.45mm, 5 Meter	450.026-5-S
0.35mm, 5 Meter	450.035-5-S
0.28mm, 5 Meter	450.028-5-S
0.30mm, 5 Meter	450.030-5-S

Kirschner Wire Fully Threaded



Length	Dia	Codes
225mm	1.5mm	452-1.5-225
225mm	2.0mm	452-2.0-225
225mm	2.5mm	452-2.5-225
150mm	1.5mm	452-1.5-150
150mm	2.0mm	452-2.0-150
150mm	2.5mm	452-2.5-150

STERILE

Length	Dia	Codes
225mm	1.5mm	452-1.5-225-S
225mm	2.0mm	452-2.0-225-S
225mm	2.5mm	452-2.5-225-S
150mm	1.5mm	452-1.5-150-S
150mm	2.0mm	452-2.0-150-S
150mm	2.5mm	452-2.5-150-S

Circlage Wire With Loop



Length 280mm	Size (mm)
457-17-280	280mm x 17 SWG (1.50mm)
457-18-280	280mm x 18 SWG (1.25mm)
457-19-280	280mm x 19 SWG (1.00mm)
457-20-280	280mm x 20 SWG (0.90mm)
457-21-280	280mm x 21 SWG (0.80mm)
457-17-600	600mm x 17 SWG (1.50mm)
457-18-600	600mm x 18 SWG (1.25mm)
457-19-600	600mm x 19 SWG (1.00mm)
457-20-600	600mm x 20 SWG (0.90mm)
457-21-600	600mm x 21 SWG (0.80mm)
457-26-600	600mm x 26 SWG (0.50mm)
457-24-600	600mm x 24 SWG (0.60mm)
457-22-600	600mm x 22 SWG (0.70mm)

STERILE

Length 280mm	Size (mm)
457-17-280-S	280mm x 17 SWG (1.50mm)
457-18-280-S	280mm x 18 SWG (1.25mm)
457-19-280-S	280mm x 19 SWG (1.00mm)
457-20-280-S	280mm x 20 SWG (0.90mm)
457-21-280-S	280mm x 21 SWG (0.80mm)
457-17-600-S	600mm x 17 SWG (1.50mm)
457-18-600-S	600mm x 18 SWG (1.25mm)
457-19-600-S	600mm x 19 SWG (1.00mm)
457-20-600-S	600mm x 20 SWG (0.90mm)
457-21-600-S	600mm x 21 SWG (0.80mm)
457-26-600-S	600mm x 26 SWG (0.50mm)
457-24-600-S	600mm x 24 SWG (0.60mm)
457-22-600-S	600mm x 22 SWG (0.70mm)



Kirschner Wire with Trocar Tip Both End



Size (mm)	Codes
Ø0.8mm x 70mm	451-0.8-70
Ø0.9mm x 70mm	451-0.9-70
Ø1.0mm x 70mm	451-1.0-70
Ø1.1mm x 70mm	451-1.1-70
Ø1.2mm x 70mm	451-1.2-70
Ø1.4mm x 70mm	451-1.4-70
Ø1.5mm x 70mm	451-1.5-70
Ø1.6mm x 70mm	451-1.6-70
Ø1.8mm x 70mm	451-1.8-70
Ø2.0mm x 70mm	451-2.0-70
Ø2.2mm x 70mm	451-2.2-70
Ø2.3mm x 70mm	451-2.3-70
Ø2.5mm x 70mm	451-2.5-70
Ø3.0mm x 70mm	451-3.0-70
Ø0.8mm x 100mm	451-0.8-100
Ø0.9mm x 100mm	451-0.9-100
Ø1.0mm x 100mm	451-1.0-100
Ø1.1mm x 100mm	451-1.1-100
Ø1.2mm x 100mm	451-1.2-100
Ø1.4mm x 100mm	451-1.4-100
Ø1.5mm x 100mm	451-1.5-100
Ø1.6mm x 100mm	451-1.6-100
Ø1.8mm x 100mm	451-1.8-100
Ø2.0mm x 100mm	451-2.0-100
Ø2.2mm x 100mm	451-2.2-100
Ø2.3mm x 100mm	451-2.3-100
Ø2.5mm x 100mm	451-2.5-100
Ø3.0mm x 100mm	451-3.0-100
Ø0.8mm x 150mm	451-0.8-150
Ø0.9mm x 150mm	451-0.9-150
Ø1.0mm x 150mm	451-1.0-150
Ø1.1mm x 150mm	451-1.1-150
Ø1.2mm x 230mm	451-1.2-230
Ø1.2mm x 150mm	451-1.2-150
Ø1.4mm x 150mm	451-1.4-150
Ø1.5mm x 230mm	451-1.5-230
Ø1.5mm x 150mm	451-1.5-150
Ø1.6mm x 150mm	451-1.6-150
Ø1.8mm x 150mm	451-1.8-150

Size (mm)	Codes
Ø2.2mm x 150mm	451-2.2-150
Ø2.3mm x 150mm	451-2.3-150
Ø2.5mm x 150mm	451-2.5-150
Ø3.0mm x 150mm	451-3.0-150
Ø0.8mm x 225mm	451-0.8-225
Ø0.9mm x 225mm	451-0.9-225
Ø1.0mm x 225mm	451-1.0-225
Ø1.1mm x 225mm	451-1.1-225
Ø1.2mm x 225mm	451-1.2-225
Ø1.4mm x 225mm	451-1.4-225
Ø1.5mm x 225mm	451-1.5-225
Ø1.6mm x 225mm	451-1.6-225
Ø1.8mm x 225mm	451-1.8-225
Ø2.0mm x 225mm	451-2.0-225
Ø2.2mm x 225mm	451-2.2-225
Ø2.3mm x 225mm	451-2.3-225
Ø2.5mm x 225mm	451-2.5-225
Ø3.0mm x 225mm	451-3.0-225
Ø1.0mm x 230mm	451-1.0-230
Ø1.5mm x 250mm	451-1.5-250
Ø2.0mm x 230mm	451-2.0-230
Ø2.5mm x 230mm	451-2.5-230
Ø0.8mm x 300mm	451-0.8-300
Ø0.9mm x 300mm	451-0.9-300
Ø1.0mm x 300mm	451-1.0-300
Ø1.1mm x 300mm	451-1.1-300
Ø1.2mm x 300mm	451-1.2-300
Ø1.4mm x 300mm	451-1.4-300
Ø1.5mm x 300mm	451-1.5-300
Ø1.6mm x 300mm	451-1.6-300
Ø1.8mm x 300mm	451-1.8-300
Ø2.0mm x 300mm	451-2.0-300
Ø2.2mm x 300mm	451-2.2-300
Ø2.3mm x 300mm	451-2.3-300
Ø2.5mm x 300mm	451-2.5-300
Ø2.8mm x 300mm	451-2.8-300
Ø3.0mm x 300mm	451-3.0-300
Ø0.8mm x 310mm	9-016-001
Ø1.0mm x 310mm	9-016-002
Ø1.2mm x 310mm	9-016-003
Ø1.4mm x 310mm	9-016-004
Ø1.6mm x 310mm	9-016-005
Ø1.8mm x 310mm	9-016-006



Kirschner Wire with Trocar Tip Both End

STERILE

Size (mm)	Codes
Ø0.8mm x 70mm	451-0.8-70-S
Ø0.9mm x 70mm	451-0.9-70-S
Ø1.0mm x 70mm	451-1.0-70-S
Ø1.1mm x 70mm	451-1.1-70-S
Ø1.2mm x 70mm	451-1.2-70-S
Ø1.4mm x 70mm	451-1.4-70-S
Ø1.5mm x 70mm	451-1.5-70-S
Ø1.6mm x 70mm	451-1.6-70-S
Ø1.8mm x 70mm	451-1.8-70-S
Ø2.0mm x 70mm	451-2.0-70-S
Ø2.2mm x 70mm	451-2.2-70-S
Ø2.3mm x 70mm	451-2.3-70-S
Ø2.5mm x 70mm	451-2.5-70-S
Ø3.0mm x 70mm	451-3.0-70-S
Ø0.8mm x 100mm	451-0.8-100-S
Ø0.9mm x 100mm	451-0.9-100-S
Ø1.0mm x 100mm	451-1.0-100-S
Ø1.1mm x 100mm	451-1.1-100-S
Ø1.2mm x 100mm	451-1.2-100-S
Ø1.4mm x 100mm	451-1.4-100-S
Ø1.5mm x 100mm	451-1.5-100-S
Ø1.6mm x 100mm	451-1.6-100-S
Ø1.8mm x 100mm	451-1.8-100-S
Ø2.0mm x 100mm	451-2.0-100-S
Ø2.2mm x 100mm	451-2.2-100-S
Ø2.3mm x 100mm	451-2.3-100-S
Ø2.5mm x 100mm	451-2.5-100-S
Ø3.0mm x 100mm	451-3.0-100-S
Ø0.8mm x 150mm	451-0.8-150-S
Ø0.9mm x 150mm	451-0.9-150-S
Ø1.0mm x 150mm	451-1.0-150-S
Ø1.1mm x 150mm	451-1.1-150-S
Ø1.2mm x 230mm	451-1.2-230-S
Ø1.2mm x 150mm	451-1.2-150-S
Ø1.4mm x 150mm	451-1.4-150-S
Ø1.5mm x 230mm	451-1.5-230-S
Ø1.5mm x 150mm	451-1.5-150-S
Ø1.6mm x 150mm	451-1.6-150-S
Ø1.8mm x 150mm	451-1.8-150-S

Size (mm)	Codes
Ø2.2mm x 150mm	451-2.2-150-S
Ø2.3mm x 150mm	451-2.3-150-S
Ø2.5mm x 150mm	451-2.5-150-S
Ø3.0mm x 150mm	451-3.0-150-S
Ø0.8mm x 225mm	451-0.8-225-S
Ø0.9mm x 225mm	451-0.9-225-S
Ø1.0mm x 225mm	451-1.0-225-S
Ø1.1mm x 225mm	451-1.1-225-S
Ø1.2mm x 225mm	451-1.2-225-S
Ø1.4mm x 225mm	451-1.4-225-S
Ø1.5mm x 225mm	451-1.5-225-S
Ø1.6mm x 225mm	451-1.6-225-S
Ø1.8mm x 225mm	451-1.8-225-S
Ø2.0mm x 225mm	451-2.0-225-S
Ø2.2mm x 225mm	451-2.2-225-S
Ø2.3mm x 225mm	451-2.3-225-S
Ø2.5mm x 225mm	451-2.5-225-S
Ø3.0mm x 225mm	451-3.0-225-S
Ø1.0mm x 230mm	451-1.0-230-S
Ø1.5mm x 250mm	451-1.5-250-S
Ø2.0mm x 230mm	451-2.0-230-S
Ø2.5mm x 230mm	451-2.5-230-S
Ø0.8mm x 300mm	451-0.8-300-S
Ø0.9mm x 300mm	451-0.9-300-S
Ø1.0mm x 300mm	451-1.0-300-S
Ø1.1mm x 300mm	451-1.1-300-S
Ø1.2mm x 300mm	451-1.2-300-S
Ø1.4mm x 300mm	451-1.4-300-S
Ø1.5mm x 300mm	451-1.5-300-S
Ø1.6mm x 300mm	451-1.6-300-S
Ø1.8mm x 300mm	451-1.8-300-S
Ø2.0mm x 300mm	451-2.0-300-S
Ø2.2mm x 300mm	451-2.2-300-S
Ø2.3mm x 300mm	451-2.3-300-S
Ø2.5mm x 300mm	451-2.5-300-S
Ø2.8mm x 300mm	451-2.8-300-S
Ø3.0mm x 300mm	451-3.0-300-S
Ø0.8mm x 310mm	9-016-001-S
Ø1.0mm x 310mm	9-016-002-S
Ø1.2mm x 310mm	9-016-003-S
Ø1.4mm x 310mm	9-016-004-S
Ø1.6mm x 310mm	9-016-005-S
Ø1.8mm x 310mm	9-016-006-S



Kirschner Wire with One Side Trocar Tip



Size (mm)	Codes
Ø0.8mm x 150mm	451-0.8-150E
Ø0.9mm x 150mm	451-0.9-150E
Ø1.0mm x 150mm	451-1.0-150E
Ø1.4mm x 150mm	451-1.4-150E
Ø1.6mm x 150mm	451-1.6-150E
Ø1.8mm x 150mm	451-1.8-150E
Ø1.2mm x 300mm	451-1.2-300E
Ø1.4mm x 300mm	451-1.4-300E
Ø1.5mm x 300mm	451-1.5-300E
Ø1.6mm x 300mm	451-1.6-300E
Ø1.8mm x 300mm	451-1.8-300E
Ø2.0mm x 300mm	451-2.0-300E
Ø2.2mm x 300mm	451-2.2-300E
Ø2.5mm x 300mm	451-2.5-300E
Ø1.5mm x 310mm	451-1.5-310E
Ø2.0mm x 310mm	451-2.0-310E
Ø2.2mm x 310mm	451-2.2-310E
Ø2.5mm x 310mm	451-2.5-310E
Ø1.8mm x 310mm	451-1.8-310E
Ø3.0mm x 310mm	451-3.0-310E
Ø2.2mm x 225mm	451-2.2-225E



Size (mm)	Codes
Ø0.8mm x 150mm	451-0.8-150E-S
Ø0.9mm x 150mm	451-0.9-150E-S
Ø1.0mm x 150mm	451-1.0-150E-S
Ø1.4mm x 150mm	451-1.4-150E-S
Ø1.6mm x 150mm	451-1.6-150E-S
Ø1.8mm x 150mm	451-1.8-150E-S
Ø1.2mm x 300mm	451-1.2-300E-S
Ø1.4mm x 300mm	451-1.4-300E-S
Ø1.5mm x 300mm	451-1.5-300E-S
Ø1.6mm x 300mm	451-1.6-300E-S
Ø1.8mm x 300mm	451-1.8-300E-S
Ø2.0mm x 300mm	451-2.0-300E-S
Ø2.2mm x 300mm	451-2.2-300E-S
Ø2.5mm x 300mm	451-2.5-300E-S
Ø1.5mm x 310mm	451-1.5-310E-S
Ø2.0mm x 310mm	451-2.0-310E-S
Ø2.2mm x 310mm	451-2.2-310E-S
Ø2.5mm x 310mm	451-2.5-310E-S
Ø1.8mm x 310mm	451-1.8-310E-S
Ø3.0mm x 310mm	451-3.0-310E-S
Ø2.2mm x 225mm	451-2.2-225E-S



TRAUMATOLOGY

The content of this catalogue is intended for healthcare professionals only. All product information, illustrations, descriptions, and specifications provided here are for general informational purposes and are subject to change without notice. Auxein Medical Pvt. Ltd. reserves the right to improve or discontinue products at any time without prior notification.

This document does not constitute medical advice. Surgeons and medical practitioners must exercise their own professional judgment when using any medical devices or instruments, including those provided by Auxein.

CLINICAL USE

The implants and instruments in this catalogue are intended for use by qualified orthopedic trauma surgeons trained in fracture fixation and trauma management procedures. Proper implant selection and surgical technique are essential for effective use. Improper selection, misuse, or insufficient surgical experience may result in complications.

REGULATORY COMPLIANCE CERTIFIED

CE Approved (EUMDR-2017/745)

USFDA 510(k)

ISO 13485:2016


CDSCO

MDSAP

DISCLAIMER: Not all products in this catalogue are covered by European CE, US FDA, or other certifications. Please contact our Sales Team for product-specific certification details.

INTELLECTUAL PROPERTY

All trademarks, trade names, brand names, product names, and logos mentioned in this catalogue are the property of Auxein Medical Pvt. Ltd. Unauthorized use or reproduction is strictly prohibited.



SCAN TO KNOW MORE

INDIA
AUXEIN MEDICAL PVT.LTD.

PLOT NO. 168-169-170, PHASE-4,
KUNDLI INDUSTRIAL AREA,
HSIIDC, SECTOR-57,
SONEPAT-131028, HARYANA
TEL: +91 99106 43638
FAX: +91 86077 70197
EMAIL: INFO@AUXEIN.COM

MÉXICO
AUXEIN MÉXICO, S.A DE CV

TEPIC 39 INT 801, COLONIA
ROMA SUR,
ALCADÍA CUAUHTÉMOC, CDMX
MÉXICO, C.P 06760
TEL: +521 55 7261 0318
EMAIL: INFO@AUXEIN.MX

USA
AUXEIN INC.

1500 NW 89TH COURT,
SUITE 107-108
DORAL, FLORIDA 33172
TEL. +1 305 395 6062
E FAX: +1 305 395 62 62
EMAIL: USOFFICE@AUXEIN.COM

DUBAI UAE
AUXEIN MEDICAL

TRADING FZCO OUD METHA
OFFICES BUILDING B BLOCK,
1ST FLOOR, OFFICE NO.26,
NEAR HEALTHCARE CITY,
PO BOX-25572
TEL. +971 438 71333
EMAIL: INFO@AUXEIN.COM

WWW.AUXEIN.COM