






BIO BRAID





AUXILOCK® BIOBRAID SUTURES

MO-6
Needle



Sutures

AUXILOCK® BioBraid Sutures are braided sterile sutures prepared from Ultra High Molecular Weight Polyethylene (UHMWPE). They can be used in soft tissue approximations along with allograft tissues in arthroscopy procedures. BioBraid sutures are non-absorbable and do not impose any significant changes in tensile strength retention known to occur in vivo. BioBraid sutures are available in a variety of colour combinations for easy suture management in complicated repairs.

ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION	SUTURE COLOUR	NEEDLE
6-011-03	AUXILOCK® Two #2 Bio-Braid: White/Blue & White/Black, 36in Total Length		
6-011-04	AUXILOCK® #2 BioBraid: Blue, With Needles: MO-6 & CE, 36in Total Length		
6-011-05	AUXILOCK® #2 BioBraid: White, With Needle: MO-6, 36in Total Length		
6-011-06	AUXILOCK® #2 BioBraid: White/Black, With Needle: MO-6, 36in Total Length		
6-011-08	AUXILOCK® #2 BioBraid: White/Blue, With Needle: MO-6, 36in Total Length		

MO-6 : - 26 mm length, 1/2 circle shape, and taper point
MO-7: - 22mm 1/2 Circle Taper Point needle
CE : - 1/2 Circle Reverse Cutting
CCS : - 1/2 Circle Trocar point


ORDERING INFORMATION

Note: Each box consist of 10 Pieces

CODE	PRODUCT DESCRIPTION	SUTURE COLOUR	NEEDLE
6-011-09	AUXILOCK® #2 BioBraid: White/Black, 36in Total Length		
6-011-10	AUXILOCK® #2 BioBraid: White/Green, With Needle: MO-6, 36in Total Length		
6-011-11	AUXILOCK® #2 BioBraid: White/Green, 36in Total Length		
6-011-12	AUXILOCK® #2 BioBraid: Blue, 36in Total Length		
6-011-13	AUXILOCK® #2 BioBraid: White, 36in Total Length		
6-011-17	AUXILOCK® #2 BioBraid: White/Blue, 36in Total Length		
6-011-32	AUXILOCK® #2 BioBraid: White, With Needle: CCS, 36in Total Length		
6-011-07	AUXILOCK® #5 BioBraid: White, With Needle: CCS, 36in Total Length		
6-011-33	AUXILOCK® #2-0 BioBraid: White/Blue, 36in Total Length		



AUXILOCK® BIOBRAID INFINITY LOOP NEEDLE



AUXILOCK® BioBraid Infinity Loop is a continuous loop of BioBraid with a straight needle. BioBraid Infinity Loop utilises the whipstitch technique on graft preparation and reduces the time spent on graft preparation. It also compresses the graft uniformly and improves strength.


ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION	SUTURE COLOUR	NEEDLE
6-012-01	AUXILOCK® #2 BioBraid Infinity Loop: White, With Needle: Straight, 20in Loop Length, 24in Total Length		
6-012-02	AUXILOCK® #2 BioBraid Infinity Loop: White/Blue, With Needle: Straight, 20in Loop Length, 24in Total Length		
6-012-03	AUXILOCK® #5 BioBraid Infinity Loop: White/Blue, With Needle: Straight, 20in Loop Length, 24in Total Length		

Note: Each box consist of 10 Pieces.



AUXILOCK® BIOBRAID SUTURE TAPE



Suture Tape

AUXILOCK® BioBraid Suture Tape made from UHMWPE The suture tape design accommodates the tape in the centre with a suture on the either sides of the tape.

It can be used in Knee, Shoulder, and Foot & Ankle surgeries including AC joint reconstruction, knotless double-row rotator cuff repair, patella fracture repair, MPFL Reconstruction etc. The broad footprint of the BioBraid Suture Tape is appropriate for repairs in degenerative cuff tissue where tissue pull-through may be a concern.

ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION	SUTURE COLOUR
6-013-01	AUXILOCK® 1.8mm BioBraid Suture Tape: White/Blue, 39in Total length	
6-013-02	AUXILOCK® 1.4mm BioBraid Suture Tape: White/Blue, 39in Total length	
6-013-03	AUXILOCK® 1.4mm BioBraid Suture Tape: White/Black, 39in Total length	
6-013-05	AUXILOCK® 2.0mm BioBraid Suture Tape: White/Blue, 39in Total length	
6-013-06	AUXILOCK® 1.0mm BioBraid Suture Tape: White/Blue, 39in Total length	
6-013-07	AUXILOCK® 1.4mm BioBraid Suture Tape: Solid Blue, 39in Total length	

Note: Each box consist of 10 Pieces.