




HIP IMPLANTS



AUXILOCK® KNOTLESS PEEK CF PUSH-IN SUTURE ANCHOR



The AUXILOCK® HIP-Length Knotless PEEK CF Push-In Suture Anchor is designed for labral repair in hip arthroscopy. Engineered for deep acetabular access, the system facilitates accurate anchor placement along the acetabular rim while ensuring reliable soft tissue fixation. The extended hip-length driver enables improved reach and control within the constrained geometry of the hip joint. The knotless push-in mechanism allows controlled suture tensioning without the need for arthroscopic knot tying, supporting procedural efficiency and reproducible fixation outcomes.

ORDERING INFORMATION


CODE	PRODUCT DESCRIPTION
6-038-02	AUXILOCK® 2.8mm Knotless PEEK CF Push-In Suture Anchor With Hip Length Driver, With Suture Passer (Hip)

RECOMMENDED INSTRUMENTATION

CODE	PRODUCT DESCRIPTION	COMPATIBLE IMPLANTS
7-138-01	8 Point Drill Guide, Long	6-038-02
7-138-02	Ø2.5mm X 360mm Knotless Drill Bit, Long	6-038-02
7-138-03	Trocar Tip Obturator, Long	6-038-02
7-138-04	Tissue Grasper, Long	6-038-02



AUXILOCK® KNOTLESS PEEK OPTIMA PUSH-IN SUTURE ANCHOR




The AUXILOCK® HIP-Length Knotless PEEK OPTIMA Push-In Suture Anchor is designed for labral repair in hip arthroscopy. Engineered for deep acetabular access, the system facilitates accurate anchor placement along the acetabular rim while ensuring reliable soft tissue fixation. The extended hip-length driver enables improved reach and control within the constrained geometry of the hip joint. The knotless push-in mechanism allows controlled suture tensioning without the need for arthroscopic knot tying, supporting procedural efficiency and reproducible fixation outcomes.

ORDERING INFORMATION


CODE	PRODUCT DESCRIPTION
6-039-02	AUXILOCK® 2.8mm Knotless PEEK OPTIMA Push-In Suture Anchor With Hip Length Driver, With Suture Passer (Hip)

RECOMMENDED INSTRUMENTATION

CODE	PRODUCT DESCRIPTION	COMPATIBLE IMPLANTS
7-138-01	8 Point Drill Guide, Long	6-039-02
7-138-02	Ø2.5mm X 360mm Knotless Drill Bit, Long	6-039-02
7-138-03	Trocar Tip Obturator, Long	6-039-02
7-138-04	Tissue Grasper, Long	6-039-02



INSTRUMENTS





7-138-01

8 Point Drill Guide, Long



7-138-02

Ø2.5mm X 360mm Knotless Drill Bit, Long (for 2.8mm Knotless Anchor)



7-138-03

Trocar Tip Obturator, Long



OPTIONAL INSTRUMENTATION

7-138-04

Tissue Grasper Long

